

COVIFIND™

By Meril

INSTRUCTIONS FOR USE

■ INTENDED USE

The COVIFIND COVID-19 Antigen Self-Test (home use) is an *in-vitro* diagnostic test kit that uses a lateral flow immunoassay for the qualitative detection of SARS-CoV-2 specific antigens in nasal swab specimens from symptomatic individuals. This test is authorized for the collection of nasal specimens from individuals aged 18 years or older who have experienced COVID-like symptoms within the last seven days. It is intended for use in the home, workplace, or other settings. Individuals under 18 years will require adult supervision or assistance. This kit is designed for self-testing by laypersons in home or non-laboratory environments.


■ PRINCIPLE

COVIFIND COVID-19 Antigen Self-Test (home use) is an immunochromatographic rapid assay kit for qualitative determination of SARS-CoV-2 infection from Nasal swab specimens. Monoclonal anti-SARS-CoV-2 nucleocapsid antibody is coated on the test line region. Antigens of SARS-CoV-2 in the specimens react with the anti-SARS-CoV-2 monoclonal nucleocapsid antibody coupled with gold conjugate and form an antigen-antibody complex followed by reaction with anti-SARS-CoV-2 monoclonal nucleocapsid antibodies immobilized in the test line. This complex migrates on the membrane, where it will be captured by the monoclonal anti-SARS-CoV-2 antibody. A coloured test line would be visible in the result window if SARS-CoV-2 antigens are present in the specimen. The intensity of coloured test line will vary depending upon the amount of SARS-CoV-2 antigen present in the specimen. If SARS-CoV-2 antigens are not present in the specimen, then no line appears in the test line. The control band is used for procedural control and should always appear if the test procedure is performed correctly.

■ WHAT ARE THE SYMPTOMS OF COVID-19

Use COVIFIND COVID-19 Antigen Self-Test if you have one or more symptoms of COVID-19. COVID-19 symptoms can vary widely, but common ones include: Fever or chills, Cough, Shortness of breath or difficulty breathing, Fatigue, Muscle or body aches, Headache, New loss of taste or smell, Sore throat, Congestion or runny nose, Nausea or vomiting, Diarrhea. Symptoms typically appear 2-14 days after exposure.

KIT CONTENTS



TEST PREPARATION

Wash your hands
Mark sure they are dry before starting the test




Remove the test components from the box and place them on a flat surface.

Remove the test device from the pouch.




Peel the Aluminium foil from the pre-filled buffer tube.

Push the buffer tube in to the **perforated circle**.



TEST PREPARATION

Remove the sterile swab from the tail end.



DO NOT TOUCH THE SWAB HEAD

SAMPLE COLLECTION: STEP 1



Carefully insert the swab into one nostril 2–4cm or until resistance is met



5 x


Slowly rotate the swab, gently rubbing it along the insides of your nasal passage at least five (5) times.



5 x

Using the same swab, repeat in the other nostril.

SAMPLE COLLECTION: STEP 2




8-10 x

Insert the swab in the buffer tube and **swirl** it 8-10 times in the fluid.


Pinch the Swab tube through the tube to remove any remaining fluid



Hold the buffer tube firmly, lift the swab till **breakpoint** and **snap** the swab handle at breakpoint




Leave the swab inside the tube and discard the stick



Place the cap on top of the buffer tube

RUN TEST




Wait 15-20 minutes

Add four (4) drops (~100µl) of the solutions into the sample well

READ RESULTS

POSITIVE



■ WHAT SHOULD I DO IF THE RESULT SHOWS POSITIVE?

Refer to your relevant health authority for advice on whether a PCR test is required to confirm your result. You should also inform the immediate contacts you have had in past 24 hours so they can take any appropriate precautions.


CONSIDER THE LINK FOR MORE DETAILS

<https://www.health.gov.au/topics/covid-19/testing-positive>

If **Control (C)** and **Test (T)** bands both appear, the test indicates for the presence of antigens to SARS-CoV-2 in the sample.


Look very closely! Even a faint line at the bottom indicates a positive result.

Different possibilities of positive result



■ WHAT SHOULD I DO IF THE RESULT SHOWS NEGATIVE?

Negative results may require additional testing to confirm your results if you are symptomatic. If symptomatic, continue antigen testing every 24 hours for three (3) days or take a laboratory PCR test. If asymptomatic, it is likely that you were not infectious at the time the test was taken. A negative test result, however, is not a guarantee that you do not have coronavirus, take any appropriate If only the Control (C) band appears, the test indicates that the result is **Negative**.



NEGATIVE



INVALID

If no Control (C) band appears, then the assay is invalid regardless of colour development on Test (T) band as indicated below. Repeat the assay with a new device.

DISPOSAL OF TEST KIT



Place the swab, tube and test device into the disposal bag



Seal the disposal bag and dispose the bag in a non-recycle waste bin



Safe Interact Pty Ltd
26/155 Cooper Road, Yagoona 2199NSW Australia

PRODUCT INFORMATION: www.covifind.au

For assistance regarding the use of the product and interpretation of test results call 1800312940. This service is available between 9 am and 7 pm (AEST), or 9 am and 8 pm (AEDT), 7 days per week.



Scan the QR code for more information about Covifind.

■ **WARNINGS AND PRECAUTIONS**

- For in vitro diagnostic use only.
- Read instructions prior to performing this test. Follow all instructions to achieve valid results.
- Do not use the kit contents beyond the expiry date.
- Do not store the test kit in direct sunlight.
- Do not freeze the kit or expose the kit over 30°C.
- Do not eat or smoke while handling specimens.
- Remove any piercings from the nose before testing.
- Wash your hands thoroughly before and after the test is completed.
- Clean up spills thoroughly using an appropriate disinfectant.
- Place all of the items in the bag provided and dispose in a non-recyclable rubbish bin.
- Use only the ingredients provided in the kit.
- Keep the test kit out of reach of children.
- Do not touch the swab head as touching will cause the test to be incorrect.
- The provided Swab should be used only for nasal specimen collection.
- Avoid contact of any liquids with eyes and skin.
- If the swab stick breaks during specimen collection, dispose the kit and recommence with a new Test Kit.
- Each single Test Device, Swab, Tube, Cap, Bottle and Bag are single use. Do not reuse individual components.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- The swab specimen should be tested immediately after collection.
- Do not re-use any of the items in the Test kit for any other purpose.

■ **KIT COMPONENTS**

- Prefilled buffer tube with cap
- Sterile Nasal Swab
- Test Device
- Disposal Bag
- Instructions for Use

■ **MATERIALS REQUIRED BUT NOT PROVIDED**

- - Timer, stopwatch, or clock

■ **STORAGE AND STABILITY**

The kit should be stored between 2° to 30°C.

Do not store in the direct sunlight.

In the event that the desiccant pouch has changed color from blue to light pink or colorless, the device should not be used.

■ **CAN THE COVIFIND COVID-19 ANTIGEN TEST DETECT VARIOUS VARIANTS OF COVID-19?**

Yes, the COVIFIND COVID-19 Antigen Self-Test (home use) can detect Alpha, Beta, Gamma, Delta and Omicron COVID-19 mutants based on the studies conducted so far.

■ **LIMITATIONS**

- The COVIFIND COVID-19 Antigen Self-Test (home use) is designed for the primary test of SARS-CoV-2.
- The results obtained with this test should be interpreted as a presumptive test only and requires confirmation by a PCR test as required in your geographic location.
- A negative or non-reactive test result does not preclude the possibility of exposure to or infection with SARS-CoV-2 virus at any time.
- Failure to follow the test procedure and interpretation of test results may adversely affect the test performance and/or produce invalid results.
- Positive test results do not exclude co-infections with other pathogens.
- The accuracy of the test depends on the sample collection process. Improper sample collection, improper sample storage may affect the test result.
- If the test is not performed within seven (7) days of symptom onset, false negatives may occur.
- A positive result does not guarantee infection.
- The test is less reliable in the later phase of infection and in asymptomatic individuals.
- Repeat testing within 1–3 days is recommended in occupational risk, high risk settings or if there is an ongoing suspicion of infection.
- Reading the results later than 20 minutes will give incorrect results.
- Negative test results are not intended to rule in other non-SARS viral or bacterial infections.

■ **TEST LIMITATIONS**

- The contents of this kit are to be used for the qualitative detection of SARS-CoV-2 antigen from nasal swab.
- Failure to follow these instructions and interpretation of test results may adversely affect test performance and/or produce invalid results.
- A negative test result may occur if the specimen was collected, extracted or transported improperly.
- If symptoms continue, you should repeat the test after 1–2 days, as the coronavirus may not be detectable in the very early phases of infection.
- You are also advised to continue following local guidelines for self-isolation and consult your doctor.
- Positive test results do not rule out co-infections with other pathogens.
- Reading the test results earlier than 15 minutes or later than 19 minutes may give incorrect results.
- The COVIFIND COVID-19 Antigen Self-Test (home use) is not intended to detect from defective (non-infectious) virus during the later stages of viral shedding that might be detected by PCR molecular tests. Therefore, a positive result cannot necessarily determine whether a person is infectious.
- Due to cross-reactivity with high concentrations of SARS-CoV, a false positive result may occur in the case of infection with SARS-CoV.
- Wait four (4) hours before repeating the test following an invalid result.
- The test is less reliable in the later phase of infection and asymptomatic individuals.
- Failure to follow the test procedure and interpretation of test results may adversely affect the test performance and/or produce invalid results.

■ **USABILITY REPORT**

A usability study was conducted with a pool of 106 lay persons in the self-testing environment. The sensitivity was found to be >90% and the specificity was confirmed to be 100% in the hand of the lay person, comparing with a professional PCR test. Therefore the Summative Evaluation has proven that the usability of the COVIFIND COVID-19 Antigen Self Test ensures a safe and proper use of the device.

■ **PERFORMANCE EVALUATION OF COVIFIND COVID-19 ANTIGEN SELF TEST:**

A. Diagnostic Sensitivity:

Overall Diagnostic Sensitivity of COVIFIND COVID-19 Antigen Self Test:
Total 120 SARS-CoV-2 positive nasal swab specimens were tested with COVIFIND COVID-19 Antigen Self Test and 114 out of 120 specimens were detected as positive.
So, diagnostic sensitivity of COVIFIND COVID-19 Antigen Self Test is calculated as 95.00% (95% CI : 89.43 % to 98.14%).

B. Diagnostic Specificity:

Overall Diagnostic Specificity of COVIFIND COVID-19 Antigen Self Test:
Total 405 SARS-CoV-2 negative specimens were tested with COVIFIND COVID-19 Antigen Self Test. All specimens were identified as negative when tested with COVIFIND COVID-19 Antigen Self Test. So, overall Diagnostic specificity of COVIFIND COVID-19 Antigen Self Test is calculated as 100.00% (95% CI: 99.09% to 100%).

C. Analytical Sensitivity (Limit of Detection):

Limit of detection for COVIFIND COVID-19 Antigen Self Test is 933TCID₅₀/ml.

D. Analytical Specificity (Cross Reactivity)

The following Cross reactants and microorganisms had no impact on the performance of COVIFIND COVID-19 Antigen Self Test:

Microbial Organisms	
Adenovirus	Mycoplasma pneumoniae
Human Metapneumovirus (hMPV)	Moraxella Catarrhalis
Parainfluenza virus -1	Chlamydia pneumoniae
Parainfluenza virus -4	Legionella pneumophila
Influenza A	Staphylococcus aureus
Influenza B	Staphylococcus epidermidis
Enterovirus	Mycobacterium tuberculosis
Respiratory syncytial virus	Human coronavirus 229E
Rhinovirus	Human coronavirus OC43
Haemophilus influenzae	Human coronavirus NL63
Streptococcus pneumoniae	MERS-coronavirus
Streptococcus pyogenes	Parainfluenza 2
Candida albicans	Parainfluenza 3
Pooled human nasal wash	
Bordetella pertussis	

■ **INTERFERING SUBSTANCES:**

The following compounds have been tested using the COVID-19 Antigen Test and no interference was observed with Whole Blood, Mucin, Mupirocin, Oxymetazoline, Dexamethasone, Flunisolide, Budesonide Nasal Spray, Phenylephrine, Rebekol, Relenza, Tamiflu, and Tobraycin.

■ **SYMBOLS USED ON LABELS:**

	Catalogue No.		Caution
	Manufacturer		Consult instruction for use
	Manufacturing date		For single use only do not reuse
	Storage temperature		Keep away from direct sunlight
	In Vitro diagnostics		Do not use if box open or damaged
	Batch No.		European health and safety product label
	Expiry No.		Authorized European Representative in European Community
	Keep dry		This CE mark concerns sterile nasal swab
	Sufficient for		

■ **PRODUCT INFORMATION**

- **Website** – www.covifind.au

■ **GENERAL INFORMATION**

- **Manufacturer**
Meril Diagnostics Pvt Ltd
Second floor, D1-D3, Meril Park, Survey No. 135/2/B & 174/2, Muktanand Marg Chala, Vapi-396191 India
- **Sponsor Details**
Safe Interact Pty Ltd
ACN 653 504 708
26/155 Cooper Road, Yagoona 2199 NSW Australia
PH: 1800312940 +61
Email: covifind@safeinteract.com.au/covifind@merillife.com
Web: www.covifind.au

■ **IMPORTANT CONTACTS**

- **Australian Capital Territory Department of Health**
General enquiries: 02 5124 9213
Coronavirus helpline (8am to 8pm daily): 02 6207 7244
Website: <https://health.act.gov.au/>
- **New South Wales Department of Health**
General enquiries: 1300 066 055
Coronavirus hotline (Service NSW, 24/7): 137 788
Website: <https://www.health.nsw.gov.au/>
- **Northern Territory Department of Health**
General enquiries: 08 8922 8044
Coronavirus hotline (National helpline): 1800 020 080
Website: <https://health.nt.gov.au/>
- **Queensland Department of Health**
General enquiries: 13HEALTH or 13 432 584
Coronavirus hotline: 134COVID or 134 268
Website: <https://www.health.qld.gov.au/>
- **South Australian Department of Health**
General enquiries: 1300 232 272
Coronavirus hotline (9am to 5pm daily): 1800 253 787
Website: <https://www.sahealth.sa.gov.au/>
- **Tasmanian Department of Health**
General enquiries: 1300 135 513
Public Health Hotline (coronavirus): 1800 671 738
Website: <https://www.health.tas.gov.au/>
- **Victorian Department of Health**
Department of Health and Human Services
Phone: 1300 650 172
Victorian coronavirus hotline (24/7) : 1800 675 398
Website: <https://www.dhhs.vic.gov.au/>
- **Western Australian Department of Health**
General enquiries 08 9222 4222
Coronavirus hotline: 13COVID (8am to 6pm, Mon–Fri)
Website: <https://www.health.wa.gov.au/>

For assistance regarding the use of the product and interpretation of test results call 1800 312 940 +61. This service is available between 9 am and 7 pm (AEST), or 9 am and 8 pm (AEDT), 7 days per week.

In the event you are experiencing problems with the test, please contact Safe Interact.

Additionally, you may wish to report poor performance or usability issues directly to the Therapeutic Goods Administration (TGA) via the Medical Device Incident Reporting scheme, email iris@tga.gov.au or call 1800 809 361.