



Australian Government

Department of Health, Disability and Ageing
Therapeutic Goods Administration

Prohibited ingredients in vaping products by Gas Chromatography-Mass Spectrometry

June 2025

Copyright

© Commonwealth of Australia 2025

This work is copyright. You may reproduce the whole or part of this work in unaltered form for your own personal use or, if you are part of an organisation, for internal use within your organisation, but only if you or your organisation do not use the reproduction for any commercial purpose and retain this copyright notice and all disclaimer notices as part of that reproduction. Apart from rights to use as permitted by the *Copyright Act 1968* or allowed by this copyright notice, all other rights are reserved and you are not allowed to reproduce the whole or any part of this work in any way (electronic or otherwise) without first being given specific written permission from the Commonwealth to do so. Requests and inquiries concerning reproduction and rights are to be sent to the TGA Copyright Officer, Therapeutic Goods Administration, PO Box 100, Woden ACT 2606 or emailed to <tga.copyright@tga.gov.au>.

Contents

Instrument set up	4
Chromatographic system	4
MRM Mode: Mass Spectrometer System	4
Limits of detection	5
Solutions	5
Prohibited ingredients stock solution	5
Working Standard	5
Quality control sample	6
Limit standard	6
Concentrated standard	6
Sample solution	6
Injection sequence	6
System suitability criteria for the Limit solution	7
Identification criteria	7
Prohibited Ingredients: -----	7

Prohibited ingredients in vaping products by GCMS

This method describes the identification of the following prohibited ingredients in vaping products by Gas Chromatography-Mass Spectrometry (GCMS): Diacetyl; 2,3-Pentanedione; Acetoin; Ethylene glycol; Benzaldehyde; Diethylene glycol and trans-Cinnamaldehyde.

Instrument set up

Chromatographic system

Table 1: Chromatographic system set up

Instrument:	GCMS		
Column:	SH-624, 0.32 mm ID, 1.8 µm DF and 30 m L.		
Oven Temperature Program:	Rate	Temperature	Hold Time
	n/a	40 °C	3.00 min
	10 °C/min	210 °C	0.00 min
	25 °C/min	300 °C	1.40 min
Injection Temperature	250 °C		
Injection Mode	Split		
Split Ratio	20.0		
Column Flow	2.00 mL/min		

MRM Mode: Mass Spectrometer System

Table 2: Ion source and interface temperatures

Ion Source Temperature	200 °C
Interface Temperature	250 °C

Table 3: MRM event table for method: Quantifier transitions

Peak	Start time	End time	Precursor → Product	CE
Diacetyl	3.70	5.50	86→43	6
2,3-pentanedione	5.50	6.84	100→57	6
Acetoin	6.84	7.52	45→27	12
Ethylene glycol	7.52	10.00	62→33	3
Benzaldehyde	10.00	12.59	106→105	6
Diethylene glycol	12.59	15.80	75→45	6
Trans-Cinnamaldehyde	17.38	18.18	131→77	27

Table 4: MRM event table for method: Qualifier transitions

Peak	Precursor → Product	CE	Ion ratio ¹	Precursor → Product	CE	Ion ratio ¹
Diacetyl	86→86	0	3.0 ± 1.5	n/a	n/a	n/a
2,3-pentanedione	100→43	12	9.4 ± 2.8	n/a	n/a	n/a
Acetoin	45→29	9	106 ± 32	n/a	n/a	n/a
Ethylene glycol	62→31	6	55 ± 16	n/a	n/a	n/a

¹ The ion ratios in this table are indicative only, these should be updated using the Limit Standard or the Concentrated standard. Specifications: Ion ratio ≤ 5 (± 50% relative); Ion ratio > 5 (± 30% relative).

Peak	Precursor → Product	CE	Ion ratio ¹	Precursor → Product	CE	Ion ratio ¹
Benzaldehyde	106→77	16.5	50 ± 15	106→78	15	4.9 ± 2.4
Diethylene glycol	45→27	12	18 ± 6	45→29	9	15 ± 4
Trans-Cinnamaldehyde	131→103	9	45 ± 14	131→131	0	18 ± 5

Limits of detection

Table 5: Analyte Limits of detection

Analyte	Limit of Detection (µg/mL)
Diacetyl	0.3
2,3-pentanedione	0.3
Acetoin	0.2
Ethylene glycol	0.8
Benzaldehyde	0.2
Diethylene glycol	0.8
Trans-Cinnamaldehyde	0.2

Solutions

Prohibited ingredients stock solution

Prepare individual solutions in acetonitrile of the following analytes with the concentrations described in the Table 6. These solutions can be stored at 4 °C in an amber flask for up to 40 days.

Table 6: Analyte stock solution concentrations in acetonitrile

Analyte	Concentration (µg/mL)
Diacetyl	1000
2,3-Pentanedione	1000
Acetoin	1000
Benzaldehyde	1000
trans-Cinnamaldehyde	1000
Ethylene glycol	2500
Diethylene glycol	2500

Working Standard

Prepare a mixed standard using the stock solutions and methanol as the diluent, adjusting each individual aliquot to obtain the concentrations in Table 7. These solutions should be made fresh prior to testing.

Table 7: Analyte stock solution concentrations in acetonitrile

Analyte	Concentration (µg/mL)
Diacetyl	10
2,3-Pentanedione	10
Acetoin	10
Benzaldehyde	10
trans-Cinnamaldehyde	10
Ethylene glycol	25

Analyte	Concentration (µg/mL)
Diethylene glycol	25

Quality control sample

Prepare a mixture of glycerol and propylene glycol by measuring 7 mL of glycerol and 3 mL of propylene glycol and mixing the two until a homogeneous mixture is obtained (QC sample).

Using a suitable positive displacement piston operated volumetric apparatus (POVA), transfer 80 µL of the resulting mixture to an HPLC vial and then transfer 1520 µL of methanol (QC sample blank).

Limit standard

Using a suitable positive displacement POVA, transfer 80 µL of the QC sample to an HPLC vial and then transfer 80 µL of the Working Standard followed by addition of 1440 µL of methanol (QC sample – limit test solution).

Concentrated standard

Using a suitable positive displacement POVA, transfer 80 µL of the QC sample to an HPLC vial and then transfer 1520 µL of the Working Standard (Second identification standard).

Sample solution

Prepare a fresh 1 in 20 dilution of the sample in methanol by transferring 80 µL of the sample to an HPLC vial followed by 1520 µL of methanol using a suitable positive displacement POVA.

For samples that contain a large amount of acetoin which is not identified due to the ion ratios being saturated, a further dilution can be prepared to improve the ion ratio match.

Injection sequence

Inject in this order the following solutions:

Methanol	2 injections
QC Sample blank	2 injections
Limit Standard	2 injections
Concentrated standard	2 injections
Samples (up to 6 samples)	1 injection
Methanol	1 injection
QC Sample blank	1 injection
Limit Standard	1 injection
Concentrated standard	1 injection
Samples (up to 6 samples)	1 injection
Methanol	1 injection
QC Sample blank	1 injection
Limit Standard	1 injection

Concentrated standard

1 injection

End

System suitability criteria for the Limit solution

The %RSD of the area of each analyte for the Limit solution (injected at least 3 times) should be NMT 20%.

Identification criteria

Prohibited Ingredients:

Disregard any peaks with areas below the areas observed in the Limit Standard.

Disregard any peaks with retention time difference (absolute) more than 0.10 min compared with the same analyte in the Limit Standard.

Disregard any peaks where the ion ratio test fails, taking into consideration the ion ratios set up either by the Limit Standard or the Concentrated standard for identification.

A prohibited ingredient is present when the area of the peak is greater than the corresponding peak in the limit standard, and the retention time difference and ion ratio requirements are met.

Therapeutic Goods Administration

PO Box 100 Woden ACT 2606 Australia

Email: info@tga.gov.au Phone: 1800 020 653 Fax: 02 6203 1605

Web: tga.gov.au