



SARS-CoV-2 Antigen Rapid Test (Nasal Swab)

For Self-Testing REF: INCP-2020 | ENGBR

A rapid test for the qualitative detection of SARS-CoV-2 Nucleocapsid Protein antigens present in nasal swab specimen. For self-testing in vitro diagnostic use.

Read the instructions carefully before taking the test.


Australian Distributors: Touch Biotechnology Pty Ltd

Customer Support Number: 1300 166 282
Hours: 9am-7pm (AEST), or 9am-8pm (AEDT), 7 days per week
Website: www.touchbiotech.com.au
Email: touch@touchbiotech.com.au
Address: 4 Tabor Road, Macquarie Park, NSW 2113

Scan and Read the "How to Use" instructions on how to use the test.



COMPONENTS PROVIDED



Other kits has Tube holder on the box.

TEST PROCEDURE

STEP-1 Wash your hands

Wash your hands with soap and water for at least 20 seconds before and after test. If soap and water are not available, use hand sanitizer with at least 60% alcohol. Make sure they are dry before starting.

STEP-2 Read Instructions for use

Read instructions for use carefully before using the test.

STEP-3 Prepare for the test

- Check the expiration date on the box. Do not use the kit if it has been damaged or has expired.
- Ensure kit is at room temperature for at least 30 minutes prior to use. Open the box carefully as it will be used in a later step.
- Do not open individual components until instructed.

NOTE: A timing device (clock, timer, phone etc.) is required, but not provided.

STEP-4 Preparation

Remove the cover of the tube with Extraction Buffer and place the tube in the tube holder.

Remove the sterile swab from the pouch.

STEP-5 Specimen Collection

Insert the swab into your nostril until you feel slight resistance. Slowly twist the swab, rubbing it along the insides of your nostril for 5-10 times against the nasal wall.

NOTE: This may feel uncomfortable. Do not insert the swab any deeper if you feel strong resistance or pain. When the nasal mucosa is damaged or bleeding, avoid swab collection is advised.

If you are wearing others, please wear a face mask. With children, you may not need to insert the swab as far into the nostril. For very young persons to steady the child's head while swabbing.

Identify remove the swab.

Using the same swab, repeat step 2 in your other nostril.

Withdraw the swab.

STEP-6 Specimen Preparation

Place the swab into the Extraction tube, ensure it is touching the bottom and stir the swab to mix well. Press the swab and against the tube and rotate the swab for 10-15 seconds.

STEP-7 Take out the swab

Remove the swab by separating the swab from the inside of the Extraction tube.

Place the swab in a biosafety bag.

Close the cap or fit the tube tip onto the tube.


STEP-8 Take out the cassette

Remove the test cassette from the sealed foil pouch and use it within one hour. Best results will be obtained if the test is performed immediately after opening the foil pouch.

Place it on a flat and clean surface.

STEP-9 Test Operation

Invert the specimen extraction tube and add 3 drops of extracted specimen to the sample well (S) of the test cassette and start the timer.



STEP-10 Wait for result



Read the result at 15 minutes.



Do not read the result after 20 minutes.



STEP-11 Read your results

To read your test results, please go to the Interpretation of the results section provided below.

STEP-12 Disposal

Please dispose all parts of the test kits and place them in the biosafety bag that can be disposed in the bin. If there are local regulations, please follow them.

STEP-13 Wash your hands

Wash your hands thoroughly after test completion.


For Customer Support Helpline:
Call (03) 5986 5465
9am-7pm (AEST), 7 days per week
for information on the correct use of this test and for interpretation of the test results.

INTERPRETATION OF THE RESULTS



Positive

SARS-CoV-2 Positive



POSITIVE: * TWO DISTINCT COLORED LINES APPEAR.

ONE COLORED LINE SHOULD BE IN THE CONTROL REGION (C) AND ANOTHER COLORED LINE SHOULD BE IN THE TEST REGION (T).


*NOTE: The intensity of the color in the test line region (T) will vary based on the amount of SARS-CoV-2 antigen present in the sample. So any faint colored lines in the test region (T) should be considered positive.

A positive result means it is very likely you have COVID-19, but the positive samples should be confirmed to reflect this. If you have COVID-19 positive result, staying at home protects the people in your community and you should not visit high-risk settings like hospitals and aged care facilities. If you feel unwell or need COVID-19 advice for someone in your care, ask your health provider, or speak to a nurse by calling the medical helpline on 1800 022 222. If you develop symptoms such as severe aches, muscle pain, call triple zero (000) immediately.



Invalid

SARS-CoV-2 Invalid




CONTROL LINE FAILS TO APPEAR.

Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test or contact with your doctor or a COVID-19 test center.



Negative

SARS-CoV-2 Negative



ONE COLORED LINE APPEARS IN THE CONTROL REGION (C).

NO COLORED LINE APPEARS IN THE TEST LINE REGION (T).

You are unlikely to have COVID-19. However, it is possible for this test to give a negative result that is incorrect (a false negative) in some people with COVID-19. This means you could possibly still have COVID-19 even though the test is negative.

If you experience symptoms such as headaches, migraines, fever, loss of sense of smell or taste, contact the nearest Covid test centre according to the rules of your local authority. In addition, you can repeat the test with a new test kit. In case of suspicion, repeat the test after 1-2 days, as the coronavirus cannot be precisely detected in all phases of an infection. Even with a negative test result, distance and hygiene rules must be observed.

