



Before testing, scan the QR code to watch the how to use video, or Visit www.mycovidtest.com.au/how-it-works/
For additional language instructions please visit www.mycovidtest.com.au/product-resources

COVID-19 Antigen Rapid Test (Oral Fluid) Instruction Guide

BEFORE STARTING


Before collecting oral fluid, nothing is to be placed in the mouth including food, drink, gum, or tobacco products for at least 10 minutes prior to collection.



Was or sanitise your hands.

Make sure they are dry before starting.

MATERIALS PROVIDED



1 PREPARE FOR THE TEST

- 1A. Check the expiration date on the box. Do not use if the kit has expired.
- 1B. Ensure the kit is at room temperature for at least 30 minutes prior to use.

Open the box carefully as it will be used in a later step.

Do not open individual components until instructed.

Note: A time device (clock, timer, etc) is required, but not provided.


2 COLLECT SAMPLE

- 2A. Deeply cough 3 – 5 times


Note: Wear a face mask or cover your mouth and nose with a tissue when you are coughing and keep distance with other people.




2B.



3 ADD BUFFER



4 ADD DROPS



5 WAIT FOR RESULTS

Do not touch the Test Device during this period.

Read the result at 15 minutes.

Do not read the result earlier than 15 minutes or after 20 minutes.



Note: A control (C) line may appear in the result window within a few minutes, but a test (T) line may take as long as 15 minutes to appear.

Note: After 20 minutes the result might become inaccurate.

6 READ TEST RESULT

Please share your test result with your healthcare provider and carefully follow your local COVID guidelines/requirements

INVALID RESULT (Test did not work)

Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test or contact our COVID-19 test centre.



POSITIVE RESULT

One coloured line should be in the control region (C) and another coloured line should be in the Test region (T).

**Note: The intensity of the colour in the test line region (T) will vary based on the amount of SARS-(CoV-2) antigen present in the sample. So any shade of colour in the test region (T) should be considered positive.* A positive result means it is very likely you have COVID-19, but the positive samples should be confirmed to reflect this. Immediately go into self-isolation in accordance with the local guidelines and immediately contact your general practitioner/doctor or the local health department in accordance with the instructions of your local authorities.

Your test result will be checked by a PCR confirmation test and you will be explained the next steps.



For ALL Positive results a Confirmatory PCR test by a laboratory is required. Please contact your local COVID Help Line on the reverse side of these instructions.

NEGATIVE RESULT

No apparent coloured line appears in the test line region (T) You are unlikely to have COVID-19. However, it is possible for this test to give a negative result that is incorrect (a false negative) in some people with COVID-19. This means you could possibly still have COVID-19 even though the test is negative.

If you experience symptoms such as headaches, migraines, fever, loss of sense of smell or taste, contact the nearest medical facility according to the rules of our local authority. In addition, you can repeat the test with a new test kit. In case of suspicion repeat the test after 1-2 days, as the corona virus cannot be precisely detected in all phases of an infection. Even with a negative test result, distance yourself and hygiene rules must be observed, migration/traveling, attending events, etc, you should follow your local COVID guidelines/requirements.



7 DISPOSE THE TEST KIT

7A. After the test is complete, place all the components into the plastic Bio-Safety bag (supplied)

7B. Dispose according to local regulations

For Customer Support Call 1800 472 743 for information on the correct use of this test and for interpretation of the test results. Customer service hours are: 24 Hours , 7 Days