

# Bionik Resurfacing Femoral Head when used in conjunction with the Bionik Acetabular Component

Sponsor: **ESKA Australia**

Manufacturer: **ESKA**

NJRR Data:

# Implanted	# Revised	Revision Rate (%)
175	9	5.1
<i>13307</i>	<i>548</i>	<i>4.1</i>
Observed Comp Yrs	Revns/100 Comp Yrs	CL on revs/100 c.yrs
327	2.8	1.26 - 5.23
<i>55420</i>	<i>1.0</i>	<i>0.90 - 1.06</i>

*The numbers in shaded italics are the comparison figures for the same type of implant. In most cases this is the numbers for all implants of the same type received by the NJRR*

Number of implanting hospitals: 25

Number of hospitals where revisions occurred: 5

Reason for Revision	N	%
Loosening/Lysis	5	55.6
Fracture (Bone)	3	33.3
Pain	1	11.1
<b>Total</b>	<b>9</b>	<b>100%</b>

Type of Revision	N	%
Femoral and Acetabular	6	66.7
Acetabular Only	3	33.3
<b>Total</b>	<b>9</b>	<b>100%</b>

## **TGA Observations on NJRR Data**

Loosening/lysis and fracture are the main reasons for revision of the Bionik Implant, Loosening/lysis and pain are over-represented. That is the proportion of Bionik implants that are revised for these reasons is greater than the proportion of implants revised for this reasons in all other implants of the same type. Femoral neck fracture is a common cause for revision with this type of implant. Femoral neck fracture is considered to be related to the implant because preservation of the femoral neck is part of the design philosophy of this type of implant (The cumulative revision rate of the implant appears to be increasing and diverging from the revision rate curve for all other implants, but once again the certainty of this trend is difficult to establish due to the relatively low number of observed years.

## **TGA Observations on Manufacturers Reply**

There are 9 revisions reported by the NJRR in this series. The sponsor states that 1 was a ceramic on ceramic implant and the remaining 8 were metal on metal. The Sponsor dismisses them all as not being implant related:

The revision of the ceramic on ceramic implant was due to a femur fracture due to AVM.

Of the implants that the sponsor claims were metal on metal hips, 3 revisions were due to neck of femur fractures (2 of which are claimed to be due to trauma), 3 revisions are due to malposition of the acetabular cup and in two cases the femoral heads were not in the varus position.

Note that this does not completely account for all revisions and is not consistent with the NJRR data, which cites 3 fractures, not 4.

ESKA have also supplied some published papers and abstracts about the performance of the ESKA implant. The following are some observations about this literature.

In a general article about hip joint surface replacement Rudert et al report a case series of 20 Bionik surface replacement prostheses inserted between 2003 and 2005. At an average follow up period of 18 months there were no infections or aseptic loosening, but there was one femoral neck fracture and one dislocation, placing the revision rate at approximately 10% (or  $(100 \times 2 \text{ revisions} / 1.5 \text{ years} \times 20 \text{ implants}) = 6.6 \text{ revisions}/100 \text{ component years}$ ).

ESKA has provided a "Data Summary" on a series involving 248 patients (number of implants not stated). Enrolments began in February 2003 and patients were followed until February 2006 (estimated average follow up of 1.5 years). During that time ESKA reports that there were 7 revisions for various reasons - mostly femoral neck fracture. The revision rate is not calculated, but based on the information provided above, an estimate would be  $100 \times 7 \text{ revisions} / 1.5 \text{ years} \times 248 \text{ implants} = 1.88 \text{ revisions}/100 \text{ component years}$ .


Beaulé et al report a retrospective review of 94 cases for which the mean follow up was 4.2 years. 13 patients are reported to have had a bad outcome. A bad outcome is defined as conversion to THR, radiolucency of greater than 1mm or narrowing of the femoral neck by greater than 10%. It is not clear whether all 13 required revision, but if they did then the revision rate was 3.29 revisions/100 component years. To achieve the same revision rate as the average revision rate of similar implants in Australia, the number of implants that were revised should can be no greater than 4. However, Beaulé et al made an important observation: The number of failures is related to a neck shaft angle  $< 130^\circ$ . The relative risk of problems with Bionik hips where the neck shaft angle is  $< 130^\circ$  is 6 times that where the neck shaft angle is  $> 130^\circ$ . It is not clear whether this piece of information is conveyed to surgeons through product literature or training.

The paper by Gerdesmeyer et al on minimally invasive surgery reports that in a series of 31 patients using a minimally invasive approach, no instances of loosening or dislocation or other sequelae were observed after 12 months of follow up.

Two abstracts to papers are submitted with no accompanying citation... and are not discussed further here.

In summary the Sponsor asserts that none of the revisions reported in the NJRR against the Bionik implant are related to the design of the implant. ESKA has also provided papers and citations as evidence of implant performance elsewhere in the world, but the revision rates reported in the literature provided appears to be higher – sometimes much higher than the revision rate reported by the NJRR for this implant.

One author noted the importance of neck shaft angle - The TGA seeks the advice of the OEWG as to whether this affects all similar resurfacing implants and whether this is a commonly known in the orthopaedics field.



## Adaptor (cementless) Total Conventional Hip Femoral Prosthesis

This analysis compares the Adaptor (cementless) Total Conventional Hip Femoral Prosthesis with all Other Total Conventional Hip prostheses. This Prosthesis has been identified as having a significantly higher revision rate.

For a detailed explanation of the process used by the Registry that results in identification of prostheses that have a higher than anticipated rate of revision please refer to the 'Prostheses with Higher than Anticipated Rates of Revision' chapter of the most recent AOANJRR Annual Report, <http://www.dmac.adelaide.edu.au/aoanjrr/publications.jsp>.

**TABLE 1**

### *Revision Rate of Primary Total Conventional Hip Replacement*

The **Revision Rate** of the Adaptor (cementless) Total Conventional Hip Femoral Prosthesis is compared to all Other Total Conventional Hip prostheses.

**Table 1: Revision Rates of Primary Total Conventional Hip Replacement**

Component	N Revised	N Total	Obs. Years	Revisions/100 Obs. Yrs (95% CI)
Adaptor (cementless)	23	567	1092	2.11 (1.34, 3.16)
Other Total Conventional Hip	5054	170537	651226	0.78 (0.75, 0.80)
<b>TOTAL</b>	<b>5077</b>	<b>171104</b>	<b>652317</b>	<b>0.78 (0.76, 0.80)</b>

**TABLE 2**

### *Yearly Cumulative Percent Revision of Primary Total Conventional Hip Replacement*

The **Yearly Cumulative Percent Revision** of the Adaptor (cementless) Total Conventional Hip Femoral Prosthesis is compared to all Other Total Conventional Hip prostheses.

**Table 2: Yearly Cumulative Percent Revision of Primary Total Conventional Hip Replacement**

CPR	1 Yr	3 Yrs	5 Yrs	7 Yrs	9 Yrs
Adaptor (cementless)	2.7 (1.6, 4.5)	5.4 (3.5, 8.3)			
Other Total Conventional Hip	1.5 (1.5, 1.6)	2.7 (2.6, 2.8)	3.5 (3.4, 3.6)	4.4 (4.3, 4.5)	5.4 (5.1, 5.7)


**FIGURE 1**

**Yearly Cumulative Percent Revision of Primary Total Conventional Hip Replacement**

The **Yearly Cumulative Percent Revision** of the Adaptor (cementless) Total Conventional Hip Femoral Prosthesis is compared to all Other Total Conventional Hip prostheses. In addition, Hazard Ratios are also reported.

Hazard Ratios are reported for specific time periods during which the Hazard Ratio is constant. This is done to enable more specific and valid comparisons of the risk of revision over time. The pattern of variation in risk has important implications with respect to the underlying reasons for any difference.

**Figure 1: Cumulative Percent Revision of Primary Total Conventional Hip Replacement**



Number at Risk	0 Yr	1 Yrs	2 Yrs	3 Yrs	4 Yrs	5 Yrs	6 Yrs	7 Yrs	8 Yrs	9 Yrs
Adapter (cementless)	567	396	265	137	16	0	0	0	0	0
Other Total Conventional Hip	170537	142264	117871	96114	75956	57166	40037	24821	11582	2968

TABLE 3

**Primary Diagnosis for Revised Primary Total Conventional Hip Replacement**

This table identifies the diagnosis of the primary procedure which was subsequently revised. This information is provided as there is a variation on outcome depending on the primary diagnosis. It is therefore important when considering the reasons for a higher than anticipated rate of revision that there is identification of the primary diagnosis. This information should be compared to the primary diagnosis for the revisions of all Other Total Conventional Hip prostheses.

**Table 3: Primary Diagnosis for Revised Primary Total Conventional Hip Replacement**

Primary Diagnosis	Adapter (cementless)		Other Total Conventional Hip	
	Number	Percent	Number	Percent
Osteoarthritis	22	95.7	4262	84.3
Avascular Necrosis			253	5.0
Fractured Neck Of Femur	1	4.3	249	4.9
Rheumatoid Arthritis			96	1.9
Developmental Dysplasia			87	1.7
Failed Internal Fixation			32	0.6
Tumour			32	0.6
Other Inflammatory Arthritis			29	0.6
Fracture/Dislocation			7	0.1
Arthrodesis Takedown			5	0.1
Other			2	0.0
<b>TOTAL</b>	<b>23</b>	<b>100.0</b>	<b>5054</b>	<b>100.0</b>

TABLE 4

**Revision Rates of Adapter (cementless) Primary Total Conventional Hip Replacement by Fixation.**

This analysis is provided as some prostheses have more than one fixation option. Additionally there are prostheses where an alternative to the recommended approach to fixation was used e.g. a cementless prosthesis that has been cemented or vice-versa.

**Table 4: Revision Rates of Adapter (cementless) Primary Total Conventional Hip Replacement by Fixation**

Fixation	N Revised	N Total	Obs. Years	Revisions/100 Obs. Yrs (95% CI)
Cemented	0	1	3	0.00 (0.00, 115.0)
Cementless	23	550	1071	2.15 (1.36, 3.22)
Hybrid	0	16	18	0.00 (0.00, 21.04)
<b>TOTAL</b>	<b>23</b>	<b>567</b>	<b>1092</b>	<b>2.11 (1.34, 3.16)</b>

011

**TABLE 5**

**Revision Rates of Adaptor (cementless) Primary Total Conventional Hip Replacement by Bearing Surface.**

This analysis is provided as some prostheses are combined with a variety of different bearing surfaces. All bearing surfaces used with this Prosthesis are listed.

**Table 5: Revision Rates of Adaptor (cementless) Primary Total Conventional Hip Replacement by Bearing Surface**

Bearing Surface	N Revised	N Total	Obs. Years	Revisions/100 Obs. Yrs (95% CI)
Ceramic/Ceramic	2	38	78	2.56 (0.31, 9.26)
Ceramic/Polyethylene	1	53	99	1.01 (0.03, 5.65)
Ceramic/Modified Polyethylene	0	11	2	0.00 (0.00, 219.1)
Metal/Metal	17	356	759	2.24 (1.30, 3.59)
Metal/Polyethylene	3	105	150	2.00 (0.41, 5.84)
Metal/Modified Polyethylene	0	4	4	0.00 (0.00, 83.74)
<b>TOTAL</b>	<b>23</b>	<b>567</b>	<b>1092</b>	<b>2.11 (1.34, 3.16)</b>

**TABLE 6****Type of Revision Performed for Primary Total Conventional Hip Replacement**

This analysis identifies the components used in the revision of the Adaptor (cementless) Total Conventional Hip Femoral Prosthesis and compares it to the components used in the revision of all Other Total Conventional Hip prostheses.

The reason this analysis is undertaken is to identify whether there is one or more components which are being replaced that differ from the components replaced for revisions of all Other Total Conventional Hip prostheses i.e. is there a difference in the type of revision undertaken for the Adaptor (cementless) Total Conventional Hip Femoral Prosthesis compared to all Other Total Conventional Hip prostheses.

**Table 6: Primary Total Conventional Hip Replacement - Type of Revision**

Revision Type	Adaptor (cementless)		Other Total Conventional Hip	
	Number	Percent	Number	Percent
Femoral Only	6	26.1	1482	29.3
Acetabular Only	5	21.7	1260	24.9
THR (Femoral/Acetabular)	2	8.7	625	12.4
Cement Spacer	1	4.3	244	4.8
Removal of Prostheses			41	0.8
Reinsertion of Components			6	0.1
Bipolar Head and Femoral			2	0.0
Saddle			2	0.0
<b>N Major</b>	<b>14</b>	<b>60.9</b>	<b>3662</b>	<b>72.5</b>
Head/Insert	5	21.7	922	18.2
Head Only	1	4.3	285	5.6
Minor Components			82	1.6
Insert Only	1	4.3	70	1.4
Head/Neck	2	8.7	30	0.6
Neck Only			2	0.0
Neck/Insert			1	0.0
<b>N Minor</b>	<b>9</b>	<b>39.1</b>	<b>1392</b>	<b>27.5</b>
<b>TOTAL</b>	<b>23</b>	<b>100.0</b>	<b>5054</b>	<b>100.0</b>



**TABLE 7****Reason for Revision of Primary Total Conventional Hip Replacement**

This is reported in two ways; a percentage of all revisions and also as a percentage of all primary procedures.

This analysis includes a comparison of reasons for revision to all Other Total Conventional Hip prostheses.

This analysis is undertaken to identify if there are differences in the reasons for revision and the number of revisions performed for those reasons between the Adaptor (cementless) Total Conventional Hip Femoral Prosthesis and all Other Total Conventional Hip prostheses.

**Table 7: Primary Total Conventional Hip Replacement - Reason for Revision**

Revision Diagnosis	Adaptor (cementless)			Other Total Conventional Hip		
	Number	% Revision	% Primary	Number	% Revision	% Primary
Loosening/Lysis	5	21.7	0.9	1514	30.0	0.9
Prosthesis Dislocation	5	21.7	0.9	1395	27.6	0.8
Infection	3	13.0	0.5	843	16.7	0.5
Fracture	2	8.7	0.4	742	14.7	0.4
Pain	4	17.4	0.7	100	2.0	0.1
Other				67	1.3	0.0
Leg Length Discrepancy	1	4.3	0.2	63	1.2	0.0
Metal Sensitivity	1	4.3	0.2	61	1.2	0.0
Malposition				54	1.1	0.0
Wear Acetabulum	1	4.3	0.2	51	1.0	0.0
Implant Breakage Stem				35	0.7	0.0
Implant Breakage Acetabular	1	4.3	0.2	32	0.6	0.0
Incorrect Sizing				29	0.6	0.0
Instability				28	0.6	0.0
Implant Breakage Head				20	0.4	0.0
Heterotopic Bone				8	0.2	0.0
Tumour				8	0.2	0.0
Avascular Necrosis				2	0.0	0.0
Synovitis				2	0.0	0.0
<b>N Revision</b>	<b>23</b>	<b>100.0</b>	<b>4.1</b>	<b>5054</b>	<b>100.0</b>	<b>3.0</b>
<b>N Primary</b>	<b>567</b>			<b>170537</b>		

**FIGURE 2****Revision Diagnosis Cumulative Incidence by Time to Revision for Primary Total Conventional Hip Replacement**

This figure details the cumulative incidence of the most common reasons for revision.

The five most common reasons for revision are included as long as each of these reasons account for more than 10 procedures or at least 5% of all revisions for the Adaptor (cementless) Total Conventional Hip Prosthesis. A comparative graph is provided of the cumulative incidence for the same reasons for revisions for all Other Total Conventional Hip prostheses.


**Figure 2: Cumulative Incidence Revision Diagnosis for Primary Total Conventional Hip Replacement**

TABLE 8

**Revision Rates of Primary Total Conventional Hip Replacement by State**

This enables a state by state variation to be identified for the Adaptor (cementless) Total Conventional Hip Femoral Prosthesis and provides the comparative data for each of the states for all Other Total Conventional Hip prostheses.

This analysis is undertaken for similar reasons as those outlined above for Table 8.

**Table 8: Revision Rates of Primary Total Conventional Hip Replacement by State**

Component	State	N Revised	N Total	Obs. Years	Revisions/100 Obs. Yrs (95% CI)
Adaptor (cementless)	NSW	14	308	575	2.44 (1.33, 4.09)
Adaptor (cementless)	VIC	3	74	141	2.13 (0.44, 6.22)
Adaptor (cementless)	QLD	1	19	55	1.83 (0.05, 10.17)
Adaptor (cementless)	WA	0	32	20	0.00 (0.00, 18.44)
Adaptor (cementless)	SA	0	5	5	0.00 (0.00, 69.02)
Adaptor (cementless)	TAS	5	128	296	1.69 (0.55, 3.95)
Adaptor (cementless)	ACT/NT	0	1	0	0.00 (0.00, 1773)
Other Total Conventional Hip	NSW	1405	50347	180311	0.78 (0.74, 0.82)
Other Total Conventional Hip	VIC	1378	46239	176185	0.78 (0.74, 0.82)
Other Total Conventional Hip	QLD	867	27407	105373	0.82 (0.77, 0.88)
Other Total Conventional Hip	WA	647	18985	76117	0.85 (0.79, 0.92)
Other Total Conventional Hip	SA	462	17571	74846	0.62 (0.56, 0.68)
Other Total Conventional Hip	TAS	163	6120	24233	0.67 (0.57, 0.78)
Other Total Conventional Hip	ACT/NT	132	3868	14161	0.93 (0.78, 1.11)
<b>TOTAL</b>		<b>5077</b>	<b>171104</b>	<b>652317</b>	<b>0.78 (0.76, 0.80)</b>

TABLE 9

**Number of Revisions of Adaptor (cementless) Primary Total Conventional Hip Replacement by Year of Implant**

This analysis details the number of prostheses reported each year to the Registry for the Adaptor (cementless) Total Conventional Hip Prosthesis. It also provides the subsequent number of revisions of the primaries reported in that year.

Primary procedures performed in later years have had less follow up time therefore the number revised is expected to be less than the number revised in earlier years. For example, a primary procedure performed in 2009 has a maximum of one year to be revised, whereas a primary performed in 2007 has a maximum of three years to be revised.

**Table 9: Number of Revisions of Adaptor (cementless) Primary Total Conventional Hip Replacement by Year of Implant**

Year of Implant	Number Revised	Total Number
2005	1	19
2006	9	140
2007	6	131
2008	4	121
2009	3	156
<b>TOTAL</b>	<b>23</b>	<b>567</b>

TABLE 10

**Revision rates of Adaptor (cementless) Primary Total Conventional Hip Replacement by Catalogue Number Range**

Many prostheses have a number of catalogue ranges. The catalogue range is specific to particular design features; more than one catalogue range usually indicates a minor difference in design in a particular Adaptor (cementless) femoral prosthesis.

This analysis has been undertaken to determine if the revision rate varies according to the catalogue number range.

**Table 10: Revision Rates of Adaptor (cementless) Primary Total Conventional Hip Replacement by Catalogue Number Range**

Catalogue Range		Catalogue Description		
Femoral Stem				
Adapter (cementless)		11090302-11090373	ADAPTER CEMENTLESS COLLARLESS	
Adapter (cementless)		11099012-11099044	HIP STEM CEMENTLESS TINB-C COLLARLESS	
Femoral Stem Range	N Revised	N Total	Obs. Years	Revisions/100 Obs. Yrs (95% CI)
11090302-11090373	22	522	1061	2.07 (1.30, 3.14)
11099012-11099044	1	45	31	3.27 (0.08, 18.23)
TOTAL	23	567	1092	2.11 (1.34, 3.16)

TABLE 11

**Revision rates of Adaptor (cementless) Primary Total Conventional Hip Replacement by Component**

A prosthesis may be combined with multiple components.

This analysis has been undertaken to determine if the revision rate varies according to the component with which it is combined.

**Table 11: Revision Rates of Adaptor (cementless) Primary Total Conventional Hip Replacement by Acetabular Component**

Acetabular Component	N Revised	N Total	Obs. Years	Revisions/100 Obs. Yrs (95% CI)
2000	0	2	6	0.00 (0.00, 61.89)
2000 Plus	2	55	121	1.65 (0.20, 5.97)
Allofit	0	13	2	0.00 (0.00, 184.6)
Bionik	21	489	948	2.22 (1.37, 3.39)
Eska Cup	0	1	3	0.00 (0.00, 121.1)
R3	0	1	1	0.00 (0.00, 561.4)
Trabecular Metal Shell	0	5	10	0.00 (0.00, 37.95)
Trilogy	0	1	2	0.00 (0.00, 201.1)
<b>TOTAL</b>	<b>23</b>	<b>567</b>	<b>1092</b>	<b>2.11 (1.34, 3.16)</b>

## Bionik/Bionik Total Resurfacing Hip Investigation

This analysis compares the Bionik/Bionik Total Resurfacing Hip Combination with all Other Total Resurfacing Hip prostheses. This Combination has been identified as having a significantly higher revision rate.

For a detailed explanation of the process used by the Registry that results in identification of prostheses that have a higher than anticipated rate of revision please refer to the 'Prostheses with Higher than Anticipated Rates of Revision' chapter of the most recent AOANJRR Annual Report, <http://www.dmac.adelaide.edu.au/aoanjrr/publications.jsp>.

TABLE 1

### Revision Rate of Primary Total Resurfacing Hip Replacement

The **Revision Rate** of the Bionik/Bionik Total Resurfacing Hip Combination is compared to all Other Total Resurfacing Hip prostheses.

Table 1: Revision Rates of Primary Total Resurfacing Hip Replacement

Component	N Revised	N Total	Obs. Years	Revisions/100 Obs. Yrs (95% CI)
Bionik/Bionik	9	175	327	2.75 (1.26, 5.23)
Total Resurfacing Hip	539	13132	55093	0.98 (0.90, 1.06)
TOTAL	548	13307	55420	0.99 (0.91, 1.08)

TABLE 2

### Yearly Cumulative Percent Revision of Primary Total Resurfacing Hip Replacement

The **Yearly Cumulative Percent Revision** of the Bionik/Bionik Total Resurfacing Hip Combination is compared to all Other Total Resurfacing Hip prostheses.

Table 2: Yearly Cumulative Percent Revision of Primary Total Resurfacing Hip Replacement

CPR	1 Yr	3 Yrs	5 Yrs	7 Yrs	9 Yrs
Bionik/Bionik	3.8 (1.7, 8.4)	6.6 (3.2, 13.6)			
Total Resurfacing Hip	1.8 (1.6, 2.1)	3.2 (2.9, 3.5)	4.5 (4.1, 4.9)	6.1 (5.5, 6.7)	7.4 (6.4, 8.6)

**FIGURE 1**

**Yearly Cumulative Percent Revision of Primary Total Resurfacing Hip Replacement**

The **Yearly Cumulative Percent Revision** of the Bionik/Bionik Total Resurfacing Hip Combination is compared to all Other Total Resurfacing Hip prostheses. In addition, Hazard Ratios are also reported.

Hazard Ratios are reported for specific time periods during which the Hazard Ratio is constant. This is done to enable more specific and valid comparisons of the risk of revision over time. The pattern of variation in risk has important implications with respect to the underlying reasons for any difference.

**Figure 1: Cumulative Percent Revision of Primary Total Resurfacing Hip Replacement**




TABLE 3

**Primary Diagnosis for Revised Primary Total Resurfacing Hip Replacement**

This table identifies the diagnosis of the primary procedure which was subsequently revised. This information is provided as there is a variation on outcome depending on the primary diagnosis. It is therefore important when considering the reasons for a higher than anticipated rate of revision that there is identification of the primary diagnosis. This information should be compared to the primary diagnosis for the revisions of all Other Total Resurfacing Hip prostheses.

Table 3: Primary Diagnosis for Revised Primary Total Resurfacing Hip Replacement

Primary Diagnosis	Bionik/Bionik		Other Total Resurfacing Hip	
	Number	Percent	Number	Percent
Osteoarthritis	9	100.0	481	89.2
Developmental Dysplasia			35	6.5
Avascular Necrosis			13	2.4
Other Inflammatory Arthritis			6	1.1
Rheumatoid Arthritis			4	0.7
<b>TOTAL</b>	<b>9</b>	<b>100.0</b>	<b>539</b>	<b>100.0</b>

TABLE 4

**Revision Rates of Bionik/Bionik Primary Total Resurfacing Hip Replacement by Fixation.**

This analysis is provided as some prostheses have more than one fixation option. Additionally there are prostheses where an alternative to the recommended approach to fixation was used e.g. a cementless prosthesis that has been cemented or vice-versa.

Table 4: Revision Rates of Bionik/Bionik Primary Total Resurfacing Hip Replacement by Fixation

Fixation	N Revised	N Total	Obs. Years	Revisions/100 Obs. Yrs (95% CI)
Cemented	0	1	1	0.00 (0.00, 349.1)
Cementless	1	15	44	2.27 (0.06, 12.67)
Hybrid	8	159	282	2.84 (1.23, 5.59)
<b>TOTAL</b>	<b>9</b>	<b>175</b>	<b>327</b>	<b>2.75 (1.26, 5.23)</b>

**TABLE 5*****Type of Revision Performed for Primary Total Resurfacing Hip Replacement***

This analysis identifies the components used in the revision of the Bionik/Bionik Total Resurfacing Hip Combination and compares it to the components used in the revision of all Other Total Resurfacing Hip prostheses.

The reason this analysis is undertaken is to identify whether there is one or more components which are being replaced that differ from the components replaced for revisions of all Other Total Resurfacing Hip prostheses i.e. is there a difference in the type of revision undertaken for the Bionik/Bionik Total Resurfacing Hip Combination compared to all Other Total Resurfacing Hip prostheses.

**Table 5: Type of Revision for Primary Total Resurfacing Hip Replacement**

Revision Type	Bionik/Bionik		Other Total Resurfacing Hip	
	Number	Percent	Number	Percent
Femoral Only	3	33.3	284	52.7
THR (Femoral/Acetabular)	6	66.7	195	36.2
Acetabular Only			41	7.6
Cement Spacer			15	2.8
Removal of Prostheses			4	0.7
<b>N Major</b>	<b>9</b>	<b>100.0</b>	<b>539</b>	<b>100.0</b>
<b>TOTAL</b>	<b>9</b>	<b>100.0</b>	<b>539</b>	<b>100.0</b>



TABLE 6

**Reason for Revision of Primary Total Resurfacing Hip Replacement**

This is reported in two ways; a percentage of all revisions and also as a percentage of all primary procedures.

This analysis includes a comparison of reasons for revision to all Other Total Resurfacing Hip prostheses.

This analysis is undertaken to identify if there are differences in the reasons for revision and the number of revisions performed for those reasons between the Bionik/Bionik Total Resurfacing Hip Combination and all Other Total Resurfacing Hip prostheses.

**Table 6: Reason for Revision of Primary Total Resurfacing Hip Replacement**

Revision Diagnosis	Bionik/Bionik			Other Total Resurfacing Hip		
	Number	% Revision	% Primary	Number	% Revision	% Primary
Fracture	3	33.3	1.7	192	35.6	1.5
Loosening/Lysis	5	55.6	2.9	178	33.0	1.4
Infection				45	8.3	0.3
Metal Sensitivity				39	7.2	0.3
Pain	1	11.1	0.6	28	5.2	0.2
Avascular Necrosis				17	3.2	0.1
Prosthesis Dislocation				15	2.8	0.1
Malposition				12	2.2	0.1
Other				6	1.1	0.0
Implant Breakage Head				3	0.6	0.0
Instability				1	0.2	0.0
Leg Length Discrepancy				1	0.2	0.0
Synovitis				1	0.2	0.0
Tumour				1	0.2	0.0
<b>N Revision</b>	<b>9</b>	<b>100.0</b>	<b>5.1</b>	<b>539</b>	<b>100.0</b>	<b>4.1</b>
<b>N Primary</b>	<b>175</b>			<b>13132</b>		

**FIGURE 2**

**Revision Diagnosis Cumulative Incidence by Time to Revision for Primary Total Resurfacing Hip Replacement**

This figure details the cumulative incidence of the most common reasons for revision.

The five most common reasons for revision are included as long as each of these reasons account for more than 10 procedures or at least 5% of all revisions for the Bionik/Bionik Total Resurfacing Hip Combination. A comparative graph is provided of the cumulative incidence for the same reasons for revisions for all Other Total Resurfacing Hip prostheses.

**Figure 2: Revision Diagnosis Cumulative Incidence by Time to Revision for Primary Total Resurfacing Hip Replacement**




TABLE 7

**Revision Rates of Bionik/Bionik Primary Total Resurfacing Hip Replacement by Hospital**

This table details the rates of revision in each of the individual hospitals in which the Bionik/Bionik Total Resurfacing Hip Combination was used. The hospitals are identified by number only.

The purpose of this analysis is to determine if the higher than anticipated rate of revision has widespread distribution between hospitals. If there is widespread distribution then the reason for the higher than anticipated rate of revision is unlikely to be surgeon specific. If the prosthesis has been used in only a small number of hospitals it is not possible to distinguish if the higher than anticipated rate of revision is related to the prosthesis, surgeon, technique or patient.

**Table 7: Revision Rates of Bionik/Bionik Primary Total Resurfacing Hip Replacement by Hospital**

Hospital Number	N Revised	N Total	Obs. Years	Revisions/100 Obs. Yrs (95% CI)
1	0	2	7	0.00 (0.00, 52.96)
2	0	21	16	0.00 (0.00, 23.13)
3	0	4	11	0.00 (0.00, 34.82)
4	2	6	15	12.92 (1.56, 46.66)
5	0	3	4	0.00 (0.00, 85.44)
6	0	1	1	0.00 (0.00, 528.4)
7	0	12	35	0.00 (0.00, 10.60)
8	0	2	3	0.00 (0.00, 141.5)
9	0	3	3	0.00 (0.00, 112.6)
10	2	10	11	17.83 (2.16, 64.41)
11	0	4	5	0.00 (0.00, 67.30)
12	0	4	6	0.00 (0.00, 60.88)
13	0	1	2	0.00 (0.00, 168.2)
14	0	22	51	0.00 (0.00, 7.17)
15	0	1	0	0.00 (0.00, 8982)
16	0	3	3	0.00 (0.00, 108.9)
17	0	2	4	0.00 (0.00, 98.93)
18	0	2	5	0.00 (0.00, 75.40)
19	1	25	53	1.90 (0.05, 10.61)
20	0	1	0	0.00 (0.00, 2750)
21	1	15	11	9.04 (0.23, 50.36)
22	3	23	68	4.39 (0.91, 12.83)
23	0	1	1	0.00 (0.00, 344.6)
24	0	5	6	0.00 (0.00, 63.17)
25	0	2	5	0.00 (0.00, 75.10)
<b>TOTAL</b>	<b>9</b>	<b>175</b>	<b>327</b>	<b>2.75 (1.26, 5.23)</b>

**TABLE 8****Revision Rates of Primary Total Resurfacing Hip Replacement by State**

This enables a state by state variation to be identified for the Bionik/Bionik Total Resurfacing Hip Combination and provides the comparative data for each of the states for all Other Total Resurfacing Hip prostheses.

This analysis is undertaken for similar reasons as those outlined above for Table 7.

**Table 8: Revision Rates of Primary Total Resurfacing Hip Replacement by State**

Component	State	N Revised	N Total	Obs. Years	Revisions/100 Obs. Yrs (95% CI)
Bionik/Bionik	NSW	7	118	225	3.12 (1.25, 6.42)
Bionik/Bionik	VIC	0	24	53	0.00 (0.00, 6.91)
Bionik/Bionik	QLD	0	2	7	0.00 (0.00, 52.96)
Bionik/Bionik	WA	0	21	16	0.00 (0.00, 23.13)
Bionik/Bionik	TAS	2	10	26	7.67 (0.93, 27.70)
Total Resurfacing Hip	NSW	180	3899	15823	1.14 (0.98, 1.32)
Total Resurfacing Hip	VIC	189	4747	21507	0.88 (0.76, 1.01)
Total Resurfacing Hip	QLD	80	2212	8156	0.98 (0.78, 1.22)
Total Resurfacing Hip	WA	16	370	1845	0.87 (0.50, 1.41)
Total Resurfacing Hip	SA	49	1268	5560	0.88 (0.65, 1.17)
Total Resurfacing Hip	TAS	9	83	279	3.22 (1.47, 6.12)
Total Resurfacing Hip	ACT/NT	16	553	1923	0.83 (0.48, 1.35)
<b>TOTAL</b>		<b>548</b>	<b>13307</b>	<b>55420</b>	<b>0.99 (0.91, 1.08)</b>

TABLE 9

**Number of Revisions of Bionik/Bionik Primary Total Resurfacing Hip Replacement by Year of Implant**

This analysis details the number of prostheses reported each year to the Registry for the Bionik/Bionik Total Resurfacing Hip Combination. It also provides the subsequent number of revisions of the primaries reported in that year.

Primary procedures performed in later years have had less follow up time therefore the number revised is expected to be less than the number revised in earlier years. For example, a primary procedure performed in 2009 has a maximum of one year to be revised, whereas a primary performed in 2007 has a maximum of three years to be revised.

**Table 9: Number of Revisions of Bionik/Bionik Primary Total Resurfacing Hip Replacement by Year of Implant**

Year of Implant	Number Revised	Total Number
2005	0	12
2006	3	33
2007	4	33
2008	0	43
2009	2	54
<b>TOTAL</b>	<b>9</b>	<b>175</b>

TABLE 10

**Revision rates of Bionik/Bionik Primary Total Resurfacing Hip Replacement by Catalogue Number Range**

Many prostheses have a number of catalogue ranges. The catalogue range is specific to particular design features; more than one catalogue range usually indicates a minor difference in design in a particular Bionik/Bionik prosthesis.

This analysis has been undertaken to determine if the revision rate varies according to the catalogue number range.

**Table 10: Revision Rates of Bionik/Bionik Primary Total Resurfacing Hip Replacement by Catalogue Number Range**

Catalogue Range		Catalogue Description			
Head					
Bionik	10260050-10260050	HIP RESURFACING CEMENTLESS SILVER			
Bionik	10270042-10270058	HIP RESURFING CEMENTED SILVER			
Bionik	10270238-10270258	CERAMIC HEAD HIP SURF. REPLACEMENT			
Bionik	10280046-10280050	HIP RESURFACING CEMENTLESS SILVER			
Bionik	10280142-10280156	HIP SURFACE REPLACEMENT CEMENTED BIOSURF SILVER			
Bionik	10280642-10280654	HIP RESURFACING CEMENTLESS SILVER CAP-COAT			
Bionik	10282038-10282058	FEMORAL HEAD SHELL BIOSURF CEMENTED			
Acetabular					
Bionik	10201050-10201064	METAL SHELL BS TINB COAT			
Bionik	10201150-10201164	METAL SHELL BS TINB COAT SCREW FIX			
Bionik	10201248-10201264	METAL SHELL CEMENTLESS TINB COAT SCREW FIX			
Bionik	10201346-10201366	METAL SHELL TINB CAP SCREW FIX			
Head Range	Acetab Range	N Revised	N Total	Obs. Years	Revisions/100 Obs. Yrs (95% CI)
10260050-10260050	10201150-10201164	0	1	4	0.00 (0.00, 85.01)
10270042-10270058	10201150-10201164	0	6	23	0.00 (0.00, 15.72)
10270042-10270058	10201248-10201264	0	4	16	0.00 (0.00, 23.28)
10270042-10270058	10201346-10201366	2	11	36	5.51 (0.67, 19.90)
10270238-10270258	10201150-10201164	0	1	1	0.00 (0.00, 343.7)
10270238-10270258	10201248-10201264	0	1	1	0.00 (0.00, 471.1)
10270238-10270258	10201346-10201366	1	9	7	14.32 (0.36, 79.77)
10280046-10280050	10201050-10201064	0	1	4	0.00 (0.00, 83.69)
10280046-10280050	10201150-10201164	0	1	4	0.00 (0.00, 83.90)
10280142-10280156	10201050-10201064	0	1	3	0.00 (0.00, 107.4)
10280142-10280156	10201150-10201164	0	3	10	0.00 (0.00, 38.45)
10280142-10280156	10201248-10201264	0	10	22	0.00 (0.00, 16.80)
10280142-10280156	10201346-10201366	2	46	87	2.30 (0.28, 8.31)
10280642-10280654	10201050-10201064	0	2	7	0.00 (0.00, 50.24)
10280642-10280654	10201150-10201164	1	2	7	15.23 (0.39, 84.83)
10280642-10280654	10201248-10201264	0	3	8	0.00 (0.00, 46.30)
10280642-10280654	10201346-10201366	0	7	18	0.00 (0.00, 20.84)
10282038-10282058	10201150-10201164	1	3	6	15.58 (0.39, 86.82)
10282038-10282058	10201248-10201264	1	36	31	3.18 (0.08, 17.73)
10282038-10282058	10201346-10201366	1	27	30	3.34 (0.08, 18.58)
TOTAL		9	175	327	2.75 (1.26, 5.23)