ORTHODYNAMICS

[REDACTED]
 Chief Biomaterials Scientist
 Director, Biomaterials and Engineering Section
 Office of Laboratories and Scientific Services
 Department of Health and Ageing
 Therapeutic Goods Administration
 PO Box 100
 Woden ACT 2606

Orthodynamics GmbH
 Grapengießerstr. 34
 23556 Lübeck
 Germany
Telefon:
 +49 (0) 451 89000-0
Fax:
 +49 (0) 451 89000-40
E-Mail:
 office@orthodynamics.de
Internet:
 www.orthodynamics.de

25th February 2011

Dear [REDACTED]

Re: ESKA Adapter (cementless) Femoral Stem Prosthesis and ESKA Bionik Resurfacing Femoral Head

I am writing in response to your letter to ESKA Australia of the 17th January 2011 regarding your observation of higher than expected revision rates for the ESKA Adapter (cementless) Femoral Stem Prosthesis and the ESKA Bionik Resurfacing Femoral Head.

I would firstly like to thank you for granting an extension of the reporting period to the 28th February, the additional time has provided us with the opportunity to carry out a thorough review of the issues raised in your letter of the 17th January to [REDACTED]. I would also take this opportunity to advise you that information supplied by ESKA Australia in September 2010 had not been verified or reviewed by Orthodynamics GmbH.

Whilst ESKA Australia is the distributor of the implants identified in your letters, you may be aware that the manufacturer of these devices was ESKA Implants AG, although this company ceased to exist on 31st March 2010.

The assets of ESKA Implants AG were acquired by The Summit Medical Group on the 31st March 2010, and a new entity was formed "Eska Orthodynamics GmbH". This business was subsequently rebranded as "Orthodynamics GmbH" in October 2010. Orthodynamics GmbH continues to manufacture and supply the above devices to the Australian market via our distributor ESKA Australia.

In response to your invitation to further review and expand on information previously supplied, the executive management at Orthodynamics GmbH commissioned a comprehensive review of the products currently being placed on the market in Australia. This review was conducted by Regulatory Affairs, Design and Development and Customer Services Managers at our facility in Lübeck, Germany. Additional product testing by [REDACTED] Director of Surgery and Orthopaedic Laboratories at UNSW, was also commissioned. ESKA Australia is currently in liaison with [REDACTED] in order to obtain the final version of this test report; it will be forwarded to you as soon as it becomes available. The main report attached does reference the findings of this test report, based on early communication of the results from [REDACTED].

O R T H O D Y N A M I C S

It is our view that this report provides data and reasoned arguments to demonstrate that these devices provide a safe and effective solution in the management of clinical procedures for which they are indicated and as such should remain available to surgeon and patients in the Australian market.

Should you have any questions regarding the report, or require any further information, please do not hesitate to contact me directly.

Yours sincerely,

[Redacted signature]

[Redacted name]

[Redacted title]

[Redacted address line]

c.c. [Redacted] ESKA Australia

Investigation into Australian Orthopaedic Association National Joint Replacement Registry Data for the ESKA Adapter (cementless) Femoral Stem Prosthesis and the ESKA Bionik Resurfacing Femoral Head

The Australian Orthopaedic Association National Joint Replacement Registry (NJRR) demonstrates a higher than expected revision rate for the ESKA Adapter (cementless) Femoral Stem and Bionik Resurfacing Femoral Head.

Interpretation of Adapter stem joint registry data

The Adapter stem has a revision rate of 2.11 revisions/100 obs. years compared to 0.78 revisions/100 obs. years for conventional total hip arthroplasty¹.

Detailed investigation into the Adapter stem NJRR data shows that the majority of its usage in Australia (63%) has been with a large metal-on-metal bearing combination. The metal-on-metal bearing also accounts for 17 of the 23 revisions (74%). It is well documented throughout the industry that this bearing combination has been associated with higher revision rates. This has led to a Medical Device Alert for all metal-on-metal bearings in the UK², and increased vigilance by the FDA³.

The UK National Joint Registry demonstrates a revision rate of 7.7% at five years with large head metal-on-metal bearings compared to 3.4% for cementless total hip arthroplasty⁴. This compares to an overall revision rate of 4.1% for the Adapter stem between 2005 and 2009 as reported by the NJRR (4.8% for metal-on-metal bearings, and 2.8% for cementless total hip arthroplasty).

It could be argued that it is therefore slightly misleading to compare the performance of the Adapter stem to conventional total hip arthroplasty. The UK National Joint Registry demonstrates that the Adapter stem falls within reported revision rates at five years when compared with other large metal head bearings, and cementless hip arthroplasty (both of which are demonstrated to have a higher revision rate than conventional cemented hip arthroplasty).

However, there is also reason to believe that there may be design factors associated with the Bionik femoral head which could account for a higher than expected revision rate. In July 2007, the external geometry of the Bionik femoral head was modified upon feedback from the Australian surgeons to avoid irritation of the iliopsoas tendon⁵. This was followed by a further design change of the Biosurf surface in September 2008 to improve wear results (this is discussed further in the Bionik resurfacing head section). A review of the NJRR data for cementless adaptor stems and large Bionik femoral heads⁶ shows a strong correlation between stem and head revisions, and a reduction in revision rate since the implementation of the latest design in 2008⁷.

Reviewing the reasons for revision in the NJRR demonstrates that the Adapter stem's performance is directly comparable to other total conventional hip prostheses in Australia for loosening, dislocation, infection and fracture (the four most common causes of revision). The incidence of revision due to pain is higher with the Adapter stem, but that could be attributed to many factors including wear particles from the metal-on-metal articulation rather than purely complications with the stem (given that three quarters of the revisions are on joint replacements with a metal-on-metal articulation).

This argument is strengthened when investigating the types of revision. This demonstrates that there were eight femoral stem revisions (excluding one case of infection which cannot be attributed to the implant). The remaining fourteen revisions all involved replacement of bearing components, rather than the stem. This gives the Adapter stem a revision rate of 0.73%/100 obs. years which is below the 0.78% reported for conventional total hip arthroplasty in the NJRR.

There was also a significant design change to the liner which was implemented into the Australian market in February 2009⁸. This design change introduced a locating pin to the liner to ensure that the liner is centrally positioned within the metal shell. This design change is also likely to have contributed to the decline in revision rates reported by the NJRR in 2009.

The Adapter stem has been widely used in the German market, in far larger numbers than those reported in Australia, with very low incidence of complaints and adverse incidents⁹. However, the use of the Adapter stem with large metal heads is not common in the German market, which may account for the discrepancy of results when compared to those witnessed in the Australian market and described in the NJRR.

Adapter stem clinical literature

Published clinical literature to support the Adapter stem is limited, although one abstract presented at the 56th Annual Meeting of the Association of South German Orthopaedic Surgeons in Baden-Baden May 2008 demonstrated good results¹⁰. This followed up 42 surgeries to a period of up to four years with excellent clinical outcomes, as measured by the Harris Hip Score of 98, and 100% survivorship.

The same design of taper cone adaptor is also used on the G2 Adapter stem marketed in Europe. The G2 Adapter stem is very similar in design to the GHE Adapter stem, and can be considered equivalent when looking at the stem/modular taper interface¹¹. Zurstege¹² followed up 100 stems over four years and concluded that the taper/stem connection was effective with no evidence of fractures or loosening.

As well as the Adapter stem, Orthodynamics GmbH manufactures a monoblock variant which is identical in geometry to the Adapter stem, but without the modular neck taper. This stem has clinical publications of larger cohorts which demonstrate that the stem design is clinically safe and effective.

Kinner et al¹³ followed up 200 surgeries up to four years and showed three revisions, none of which were directly attributable to the stem. Patient outcome was good with a high Harris Hip Score of 97, and the three revisions account for a revision rate of 1.5%. Ogawa et al¹⁴ followed up 48 stems for a period of just over six years. They witnessed no revisions and an average Merle d'Aubigne score of 16.1. Similarly, Kusaba et al¹⁵ also followed up a cohort of 166 patients for a period of just over six years. They also witnessed no revisions, and good clinical outcomes with a Harris Hip Score of 91.

Overall, these three studies of the monoblock GHE stem followed up over 400 patients and witnessed just three revisions at a minimum of four years follow-up (none of which could be attributed to the stem). This amounts to a low revision rate of 0.7% for the monoblock variant of the Adapter stem.

As well as the published clinical data in support of the Adapter stem and its variants, Orthodynamics GmbH has also obtained surgical testimonies from experienced Adapter stem users in our domestic market^{16,17}. This is further evidence that the Adapter stem is achieving satisfactory clinical results, although again it should be stressed that these users are not using the stem in combination with metal-on-metal bearings as is often the case in Australia.

Adapter stem mechanical testing

A review of the pre-clinical testing for Adapter stems showed that whilst the worst case stem, taper, head combination had been fatigue tested to ISO 7206-6, the testing had not been conducted in saline. Upon reviewing the NJRR data, Orthodynamics GmbH commissioned a further fatigue test in saline to determine if fretting corrosion could lead to mechanical failures of the Adapter stem. This test¹⁸ was performed by an eminent, independent party in Australia, Professor Bill Walsh of the University of New South Wales, to 3.5kN and 10m cycles (beyond the ISO 7206-6 requirements).

The results showed that the stem exceeded the fatigue requirements of ISO 7206-6. Analysis of the taper adapter after testing showed that there was evidence of fretting, particularly on the medial side of the taper adapter, and the corresponding area of the female femoral stem connection. However, fretting is observed in all modular taper connections after fatigue testing. The amount of fretting documented varies with Kretzer et al¹⁹ observing low amounts of Titanium release (in the region of μg) and Traina²⁰ et al witnessing Titanium release in the region of 0.6mg/million cycles. Both studies conclude that the amount of wear debris and associated risk of fretting is acceptable, given the clinical benefits of modularity. Grupp et al²¹ demonstrate that a Cobalt chrome modular neck adapter, as utilised with the Adapter stem, offers significant reductions in fretting when compared to Titanium neck adapters. A subjective comparison of fretting witnessed with these accepted levels in the published literature shows no evidence of increased risk associated with the Adapter stem. The level of wear particles is certainly low when compared to those witnessed from metal-on-metal bearings²² (1mm^3 is equivalent to 8.28mg).

Adapter stem conclusion

A review of the joint registry data shows there to be identifiable factors and reasons to explain the higher than expected rate of revision for the Adapter stem when compared to conventional hip arthroplasty.

Firstly, the high revision rate is believed to be linked to the comparatively high metal-on-metal usage in Australia, and problems with earlier designs of large metal heads irritating the iliopsoas tendon and potential high wear rates. However, the five year results are favourable when compared to the UK National Joint Registry data where there has been large metal-on-metal implant usage, and there is a clear downward trend in revisions being recorded by the NJRR as a consequence of the design changes to the Bionik femoral head and liner.

Aside from pain, which is believed to be linked to metal-on-metal wear particles from the metal-on-metal articulation, the Adapter stem has a comparable record with the conventional hip replacement results reported in the NJRR. Breaking down the data to femoral stem revisions, shows a revision rate of 0.73 revisions/100 obs. years (lower than the standard for conventional hip arthroplasty of 0.78 revisions/100 obs. years).

The successful clinical performance of the stem is supported by good clinical results and published data in Europe, with lower revision rates than those witnessed in Australia. Again, this is believed to be attributable to the fact that the stem is not used with a large metal-on-metal bearing in Europe.

Whilst the stem/adapter stem interface has shown evidence of fretting during mechanical testing, this is believed to be comparable to similar products and acceptable for a modular implant. Wear rates associated with taper junctions are very low for cobalt chrome taper adapters when compared to the wear rates seen in metal-on-metal articulations.

The NJRR and UK NJR analysis, the published literature, the surgeon testimonies and the fatigue test results all demonstrate that the Adapter stem is clinically safe and effective when used for its intended purpose of primary hip arthroplasty. The NJRR data which demonstrates a higher than expected revision rate for the Adapter stem has been distorted by metal-on-metal failures which can be attributed to a former design of the femoral head component and liner.

Given the low risk associated with the Adapter stem, and the clinical benefits of the modular neck system, it is Orthodynamics GmbH's belief that there is sufficient data to enable continued usage of the device in the Australian market.

Interpretation of Bionik resurfacing head registry data

The Bionik resurfacing head has a revision rate of 2.75 revisions/100 obs. years compared to 0.99 revisions/100 obs. years for all hip resurfacing²³.

The Bionik resurfacing system data supplied by the NJRR covers four different types of femoral component²⁴: a conventional hip resurfacing component, a conventional hip resurfacing component with "BioSurf", a femoral head shell "Onlay" with "Biosurf", and a Cera-Metal device (no longer supplied to the Australian market).

The Cera-Metal device was a ceramic coated metal head, but it has never been manufactured or supplied by Orthodynamics GmbH and is no longer available on the market worldwide. There was one Cera-Metal failure in the market, but the revision rate was 11.11 revisions/100 obs. years. This figure obviously distorts the remaining Bionik resurfacing data and should be discounted. The revision rate for the other Bionik resurfacing variants is 2.52 revisions/100 obs. years.

The data is also potentially skewed by a 20% revision rate in Tasmania (2 revisions from 10 implantations, at 7.67 revisions/100 obs. years). Although it is clear that the revision rate is higher than expected when compared to other resurfacing data from the NJRR, there have only been nine revisions of the Bionik resurfacing head. New South Wales and Tasmania account for 131 of 178 implantations (73%), and have witnessed all of the revisions. Even discounting the Bionik data, these

centres have the highest resurfacing revision rates of 1.11 revisions/100 obs. years and 2.77 revisions/100 obs. years respectively.

The UK National Joint Registry reports a resurfacing revision rate of 6.2% at five years⁴. The NJRR data shows that the Bionik resurfacing head has a revision rate of 5.1%.

The NJRR data also demonstrates that all of the revisions involved the femoral component, and that conventional resurfacing Biosurf and femoral shell Onlay Biosurf components account for 75% of the revisions excluding Cera-Metal (3 conventional revisions and 3 Onlay revisions). The dimpled Biosurf surface is identical on the conventional resurfacing and Onlay devices, and is also used on the Bionik femoral head. It is believed that earlier design iterations of the Biosurf components could have led to higher than expected revision rates due to high wear (see mechanical testing section below for further information). The earlier design of liners without the centralising peg could also be a factor causing slight elevations in wear. These factors were also discussed as reasons for the higher than expected revision rates of the Adapter stem.

The case for metal wear particles being a cause for revision is strengthened by the fact that 5 of the 9 revisions (56%) are due to loosening or lysis. The other predominant failure rate is fracture; this was the cause of 3 of the 9 revisions (33%). However, the incidence of fractures is in line with those witnessed in the NJRR for resurfacing at around 1.5%. The cause of fractures may be more related to the general resurfacing surgical technique rather than specific implant features.

The resurfacing components have been widely sold into the German market, but currently the revision rate appears to be much lower than that reported in Australia²⁵. 9 of 181 stems sold into Australia have been revised (4.97%) compared to 17 of 1973 worldwide (0.86%).

Bionik resurfacing head clinical literature

The clinical literature for the conventional resurfacing and Onlay resurfacing prostheses is good. Gerdesmeyer's most extensive follow-up is of 104 surgeries at 3 years²⁶. The follow-up is conducted on an Onlay device with Biosurf and shows excellent clinical outcomes with a good Harris Hip Score of 98 and a revision rate of 0.96%. The one revision is due to femoral neck fracture, and therefore not strictly implant related. This relates favourably to the revision rate reported in the NJRR data.

Gerdesmeyer et al²⁷ performed a second study that investigates the same device but with a minimally invasive approach. This time 31 surgeries were followed up at one year. Again, the clinical outcome was very good with an excellent Harris Hip Score of 97 and no revisions.

Rudert et al²⁸ followed up 20 patients with the conventional hip resurfacing and Biosurf to 1.5 years. The clinical outcomes were satisfactory with Harris Hip Scores of 92, but a high revision rate of 10% was reported. It is believed that this data is skewed by the low numbers in the study. Of the two revisions, one was due to dislocation and the other was due to femoral neck fracture. There were no implant failures due to loosening as reported in the NJRR.

Literature suggests that implant positioning and alignment is of increased importance during hip resurfacing procedures, and this can significantly affect the long-term results²⁹. This may also be a factor in the differences in revision rates reported in Germany compared to Australia.

Bionik resurfacing head mechanical testing

The Bionik femoral and resurfacing heads were brought to market on the basis of a successful hip simulator wear test to ISO 14242 in 2003³⁰. Conventional resurfacing parts were used, one with Biosurf and one without. The results of the test were 2.1mg/million cycles for the Biosurf and 5.5mg/million cycles for the non-dimpled surface. These were converted into volumetric wear rates of 0.25mm³/year and 0.66mm³/year, and shown to be comparable with the lower end of reported wear rates in the literature²¹. Based upon the successful wear test, established cobalt chrome bearing materials, and established manufacturing tolerances³¹, the Bionik femoral heads were CE marked.

Following design changes to the external geometry of the Biosurf surface, a further wear test was commissioned in 2007³². This test was performed using three conventional femoral heads. The results identified a higher mean wear rate of 10.83mg/million cycles which equates to 1.31mm³/year. These wear results are in the middle range of the reported values in the literature. The results were also inconsistent with a large standard deviation.

Concern about these results led to an improvement in polishing of the Biosurf components and a further test of four components³³. This testing demonstrated a significantly higher wear rate again, and a wider range of results. The wear rate varied from 2.2mg/million cycles to 52.2mg/million cycles. Conversion of these results gives a mean volumetric wear rate of 3.29mm³/year which is beyond the levels reported in the published literature²².

Due to the wide range of results witnessed, it is not clear whether this problem related to all parts that had been previously manufactured or just the latest design iteration. The devices were left on the market by ESKA Implants AG until a design change was approved in September 2008. This remains the latest design iteration, and it was approved on the basis of achieving three consistently low results from the wear simulator testing for both surfaces^{34,35}. These latest wear results of 1.51mg/million cycles for the Biosurf and 1.81mg/million cycles (0.18mm³/year and 0.39mm³/year) for the conventional surface were lower than those seen in the original testing, and are at the low end of reported results in the literature.

The design was improved by moving the region of the dimples closer to the equator and reducing the number and depth of the dimples³⁶. Following the successful wear results, the latest design of device replaced all previous issue stock in the ESKA Implants warehouse and all previous issue consigned European stock. It is understood that Australian stock was not replaced by the then owners of ESKA Implants AG and/or ESKA Australia. It should be stressed that this action was performed long before the change of ownership, and is in no way condoned by Orthodynamics GmbH. A review of stock currently held by the Australian distributor, ESKA Australia, and their customers, confirmed that there were still previous issue femoral heads in stock at their facility. Arrangements have now been put in place to ensure that this previous issue stock is returned to Orthodynamics GmbH to prevent its usage in the Australian market.

Bionik resurfacing head conclusion

The NJRR data demonstrated that even allowing for the small numbers, there is a higher than expected revision rate for the Bionik revision femoral head. However, given the low numbers and possible reasons for a slightly higher rate such as Cera-Metal and poor results at one centre, the revision rate of 5.1% is not excessive for resurfacing devices, as demonstrated by the UK National Joint Registry data.

The NJRR data indicates that there may be a problem with femoral heads, with them being revised in all 9 cases, and with Biosurf, as that product was involved in six of the eight metal revisions. The fact that loosening and lysis is higher than expected also suggests that metal wear particles may be a contributory factor.

The clinical literature does not match the conclusions of the NJRR with predominately positive clinical results published. The sales and complaint data from the rest of the world also indicate more positive results than the high revision numbers from the NJRR.

The mechanical testing results indicate a possible problem with potentially high and unreliable wear rates for earlier designs of the Biosurf femoral head. It is not clear whether this will affect all previous Bionik femoral head designs, but there is certainly evidence that heads manufactured in 2007 and before September 2008 could have a higher than normal wear rate. It is believed that this could be a contributory factor to the higher than expected revision levels of the Bionik resurfacing head and the cementless Adapter stem when used with a metal-on-metal bearing.

However, the problem was rectified by a design change in September 2008, and the NJRR is reporting lower revision rates, certainly for the Adapter stem and conventional Bionik femoral head, since the design change. It can be confirmed that our distributor has advised there is no stock available for implantation in Australia to the previous issue design, so Orthodynamics GmbH would expect to witness an increase in survivorship of all of its metal-on-metal bearing combinations over the next few years.

Even assuming the potential high wear rate, the clinical data gathered by Professor Dr. Gerdesmeyer et al and Dr. Rudert demonstrates that good clinical results are obtainable with this device when it used correctly and for suitable indications.

Given that the higher than expected revision rate can be attributed to a former design of femoral head, and the successful clinical and mechanical test results outside of the NJRR data, it is believed that the Bionik resurfacing head is clinically safe and effective for its intended use, and this information is sufficient to enable the continued usage of the device in the Australian market.

Report completed by:

Quality Management and Regulatory Affairs Manager - Implants

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Ref: MDA/2010/033 Issued: 22 April 2010 at 14:00

Device

All metal-on-metal (MoM) hip replacements.

Problem

The MHRA has received reports of revisions of MoM hip replacements involving soft tissue reactions. These reactions may be associated with unexplained hip pain.

Action

Put systems in place for the follow-up of patients implanted with MoM hip replacements including, where appropriate, blood metal ion measurements and cross sectional imaging.

Action by

- Medical directors
- Orthopaedic departments
- Orthopaedic surgeons
- Staff involved in the management of patients with joint replacement implants.

CAS deadlines

Action underway: 21 May 2010
Action complete: 17 June 2010

Note: These deadlines are for systems to be in place to take actions and not for the completion of patient follow-up and testing.

Problem

The majority of patients implanted with MoM hip replacements have well functioning hips and are thought to be at a low risk of developing serious problems.

A small number of patients implanted with these hips may, however, develop progressive soft tissue reactions to the wear debris associated with MoM articulations. The debris can cause soft tissue necrosis and adversely affect the results of revision surgery. Early revision of poorly performing MoM hip replacements should give a better revision outcome.

Following extensive consultation with orthopaedic experts and using information from the National Joint Registry for England and Wales, the MHRA is issuing this interim advice to healthcare professionals involved in the management of patients implanted with MoM hip replacements.

The MHRA is continuing to monitor the situation in consultation with orthopaedic experts and may issue further advice.

Action

For patients implanted with MoM hip replacements:

- follow up patients at least annually for five years postoperatively and more frequently in the presence of symptoms. Beyond five years, follow up in accordance with locally agreed protocols
- investigate patients with painful MoM hip replacements. Specific tests should include evaluation of cobalt and chromium ion levels in the patient's blood and cross sectional imaging including MRI or ultrasound scan
- consider measuring cobalt and chromium ion levels in the blood and/or cross sectional imaging for the following patient groups:
 - > patients with radiological features associated with adverse outcomes including component position
 - > patients with small component size (hip resurfacing arthroplasty only)
 - > cases where the patient or surgeon is concerned about the MoM hip replacement
 - > cohorts of patients where there is concern about higher than expected rates of failure
- if either cobalt or chromium ion levels are elevated above seven parts per billion (ppb), then a second test should be performed three months after the first in order to identify patients who require closer surveillance, which may include cross sectional imaging
- if imaging reveals soft tissue reactions, fluid collections or tissue masses then consider revision surgery.

Note: Measurements of cobalt or chromium ions should be carried out by laboratories participating in the Trace Elements External Quality Assessment Scheme (TEQAS) (<http://www.sas-centre.org/home.html>).

Distribution

This MDA has been distributed to:

- NHS trusts in England (Chief Executives)
- HSC trusts in Northern Ireland (Chief Executives)
- NHS boards in Scotland (Chief Executives)
- NHS boards and trusts in Wales (Chief Executives)
- Care Quality Commission (Headquarters)
- Primary care trusts in England (Chief Executives)

Issued: 22 April 2010 at 14:00

105
Ref: MDA/2010/033

Onward distribution

Please bring this notice to the attention of all who need to know or be aware of it. This may include distribution by:

Trusts to:

CAS and SABs (NI) liaison officers for onward distribution to all relevant staff including:

- Clinical governance leads
- Medical directors
- Nursing executive directors
- Orthopaedic departments
- Orthopaedic surgeons
- Outpatient clinics
- Outpatient theatre managers
- Outpatient theatre nurses
- Pathologists
- Radiology departments
- Radiology directors
- Risk managers
- Theatre managers

Care Quality Commission (CQC) (England only) to:

The MHRA considers this information to be important to:

- Hospitals in the independent sector
- Independent treatment centres
- Private medical practitioners

Primary care trusts to:

- Directors of public health
- General practitioners
- NHS walk-in centres

England

If you are in England, please send enquiries about this notice to the MHRA, quoting reference number **MDA/2010/033 or 2010/004/019/291/007**

Technical aspects

Ms Feza Haque or Dr Crina Cacou or Dr Khalid Razak
Medicines & Healthcare products Regulatory Agency
Market Towers
1 Nine Elms Lane
London SW8 5NQ

Tel: 020 7084 /3066/ 3338/ 3200

Fax: 020 7084 3106

Email: feza.haque@mhra.gsi.gov.uk

crina.cacou@mhra.gsi.gov.uk

khalid.razak@mhra.gsi.gov.uk

Clinical aspects

Dr Susanne Ludgate
Medicines & Healthcare products Regulatory Agency
Market Towers
1 Nine Elms Lane
London SW8 5NQ

Tel: 020 7084 3123

Fax: 020 7084 3111

Email: susanne.ludgate@mhra.gsi.gov.uk

How to report adverse incidents

Please report via our website <http://www.mhra.gov.uk>

Further information about CAS can be found at <https://www.cas.dh.gov.uk/Home.aspx>

Northern Ireland

Alerts in Northern Ireland will continue to be distributed via the NI SABS system.

Enquiries and adverse incident reports in Northern Ireland should be addressed to:

Northern Ireland Adverse Incident Centre

Health Estates Investment Group

Room 17

Annex 6

Castle Buildings

Stormont Estate

Dundonald BT4 3SQ

Tel: 02890 523 704

Fax: 02890 523 900

Email: NIAIC@dhsspsni.gov.uk

<http://www.dhsspsni.gov.uk/index/hea/niaic.htm>

How to report adverse incidents in Northern Ireland

Please report directly to NIAIC, further information can be found on our website <http://www.dhsspsni.gov.uk/niaic>

Further information about **SABS** can be found at <http://sabs.dhsspsni.gov.uk/>

Scotland

Enquiries and adverse incident reports in Scotland should be addressed to:

Incident Reporting and Investigation Centre

Health Facilities Scotland

NHS National Services Scotland

Gyle Square

1 South Gyle Crescent

Edinburgh EH12 9EB

Tel: 0131 275 7575

Fax: 0131 314 0722

Email: nss.irc@nhs.net

<http://www.hfs.scot.nhs.uk/online-services/incident-reporting-and-investigation-centre-irc/>

Wales

Enquiries in Wales should be addressed to:

Dr Sara Hayes

Senior Medical Officer

Medical Device Alerts

Welsh Assembly Government

Cathays Park

Cardiff CF10 3NQ

Tel: 029 2082 3922

Email: Haz-Aic@wales.gsi.gov.uk

MHRA is an executive agency of the Department of Health

© Crown Copyright 2010

Addressees may take copies for distribution within their own organisations

FDA U.S. Food and Drug Administration


[Home](#) > [Medical Devices](#) > [Products and Medical Procedures](#) > [Implants and Prosthetics](#)

Medical Devices**Metal-on-Metal Hip Implants**


Hip joint deterioration can lead to symptoms such as pain, stiffness or difficulty walking. When symptoms do not respond to conservative treatment, patients may be advised to undergo total hip replacement or hip resurfacing. Patients may receive a "metal-on-metal" hip implant in which the "ball and socket" of the device are both made from metal.

Metal-on-Metal Hip Implant Systems

Total Hip Replacement

Hip Resurfacing




*ADAM

In this website, the FDA describes hip implants, how metal-on-metal implants differ from other hip implants and gives information and recommendations for patients and physicians about the benefits and risks of these products.

Links on this page:

Table 3.2 Revision at three and five years for primary hip replacement procedures, undertaken between 1st April 2003 and 31st December 2009, which were linked to a HES/PEDW episode.

Category	Number of patients	Revision rate at three years ²³ (95% CI)	Revision rate at five years ²³ (95% CI)	Hazard ratio ²⁴ (95% CI)
Prosthesis type				
Males <55 years				
Total replacement using cement	1,669	2.7% (1.9% to 3.8%)	4.4% (3.2% to 6.0%)	1
Total replacement not using cement	3,858	2.9% (2.3% to 3.6%)	4.0% (3.1% to 5.1%)	1.10 (0.77 to 1.57)
Hybrid total replacement	1,498	1.8% (1.2% to 2.7%)	2.6% (1.7% to 3.8%)	0.68 (0.42 to 1.09)
Hip resurfacing	4,215	3.6% (3.0% to 4.3%)	5.6% (4.5% to 6.9%)	1.37 (0.98 to 1.92)
LHMoM THR	1,564	4.5% (3.3% to 6.1%)	6.4% (4.5% to 8.9%)	1.52 (1.01 to 2.28)
Males 55 - 64 years				
Total replacement using cement	5,560	1.7% (1.3% to 2.1%)	2.6% (2.1% to 3.2%)	1
Total replacement not using cement	8,035	2.8% (2.3% to 3.2%)	3.5% (2.9% to 4.2%)	1.56 (1.22 to 2.01)
Hybrid total replacement	2,990	1.8% (1.3% to 2.4%)	2.7% (2.0% to 3.7%)	1.09 (0.78 to 1.52)
Hip resurfacing	3,852	3.3% (2.7% to 3.9%)	4.5% (3.6% to 5.6%)	2.14 (1.63 to 2.80)
LHMoM THR	1,857	4.0% (2.9% to 5.4%)	5.5% (3.9% to 7.8%)	2.19 (1.56 to 3.07)
Males 65+ years				
Total replacement using cement	27,694	1.5% (1.3% to 1.7%)	2.0% (1.8% to 2.3%)	1
Total replacement not using cement	14,499	2.5% (2.2% to 2.8%)	3.1% (2.7% to 3.6%)	1.73 (1.48 to 2.02)
Hybrid total replacement	7,748	2.2% (1.9% to 2.6%)	3.0% (2.4% to 3.7%)	1.48 (1.22 to 1.80)
Hip resurfacing	980	5.4% (4.1% to 7.3%)	6.0% (4.4% to 8.2%)	3.85 (2.83 to 5.23)
LHMoM THR	1,444	3.4% (2.2% to 5.1%)	6.4% (3.7% to 11.0%)	2.21 (1.55 to 3.16)
Females <55				
Total replacement using cement	2,174	2.5% (1.8% to 3.4%)	3.6% (2.7% to 4.8%)	1
Total replacement not using cement	5,082	2.7% (2.2% to 3.3%)	4.3% (3.4% to 5.5%)	1.22 (0.89 to 1.71)
Hybrid total replacement	1,840	2.2% (1.5% to 3.1%)	4.0% (2.8% to 5.8%)	0.97 (0.64 to 1.48)
Hip resurfacing	2,647	5.4% (4.5% to 6.5%)	8.3% (6.8% to 10%)	2.32 (1.67 to 3.21)
LHMoM THR	1,037	6.3% (4.6% to 8.6%)	9.2% (5.9% to 14.1%)	2.22 (1.48 to 3.32)
Females 55-64				
Total replacement using cement	8,239	1.6% (1.3% to 2.0%)	2.5% (2.1% to 3.0%)	1
Total replacement not using cement	10,735	2.5% (2.2% to 2.9%)	3.7% (3.1% to 4.4%)	1.54 (1.24 to 1.91)
Hybrid total replacement	4,474	1.6% (1.2% to 2.1%)	2.9% (2.2% to 3.8%)	1.06 (0.79 to 1.43)
Hip resurfacing	1,854	5.4% (4.4% to 6.6%)	7.9% (6.3% to 9.8%)	3.61 (2.78 to 4.68)
LHMoM THR	1,526	6.7% (5.2% to 8.7%)	9.5% (7.4% to 12.3%)	3.73 (2.80 to 4.96)
Females 65+				
Total replacement using cement	54,023	1.1% (1.0% to 1.2%)	1.6% (1.5% to 1.8%)	1
Total replacement not using cement	20,728	2.3% (2.0% to 2.5%)	3.0% (2.7% to 3.4%)	2.21 (1.94 to 2.53)
Hybrid total replacement	13,112	1.5% (1.3% to 1.8%)	2.0% (1.6% to 2.4%)	1.43 (1.20 to 2.53)
Hip resurfacing	305	7.1% (4.5% to 11.1%)	8.8% (5.2% to 14.7%)	6.59 (4.17 to 10.41)
LHMoM THR	1,454	3.6% (2.5% to 5.2%)	10.5% (6.1% to 17.8%)	3.65 (2.67 to 5.00)
Patient physical status				
P1 – fit and healthy	37,244	2.3% (2.2% to 2.5%)	3.4% (3.1% to 3.7%)	1
P2 – mild disease not incapacitating	138,701	1.9% (1.9% to 2.0%)	2.8% (2.6% to 2.9%)	1.06 (0.98 to 1.16)
P3+ – incapacitating systemic disease or worse	40,748	2.1% (2.0% to 2.3%)	2.8% (2.6% to 3.1%)	1.29 (1.16 to 1.43)
Type of provider				
NHS hospital	180,937	2.1% (2.0% to 2.2%)	2.9% (2.8% to 3.0%)	1
NHS treatment centre	13,554	1.8% (1.5% to 2.1%)	2.9% (2.4% to 3.4%)	0.83 (0.73 to 0.96)
Independent hospital	13,378	2.0% (1.7% to 2.4%)	2.6% (2.2% to 3.1%)	0.85 (0.73 to 1.00)
ISTC	8,824	2.3% (1.8% to 2.8%)	2.9% (2.2% to 3.9%)	1.09 (0.91 to 1.31)

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ARTNR	ETIBEZ1	ETIBEZ2	GROESSE	SEITE
10252438	METALLKOPF SILVER "BS"	K12/14	AD=38MM	X-LONG
10252444	METALLKOPF SILVER "BS"	K12/14	AD=44MM	X-LONG
10252448	METALLKOPF SILVER "BS"	K12/14	AD=48MM	X-LONG
10250438	METALLKOPF BIOSURF SILVER "BS"	K12/14	AD=38MM	X-LONG
10250444	METALLKOPF BIOSURF SILVER "BS"	K12/14	AD=44MM	X-LONG
10250448	METALLKOPF BIOSURF SILVER "BS"	K12/14	AD=48MM	X-LONG

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


16107

037

Bionik Total Conventional Hip Acetabular Prosthesis

This analysis compares the Bionik Total Conventional Hip Acetabular Prosthesis with all Other Total Conventional Hip prostheses. This Acetabular Prosthesis has been identified as having a significantly higher revision rate.

For a detailed explanation of the process used by the Registry that results in identification of prostheses that have a higher than anticipated rate of revision please refer to the 'Prostheses with Higher than Anticipated Rates of Revision' chapter of the most recent AOANJRR Annual Report, <http://www.dmac.adelaide.edu.au/aoanjrr/publications.jsp>.

TABLE 1

Revision Rate of Primary Total Conventional Hip Replacement

The **Revision Rate** of the Bionik Total Conventional Hip Acetabular Prosthesis is compared to all Other Total Conventional Hip prostheses.

Table 1: Revision Rates of Primary Total Conventional Hip Replacement

Component	N Revised	N Total	Obs. Years	Revisions/100 Obs. Yrs (95% CI)
Bionik	28	565	1103	2.54 (1.69, 3.67)
Other Total Conventional Hip	5049	170539	651215	0.78 (0.75, 0.80)
TOTAL	5077	171104	652317	0.78 (0.76, 0.80)

TABLE 2

Yearly Cumulative Percent Revision of Primary Total Conventional Hip Replacement

The **Yearly Cumulative Percent Revision** of the Bionik Total Conventional Hip Acetabular Prosthesis is compared to all Other Total Conventional Hip prostheses.

Table 2: Yearly Cumulative Percent Revision of Primary Total Conventional Hip Replacement

CPR	1 Yr	3 Yrs	5 Yrs	7 Yrs	9 Yrs
Bionik	3.2 (2.0, 5.2)	6.9 (4.6, 10.2)			
Other Total Conventional Hip	1.5 (1.5, 1.6)	2.7 (2.6, 2.8)	3.5 (3.4, 3.6)	4.4 (4.3, 4.5)	5.4 (5.1, 5.7)


FIGURE 1

Yearly Cumulative Percent Revision of Primary Total Conventional Hip Replacement

The **Yearly Cumulative Percent Revision** of the Bionik Total Conventional Hip Acetabular Prosthesis is compared to all Other Total Conventional Hip prostheses. In addition, Hazard Ratios are also reported.

Hazard Ratios are reported for specific time periods during which the Hazard Ratio is constant. This is done to enable more specific and valid comparisons of the risk of revision over time. The pattern of variation in risk has important implications with respect to the underlying reasons for any difference.

Figure 1: Cumulative Percent Revision of Primary Total Conventional Hip Replacement

Number at Risk	0 Yr	1 Yrs	2 Yrs	3 Yrs	4 Yrs	5 Yrs	6 Yrs	7 Yrs	8 Yrs	9 Yrs
Bionik	565	412	267	132	10	0	0	0	0	0
Other Total Conventional Hip	170539	142248	117869	96119	75962	57166	40037	24821	11582	2968

TABLE 3**Primary Diagnosis for Revised Primary Total Conventional Hip Replacement**

This table identifies the diagnosis of the primary procedure which was subsequently revised. This information is provided as there is a variation on outcome depending on the primary diagnosis. It is therefore important when considering the reasons for a higher than anticipated rate of revision that there is identification of the primary diagnosis. This information should be compared to the primary diagnosis for the revisions of all Other Total Conventional Hip prostheses.

Table 3: Primary Diagnosis for Revised Primary Total Conventional Hip Replacement

Primary Diagnosis	Bionik		Other Total Conventional Hip	
	Number	Percent	Number	Percent
Osteoarthritis	27	96.4	4257	84.3
Avascular Necrosis			253	5.0
Fractured Neck Of Femur	1	3.6	249	4.9
Rheumatoid Arthritis			96	1.9
Developmental Dysplasia			87	1.7
Failed Internal Fixation			32	0.6
Tumour			32	0.6
Other Inflammatory Arthritis			29	0.6
Fracture/Dislocation			7	0.1
Arthrodesis Takedown			5	0.1
Other			2	0.0
TOTAL	28	100.0	5049	100.0

TABLE 4**Revision Rates of Bionik Primary Total Conventional Hip Replacement by Fixation.**

This analysis is provided as some prostheses have more than one fixation option. Additionally there are prostheses where an alternative to the recommended approach to fixation was used e.g. a cementless prosthesis that has been cemented or vice-versa.

Table 4: Revision Rates of Bionik Primary Total Conventional Hip Replacement by Fixation

Fixation	N Revised	N Total	Obs. Years	Revisions/100 Obs. Yrs (95% CI)
Cemented	0	2	5	0.00 (0.00, 76.12)
Cementless	22	481	942	2.34 (1.46, 3.54)
Hybrid	6	82	156	3.85 (1.41, 8.37)
TOTAL	28	565	1103	2.54 (1.69, 3.67)

TABLE 5

Revision Rates of Bionik Primary Total Conventional Hip Replacement by Bearing Surface.

This analysis is provided as some prostheses are combined with a variety of different bearing surfaces. All bearing surfaces used with this Acetabular Prosthesis are listed.

Table 5: Revision Rates of Bionik Primary Total Conventional Hip Replacement by Bearing Surface

Bearing Surface	N Revised	N Total	Obs. Years	Revisions/100 Obs. Yrs (95% CI)
Ceramic/Ceramic	1	8	4	25.87 (0.65, 144.1)
Ceramic/Polyethylene	1	43	75	1.33 (0.03, 7.38)
Metal/Metal	22	418	885	2.49 (1.56, 3.77)
Metal/Polyethylene	4	96	139	2.88 (0.78, 7.37)
TOTAL	28	565	1103	2.54 (1.69, 3.67)

TABLE 6

Type of Revision Performed for Primary Total Conventional Hip Replacement

This analysis identifies the components used in the revision of the Bionik Total Conventional Hip Acetabular Prosthesis and compares it to the components used in the revision of all Other Total Conventional Hip prostheses.

The reason this analysis is undertaken is to identify whether there is one or more components which are being replaced that differ from the components replaced for revisions of all Other Total Conventional Hip prostheses i.e. is there a difference in the type of revision undertaken for the Bionik Total Conventional Hip Acetabular Prosthesis compared to all Other Total Conventional Hip prostheses.

Table 6: Primary Total Conventional Hip Replacement - Type of Revision

Revision Type	Bionik		Other Total Conventional Hip	
	Number	Percent	Number	Percent
Femoral Only	7	25.0	1481	29.3
Acetabular Only	5	17.9	1260	25.0
THR (Femoral/Acetabular)	4	14.3	623	12.3
Cement Spacer	1	3.6	244	4.8
Removal of Prostheses			41	0.8
Reinsertion of Components			6	0.1
Bipolar Head and Femoral Saddle			2	0.0
			2	0.0
N Major	17	60.7	3659	72.5
Head/Insert	7	25.0	920	18.2
Head Only	1	3.6	285	5.6
Minor Components			82	1.6
Insert Only	1	3.6	70	1.4
Head/Neck	2	7.1	30	0.6
Neck Only			2	0.0
Neck/Insert			1	0.0
N Minor	11	39.3	1390	27.5
TOTAL	28	100.0	5049	100.0

TABLE 7**Reason for Revision of Primary Total Conventional Hip Replacement**

This is reported in two ways; a percentage of all revisions and also as a percentage of all primary procedures.

This analysis includes a comparison of reasons for revision to all Other Total Conventional Hip prostheses.

This analysis is undertaken to identify if there are differences in the reasons for revision and the number of revisions performed for those reasons between the Bionik Total Conventional Hip Acetabular Prosthesis and all Other Total Conventional Hip prostheses.

Table 7: Primary Total Conventional Hip Replacement - Reason for Revision

Revision Diagnosis	Number	Bionik		Other Total Conventional Hip		
		% Revision	% Primary	Number	% Revision	% Primary
Loosening/Lysis	7	25.0	1.2	1512	29.9	0.9
Prosthesis Dislocation	7	25.0	1.2	1393	27.6	0.8
Infection	3	10.7	0.5	843	16.7	0.5
Fracture	2	7.1	0.4	742	14.7	0.4
Pain	5	17.9	0.9	99	2.0	0.1
Other				67	1.3	0.0
Leg Length Discrepancy	1	3.6	0.2	63	1.2	0.0
Metal Sensitivity	1	3.6	0.2	61	1.2	0.0
Malposition				54	1.1	0.0
Wear Acetabulum	1	3.6	0.2	51	1.0	0.0
Implant Breakage Stem				35	0.7	0.0
Implant Breakage Acetabular	1	3.6	0.2	32	0.6	0.0
Incorrect Sizing				29	0.6	0.0
Instability				28	0.6	0.0
Implant Breakage Head				20	0.4	0.0
Heterotopic Bone				8	0.2	0.0
Tumour				8	0.2	0.0
Avascular Necrosis				2	0.0	0.0
Synovitis				2	0.0	0.0
N Revision	28	100.0	5.0	5049	100.0	3.0
N Primary	565			170539		

FIGURE 2

Revision Diagnosis Cumulative Incidence by Time to Revision for Primary Total Conventional Hip Replacement

This figure details the cumulative incidence of the most common reasons for revision.

The five most common reasons for revision are included as long as each of these reasons account for more than 10 procedures or at least 5% of all revisions for the Bionik Total Conventional Hip Acetabular Prosthesis. A comparative graph is provided of the cumulative incidence for the same reasons for revisions for all Other Total Conventional Hip prostheses.

Figure 2: Cumulative Incidence Revision Diagnosis for Primary Total Conventional Hip Replacement


TABLE 8

Revision Rates of Primary Total Conventional Hip Replacement by State

This enables a state by state variation to be identified for the Bionik Total Conventional Hip Acetabular Prosthesis and provides the comparative data for each of the states for all Other Total Conventional Hip prostheses.

This analysis is undertaken for similar reasons as those outlined above for Table 8.

Table 8: Revision Rates of Primary Total Conventional Hip Replacement by State

Component	State	N Revised	N Total	Obs. Years	Revisions/100 Obs. Yrs (95% CI)
Bionik	NSW	18	306	610	2.95 (1.75, 4.67)
Bionik	VIC	3	75	145	2.07 (0.43, 6.04)
Bionik	QLD	1	20	55	1.81 (0.05, 10.08)
Bionik	WA	0	22	11	0.00 (0.00, 32.84)
Bionik	SA	0	5	5	0.00 (0.00, 69.02)
Bionik	TAS	6	136	276	2.18 (0.80, 4.74)
Bionik	ACT/NT	0	1	0	0.00 (0.00, 1773)
Other Total Conventional Hip	NSW	1401	50349	180275	0.78 (0.74, 0.82)
Other Total Conventional Hip	VIC	1378	46238	176180	0.78 (0.74, 0.82)
Other Total Conventional Hip	QLD	867	27406	105373	0.82 (0.77, 0.88)
Other Total Conventional Hip	WA	647	18995	76125	0.85 (0.79, 0.92)
Other Total Conventional Hip	SA	462	17571	74846	0.62 (0.56, 0.68)
Other Total Conventional Hip	TAS	162	6112	24253	0.67 (0.57, 0.78)
Other Total Conventional Hip	ACT/NT	132	3868	14161	0.93 (0.78, 1.11)
TOTAL		5077	171104	652317	0.78 (0.76, 0.80)

TABLE 9

Number of Revisions of Bionik Primary Total Conventional Hip Replacement by Year of Implant

This analysis details the number of prostheses reported each year to the Registry for the Bionik Total Conventional Hip Acetabular Prosthesis. It also provides the subsequent number of revisions of the primaries reported in that year.

Primary procedures performed in later years have had less follow up time therefore the number revised is expected to be less than the number revised in earlier years. For example, a primary procedure performed in 2009 has a maximum of one year to be revised, whereas a primary performed in 2007 has a maximum of three years to be revised.

Table 9: Number of Revisions of Bionik Primary Total Conventional Hip Replacement by Year of Implant

Year of Implant	Number Revised	Total Number
2005	0	11
2006	13	147
2007	6	136
2008	5	137
2009	4	134
TOTAL	28	565

TABLE 10

Revision rates of Bionik Primary Total Conventional Hip Replacement by Catalogue Number Range

Many prostheses have a number of catalogue ranges. The catalogue range is specific to particular design features; more than one catalogue range usually indicates a minor difference in design in a particular Bionik prosthesis.

This analysis has been undertaken to determine if the revision rate varies according to the catalogue number range.

Table 10: Revision Rates of Bionik Primary Total Conventional Hip Replacement by Catalogue Number Range

Catalogue Range		Catalogue Description		
Acetabular				
Bionik	10201050-10201064	METAL SHELL BS TINB COAT		
Bionik	10201150-10201164	METAL SHELL BS TINB COAT SCREW FIX		
Bionik	10201248-10201264	METAL SHELL CEMENTLESS TINB COAT SCREW FIX		
Bionik	10201346-10201366	METAL SHELL TINB CAP SCREW FIX		
Bionik	10202154-10202164	METAL SHELL CEMENTLESS DYSPL. TINB C.: SCREW FIX		
Acetabular Range	N Revised	N Total	Obs. Years	Revisions/100 Obs. Yrs (95% CI)
10201050-10201064	1	23	54	1.84 (0.05, 10.23)
10201150-10201164	6	126	271	2.22 (0.81, 4.83)
10201248-10201264	4	52	101	3.95 (1.08, 10.12)
10201346-10201366	17	361	672	2.53 (1.47, 4.05)
10202154-10202164	0	3	4	0.00 (0.00, 87.78)
TOTAL	28	565	1103	2.54 (1.69, 3.67)

TABLE 11***Revision rates of Bionik Primary Total Conventional Hip Replacement by Component***


A prosthesis may be combined with multiple components.

This analysis has been undertaken to determine if the revision rate varies according to the component with which it is combined.

Table 11: Revision Rates of Bionik Primary Total Conventional Hip Replacement by Femoral Stem Component

Femoral Stem Component	N Revised	N Total	Obs. Years	Revisions/100 Obs. Yrs (95% CI)
Adapter	27	560	1097	2.46 (1.62, 3.58)
Corail	0	1	0	0.00 (0.00, 806.8)
Custom Made (Eka)	0	1	1	0.00 (0.00, 246.8)
Edinburgh	0	1	3	0.00 (0.00, 128.1)
K2	1	1	0	1304.46 (33.03, 7268)
S-Rom	0	1	1	0.00 (0.00, 687.4)
TOTAL	28	565	1103	2.54 (1.69, 3.67)

Graph comparing revisions of Adapter stems and Bionik femoral heads

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

09106

086

Orthodynamics GmbH Sales and Complaints Data - Cementless Adapter Stem

Year	Product Name	# supplied in Australia	# supplied World Wide	# of complaints in Australia	# of complaints World Wide (incl. Australia)	# of Adverse Events in Australia	# of Adverse Events World Wide (incl. Australia)
2005	Adapter (Cementless) Femoral Stem	20	475	0	0	1	1
2006	Adapter (Cementless) Femoral Stem	143	338	0	0	9	9
2007	Adapter (Cementless) Femoral Stem	136	398	0	0	6	6
2008	Adapter (Cementless) Femoral Stem	129	351	0	0	4	4
2009	Adapter (Cementless) Femoral Stem	169	427	0	0	3	3

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Experience with a Hydroxyapatite-Coated and Macroporous-surfaced Hip Prosthesis

Introduction: Cementless hip arthroplasty is not completely satisfactory – even with macroporous structured surface. Medium term results are sometimes disappointing because of insufficient secondary stability. Thus our aim was to improve fixation by additional coating with hydroxyapatite.

Methods: 200 consecutive patients who had hydroxyapatite-coated cementless hip replacement with a macroporous hip prosthesis (09/92–03/96) were studied prospectively. All patients were included in a prospective follow up schedule according to the criteria of Johnston et al. As measure of clinical outcome we calculated Harris Hip Score as well as Enghs Score to assess radiologic fixation and stability.

Results: 91% could be followed up regularly. No revision because of aseptic loosening had to be done. Analysis of clinical results showed almost painless patients from early after the operation – especially no thigh pain. Average HHS after 2 years was 97.

Radiological evaluation showed early and complete osteointegration of all components. According to Enghs Score they are stable and well fixed.

Conclusion: HA coated macroporous implants provide some striking advantages, which encouraged us to continue with this system. By early and secure bony ingrowth a fibrous interface is avoided and thereby also long-lasting thigh pain.

Key words: Hip arthroplasty – surface structure – hydroxyapatite

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Osteoblastic response to osteoarthritis of the hip does not predict outcome of cementless cup fixation

79 patients followed for 5-11 years

Takashi Nishii, Nobuhiko Sugano, Takashi Sakai, Keiji Haraguchi, Kenji Ohzono and Hideki Yoshikawa

Department of Orthopaedic Surgery, Osaka University Medical School E3, 2-2 Yamadaoka, Suita, Osaka 565-0871, Japan
E-mail: nishii@ort.med.osaka-u.ac.jp
Submitted 01. Accepted 01


ABSTRACT - We evaluated the influence of osteoblastic response to osteoarthritis of the hip on the outcome of cementless acetabular cup after 91 total hip replacements in 79 patients. Of the 91 hips, 23 were atrophic, 37 normotrophic, and 31 hypertrophic, according to Bombelli's criteria. There were no clinical or radiographic differences among the three groups at the final follow-up (average 7 (5-11) years), when stable bone growth had been achieved by all of the acetabular cups in patients with the atrophic type, 35/37 of the normotrophic type, and all the hypertrophic type. Revision of the acetabular cup was performed on 1 hip of the normotrophic type, in connection with severe polyethylene liner wear and progressive osteolysis. ■

Osteoblastic response by host bone to implant surface may be an important factor affecting longevity of total hip arthroplasty (THA). In patients with osteoarthritis (OA), the type of osteoblastic response, in terms of Bombelli's (1983) classification, has been reported to predict loosening of the cemented acetabular socket (Saito et al. 1987, Kobayashi et al. 1997). The atrophic OA type hips showed a higher incidence of unsatisfactory clinical results and radiographic loosening of the cemented cup than did the hypertrophic or normotrophic OA type hips. The influence of osteoblastic response on cementless THA has not been studied. We thought that cementless fixation may need a higher osteoblastic response to achieve good bone ingrowth or ongrowth fixation. Therefore we inves-

tigated whether there was any association between osteoblastic response and outcome in cementless THA, as previously reported in cemented THA.

Patients and methods

Between 1987 and 1993, a cementless THA, using the Spongiosa Metal Total Hip System (ESKA, Lübeck, Germany), was performed in a consecutive series of 100 hips for 87 patients with OA. All the cases were secondary to acetabular dysplasia or congenital subluxation or dislocation of the hip. No cases of primary osteoarthritis, or completely dislocated hip were included. The stem and metal socket of this system are made of cobalt chrome molybdenum alloy with the entire surface covered by a spongy metal structure similar to cancellous bone. The metal cup has a hemispherical shape with two spikes wedged into the anterior and posterior acetabular rims and one peg inserted into the ischium. The curved stem with a collar anatomically matches the medullary canal of the proximal femur. An interior polyethylene liner (ESKA, Lübeck, Germany) and a modular aluminum ceramic head, having a diameter of 28 mm (BioloX: Ceramtech, Plochingen, Germany), are used for the articulation. As a rule, the cup was placed in the anatomical position of the acetabulum, and a block bone graft from the resected femoral head was used in 14 hips and chip grafts in 23 hips to fill the superolateral bony defect in the acetabulum. Partial weight bearing was started at

Morphological types of osteoarthritis of the hip according to Bombelli's criteria: a) atrophic OA type, showing poor osteophyte formation and gradual reduction in size of the femoral head, b) normotrophic OA type, showing a moderately sized osteophyte formation, c) hypertrophic OA type, having large osteophytes on the margins and floor of the acetabulum as well as around the femoral head.

6 weeks and full weight bearing was allowed at 12 weeks after the operation. Clinical and radiographic assessments were done annually.

2 patients, accounting for 3 hips, died within 5 years after surgery due to causes unrelated to

the hip. Deep infection occurred in 1 hip 4 years postoperatively, and both the acetabular cup and stem were removed. 5 patients (5 hips) were lost to follow-up within 5 years after surgery, but they showed excellent clinical and radiographic results when last reviewed. The remaining 91 hips of 79 patients were the subjects for subsequently analyzed with a minimal follow-up of 5 years. The average age of the patients was 52 (29-67) years at the time of operation, and the average follow-up 7 (5-11) years.

The clinical assessments were based on the Merle d'Aubigné and Postel (1954) hip score, which allocates up to 6 points for each type of pain, mobility, and function, with a total of 18 points given to a normal hip. Osteoblastic response of the arthrotic hips was evaluated with Bombelli's (1983) criteria on preoperative radiographs (Figure). These radiographs were interpreted blindly by two of the authors without knowledge of the clinical and radiographic outcomes. Both observers had the same view about 78 hips, but came to an agreement about the remaining cases. 23 hips were classified as atrophic OA type (A-type), 37 as normotrophic OA type (N-type), and 31 as hypertrophic OA type (H-type). All patients were female, except 1 in the A-type group and 2 in the N-type group. The follow-up period was about the same in the three groups (mean months (SD): 85 (20) in the A-type group, 85 (19) in the N-type group, 87 (15) in the H-type group). The patients' ages at operation differed significantly: 56 years (SD 6) in the A-type group, 53 years (7) in the N-type group, 49 years (8) in the H-type group, the latter being younger than the A-type group ($p = 0.003$, Scheffe test).

At the annual postoperative follow-ups, anteroposterior and lateral radiographs were evaluated for evidence of radiolucent lines, implant fixation, and osteolysis by a zonal interface analysis of the acetabular cup with DeLee and Charnley's method (1976). Fixation of the cup was evaluated by a modified DeLee and Charnley classification (McPherson et al. 1995), according to which fixation by bone ingrowth (type 1) has three subtypes: 1A means the absence of radiolucent lines or migration of the cup, 1B has a radiolucent line in one zone, and 1C has a radiolucent line in two zones. Stable fibrous fixation (type 2) has a com-

Incidence of radiolucent lines around the acetabular cup and cup fixation at the final follow-up

Type of osteoarthrosis	Radiolucent line				Cup fixation				
	Present in one or more zones	Zone 1	Zone 2	Zone 3	1A	1B	1C	2	3
Atrophic	2/2 (CI 0-0.22)	2/2	0	0	20/21	2/2	0	0	0
Normotrophic	7/7 (CI 0.06-0.35)	1/1	4/4	5/5	27/29	4/4	1/2	0	2/2 (CI 0-0.14)
Hypertrophic	3/3 (CI 0.02-0.24)	0	0	3/3	24/28	3/3	0	0	0

Values are number of patients/hips, CI = 95 % confidence interval of incidence ratios

plete radiolucent line less than 2 mm in width. Unstable fibrous fixation (type 3) has a complete radiolucent line equal to or greater than 2 mm in width, a progressive zone 3 radiolucent line, or cup migration by more than 2 mm.

Statistical analysis was performed with StatView software (version 5.0; SAS Institute, Cary, North Carolina) and a p-value less than 0.05 was considered significant.

Results

Clinical results

The mean Merle d'Aubigné and Postel hip score of all patients improved from 7.6 (pain 1.2, mobility 3.7, function 2.7) before operation, to 17 (pain 5.9, mobility 5.3, function 5.6) 2 years post-operatively, and to 17 (pain 5.9, mobility 5.3, function 5.5) at the final follow-up. All three groups achieved good hip score improvement from preoperative examinations to the final follow-up: from 8.0 points to 17 in the A-type group, from 7.7 to 17 in the N-type group, and from 7.2 to 17 in the H-type group. There was no statistically significant difference in the hip scores among the three groups at the final follow-up (one-way analysis of variance). The results were unsatisfactory, indicated by a total score of 14 or less, in 2 N-type hips due to stem fracture or disability from osteoarthritis of the contralateral hip, and in 1 H-type hip due to disability from lumbar spondyloarthritis.

Radiographic results

Radiolucent lines in one or more zones around the acetabular cup were seen in 2 patients in the A-type group, 7 patients in the N-type group, and 3 patients in the H-type group at the final follow-up (Table). When adjusted for other factors, including age, sex, follow-up period, and performance of block bone graft, osteoblastic response type showed no significant correlation with the occurrence of a radiolucent line (logistic regression analysis). Satisfactory fixation of the acetabular cup was achieved in all 3 groups at the final follow-up. All patients in A-type and H-type groups, and 32 patients in the N-type group showed bone ingrowth fixation (grades 1A, 1B, or 1C). Osteolysis around the acetabular cup was most marked in the N-type group (6 patients), while only 1 patient in the A-type group and none in the H-type group showed osteolysis.

Two revision operations were performed during the study: 1 hip in the N-type group underwent a revision of the acetabular cup due to unstable fibrous fixation, and 1 hip in the N-type group underwent revision of the femoral stem due to stem fracture without loosening of the cup.

Discussion

The morphological type of osteoarthritis is said to affect the outcome of cemented THA (Saito et al. 1987, Kobayashi et al. 1997): atrophic type hips had a high incidence of unsatisfactory clinical results and acetabular loosening. Histological

examinations showed that bone specimens from hypertrophic type hips had many osteoblasts with few osteoclasts on the surface of subchondral trabeculae, while bone specimens from atrophic type hips had many osteoclasts with fewer osteoblasts (Saito et al. 1987).

These findings suggest that the results with morphological types of osteoarthritis may be similar to those with acetabular cups after cementless THA. Hernández-Vaquero et al. (1996) compared the outcome of threaded acetabular cups after cementless THA among the atrophic, normotrophic, and hypertrophic types of hips at a mean follow-up of 6 years, and found that periprosthetic radiolucent lines and loosening of the acetabular cups occurred significantly oftener in the atrophic hips. Our study, however, showed no significant differences in clinical results, incidence of radiolucent lines or cup fixation after 7 years among the three groups. The discrepancy between the findings of Hernández-Vaquero et al. (1996) and ours may be due to the tendency of threaded cups to fail after a few years (Snorrason and Kärrholm 1990, Yahiro et al. 1995). Adequate bone ingrowth and long-term fixation of the acetabular cup can be achieved only by stable initial fixation into the acetabular bone. If stability is insufficient, minimal bone ingrowth can be expected because of weak osteoblastic activity in atrophic type hips. Here we used a press-fitting technique with underreaming to obtain stable fixation of the acetabular cup. The results half-way through the study using this technique for cup fixation were satisfactory (Morscher 1992, Dorr et al. 1998). Linder et al. (1988) showed that osseointegration with the implant surface was usually attained even in patients with compromised bone quality, such as osteoporosis and rheumatoid arthritis, and adequate initial stability with the acetabular cup we used should have resulted in adequate subsequent bone ingrowth fixation and a successful outcome, even in hips of the atrophic type.

Our study has certain limitations: almost all patients were female, and osteoarthritis secondary to dysplasia or congenital subluxation or dislocation of the hip was evaluated exclusively because these are the main causes of osteoarthritis in our country. This may have affected the relation between the outcome of the arthroplasty and the osteoblastic types of hip. We found no significant

difference among the three groups as regards clinical and radiographic outcomes. However, there may be a true difference among the groups (type II error), because of the limited sample size (Freedman and Bernstein 1999). Further studies in more patients including more males and various causes of osteoarthritis are needed to increase the reliability of the statistical analysis and extrapolate the findings to the usual type of osteoarthritis of the hip or a population with more males.

No benefits in any form have been received or will be received from a commercial party related directly or indirectly to the subjects of this article. No funds were received to support this study.

Bombelli R. Osteoarthritis of the hip: classification and pathogenesis: the role of osteotomy as a consequent therapy. 2nd ed. Springer-Verlag, Berlin Heidelberg New York 1983.

DeLee J G, Charnley J. Radiological demarcation of cemented sockets in total hip replacement. *Clin Orthop* 1976; 121: 20-32.

Dorr L D, Wan Z, Cohen J. Hemispheric titanium porous-coated acetabular component without screw fixation. *Clin Orthop* 1998; 351: 158-68.

Freedman K B, Bernstein J. Sample size and statistical power in clinical orthopaedic research. *J Bone Joint Surg (Am)* 1999; 81 (10): 1454-60.

Hernández-Vaquero D, Suárez-Vázquez A, Fernández-Corona C, Menéndez-Viñuela G, Alegre-Mateo R, García-Sandoval M A. Loosening of threaded acetabular cups in arthroplasty of the hip: the association with different types of coxarthrosis. *International Orthopaedics* 1996; 20: 70-4.

Kobayashi S, Takaoka K, Saito N, Hira K. Factors affecting aseptic failure of fixation after primary Charnley total hip arthroplasty. *J Bone Joint Surg (Am)* 1997; 79 (11): 1618-27.

Linder L, Carlsson A, Marsal L, Bjursten L M, Brånemark P-I. Clinical aspects of osseointegration in joint replacement: a histological study of titanium implants. *J Bone Joint Surg (Br)* 1988; 70 (4): 550-5.

McPherson E J, Dorr L D, Gruen T A, Saberi M T. Hydroxyapatite-coated proximal ingrowth femoral stems: a matched pair control study. *Clin Orthop* 1995; 315: 223-30.

Merle d'Aubigné R, Postel M. Functional results of hip arthroplasty with acrylic prosthesis. *J Bone Joint Surg (Am)* 1954; 36 (3): 451-75.

Morscher E W. Current status of acetabular fixation in primary total hip arthroplasty. *Clin Orthop* 1992; 274: 172-93.

Saito M, Saito S, Ohzono K, Ono K. The osteoblastic response to osteoarthritis of the hip: its influence on the long-term results of arthroplasty. *J Bone Joint Surg (Br)* 1987; 69 (5): 746-51.

Snorrason F, Kärrholm J. Primary migration of fully-threaded acetabular prostheses: a roentgen stereophotogrammetric analysis. *J Bone Joint Surg (Br)* 1990; 72 (4): 647-52.

Yahiro M A, Gantenberg J B, Nelson R, Lu H T C, Mishra N K. Comparison of the results of cemented, porous-ingrowth, and threaded acetabular cup fixation: a meta-analysis of the orthopaedic literature. *J Arthroplasty* 1995; 10 (3): 339-50.

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Literaturverzeichnis

- Kusaba A, Kondo S, Kuroki Y: A five year follow-up of Spongiosa Metal II ceramic-ceramic total hip replacement system for dysplastic osteoarthritis. *Orthopaedic Ceramic Implant* 25: 53-6, 2005
- Kalamchi A, MacEwen GD: Avascular necrosis following treatment of congenital dislocation of the hip. *J Bone Joint Surg Am* 62:876-88, 1980
- Crowe JF, Mani VJ, Ranawat CS: Total hip replacement in congenital dislocation and dysplasia of the hip. *J Bone Joint Surg Am* 61: 15-23, 1979
- Walter WL, O'toole GC, Walter WK, et al.: Squeaking in ceramic-on-ceramic hips: the importance of acetabular component orientation. *J Arthroplasty* 22:496-503, 2007
- Callaghan JJ, Dysart SH, Savory CG: The uncemented porous-coated anatomic total hip prosthesis: two-year results of a prospective consecutive series. *J Bone and Joint Surg Am* 70: 337-46, 1988
- Lusty PJ, Watson A, Tuke MA, et al.: Orientation and wear of the acetabular component in third generation alumina-on-alumina ceramic bearings. An analysis of 33 retrievals. *J Bone Joint Surg Br* 89: 1158-64, 2007
- Santavirta S, Böhler M, Harris WH et al.: Alternative materials to improve total hip replacement tribology. *Acta Orthop Scand*: 74: 380-8, 2003
- Nam KW, Yoo JJ, Lae Kim et al.: Alumina-debris-induced osteolysis in contemporary alumina-alumina total hip arthroplasty. *J Bone Joint Surg Am* 89: 2499-503, 2007
- Bader R, Steinhauser E, Zimmermann S, et al.: Differences between the wear couples metal-on-polyethylene and ceramic-on-ceramic in the stability against dislocation of total hip replacement. *J Mater Sci Mater Med* 15: 711-8, 2004
- Garcia-Cimbrelo E, Garcia-Rey E, Murcia-Mazón, A, et al.: Alumina-on-alumina in THA: a multicenter prospective study. *Clin Orthop Relat Res* 466:309, 2008
- Bizot P, Banallec L, Sedel L, et al.: Alumina-on-alumina total hip prostheses in patients 40 years of age or younger. *Clin Orthop Relat Res* 379: 68-76, 2000
- Yoo JJ, Kim HJ, Kim YM: Damage of an alumina-on-alumina bearings surface from a difficult reduction of a total hip arthroplasty. *J Bone Joint Surg Am* 86: 376-8, 2004
- Yoo JJ, Kim YM, Yoon KS, et al.: Alumina-on-alumina total hip arthroplasty. A five-year minimum follow-up study. *J Bone Joint Surg Am* 87: 530-5, 2005
- Taylor S, Manley MT, Sutton K: The role of stripe wear in causing acoustic emissions from alumina ceramic-on-ceramic bearings. *J Arthroplasty* 22 (7 Suppl 3): 47, 2007

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Mr. S. Meakins

Offenburg Ebertplatz

Quality Management and Regulatory Affairs Manager – Implants

Offenburg St. Josefsklinik

Orthodynamics GmbH

Gengenbach

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23556 Lübeck

Traumatologisches Zentrum

Offenburg Ebertplatz

Tel. 0781/472-3301 Kliniksekretariat

traumatologie@og.ortenau-klinikum.de

Tel. 0781/472-3311 P.O. Ambulanz

www.ortenau-klinikum.de

22nd February 2011

Re: GHE Adapter stems

Dear Mr. Meakins,

I can confirm that I have been using the cemented GHE Adapter stem for primary hip replacement since 2006.

During this time, I have implanted approximately 35 Adapter stems per year, over 150 implantations in total.

Whilst there has been no official follow-up of these devices, from my experience the incidence of problems is low and the implant survivorship is good. In my opinion, the cemented Adaptor stem offers a good clinical outcome that is comparable to other stems used to treat primary hip arthroplasty indications.

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
Re: GHE Adapter stems

Dear Mr. Meakins,

I can confirm that I have been using the cementless GHE Adapter stem for primary hip replacement since 2004.

During this time, I have implanted approximately 38 Adapter stems per year, over 200 implantations in total.

Whilst there has been no official follow-up of these devices, from my experience the incidence of problems is low and the implant survivorship is good. In my opinion, the cementless GHE Adaptor stem offers a good clinical outcome that is comparable to other stems used to treat primary hip arthroplasty indications.

Yours sincerely,

Dr. Seifert

Apolda

Metal release and corrosion effects of modular neck total hip arthroplasty

J. Philippe Kretzer · Eike Jakubowitz ·
Michael Krachler · Marc Thomsen · Christian Heisel

Received: 5 December 2008 / Accepted: 20 December 2008 / Published online: 14 February 2009
© Springer-Verlag 2009

Abstract Modular neck implants are an attractive treatment tool in total hip replacement. Concerns remain about the mechanical stability and metal ion release caused by the modular connection. Five different implant designs were investigated in an experimental set-up. In vivo conditions were simulated and the long-term titanium release was measured. Finally, the modular connections were inspected for corrosion processes and signs of fretting. No mechanical failure or excessive corrosion could be identified for the implants tested. The titanium releases measured were extremely low compared to in vivo and in vitro studies and were not in a critical range.

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Introduction

In current total hip arthroplasty, modularity between the head-neck connection of the femoral stem is well established [1]. More recently, implants with an additional modular connection between the neck and stem have become more popular. The ability to adjust the head centre relative to the stem axis intraoperatively, by connecting the neck in different positions, guarantees high versatility during surgery.

The mechanical loading conditions of the neck-stem connection fundamentally differ from the well-established head-neck connection. The forces at the established head-neck connection are transmitted centrally through the head resulting in low stresses. In contrast, the eccentric load on the modular neck-stem connection leads to higher stresses at the connection.

Fretting caused by micro-motion between the modular components may lead to particle and ion release to the surrounding tissue [4, 12, 21]. Titanium-containing particles may cause the release of cytokines that are associated with periprosthetic bone loss [25]. Such particles may get trapped between the articulating surfaces of the artificial joint resulting in third-body wear [21]. Furthermore, mechanically assisted crevice corrosion (as a combination of fretting and crevice corrosion) may attack the alloy and

J. P. Kretzer (✉) · E. Jakubowitz
Laboratory of Biomechanics and Implant Research,
Department of Orthopaedics, University of Heidelberg,
Schlierbacher Landstrasse 200A,
69117 Heidelberg, Germany
e-mail: philippe.kretzer@web.de

M. Krachler
Institute of Earth Sciences, University of Heidelberg,
Heidelberg, Germany

M. Thomsen
German Red Cross Hospital,
Lilienmattstr. 5,
76530 Baden-Baden, Germany

C. Heisel
ARCUS Kliniken Pforzheim,
Department of Orthopaedics and Traumatology,
Rastatter Strasse 17-19,
75179 Pforzheim, Germany

lead to implant failure [7]. Modularity is often associated with corrosion effects and is seen as a weak link with a possible source of complications [24]. A multicentre retrieval analysis of 231 modular hip implants revealed corrosion processes at the modular head-neck connection in more than 28% of the cases [8].

The aim of this study was to analyse and compare the corrosion effects and to evaluate the particle and ion release of the modular neck hip implants currently used. The following questions were addressed: (1) Are there differences in the total titanium release between different implant designs? (2) How does the geometry of the taper influence the titanium release? (3) Are there differences in the titanium release over time?

Material and methods

Three implants were analysed (Fig. 1): Eco-Modular® (Endoplast, Marl, Germany), Varicon® (Falcon Medical, Mödling, Austria), Metha® (Aesculap, Tuttlingen, Germany) and SPS-Modular® (Symbios, Yverdon, Switzerland). All implants were made of Ti-6Al-4V alloy. Additionally, a universal modular neck adapter (Bio-Ball®, Merete, Berlin, Germany) which is connectable to standard hip stems (12/14 taper) was included. This adapter was combined with a CF30®-stem (Zimmer, Warsaw, IN, USA). The Metha® and SPS-Modular® taper shapes were oval, designed with a female stem connection. All other taper geometries were round and had a male stem connection (Fig. 1).

The stem and neck sizes were chosen on the basis of the femoral size of an average patient. Comparable biomechanical loading conditions were achieved by selecting similar


femoral anteversion and head neck diaphyseal (CCD) angles and lateral offsets for all implants (anteversion: 6.8°; CCD: 133.1°; lateral offset: 46.5 mm).

The taper angles were determined using a coordinate measuring machine (CMM) (MarVision MS 222, Mahr, Göttingen, Germany; accuracy: $\pm 2.3 \mu\text{m}$) and the surface roughness was analysed using a roughness measurement instrument (Perthometer M2, Mahr, Göttingen, Germany; accuracy: 12 nm). The stems and modular necks of each design were randomly paired. The resulting taper angle differences (cone angle differences) were calculated and the total surface roughness of each couple was determined quadratically. All geometrical measurements were repeated five times.

The assembly load and environment play a major role in the initial stability of modular connections [16]. To simulate surgical conditions the modular connections were connected wet by a hammer stroke. The stems were orientated in 10° adduction and in 9° flexion resulting in bending and torsional loading of the stem and stem-modular neck connection according to ISO 7206 (Fig. 2). A sinusoidal load of between 0.3 kN (min.) and 2.3 kN (max.) was applied for 10×10^6 cycles. The frequency was repeatedly altered between 3 Hz (100,000 cycles) followed by 15 Hz (900,000 cycles).

The particles and ions released from the stem-modular neck connection were measured using a fluid reservoir filled with 800 ml calf serum (diluted to a protein content of 30 g/l, at 37°C, pH: 7.4). Serum was chosen because proteins are reported to affect the corrosion resistance of titanium alloy [11, 14]. The serum was continuously circulated by a peristaltic pump from the reservoir to a small simulation chamber which surrounded the modular connection (Fig. 2).

Fig. 1 Implant designs investigated. The oval modular taper of the two designs on the right side are female connections


Fig. 2 Alignment of the stem and modular neck relative to spatial coordinates. The peristaltic pump continuously circulates serum between the fluid reservoir and simulation chamber

The titanium concentration of the serum was analysed using HR-ICP-MS (Element2, Thermo Fisher Scientific, Bremen, Germany) at intervals of 0 , 1×10^6 , 2×10^6 , 4.5×10^6 , 7×10^6 and 10×10^6 cycles. In order to detect not only ions but also particles in the serum, the samples were first digested with high purity nitric acid (HNO_3) and hydrogen peroxide (H_2O_2) following a standardised, well-established procedure [10]. The titanium mass was calculated based on the test medium's volume and the measured titanium concentration.

All data are presented as mean $\pm$ standard deviation (SD). An analysis of variance (ANOVA) was used to prove whether there were differences between the prostheses in total titanium release. Pearson's correlation test was performed to deal with geometrically related influences (total surface roughness, mean cone angle and delta cone angle) on total titanium release. Student's t -test was used to detect possible influences of the taper shape (round vs oval). Regression analysis was performed to demonstrate a progression in titanium release over time. All statistical

analyses were performed using SPSS® software (SPSS® for Windows 16.0.1, SPSS Inc. Chicago, IL, USA). A p value < 0.05 was considered significant.

Finally, the tapers were carefully disconnected and inspected by scanning electron microscope (SEM) (Type 440, Carl Zeiss, Cambridge, UK) for corrosion effects.

Results

The total titanium release of all systems varied between 12 and 44 μg at the end of the test (Fig. 3). No significant differences in total titanium release were found between the tested implants ($p=0.66$). Table 1 summarises the results of the geometrical measurements. Neither the mean cone angle ($p=0.45$), the cone angle difference ($p=0.30$) nor the shape of the taper ($p=0.44$) had a significant influence on the total titanium release. Based on an almost significant ($p=0.06$) negative correlation coefficient of $r=-0.86$, the total surface roughness indeed seemed to have an effect. A significant linear progression in titanium release over time was found for the Eco-Modular® implant ($p<0.01$, $r=0.98$), the SPS-Modular® implant ($p<0.01$, $r=0.99$) and the Bio-Ball® implant ($p=0.04$, $r=0.90$). No linear correlation was found for the Varicon® implant ($p=0.30$, $r=0.59$) and the Metha® implant ($p=0.19$, $r=0.70$).

SEM revealed moderate signs of fretting in all tested connections. Figure 4 shows scars caused by micro-motion which perpendicularly obliterate machine lines. No severe corrosion effects were detected in any of the implants investigated.

Discussion

The literature unequivocally demonstrates that design- and material-related parameters influence micro-motion and


Fig. 3 Titanium release of the five implant designs during the course of simulation

Table 1 Results of the cone geometry measured

	Total surface roughness (μm)		Mean cone angle ($^{\circ}$)		Cone angle difference ($^{\circ}$)	
Endoplant Eco-Modular®	1.28	± 0.31	5.67	± 0.02	0.04	± 0.04
Falcon Varicon®	3.48	± 0.17	5.80	± 0.01	0.31	± 0.03
Aesculap Metha®	0.72	± 0.02	3.69	± 0.01	0.09	± 0.04
Symbios SPS®	0.34	± 0.01	5.01	± 0.01	0.10	± 0.02
Merete Bio-Ball®	12.15	± 0.39	5.68	± 0.00	0.11	± 0.03


fretting of implant connections [8, 9, 11, 13, 15]. However, no differences in total titanium release were found for the implant designs tested here. Compared to other studies the titanium release of all systems investigated in this study is low. A similar study on a modular neck hip implant design made of a titanium alloy revealed material loss of between 280 and 1,640 μg after 5.5×10^6 loading cycles [22]. Another study on a titanium plate-screw connection resulted in a titanium release of 849 μg after 1.2×10^6 loading cycles [13]. Lower concentrations ranging from 3 to 21 μg were also reported for plate-screw connections [11]. However, this study was only performed for 0.4×10^6 loading cycles.

The surface roughness showed a negative linear correlation to the total titanium release. With increasing surface roughness the titanium release decreased. Higher surface roughness may locally increase and evenly distribute the contact pressure [2] and thus decrease micro-motions between both components. The reported effect of increased micro-motion and fretting depending on the cone angle difference [20, 23] was not seen in this study. The implants showed different patterns with respect to progression in titanium release over time. For the Symbios SPS®, Endoplant Eco-Modular® and Merete Bio-Ball® implants the titanium release proceeds in a linear manner. However, the Aesculap Metha® and Falcon Varicon® systems showed a titanium release which increased during the first 2×10^6

cycles and subsequently remained stable (Fig. 3). Similar progressions are reported for wear studies on metal-on-metal joint bearings in total hip replacement [10]. During the initial run-in phase of those implants, the surface geometries of both counterparts adapt to each other. Such a mechanism may also occur at the modular connection leading to decreased micro-motion and resulting in a more stable connection. The explanation for the measured steady-state titanium concentrations over time in two of the five tested implant designs may be related to a limited fluid transfer between the crevice of the modular connection and the surrounding fluid. If the fluid transfer is restricted, metal ions caused by fretting may not enter into the surrounding fluid. In this case mechanically assisted crevice corrosion as reported by Gilbert et al. [7] may occur. Based on fretting and subsequent repassivation of the titanium, the local concentration of free oxygen will drop and this will result in an increased concentration of free metal ions in the crevice. The excess of metal ions then attracts chloride ions to form metal chlorides. These metal chlorides will react with water to form metal hydroxide and hydrochloric acid, lowering the pH and resulting in a hydrochloric acid solution with a very low pH in the crevice. The metal then loses its passive film and becomes thermodynamically unstable. Subsequent failure of the modular connection might occur. However, in this study neither mechanically assisted crevice corrosion nor other severe corrosion effects could be identified by final SEM analysis. No gross material failure occurred.

Forces of normal gait were applied in this study [3]. For active or obese patients higher forces may occur in vivo, especially in combination with high offset necks.

Titanium is the commonly used implant material. In spite of its functional benefits and local biocompatibility, concerns exist about the long-term risks caused by the metal particles and ion release. Macrophage activation caused by titanium particles and increased cell mediator release has been observed in vitro and in vivo [18, 19]. In vitro chromosomal damage caused by high concentrations of titanium particles have also been reported [5, 6]. Nevertheless, such reactions depend strongly on the applied particle or ion concentration. Rogers et al. investigated the effect of different titanium particle concentrations on the cell mediator release and cytotoxicity [18]. They applied

**Fig. 4** Fretting scars run perpendicular to machine lines

particle concentrations of between 14 and 58 $\mu\text{g/ml}$ to human monocytes in vitro and reported an increase in cell mediator release depending on the particle concentration. However, a significant toxicity was not seen in their study. Assuming a patient activity of 2×10^6 cycles/year [10], the mean titanium release of all tested systems in our study would correspond to an annual total titanium release of about 6 μg to the whole organism. Compared to in vivo or other in vitro studies this concentration is extremely low. Clinically, reactions to titanium are rare and in vivo studies also confirmed no genotoxicity or cytotoxicity of titanium [17]. Consequently, the authors believe that the clinical use of these implants would appear to be appropriate.

The titanium concentrations exhibited quite a high variability. The measurement of small values near the threshold of a specific measuring method generates wide ranges. Thus, the variability measured for some values in this study may be related to the low values generally measured.

Conclusion

The results of this study are promising and can support limited clinical use of these implants. Very low concentrations of titanium release were found. No differences between the different designs in total titanium were seen. However, the surface roughness of the taper connection seems to affect the total titanium release of the implants tested. Different modes of titanium release over time were measured. No mechanical failure or excessive corrosion processes could be identified on the implants tested. Compared to in vivo and other in vitro studies extremely low concentrations of titanium release were seen. From the authors' point of view the measured titanium concentrations are within a clinically non-critical range. The safety of the modular neck-stem connection still needs to be proven in clinical studies.

Acknowledgments The authors thank the Ministry of Arts and Science of Baden-Württemberg (Germany) for supporting this work with research grants and Tom Bruckner (Ph.D.) from the Institute of Biometrics of the University of Heidelberg for expert advice on statistics.

References

1. Aldinger PR, Thomsen M, Mau H, Ewerbeck V, Breusch SJ (2003) Cementless Spotorno tapered titanium stems: excellent 10–15-year survival in 141 young patients. *Acta Orthop Scand* 74:253–258
2. Bachnick M, Hasenpusch M, Richter H, Boenick U (1994) Effect of hardness and surface quality of metal cones on fracture load of ceramic ball heads of hip endoprostheses (in German). *Biomed Tech (Berl)* 39:302–306
3. Bergmann G, Deuretzbacher G, Heller M, Graichen F, Rohlmann A, Strauss J, Duda GN (2001) Hip contact forces and gait patterns from routine activities. *J Biomech* 34:859–871
4. Bobyn JD, Tanzer M, Krygier JJ, Dujovne AR, Brooks CE (1994) Concerns with modularity in total hip arthroplasty. *Clin Orthop Relat Res* 298:27–36
5. Coen N, Kadhim MA, Wright EG, Case CP, Mothersill CE (2003) Particulate debris from a titanium metal prosthesis induces genomic instability in primary human fibroblast cells. *Br J Cancer* 88:548–552
6. Daley B, Doherty AT, Fairman B, Case CP (2004) Wear debris from hip or knee replacements causes chromosomal damage in human cells in tissue culture. *J Bone Joint Surg Br* 86:598–606
7. Gilbert JL, Buckley CA, Jacobs JJ (1993) In vivo corrosion of modular hip prosthesis components in mixed and similar metal combinations. The effect of crevice, stress, motion, and alloy coupling. *J Biomed Mater Res* 27:1533–1544
8. Goldberg JR, Gilbert JL, Jacobs JJ, Bauer TW, Paprosky W, Leurgans S (2002) A multicenter retrieval study of the taper interfaces of modular hip prostheses. *Clin Orthop Relat Res* 401:149–161
9. Hallab NJ, Messina C, Skipor A, Jacobs JJ (2004) Differences in the fretting corrosion of metal-metal and ceramic-metal modular junctions of total hip replacements. *J Orthop Res* 22:250–259
10. Heisel C, Streich N, Krachler M, Jakubowitz E, Kretzer JP (2008) Characterization of the running-in period in total hip resurfacing arthroplasty: an in vivo and in vitro metal ion analysis. *J Bone Joint Surg Am* 90(Suppl 3):125–133
11. Hol PJ, Molster A, Gjerdet NR (2008) Should the galvanic combination of titanium and stainless steel surgical implants be avoided? *Injury* 39:161–169
12. Jacobs JJ, Urban RM, Gilbert JL, Skipor AK, Black J, Jasty M, Galante JO (1995) Local and distant products from modularity. *Clin Orthop Relat Res* 319:94–105
13. Kawalec JS, Brown SA, Payer JH, Merriitt K (1995) Mixed-metal fretting corrosion of Ti6Al4V and wrought cobalt alloy. *J Biomed Mater Res* 29:867–873
14. Khan MA, Williams RL, Williams DF (1999) The corrosion behaviour of Ti-6Al-4V, Ti-6Al-7Nb and Ti-13Nb-13Zr in protein solutions. *Biomaterials* 20:631–637
15. Lewis G (2004) Geometric element analysis of fretting in a model of a modular femoral component of a hip implant. *Biomed Mater Eng* 14:43–51
16. Mroczkowski ML, Hertzler JS, Humphrey SM, Johnson T, Blanchard CR (2006) Effect of impact assembly on the fretting corrosion of modular hip tapers. *J Orthop Res* 24:271–279
17. Piozzi R, Ribeiro DA, Padovan LE, Filho HN, Matsumoto MA (2009) Genotoxicity and cytotoxicity in multiple organs induced by titanium miniplates in Wistar rats. *J Biomed Mater Res A* 88:342–347
18. Rogers SD, Howie DW, Graves SE, Percy MJ, Haynes DR (1997) In vitro human monocyte response to wear particles of titanium alloy containing vanadium or niobium. *J Bone Joint Surg Br* 79:311–315
19. Schwab LP, Xing Z, Hastly KA, Smith RA (2006) Titanium particles and surface-bound LPS activate different pathways in IC-21 macrophages. *J Biomed Mater Res B Appl Biomater* 79:66–73
20. Shareef N, Levine D (1996) Effect of manufacturing tolerances on the micromotion at the Morse taper interface in modular hip implants using the finite element technique. *Biomaterials* 17:623–630
21. Urban RM, Jacobs JJ, Gilbert JL, Galante JO (1994) Migration of corrosion products from modular hip prostheses. Particle micro-

- analysis and histopathological findings. *J Bone Joint Surg Am* 76:1345–1359
22. Viceconti M, Baleani M, Squarzone S, Toni A (1997) Fretting wear in a modular neck hip prosthesis. *J Biomed Mater Res* 35:207–216
23. Viceconti M, Ruggeri O, Toni A, Giunti A (1996) Design-related fretting wear in modular neck hip prosthesis. *J Biomed Mater Res* 30:181–186
24. Virtanen S, Milosev I, Gomez-Barrena E, Trebse R, Salo J, Kontinen YT (2008) Special modes of corrosion under physiological and simulated physiological conditions. *Acta Biomater* 4:468–476
25. Warashina H, Sakano S, Kitamura S, Yamauchi KI, Yamaguchi J, Ishiguro N, Hasegawa Y (2003) Biological reaction to alumina, zirconia, titanium and polyethylene particles implanted onto murine calvaria. *Biomaterials* 24:3655–3661

Scientific Exhibit at the 71st AAOS Annual Meeting

MODULAR NECK PRIMARY PROSTHESIS: EXPERIMENTAL AND CLINICAL OUTCOMES

Francesco Traina^{1,2}, MD, Massimiliano Baleani², M Eng,
Marco Viceconti², PhD, M Eng, Aldo Toni^{1,2}, MD

1. I Department of Orthopaedic Surgery, Rizzoli Orthopaedic Institute, Bologna, Italy
2. Medical Technology Laboratory, Rizzoli Orthopaedic Institute, Bologna, Italy

San Francisco, March 10-14, 2004

INTRODUCTION

Modular acetabular designs, due to their versatility, are widespread used in primary total hip replacements while primary modular stem designs are rarely employed. Nevertheless primary stem modularity could be useful to manage difficult cases, when the anatomy is overthrown, and for mini-incision approaches. Stem modularity, as a matter of fact, provides an increased adaptability without any need for a large inventory of monoblock prostheses or expensive custom made prostheses. Besides the availability of a variable geometry stem could be useful to restore the leg length or the hip offset and to correct a pathological femoral anatomy.

While revision modular stems, with different interchangeable components^{1,2,3} are every day more often employed and well documented, very little is reported on primary modular stems.

THE MODULAR PROSTHESIS

The stem employed for the pre-clinical validation and in the clinical practice is a titanium alloy stem (Ti6Al4V) anatomically shaped (ANCA-Fit, Wright-Cremascoli), grit blasted (Ra 6 μm) and coated with

80 μm high crystalline plasma sprayed hydroxyapatite (Ca/P rate of 1.67 ± 0.05 ; purity >97%; crystallinity >60%) in the proximal third. In the proximal end a double tapered housing for the modular necks is provided. The distal 2/3 of the stem are tapered and grit blasted (Ra 2 μm).


Fig. 1 a) Schematic drawing showing the possible implant solution with a straight, an 8° and a 15° retroverted neck. b) Different possible implant lengths and offset with a straight or a varus neck of the two lengths. c) Picture of the ANCA Fit hip system.

MODULAR NECKS

The modular necks are made of titanium (Ti6Al4V), present an oblong section and a conical design to limit the risk of impingement. A thin titanium oxide layer (2-5 nm thickness) cover the whole neck surface providing an exceptional corrosion resistance. The neck is connected with the modular head by a cylindrical taper fitting, while is connected to the stem by a taper fitting having a rectangular cross section with semicircular short sides. Necks present two different lengths: short (28 mm) and long (38.5 mm). For both lengths there are 6 different models: straight, varus-valgus of 8°, antverted-retroverted of 8° or 15°, the combination of 6° of varus and 4.5° of retroversion for the left and the right side. The modular heads are available in three sizes: short (-3.5 mm), medium (0 mm) and long (+3.5 mm) providing with the two necks a range of possible lengthening of 17.5 mm.


Fig. 2 Modular necks: drawings representing the possible length solution and the different geometrical designs. AR-VV1: retroverted of 6° and varus of 4.5° for the left side; AR-VV2: retroverted of 6° and varus of 4.5° for the right side.

AIM OF THE STUDY

The aim of the study was to present the pre-clinical validation and the clinical results of 1734 modular neck primary prostheses at 8 years of survival.

PRE-CLINICAL VALIDATION

The fretting-corrosion behaviour of the neck-stem coupling was investigated. A validated internal procedure^{4,5}, defined on the basis of the ISO 7206 recommendations for fatigue test of hip stem, was used. Different testing conditions were fixed to investigate the effect of the environmental condition, applied load, stem size, and test duration (that represents the *in vivo* prosthesis life) on the damage process of the coupling. The coupling surface underwent to microscope analysis to determine the morphology of the damage. Finally, the amount of particulate produced during the test was calculated by means of the weight loss of the neck.

Experimental Results

Although long load history (3300 N applied for 20 million cycles that conventionally represent 20 years of active patient use) or high load level (4200 N applied for 5.5 million cycles that would represent 5 years of jogging activity), **no mechanical fracture of the neck-stem coupling occurred during the tests.** The neck was firmly assembled to the stem at the end of the test.

The microscopical analysis of the coupling surface showed no evidence of primary corrosion also when the test was conducted in chemically aggressive environment (FeCl_3 solution). Flakes of spalled material were observed in the burnished areas suggesting the presence of a mechanical damage due to micromotion


Fig. 3 a) Schematic drawing showing the experimental set-up. b) Picture of a modular neck during the tests. c) Modular neck after 5.5 million cycles at 3300 N: red circles show the burnished areas.

(fretting). Nevertheless, the extension of the damaged area was affected by the test environment suggesting a mechanical damage chemically enhanced. The damage extension is linearly related to the load level, the stem size and the number of applied load cycles. However, **the amount of fretted material was estimated 0.6 mg/million cycles.**

CLINICAL EXPERIENCE

From August 1995 to August 2003, 1752 primary surgeries were performed with a modular neck stem in our Hospital; all the prostheses present a ceramic coupling with 28 mm head. To evaluate the influence of the learning curve with modular neck prostheses the population was divided in two cohorts. The first 864 consecutive surgeries performed from August 1995 to December 2000 were included in the first cohort, and of these patient 16 were lost to follow-up leaving under observation 848 prostheses. The following 888 consecutive surgeries performed between January 2001 and August 2003 were included in the second cohort: in this group 2 patients were lost to follow-up leaving 886 prostheses under control. Summarizing 18 patients (1.02%) out of 1752 were lost to follow-up, the remaining 1734 patients were investigated.

Considering the pre-operative diagnosis, with the conventional related anatomical deformity of the hip joint, the 1734 prostheses were divided in 4 major groups: A) Normal hip morphology; B) Mild hip deformity; C) Moderate hip deformity; and D) Severe Hip Deformity (Figure 4).

The demographic features of the 2 cohorts and of the whole population are shown in Table 1: there is not any statistically difference for age, gender and pre-operative diagnosis between the three.

Modular neck distribution was investigated in relation to the four different groups of hip morphology and in relation to the surgeons' learning curve.

The modular stem survival was estimated by the Kaplan-Meier method⁶. All the patients considered followed up had at least one follow-up control in the last 18 months; patients who were not reviewed in the last 18 months were checked by phone to assess their clinical status and the survival of the stem. Patients who were not reached by phone were considered truly lost to follow-up.


Fig. 4 Patients pre-operative diagnoses were grouped in 4 major categories considering the conventional related anatomical deformity.

	Aug 1995 - Dec 2000 848 prostheses	Jan 2001 – Aug 2003 886 prostheses	Aug 1995 – Aug 2003 1734 prostheses
Surgeries/year			
Average age	54.7 years (16-81)	56.4 years (14-86)	55.5 years (14-86)
Gender	52.6% ♀ 47.4% ♂	54.1% ♀ 45.9% ♂	53.3% ♀ 46.7% ♂
Pre-op diagnosis	 6.6% 66.7% 1.8% 24.9%	 9.3% 66% 1.5% 23.3%	 8% 66.4% 1.6% 24%

Tab. 1 Demographic features of the two cohorts and the whole population.

Clinical Results

None of the 1734 prostheses needed revision for a modular neck failure. In our experience the straight neck was implanted in 56.5% of the cases of group A, 55.4% of B, 39.3% of C, and 21.6% of D (Figure 5). This distribution shows how modularity is always more over useful increasing the difficulty of the case, but also how frequently a non-straight neck is employed in easy or standard cases (A and B). Figure 6 shows how non-straight necks were more employed in the second 886 cases because of the increased familiarity of the surgeons with modularity opportunities.

Dislocation rate of the 1734 prostheses, with a ceramic coupling without a liner lip, was 1.4% (25 patients), 6 cases of recurrent dislocation led to revision surgery (0.3%). Eighteen implants failed, 6 for recurrent dislocation, 5 for late aseptic loosening, 3 for early failure (2 cases for intra-operative cracks leading to femoral fractures, 1 for stem undersizing), and 3 for septic loosening. The survival rate of the modular neck stem, without septic loosening, is 97.5% (C.I. 94.9-100%) at 8 years (Figure 7).


Fig. 5 Necks distribution for different anatomical deformity. **Top:** non-straight necks were used mostly for difficult cases (C-D) but also in about 45% of the easier cases. **Bottom:** Neck length choice is not related to the preoperative diagnosis.


Fig. 6 Necks distribution in relation to surgeons learning curve. Non-straight necks were more employed in the second cohort when surgeons familiarity with modularity opportunities increased.


Fig. 7 Survival analysis at 8 years of the modular neck stem. No stem failure was due to the modular neck.

MODULAR NECK CLINICAL ADVANTAGES

Modular necks present several advantages, the 2 neck lengths and the 3 head lengths enable a good modulation of the implant length and offset: indeed, the former can vary between 24.5 mm (28 mm - 3.5 mm) and 42 mm (38.5 mm + 3.5 mm) considering the 2 different neck lengths matched with the available 3 head lengths. Furthermore, neck modularity allows the surgeon to correct intra-operatively the implant configuration, based on intra-operative trials of the hip range of motion, reducing the risk of impingement (Figure 8). Besides, different neck designs allow the correction of a pathological anatomy respecting the hip center of rotation or, in very difficult cases, allow to change the hip center of rotation respecting the hip offset (Figure 10). In addition, the modular prosthesis allows the surgeon to replace only the modular components of the implant; in case of revision this could be very useful, for example, for revisions due to recurrent dislocation when a very easy substitution of the neck could be the only surgery needed.


Fig. 8 In the drawing is shown how modular necks could be useful in case of impingement: a long retroverted neck prevents the impingement inducible extra-rotating the leg with a short straight neck.


Fig. 9 a) Two superimposed CT scan slices of a DDH case showing a highly pathological femoral neck anteversion. b) a straight neck does not allow an anatomic reconstruction. c) a 15° retroverted neck allows the correction of the deformity.



Fig. 10 a) X-ray of a DDH. Single leg stance. C: hip center of rotation; F: resultant force of abductor muscles; W: body weight; a_F : moment arm of the abductor muscles force; a_W : moment arm of the body weight. The condition for the rotational equilibrium of the pelvis about C is $F \cdot a_F = W \cdot a_W$. b) Ideal solution of hip replacement with a straight short neck and the cup in the ideal hip center. c) By this solution the hip offset is restored and a good leg lengthening is achieved, but an augmentation technique with a bone graft is needed for the cup stability. d) To avoid bone grafting the cup should be implanted with a higher center of rotation. e) Despite the high hip center of rotation, with a long varus modular neck, the hip offset, the leg lengthening and the moment arm of the system are respected ($a_R = a_Y$).

CASE REPORTS


Case 1. Female, 54 years old, showing a bilateral DDH.

Right hip: a) Pre-op frontal view; b) Post-op frontal view: a short 15° retroverted neck was employed; b') Post-op frog view; c) and c') 7.5 years follow-up.

Left hip: a) Pre-op frontal view; d) Post-op frontal view: a long 8° retroverted neck was employed; d') Post-op frog view; e) and e') 6.5 years follow-up; f) frontal view at the last follow-up.


SHORT VARUS

Case 2. Female, 64 years old, primary arthritis.

a) Pre-op frontal view; b) Post-op: a short 8° varus neck was employed; c) Post-op frog view; d) 8 years follow-up; e) Follow-up, frog view.

LONG STRAIGHT

Case 3. Female, 57 years old, primary arthritis.

a) Pre-op frontal view; b) Post-op: a long straight neck was employed; c) Post-op frog view; d) 7 years follow-up; e) Follow-up, frog view.

DISCUSSION AND CONCLUSION

Modular prostheses allow the surgeon to correct the length of the lower limb by choosing the offset, the length and the anti-retroversion of the neck; besides, this is possible avoiding the long-waiting times and high costs of custom-made prostheses, and assuring the reliability of mechanical tests (done on standard prostheses). Such modularity is particularly useful in cases where the femoral deformity is threatening (only in the 21.6% of group D and in the 39.3 % of group C a straight neck was employed) but also in standard cases as shown in Figure 5. Indeed in about 45% of the prostheses implanted with a normal or a mild anatomically deformed hip (groups A and B) a non-straight neck has been employed. This results show the poor predictive value of a 2D templating that usually underestimate the femoral neck anteversion and is not useful to predict the possible causes of impingement. Furthermore the impingement due to a cup malpositioning is usually corrigible, during surgery, by a modular neck. Besides, the possibility of increase the stem length enhances the implant stability, this is particularly useful with the old age patients with muscular weakness and in the cases where the acetabular anatomy is badly compromised (group D). Indeed in our series the rate of dislocation is only 1,4% and there is not any significant difference between patients over and under 65 years old and between the 4 joint morphological cohorts (A; B; C; D) contradicting what was expected considering the letterature.

In spite of these advantages, a modular prosthesis has some concerns related to the risk of corrosion, fretting^{7,8,9,10,11}, dislocation of the modular components^{12,13,14} and mostly of fracture of the modular neck¹⁵. These potential problems risen in presence of two coupled components (head-neck) could be worsened in case of a further modularity between the neck and the stem. The experimental study tried to address the causes of concern. The investigated modular necks have shown a good mechanical behaviour *in vitro*. The neck-stem coupling undergoes to a fretting damage process. However, the amount of produced fretting debris (0.6 mg per million cycles) seems negligible considering that a stable prosthesis is likely to produce more than 10 mg/year of metal debris^{4,5}.

These experimental findings are in agreement with the clinical outcomes. Besides, among the 1734 modular necks implanted in our hospital, out of more than 50.000 modular necks already implanted in Europe, no one needed a revision for mechanical failure. In the explanted cases, the cause of revision was not related to the neck-stem modularity. Furthermore, the comparative analysis between retrieved necks and those experimentally studied confirmed the absence of corrosion¹⁶ and the completely absence of periprosthetic metallosis. These good results are to be related to the corrosion resistant features of the titanium-titanium coupling and the design of the taper.

In conclusion the use of a modular prosthesis should be routinely reserved for difficult primary surgeries or revision surgeries, where the non-modular prostheses could be not satisfactory. Nevertheless in our experience modular necks are very useful even for standard procedure to achieve a good implant stability and the best implant offset.

REFERENCES

- Wirtz DC, Heller KD, Holzwarth U, Siebert C, Pitto RP, Zeiler G, Blencke BA, Forst R. A modular femoral implant for uncemented stem revision in THR. *Int Orthop*. 2000;24(3):134-8. Review.
- Berry DJ. Treatment of Vancouver B3 periprosthetic femur fractures with a fluted tapered stem. *Clin Orthop*. 2003 Dec;(417):224-31.
- Sporer SM, Paprosky WG. Revision total hip arthroplasty: the limits of fully coated stems. *Clin Orthop*. 2003 Dec;(417):203-9.
- Viceconti M, Ruggeri O, Toni A, Giunti A. Design-related fretting wear in modular neck hip prosthesis. *J Biomed Mater Res*. 1996 Feb;30(2):181-6.
- Viceconti M, Baleani M, Squarizoni S, Toni A. Fretting wear in a modular neck hip prosthesis. *J Biomed Mater Res*. 1997 May;35(2):207-16.
- Kaplan EL, Meier P. Non parametric estimation from incomplete observations. *J Am Stat Assoc*. 1958;53:457-81.
- Mathiesen EB, Lindgren JU, Blomgren GG, Reinholdt FP. Corrosion of modular hip prostheses. *J Bone Joint Surg Br*. 1991 Jul;73(4):569-75.
- Bobyn JD, Tanzer M, Krygier JJ, Dujovne AR, Brooks CE. Concerns with modularity in total hip arthroplasty. *Clin Orthop*. 1994 Jan;(298):27-36.
- Chimell MJ, Rispler D, Poss R. The impact of modularity in total hip arthroplasty. *Clin Orthop*. 1995 Oct;(319):77-84.
- Collier JP, Surprenant VA, Jensen RE, Mayor MB, Surprenant HP. Corrosion between the components of modular femoral hip prostheses. *J Bone Joint Surg Br*. 1992 Jul;74(4):511-7.
- Cook SD, Barrack RL, Baffes GC, Clemow AJ, Serekian P, Dong N, Kester MA. Wear and corrosion of modular interfaces in total hip replacements. *Clin Orthop*. 1994 Jan;(298):80-8.
- Pellicci PM, Haas SB. Disassembly of a modular femoral component during closed reduction of the dislocated femoral component. A case report. *J Bone Joint Surg Am*. 1990 Apr;72(4):619-20.
- Star MJ, Colwell CW Jr, Donaldson WF 3rd, Walker RH. Dissociation of modular hip arthroplasty components after dislocation. A report of three cases at differing dissociation levels. *Clin Orthop*. 1992 May;(278):111-5.
- Woolson ST, Pottorff GT. Disassembly of a modular femoral prosthesis after dislocation of the femoral component. A case report. *J Bone Joint Surg Am*. 1990 Apr;72(4):624-5.
- Salvati EA, Lieberman JR, Huk OL, Evans BG. Complications of femoral and acetabular modularity. *Clin Orthop*. 1995 Oct;(319):85-93.
- Toni A, Paderni S, Sudanese A, Guerra E, Traina F, Giardina F, Antonietti B, Giunti. Anatomic cementless total hip arthroplasty with ceramic bearings and modular necks: 3 to 5 years follow-up. *Hip International* 2001 Jan-Mar;11(1):1-7.

We wish to thank Luigi Lena* for preparing drawings, pictures and layout; Barbara Bordini* and Manuela Declerico* for the statistics of this scientific exhibit.

* Medical Technology Laboratory, Rizzoli Orthopaedic Institute, Bologna, Italy

RESEARCH ARTICLE

Open Access

Modular titanium alloy neck adapter failures in hip replacement - failure mode analysis and influence of implant material

Thomas M Grupp^{1*}, Thomas Weik^{1†}, Wilhelm Bloemer¹, Hanns-Peter Knaebel²

Abstract

Background: Modular neck adapters for hip arthroplasty stems allow the surgeon to modify CCD angle, offset and femoral anteversion intraoperatively. Fretting or crevice corrosion may lead to failure of such a modular device due to high loads or surface contamination inside the modular coupling. Unfortunately we have experienced such a failure of implants and now report our clinical experience with the failures in order to advance orthopaedic material research and joint replacement surgery.

The failed neck adapters were implanted between August 2004 and November 2006 a total of about 5000 devices. After this period, the titanium neck adapters were replaced by adapters out of cobalt-chromium. Until the end of 2008 in total 1.4% (n = 68) of the implanted titanium alloy neck adapters failed with an average time of 2.0 years (0.7 to 4.0 years) postoperatively. All, but one, patients were male, their average age being 57.4 years (36 to 75 years) and the average weight 102.3 kg (75 to 130 kg). The failures of neck adapters were divided into 66% with small CCD of 130° and 60% with head lengths of L or larger. Assuming an average time to failure of 2.8 years, the cumulative failure rate was calculated with 2.4%.

Methods: A series of adapter failures of titanium alloy modular neck adapters in combination with a titanium alloy modular short hip stem was investigated. For patients having received this particular implant combination risk factors were identified which were associated with the occurrence of implant failure. A Kaplan-Meier survival-failure-analysis was conducted. The retrieved implants were analysed using microscopic and chemical methods. Modes of failure were simulated in biomechanical tests. Comparative tests included modular neck adapters made of titanium alloy and cobalt chrome alloy material.

Results: Retrieval examinations and biomechanical simulation revealed that primary micromotions initiated fretting within the modular tapered neck connection. A continuous abrasion and repassivation process with a subsequent cold welding at the titanium alloy modular interface. Surface layers of 10 - 30 µm titanium oxide were observed. Surface cracks caused by fretting or fretting corrosion finally lead to fatigue fracture of the titanium alloy modular neck adapters. Neck adapters made of cobalt chrome alloy show significantly reduced micromotions especially in case of contaminated cone connection. With a cobalt-chromium neck the micromotions can be reduced by a factor of 3 compared to the titanium neck. The incidence of fretting corrosion was also substantially lower with the cobalt-chromium neck configuration.

Conclusions: Failure of modular titanium alloy neck adapters can be initiated by surface micromotions due to surface contamination or highly loaded implant components. In the present study, the patients at risk were men with an average weight over 100 kg. Modular cobalt chrome neck adapters provide higher safety compared to titanium alloy material.

* Correspondence: thomas.grupp@aesculap.de

† Contributed equally

¹Aesculap AG Research & Development, Tuttlingen, Germany

Background

Total Hip Arthroplasty (THA) has become a successful clinical treatment to restore the function of the joint, with a positive impact on the patient's quality of life. Modular connections for hip prostheses have been used since the early 70ies for heads with different neck sizes or diameters. Later in the 90ies, modular neck adapters have been introduced [1,2] for intraoperative adjustment of collum-caput-diaphysis (CCD) angle and femoral anteversion to optimise offset and leg length, irrespective of the hip stem implanted. These solutions have been proven their relevance in total hip arthroplasty [3,4]. Several failures of modular connections in hip replacement including primary and revision stem components [5-8] were reported. Long-term experiences with modular heads made of cobalt-chromium alloy (CoCr29Mo6) in combination with the cone of the stem out of titanium alloy (TiAl6V4) have revealed, apart from traces of fretting corrosion, no adverse events such as fractures of the tapers in clinical use [9-17]. Failures of modular neck adapters have been rarely documented [5-7]. In 2007, we reported three cases of a failure of a modular short hip stem [8]. The purpose of this paper is to present the state-of-the-art research and the recent findings of the failure analysis.

Case history

The failed neck adapters were implanted between August 2004 and November 2006. After the third incident, the titanium neck adapters were replaced by adapters out of cobalt-chromium [8]. A typical x-ray of a failed titanium neck adapter is shown in Figure 1.

The failure of the neck did not cause further damage on the acetabulum side. Subsequently, there was no need to revise the cups in all of these cases. Modular cup inserts were replaced in some cases due to visible signs of damage originated by the neck failure or for reasons of precaution. About 5000 hip joints have been implanted with the titanium stem and titanium neck adapter material combination. Until the end of 2008 1.4% ($n = 68$) of the titanium alloy neck adapters failed after an average time in vivo of 24 months (Figure 2). All, but one, patients were male and most patients (59%) had a weight above 100 kg and an average body mass index (BMI) of 31.6 (24 to 42) (Figure 3).

An overview about the average age, weight and post-operative time until failure of the titanium alloy neck adapter in 68 patients and the distribution of neck adapter geometry, head size and length is given in Table 1. The implant size did not have a detectable influence on the occurrence of neck failure. The neck adapter failures were divided into 66% with a CCD angle of 130°, 34% with a CCD angle 135° and none with 140°.

A Kaplan-Meier survival analysis (Statistica 7, StatSoft Europe GmbH, Hamburg, Germany) was undertaken


Figure 1 X-ray of a hip with a failed neck.

with revision due to failed neck adapter as an endpoint (confidence interval $\pm 95\%$). A cumulative survival rate of 98.6% was calculated with fractured neck adapter as the causative factor of failure (Figure 4).

There was no correlation between implant failures and a specific clinic or surgeon.

Methods

Implant component description

The modular neck adapters and stems used in our study were in all cases the Metha Short Hip Stem Prosthesis (Aesculap AG, Tuttlingen, Germany). Adapters and stems were made of titanium alloy. With its circumferential pure titanium porous coating and additive thin layer of dicalcium phosphate dehydrate (Plasmapore® μCaP , Aesculap, Tuttlingen, Germany), the stem is designed for cementless anchorage (Figure 5). The modular neck adapters were used in three different CCD-angles (130°, 135° and 140°) in combination with neutral, anteverted and retroverted ($\pm 7.5^\circ$) versions to adjust the CCD-angle intraoperatively [4].

Retrieval analysis

In the retrieval analysis, the failed neck adapters were subjected to a detailed examination to which most


Figure 2 Occurrence of clinical neck failures in a postoperative time period in months.


Figure 3 Correlation weight of patients and failures.

Table 1 Overview about the average age, weight and time in vivo until failure of the titanium alloy neck adapter in 68 patients and distribution of neck adapter and combined head geometry

age	Ø 57.2 years (36 - 75 years)					
weight	Ø 102.3 kg (75 to 130 kg)					
time in vivo	Ø 24 months (8 to 48 months)					
neck adapter (CCD-angle)	130°/0°	130° ± 7.5°	135°/0°	135° ± 7.5°	140°/0°	140° ± 7.5°
	29	16	18	5	-	-
head length	S	M	L	XL	XXL	unknown
	4	22	26	7	6	3
head diameter	28 mm		32 mm		36 mm	
	11		50		7	
head material	CoCrMo		BioloX forte		BioloX delta	
	13		53		2	


Figure 4 Kaplan-Meier survival analysis for failed neck adapter as reason of revision ($\pm$ 95% confidence interval).

patients consented. The analysis is in every case based on the failed neck adapter and the hip stem including the remaining distal part of the neck adapter. The revised modular heads and cup inserts were in most of the cases not available and not subject of this examination. The present report is based on the retrieval analysis of 47 devices, the complete investigation data being compiled at the end of 2008. The fracture and modular

connection surfaces were evaluated using light microscope and scanning electron microscope (SEM) (Zeiss EVO 50, Carl Zeiss NTS GmbH Oberkochen, Germany).

The fragment of the broken adapters which remained inside the shaft was cut through axially. One half of the adapters was cautiously removed from the shaft to examine the contact zone between the two components. The other half was left inside the shaft to enable the metallographic investigation of the interface and the microstructure of the materials. To enable the metallographic investigation of the neck adapter/stem interface the samples were prepared by mechanical cutting using a cut off-wheel (Secotom 10 Struers A/S, Denmark). After embedding into epoxy resin (Polyfast Struers A/S, Denmark) the cross section was prepared by grinding, polishing and etching according to the method by Kroll (composition 100 ml distilled H_2O , 1 ml HF , 3 ml HNO_3). The investigation was performed by light microscope (Wilde M3Z Herrenbrugg, Switzerland) with a magnification up to 500-times.

The chemical composition of the abrasive substances found in the crevice of the cone connection was determined using energy dispersive analysis (EDX) (Oxford EDX-Analysensystem INCA 7059 Oberkochen, Germany) and wet chemical digestion. The surface of the cone was rinsed cautiously using dilute acid. After hot digestion and dilution the analysis of the elements was performed using inductively coupled plasma optical emission spectrometry (ICP-OES Horiba Ultima II Jobin Yvon Longjumeau, France).

Micromotion analysis of the modular neck stem interface
To quantify the relative motion in the interface of the modular hip stem out of titanium alloy a total of eight


Figure 5 Metha Short Hip Stem System.


Figure 6 Test setup for measurement of micromotions in the modular cone connection (left) and particle-contaminated joining area (right) [18,19].

specimens with neck adapters out of titanium alloy and cobalt-chromium alloy were examined. As it appears that intraoperative contamination of the cone connection with bone particles has an considerable impact on the magnitude of fretting in the interface due to increased micromotions, a specific test setup was developed at the Biomechanics Section of the Hamburg University of Technology (Prof. Michael M. Morlock) and each device being tested with a clean and with a particle contaminated joining area (small porcine bone grafts approximately 1 mm in diameter) [18,19]. To determine the relative motion between neck adapter and stem a contactless measurement system (Micro-Epsilon Type U05(78) Messtechnik Ortenburg, Germany) was used. Two sensors with a range of 500 μm and a sensitivity of 0.025 μm were fixed on an aluminium plate and mounted on the neck adapter. Two steel plates fixed on the resection plane of the stem served as targets for the sensors (Figure 6).

The stem was combined with neck adapters with 130° CCD angle, embedded in bone cement (Palacos R, Heraeus Medical GmbH Wehrheim, Germany) and tested on a servohydraulic testing machine (MTS 850.2, MTS Systems Corporation Eden Prairie MN, USA). A sinusoidal axial force between 50 and 2500 N was applied via a ceramic head with neck length L at a frequency of 1 Hz for 2000 cycles to measure the relative displacement between neck adapter and stem with regard to irreversible settling and micromotions. A statistical analysis was performed to distinguish between independent groups (clean and particle contaminated joining area) (paired Student's t test) for both neck adapter materials (SPSS 15.0).

Pre-clinical fatigue testing

To determine the endurance properties of the neck region comparative tests were performed according to ASTM F 2068-03 and ISO 7206-6:1992(E) (MTS Mini Bionix II, MTS Systems GmbH Berlin, Germany). The hip contact force was set at 5340 N with a sinusoidal loading mode (ratio $F_{\min}/F_{\max} = 0.1$). In a saline

medium (0.9%, pH 7) the number of cycles was set at 10 million load alterations at a frequency of 15 Hz. In a worst case scenario using a 32 mm diameter XL head, neck fatigue was tested on the smallest modular stem in combination with neck components with 130° CCD angle and 7.5° retrotorsion out of titanium and cobalt-chromium alloy, respectively. A paired Student's t test was used to differentiate the fatigue behaviour of the two neck adapter materials (Statistica 7, StatSoft Europe GmbH, Hamburg, Germany).

To evaluate the coaction of the stem/neck modularity with the hip stem, the modular stems were tilted both in the frontal plane ($\alpha = 10^\circ$) and in the sagittal plane ($\beta = 9^\circ$). Based on CT scans on several human femurs and the specific stem design the level of embedding was set to 45 mm below the centre of the head to simulate a deficiency of proximal support due to bone loss in the neck and trochanter region (Figure 7).

For these customised test, the hip contact force was set at 2300 N for 5 million cycles according ISO 7206-8:1995(E) followed by a stepwise load increase (Locati method 500 N, 1 million cycles) until failure. These customised test setup was introduced to analyse if the mode of fatigue failure occurs in the modularity or in the stem region outside of the neck/stem connection.

To simulate the clinical failure modes in vitro, parameters for the titanium neck adapters were varied as shown in Table 2.

The resultant hip contact force was set at 3800 N to simulate a higher strain as in overweight or active patients. The test frequency was decreased from 15 Hz to 1 Hz to detect a possible influence on fretting corrosion in the modular neck interface. To simulate the patient's situation in daily life more realistically frequent stress phases alternated with rest phases. Taking into account the clinical conditions of hip replacement assembly, the connection was contaminated with blood, serum and cortical bone particles (approximately 1 mm in diameter) to provoke a mechanical disturbance in the interface. Additionally, the surrounding chemical


Figure 7 Test setup for neck test (left) and stem test (right) with reference electrode to measure corrosion potential.

conditions were altered by modifying pH values and enhancing chloride concentrations with the use of NaCl, CaCl₂ and FeCl₃ to accelerate any eventual tribochemical processes. To determine the continuous abrasion and repassivation of the titanium alloy surfaces in the neck/stem connection the redox potential was measured (InLab® 301 Mettler Toledo Balingen, Germany).

Results

Findings of the retrieval examination in clinical failures

All retrieved modular adapters showed similar breakage cross-sections. The fatigue fracture starts in the antero-lateral area at the upper part of the conical connection, the area subjected to maximum biomechanical stress. The fatigue fracture is followed by the residual forced fracture (Figure 8).

Metallographic analyses showed that microcracks developed on the surface of the cone in the clamping range (Figure 9). It appears that these microcracks induced the fatigue crack which finally led to the implant failure. Microcracks were analyzed using a scanning electron microscope. Figure 9 shows a potential micro crack in an area where fretting marks can be seen on the surface of the cone.

The examination did not indicate any product deviation, manufacturing failure or batch correlation, thus excluding any manufacturing influence as cause for implant breakage. In the majority of cases, the surface of the cones was significantly modified and showed signs of fretting corrosion. Among the forty-seven of the revised cones available for examination, a total of 87%

presented substantial fretting overlaid by laminar corrosion of different intensity (Figure 10).

The findings of the corrosive alteration apart of fretting have not been corroborated by biomechanical tests so far [8,14]. Metallographic cross sections of the cone show a 20 µm up to 30 µm thick layer developing in the modular joining area of most of the revised implants (Figure 11), a layer composed of different abrasive and corrosive substances including a high amount of titanium dioxide. The layer was formed due to abrasion of the passive oxide layer protecting the titanium alloy against corrosion. Continued abrasion and repeated repassivation of the titanium alloy surface produced the oxide layer's characteristic thickness.

Additionally in some cases this process accelerated through crevice corrosion as the reduction of the pH to values of about 2-3 indicates.

In one patient weighing 100 kg, a revision had to be carried out after 26 months to improve the offset due to cup loosening, thus providing the possibility to examine a faultless titanium adapter. (Figure 12).

The surface of the cone adapter which was not cleaned after revision was examined using a scanning electron

Table 2 In vitro model to simulate the clinical failure modes with various parameters

Hip contact force	3800 N
Frequency	1 Hz/15 Hz
Stress and rest phases	1,000 cycles, 1.5 minutes rest, 1,000 cycles, 1.5 minutes rest repeat
Contamination	Dry, Blood, Serum, Cortical bone particles
Surrounding medium	Saline 0.9% pH 7 Saline 1.8% pH 2 CaCl ₂ solution pH 2 FeCl ₃ solution pH 0.5


Figure 8 Fatigue fracture surface of a clinically failed titanium neck adapter.


Figure 9 Metallographic analyses revealed microcracks on the cone surface (left) and a potential microcrack in a fretting area (right).

microscope. No calcium or phosphor could be detected. There was no indication of a contamination of the cone connection by bone particles. Besides mild abrasion on the surface, no signs of fretting corrosion could be found.

Influence of micromotions on the modular neck stem interface

The major part of the observed displacement was caused by elastic deformation of the neck adapter. The micromotions of the titanium neck adapter were in a clean condition $6 \pm 4 \mu\text{m}$ (medial) and $6 \pm 2 \mu\text{m}$ (lateral) and with particle contamination $18 \pm 3 \mu\text{m}$ (medial) and $15 \pm 4 \mu\text{m}$ (lateral). In each contaminated joining area, the displacement and the micromotions were larger than in a clean condition. The fully uncontaminated components showed a distinct and definitive settling within the first 20 cycles. In case of a contaminated joining area the settling process of the neck adapter did not come to an end and lasted over the entire test duration (2000 cycles) (Figure 13).

The titanium alloy neck adapter showed a significant increase of micromotions in the interface of the particle-contaminated joining ($p < 0.01$) in comparison to clean conditions (Figure 14). The micromotions of the

neck adapter out of cobalt-chromium were in a clean condition $3 \pm 1 \mu\text{m}$ (medial) and $3 \pm 1 \mu\text{m}$ (lateral) and with particle contamination $5 \pm 6 \mu\text{m}$ (medial) and $3 \pm 2 \mu\text{m}$ (lateral).

In the case of the cobalt-chromium alloy neck adapter, the contamination of the joining area did not influence the micromotions significantly ($p = 0.43$).

In vitro fatigue behaviour of the stem/neck modularity

The neck adapters made of titanium and of cobalt-chromium alloy showed significantly different ($p = 0.009$) endurance behaviour. The Ti6Al4V adapters combined with a 28 L Biolog forte head failed after 2.45 million cycles (range 0.189 to 4.43). No failure occurred with the cobalt-chromium adapters (32 XL Biolog option) up to 20 million cycles.

Simulating a deficiency in proximal stem support the fatigue behaviour of the coaction of the stem/neck modularity with the hip stem was examined in a customised test setup. The titanium alloy neck adapters fulfilled the ISO 7206-8: 1995(E) requirements (2300 N for 5 million cycles) and failed during a stepwise increase of the applied force (Locati method) at a high load level of 5300 N after 0.7 million cycles (range 0.63 to 0.75). In a




surface morphology	partial large flaking	partial light flaking	no change
			
amount of the cones	23 %	30 %	13 %
time to failure	10 – 36 month	8 – 39 month	16 – 38 month
averaged time	24 month	22 month	27 month

Figure 10 Overview of interface characteristics and time to failure of the retrievals.


Figure 11 Brittle layer in cone connection between stem and neck adapter.

direct comparison, the cobalt-chromium neck adapters did not fail in the Locati test until the stem fractured at 6800 N load after 0.3 million cycles (Figure 15). Furthermore, the cobalt-chromium neck components did not fracture at 5300 N during 10 million cycles.

Contamination with cortical bone particles in the interface was found to be the main influence factor on the fatigue behaviour of the titanium alloy neck adapter. With a hip contact force of 3800 N the contaminated titanium alloy neck adapters failed after 3 million cycles (range 2.7 to 3.3) while no failure occurred under clean interface conditions (run-out limit 10 million cycles). This is mainly due to the increase of micromotions and the resulting fretting corrosion in the modular neck interface. In clean conditions, the load has to be increased up to 5300 N to provoke a similar failure mode. The redox potential measures fretting corrosion and repassivation on metallic surfaces. It shows impressively that introducing cortical particles into the modular neck interface generates a

mechanical disturbance which triggers repassivation processes in high frequency (Figure 16).

Discussion

This study confirms that the failure of the adapter cannot be attributed to any material or processing deviation or incorrect dimensioning. A twisted assembling of the adapter can also be excluded.

Due to mechanical perturbation in the modular connection, micro-movements caused fretting on the surface. Microcracks developed in the fretting zone, ultimately leading to the dynamic fatigue fracture of the implant. In 87% of the cases, fretting was accompanied by corrosion [17,20-25] as it initiates crevice corrosion [26,27]. The combination of fretting and crevice corrosion destroys the passive layer permanently. The electrochemical reactions provoke a shift in the pH value into the acidic range. The corrosion process generates a fissured surface of the connecting cone and can also generate microcracks. The examinations revealed that corrosion does not trigger the implant failure but does accelerate the process [26]. Cone adapters free from corrosive attack last the longest period up to failure (Figure 10).

In general, the titanium alloy provides excellent corrosion resistance and high biocompatibility because it quickly develops a thin passive layer with a thickness under 1 μm. Under inadequate conditions of cone assembling and in case of micro-movements within the modular connection, the protective passive layer is abraded initiating a continued abrasion and repassivation process which depletes the oxygen inside the crevice [17]. The metallic surface in the crevice becomes anodic relative to the outer surface, thus triggering an electrolytic process. The presence of chloride ions which disturb the passivity, induces pitting corrosion in the crevice leading to a quick dissolution of the metal



Figure 12 Retrieved faultless titanium alloy neck adapter (left) and stem (right) without contamination and any signs of corrosion.

accompanied by an acidification of the electrolyte [12,20,21,28].

Inside the layer, the elements calcium and phosphor could be detected and verified by different analytical methods. The crystalline composition of the calcium phosphate compound seems to indicate that the wear debris contained bone particles. Impurities which get inside the cone during intraoperative assembly may contribute to a mechanical perturbation which manifests itself as micro-movements which initiates fretting.

The study suggests the following hypothesis as to the cause of the damage. Fretting occurs when two surfaces in contact experience small amplitude oscillary relative

motion; damage is induced on the fretting region. If the fretting fatigue strength of the material is exceeded, microcracks develop on the surface. In addition, tribochemically activated particles can discharge their content from the surface. These particles react with oxygen spontaneously, thus leading to fretting corrosion.

The severity of corrosion depends mainly from the frequency of the fretting action. For lower cyclic rates such as in a hip implant, the debris is usually oxidized if the environment is chemically active. The repeated removal of oxide films due to continued abrasion and repassivation produces thick oxide layer's. Such oxides are normally harder than the virgin titanium alloy


Figure 14 Micromotions in a clean and particle-contaminated interface (Mean $\pm$ STD) for titanium and cobalt-chromium alloy neck adapters (medial and lateral).


Figure 15 Endurance behaviour of the modular stem and different failure mechanisms for neck components made of titanium (failure of the neck adapter) and of cobalt-chromium alloy (fractured stem at a level below the embedding).

leading to greater surface damage [29-31]. These damages then accelerate the crack nucleation [32,33].

Fretting reduces the fatigue strength of the titanium significantly [34]. To a great extent, this can be avoided by using of cobalt-chromium. The criteria listed in Table 3 show the comparison of the different materials which have been verified by biomechanical investigation. At a load of 5300 N the adapter made of the cobalt-based alloy did not break. At this high load level, the titanium adapter fractured after 2.45 million cycles.

The surface damage of the titanium alloy adapters caused by the microcracks or by corrosive deterioration accelerates the propagation of cracks by the cyclic loads bringing about the dynamic fatigue failure of the adapters. Micro-movements cause fretting in the cone connection. They can be increased by contamination of the cone connection through tissue or other particles intraoperatively. To anticipate this process any contamination of the connection should be avoided and the components dried before assembling. For this purpose


Figure 16 Free redox potential with frequent repassivation processes in the particle-contaminated modular neck interface (red curve) compared with a clean joining (blue curve).

Table 3 Assessment of the cone adapter made of the different materials - (ooo = excellent oo = good o = moderate)

	TiAl6V4	CoCr29Mo6
mechanical properties		
fatigue strength	ooo	ooo
stiffness/modulus of elasticity	oo	ooo
notch sensitivity	o	ooo
crack propagation	o	oo
abrasion	o	ooo
corrosion characteristics		
passive layer	ooo	ooo
re-passivation	ooo	oo
fretting corrosion	o	ooo
crevice corrosion	ooo	oo
allergic potential	ooo	o

abrasion-resistant cleaning rods are supplied together with the implants.

Conclusions

The change of the material of the adapter from titanium alloy to a cobalt-based alloy (CoCr29Mo6) increases the safety of the cone connection significantly. The combination of the cobalt-based alloy and the titanium alloy of the shaft shows a considerably higher rigidity. The smaller micro-movements reduce abrasion. Furthermore, the highly stable passive layer of the cobalt-based alloy provides an improved resistance against fretting. Due to its structure, the cobalt alloy has a much lower notch sensitivity compared to the titanium alloy. This enhances fatigue strength.

Among patients treated with the titanium alloy neck adapter, a combination of different parameters was identified as risk factors of implant failure. The parameters are intraoperative particle contamination of the cone connection, excessive loading due to a patient weight above 100 kg or high activity level, and male gender. In addition, the risk for failure rises with CCD angles of the cone adapter of 135° and smaller.

Acknowledgements

The authors would like to thank Elisa Hoenig, M.Sc. and Michael M. Morlock, Ph.D. for the examination of micromotions in the neck adapter/stem interface, Thomas Hermle, M.Sc., MBA for the performance of the survival-failure analysis and Christoph Schilling, M.Sc. for the statistical analysis of the in vitro tests.

Author details

¹Aesculap AG Research Development, Tuttlingen, Germany. ²Aesculap AG Clinical Science, Tuttlingen, Germany.

Authors' contributions

The authors TMG, TW, WB, HPK have made substantial contributions to conception and design, or acquisition of data, or analysis and interpretation

of data; 2) TMG, TW, WB, HPK have been involved in drafting the manuscript or revising it critically for important intellectual content; and 3) TMG, TW, WB, HPK have given final approval of the version to be published.

Competing interests

The authors TMG, TW, WB, HPK are all employees of Aesculap AG Tuttlingen, Germany

Received: 27 February 2009

Accepted: 4 January 2010 Published: 4 January 2010

References

1. Toni A, Paderni S, Sudanese A, Guerra E, Traina F, Antonietti B, Giunti A: Anatomic cementless total hip arthroplasty with ceramic bearings and modular necks: 3 to 5 years follow-up. *Hip International* 2001, **11**(1):1-17.
2. Traina F, Baleani M, Viceconti M, Toni A: Modular Neck Primary Prosthesis: Experimental and Clinical Outcomes. *Scientific Exhibit at the 71st AAOS Annual Meeting, San Francisco* 2004.
3. Toni A, Sudanese A, Paderni S, Guerra E, Bianchi G, Antonietti B, Giunti A: Artroprotesi d'anca non cementata con collo modulare - Cementless hip arthroplasty with a modular neck. *Chir Organi Mov* 2001, **LXXXVI**:73-85.
4. Braun A, Lazovic D, Zigan R: Modular short stem prosthesis- Implant positioning and the impact of navigation. *Orthopedics Suppl* 2007, **30**(Suppl 10):144-147.
5. Keggi KJ: My experience with proximal modular stems. *Symposium "Cutting-edge developments on proximal modularity in THA" Annual AAHKS Meeting* 2008, 16-22.
6. Keggi KJ, Kennon RE, Keggi JM: Hip reconstruction - primary total hip arthroplasty with modular stems. *Current Orthopaedic Practice* 2008, **19**(2):124-126.
7. Cameron HU, Leslie CJ, McTighe T: Design considerations and results for a modular neck in cemented THA. *Abstract, Australian Orthopaedic Association 67th Annual Scientific Meeting* 2007.
8. Hermle T, Zeller R, Grupp TM, Blömer W: Metha short stem hip prosthesis - Examination of the modular cone connection. *Springer Medizin Verlag, Heidelberg, implant* 2007, **1/07**: 12-17.
9. Bobyn JD, Tanzer M, Krygier JJ, Uujovne AR, Brooks CE: Concerns with modularity in total hip arthroplasty. *Clin Orth Rel Res* 1994, **298**: 27-38.
10. Brown SA, Flemming CA, Kawalec JS, Placko HE, Vassaux C, Merritt K, Payer JH, Kraay MJ: Fretting corrosion accelerates crevice corrosion of modular hip tapers. *J Appl Biomater* 1995, **6**(1):19-26.
11. Lieberman JR, Rimnac CM, Garvin KL, Klein RW, Salvati EA: An analysis of the head-neck taper interface in retrieved hip prostheses. *Clin Orth Rel Res* 1994, **300**: 162-167.
12. Gilbert JL, Buckley CA, Jacobs JJ: In vivo corrosion of modular hip prosthesis components in mixed and similar metal combinations. *The effect of crevice, stress, motion, and alloy coupling. Journal of Biomedical Materials Research* 1993, **27**:1533-1544.
13. Goldberg JR, Gilbert JL, Jacobs JJ, Bauer TW, Paprosky W, Leurgans S: A multicenter retrieval study of the taper interfaces of modular hip prostheses. *Clin Orth Rel Res* 2002, **401**: 149-161.
14. Windler M: Korrosionsverhalten von modularen Verbindungen bei Hüftprothesen. *Dissertation Eidgenössisch Technische Hochschule Zürich No. 15284* 2003.
15. Collier JP, Surprenant VA, Jensen RE, Mayor MB: Corrosion at the interface of cobalt-alloy heads on titanium-alloy stems. *Clin Orth Rel Res* 1991, **271**: 305-312.
16. McKellop H, Sarmiento A, Brien W, Park SH: Interface corrosion of a modular head total hip prosthesis. *J Arthroplasty* 1992, **7**(3):291-294.
17. Jacobs JJ, Gilbert JL, Urban RM: Current concepts review - Corrosion of metal orthopaedic implants. *J Bone Joint Surg [Am]* 1998, **80-A**:268-282.
18. Hoenig E, Morlock MM: "Bestimmung der Relativbewegung an der Konusverbindung der modularen Metha-Kurzschafthoprothese". *Study to determine the relative motion between the neck adapter and the stem Biomechanics Section - Hamburg University of Technology, Aesculap test report No. V1039* 2008.
19. Hoenig E, Morlock MM: Estimation of the relative motion in the taper lock interface of a modular hip prosthesis. *Abstract Orthopaedic Research Society, Las Vegas, Nevada* 2009.

20. Kaesche H: Die Korrosion der Metalle - Physikalisch-chemische Prinzipien und aktuelle Probleme. Springer-Verlag Berlin Heidelberg New York 1979, 274-276.
21. Gilbert JL, Jacobs JJ: The mechanical and electrochemical processes associated with taper fretting crevice corrosion - a review. *Modularity of orthopedic implants, ASTM STP 1301* J ASTM International Parr JE, Mayor MB 1997, 45-59.
22. Urban RM, Gilbert JL, Jacobs JJ: Corrosion of modular titanium alloy stems in cementless hip replacement. *J ASTM International* 2005, **2**(10):1-10.
23. Rodrigues DC, Urban RM, Jacobs JJ, Gilbert JL: In vivo severe corrosion and hydrogen embrittlement of retrieved modular body titanium alloy hip-implants. *J Biomed Mater Res Part B Appl Biomater* 2009, **88B**(1):206-219.
24. Collier JP, Surprenant VA, Jensen RE, Mayor MB, Surprenant HP: Corrosion between the components of modular femoral hip prostheses. *J Bone Joint Surg [Br]* 1992, **74-B**:511-517.
25. Cook SD, Barrack RL, Baffes GC, Clemow AJT, Serekian P, Dong N, Kester MA: Wear and corrosion of modular interfaces in total hip replacements. *Clin Orth Rel Res* 1994, **298**: 80-88.
26. Goldberg JR, Buckley CA, Jacobs JJ, Gilbert JL: Corrosion testing of modular hip implants. *Modularity of orthopedic implants, ASTM STP 1301* J ASTM International Parr JE, Mayor MB 1997, 157-176.
27. Goldberg JR, Gilbert JL: In vitro corrosion testing of modular hip tapers. *J Biomed Mater Res Part B: Appl Biomater* 2003, **64B**:78-93.
28. Duisabeau L, Combrade P, Forest B: Environmental effect on fretting of metallic materials for orthopaedic implants. *Wear* 2004, **256**:805-816.
29. Hoepfner DW, Gates FL: Fretting Fatigue Consideration in Engineering Design. *Wear*. 1981, **70**:155-164.
30. Viceconti M, Baleani M, Squarizoni S, Toni A: Fretting wear in modular neck hip prosthesis. *Journal of Biomedical Materials Research* 1997, **35**:207-216.
31. Viceconti M, Ruggeri O, Toni A, Giunti A: Design-related fretting wear in modular neck hip prosthesis. *Journal of Biomedical Materials Research* 1996, **30**:181-186.
32. Hoepfner DW, Chandrasekaran V: Fretting In orthopaedic implants: a review. *Wear*. 1994, **173**:189-197.
33. Martinez SA, Sathish S, Mall S, Blodgett MP: Evolution of fretting fatigue damage and relaxation of residual stress in shot-peened Ti6Al4V. *Metallurgical and Materials Trans A* 2005, **36A**:3385-3396.
34. Forschungskuratorium Maschinenbau EV: Reibkorrosion - Vorhaben Nr. 6. *Maschinenbau-Verlag GmbH, 6000 Frankfurt/M, Forschungsheft* 1976, **56**: 110-113.

Pre-publication history

The pre-publication history for this paper can be accessed here: <http://www.biomedcentral.com/1471-2474/11/3/prepub>

doi:10.1186/1471-2474-11-3


Cite this article as: Grupp et al.: Modular titanium alloy neck adapter failures in hip replacement - failure mode analysis and influence of implant material. *BMC Musculoskeletal Disorders* 2010 **11**:3.

Submit your next manuscript to BioMed Central and take full advantage of:

- Convenient online submission
- Thorough peer review
- No space constraints or color figure charges
- Immediate publication on acceptance
- Inclusion in PubMed, CAS, Scopus and Google Scholar
- Research which is freely available for redistribution

Submit your manuscript at
www.biomedcentral.com/submit

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The numbers in shaded italics are the comparison figures for the same type of implant. In most cases this is the numbers for all implants of the same type received by the NJRR

Number of implanting hospitals: 34

Number of hospitals where revisions occurred: 10

Reason for Revision	N	%
Loosening/Lysis	5	21.7
Dislocation	5	21.7
Pain	4	17.4
Infection	3	13.0
Fracture (Bone)	2	8.7
Leg Length Discrepancy	1	4.3
Metal Sensitivity	1	4.3
Wear Acetabulum	1	4.3
Implant Brakage Acetabular	1	4.3
Total	23	100%

Type of Revision	N	%
Femoral Only	6	26.1
Acetabular Only	5	21.7
Femoral and Acetabular	2	8.7
Cement Spacer	1	4.4
Head/Insert	5	21.7
Head/Neck	2	8.7
Head Only	1	4.3
Insert Only	1	4.3
Total	23	100%

TGA Observations on NJRR Data

The revisions seem to be evenly distributed among using hospitals, which would imply that if there is a problem it is not technique or surgeon related. The cumulative revision rate appears to get steadily worse and diverging from the cumulative revision rate of all other similar implants. Pain, leg leg discrepancy, metal sensitivity, wear and implant breakage are all over-represented as reasons for revision in this set.

TGA Observations on Manufacturers Reply

In this case the sponsor approach has been to explain, on a case by case basis why the implants have been revised, and then to provide published literature that suggests that this femoral stem prosthesis has had excellent results elsewhere in the world.

The 6 femoral only revisions and the 5 acetabular only revisions are said not to be due to "the failure of the device", but rather due to the "learning stage of the surgeons". However no supporting such as the manufacturer's own investigation records etc are provided to back up the statement, and the NJRR report contains no information about the surgeons who put these implants in or when.

One of the total revisions is blamed on the "high level of chromium that the patient had" requiring a change from a metal on metal bearing to a ceramic on ceramic bearing – the sponsor considers this not to be implant related. The second total revision is said to have happened because an inexperienced surgeon chose to use too small an acetabulum, the revision was required to rectify this. (Note: this explanation was provided for one of the two femoral + acetabular revisions.

The sponsor provides an explanation for all all of the revisions in the "Type of Revision table in terms of surgical error, but no supporting evidence for the assertions (eg a cross reference to investigation reports, implant retrieval analysis etc is provided) is provided. However, the sponsor states that he has received no complaints about the Adapter stem and 23 adverse event reports. Presumably this means that the sponsor has received an adverse event report from a surgeon that matches every revision cited by the NJRR. The sponsor states that there have been no complaints or adverse events reported elsewhere in the world.

As with other sponsors, ESKA have supplied a table summarising the available literature and copies of key publications.

Closer inspection of the paper by Sielewicz et al - reveals that the studies described in that paper were not without their complications (eg loosening, intraoperative stem fracture, pain). Also the relationship between the stem used by Sielewicz in 1983-88 and the Adapter stem is not clear other than the fact that both have a spongy surface.

Gollwitzer et al reports very good results with an ESKA hip system - presumably Adapter (76 hips, mean follow up 7.9 years, 1 femoral component revision and 2 acetabular component revision and good hip scores).

Götze et al report a cumulative survival rate on $90 \pm 8\%$ at for the acetabular and $86 \pm 5\%$ for the stems at 14.9 years. 137 patients were followed in this study and the authors state that four prosthesis (a relatively large number) were revised due to implant fracture. These are not great results. The authors conclude that

"The long term results of the spongy metal cup are good, whereas the high loosening and fracture rate of fully coated stem are a source of concern especially with regard to the difficult revision scenario with frequent massive bone loss."

Götze et al state that the stem has since been modified to proximal coating only but it is not clear how similar the Adapter stem is to the stem used in this study.

Sholtz et al reviewed 165 consecutive ESKA hip arthroplasties conducted between 1983 and 1985. 32 patients had died of unrelated causes, 14 hips had been revised. The minimum follow up was 12 years the mean follow up time was 13 years. Of the 14 that were revised: 1 was due to infection – 11 (femoral only) were due to loosening – 1 was due to stem fracture and one (acetabulum only) due to loosening. The survival rate was calculated to be 88% at 15 years with good hip scores.

Powerpoint slides of the Scholtz data set were also provided.