



Device Incident Report : Medical Devices Branch - Device Vigilance and Monitoring

DIR : 41 - ID : 510535

16/10/2020

UNSIGNED

Released by [REDACTED] on 19/11/2020 14:40:19

Report #:	Records Management #:	Reporter's Reference #:	Report Type:
65968	E20-373384	RM4949	Final

ARTG: [Document Container URL](#)

Report Information Section

Report Status:	Sponsor's Reported Category:	Date of Adverse Event:	Date of In tial Report:
Active	Trend	25/09/2020	16/10/2020
Date of Final Report:	Date of In tial TGA Action:	Reviewed by Team:	Date Response Received:
30/06/2021	16/10/2020	27/10/2020	
Date Completed:	Operator at Time of Event:	If 'Other' Operator Selected:	Reporter consents to contact by sponsor:
			N/A
Source of Report:	If 'Other' Source Selected:	Type of In tial Act on:	
Sponsor		For Team Meeting	

Event Description for Website Publ cat on:

Clin cal Event Informat on:

The outcome following a surveillance aud t on Repromed SA Genet cs laboratory's cell free (non-invasive) PGT-A program that has triggered a detailed internal investigation. Whilst our investigation has not been completed, irregular ties have been found in the Val dation document that supported the implementation of this test and the accred tat on granted by NATA. Consequently, Repromed have decided that the test will cease to be offered clin cally to its patients, or any patients it performed the test for in the broader Monash IVF Group's network of clin cs. NATA have been informed accordingly, and a repeat of the val dation study performed as a matter of urgency. Monash IVF Group will work through a process of Open disclosure w th our clin cians and patients wh ch will include remediation where found to be appropriate. I have also included below a sample of the sort email that has been sent to patients, this one is for patients w th aneuploidy only embryos.

Background
Monash IVF's NI PGT program was launched following an almost four-year long process of research and Val dation trials. The program was approved by Monash IVF GMAC and had human eth cal trials approval from an accredited HREC. The research and validat on, led by [REDACTED], focussed on whether the genet c results from NI PGT systems were as accurate as those from gold standard PGT-A with trophectoderm b opsy.

The val dation data was obtained on a cohort of 121 embryos, and showed the two methods were very similar in their outcomes of detecting aneuplo dy. There was a correlation of 98% between invasive PGT and NI-PGT; consequently the test was subm tted for NATA accred tation. The NATA accred tation audit with peer review from a techn cal expert and review of val dation report occurred in early 2019. Monash IVF Group launched the NI PGT-A program across New South Wales, Northern Territory, Queensland, South Australia, Tasmania and V ctoria in May 2019. The test has also been performed on behalf of [REDACTED] in [REDACTED].

Findings
As part of our routine surveillance program, a review of the NI-PGT outcomes was undertaken in June 2020. This review cons dered not only how accurately the two alternatives performed in terms of detecting abnormality, but also in clinical pregnancy rates. In add tion, the post launch surveillance sample size was considerably larger (n = 805) than the val dation studies (n = 121). This review highlighted some variat ons of the performance of the test in clinical use, as compared to the validat on results:

- Failed DNA amplif cation rates were pleasingly lower than the validat on study (2.6% Surveillance vs 5.0% Val dation).
- There was an increase in inconclusive rates compared to val dation study (6.3% Surveillance vs 1.6% Validat on); however both outcomes were within parameters experienced in routine clinical practice.
- The significant unexpected finding, was that there was a signif cant increase in aneuplo dy rates in the NI PGT tested embryos when compared to the current invasive PGT tested embryos. This was evident across all ages and for non-delayed and delayed embryos. The increase in aneuploidy was 20-30% higher in the NI PGT group.
- This implies a higher false positive rate compared w th PGT-A w th b opsy, indicating that more embryos may have been called abnormal when in fact they may be normal, compared w th b opsy PGT-A.

The discordant results between Biopsy PGT and NI PGT prompted a full interrogation of the validation study data files to try and better understand these unexpected outcomes. This review revealed some discrepancies in the val dation data that are not yet fully understood, but bring into quest on its scientific and clinical val d ty.

Next Steps
Monash IVF Group are working with clinicians to prior tise patient communications and support programs, as we notify patients that the test is being suspended and the impact that our findings have on any stored aneuplo dy embryos. A further val dation study is underway, under the supervision of a multidisciplinary Steering Committee, to assess the possibil ty of continuing to prov de the test. NATA have been advised of our decision to suspend the test and we will work with them in the event that we are in a position to repeat the accred tation process in the coming months. If you have any quest ons regarding any of the informat on or would like further information, please do not hes tate to call me.

Number of Inc dents in Report:	Contact:	Alternative Person Title:	Alternative Person First Name:
1			
Alternative Person Surname:	Alternative Person Phone:	Alternative Person Fax:	Alternative Person Email:

Recorded Problems Observed

Recorded Problems Observed:

Output Problem -> Incorrect, Inadequate or Imprecise Result or Readings -> False Positive Result

Clin cal Signs, Symptoms and Condit ons

Recorded Clinical Signs, Symptoms and Conditions:

Others -> Insufficient Information ->

Health Impact

Recorded Health Impacts:

Delay to Treatment/ Therapy -> ->

Patient Information

Sex:

Weight:

Age:

Patient Focused Corrective Action Taken:

Patient History:

Patient Outcome/Consequences:

Additional Event Description:

Describe any test (Lab, xray, etc.):

Injured - Extent of Injury:

Other medical devices currently using/implanted:

Medical Problem Device Used For:

Additional Patients Added:

0

Submitting Reporter Section

Search Reporter By Surname:

Reporter #:

Preferred Contact Method:

Reporter Title:

First Name:

Surname:

Position:

Company/Institution:

Address 1:

Address 2:

Monash IVF Group

Country:

Postcode:

Town/Suburb:

State:

Australia

3121

Richmond

VIC

Mobile:

Email:

Phone:

Fax:

Last External Submission By:

104241_65620 - 16/10/2020 17:03

Initial Reporter Section

As Above?:

If No, fill out the following:

Initial Reporter Confidential:

No

Yes

Search Reporter By Surname:

Initial Reporter #:

Preferred Contact Method:

Title:

First Name:

Surname:

Position:

Company/Institution:

Address 1:

Address 2:

79 631 193 489

Postcode:

Country:

Town/Suburb:

State:

Mobile:

Email:

Phone:

Fax:

Allow the device company to contact you about the incident:

☐

Device Information Section

Product Exempt (Note: If not exempt, enter ARTG No):

No

Product Licence Category:

Brand Name:

In-house IVD - cell free PGT-A

Model #:

Purchase Date:

Date of Initial Procedure:

Access Contact First Name:

Access Contact Email:

Search Device ARTG:

DV-2017-IVI-10175-1

Device Class:

Initial Device Description:

In-house IVD - cell free PGT-A

Serial #:

Expiry Date:

Place of Implantation:

Access Contact Surname:

Licence Status:

Device ARTG #:

GMDN / UMDN Code:

Usage of Device:

Batch #:

Date of Implant:

Reported Device Location:

Access Contact Phone:

Status Effective Date:

Therapeutic Licence Type:

GMDN / UMDN Text:

Software Version:

Lot #:

Date of Explant:

Access Contact Title:

Access Contact Fax:

Additional Devices Added:

0

Manufacturer Information Section

Manufacturer Name:

Monash IVF Group

Address 2:

Postcode:

Manufacturer Informed:

Yes

Contact Surname:

Town/Suburb:

Phone:

Date Aware of Adverse Event:

25/09/2020

Manufacturer Client Id:

State/Province:

Fax:

Contact Title:

Address 1:

Country:

Email:

Contact First Name:

Supplier Information Section

Supplier Name:

Town/Suburb:

Phone:

Supplier Informed:

Contact Surname:

State:

Fax:

Date of Supplier Contact:

Contact Phone:

Address 1:

Country:

Email:

Contact Title:

Contact Fax:

Address 2:

Postcode:

Website:

Contact First Name:

Contact Email:

Report Status

For website publication:

Ready for Publication:

No

Investigated:

Yes

Investigation Reason:

More information is required

Team Assignment:

Team C (IIa, I and IVD 1-4)

Team Priority:

Routine

Team Review

Reviewed by Team: 27/10/2020	Reason Sent To Meeting: Advice from IVD section given on requirement of more information	Outcome from team meeting: Investigation (within DIR) is recommended
---------------------------------	---	---

Notes for Team meeting:

Exempt in-house IVD - still required to comply with post-market reporting. Request investigation report - possibly request updates (depending on expected time-frame of final report). As sponsor had not thought it necessary to notify TGA of the incident, they may be unaware of their requirement to submit a final report. Contact sponsor and advise them they are required to submit a final report within 100 days.

Outcomes from Team Meeting:

Investigation to enquire: Rate information, Investigation report. What standards is this IVD done under? Provide evidence that testing was conducted according to these standards. To clarify with SW the role and scope of DPMM in regulating this IVD.

DPRC Review

Reviewed by DPRC:	DPRC Reason Sent To Meeting:	Outcome from DPRC Meeting:
-------------------	------------------------------	----------------------------

Meeting Notes:

Initial Risk Analysis

Background Information	Risk Assessment - Section A	Risk Assessment - Section B	Risk Assessment - Section C	Risk Assessment - Section D
Date: 16/10/2020	Severity: 1 - No harm has been reported to occur	Incidents in the last 12 months: 	Manufacturer analysis: Yes	Assessor: [REDACTED]
Incidents in last 24 months: 	Manufacturer act on: Yes	ESTIMATED LEVEL OF INVESTIGATION: Level 3 Investigation (for multiple DIRs/a DII)	FINAL LEVEL OF INVESTIGATION: 	Injured Party: Patient
Incidents in last 36 months: 	IVD status: A point of care IVD	EXCEPTION TO INVESTIGATION LEVEL: 	Found Prior To Use: No	Is AE covered by current recall: No
Incidents Worldwide: 	Number of potential contributing factors: Yes - some potential factors (up to 3)		Reusable: Yes	Similar events (past 6 months): 0 incidents
Products supplied the last 12 months: 	Specific factors identified: Compatibility of device - patient characteristics, Use of device - fit for purpose	ESTIMATED LEVEL OF PRIORITY: Routine	FINAL LEVEL OF PRIORITY: 	3 or more events - batch/model:
Products supplied last 24 months: 	Number of potential sensitivities: Yes - multiple potential sensitivities (more than 3)	EXCEPTION TO PRIORITY LEVEL: 		3 or more events - health district:
Products supplied last 36 months: 	Specific sensitivities identified: Device uses new technology, New evidence about device efficacy, Public interest/requests for information, Media interest/requests for information, Industry interest/requests for information, Investigations/requests by government, Enquiries/investigations/actions by the TGA, Potential for current incident to become sensitive			3 or more events - organisation:
Products supplied Worldwide: 	Consultations during risk assessment: I discussed issues with one of my peers	Final Risk Assessment: No		

Additional Risk Analysis

Click 'N' to start a new risk analysis


Analysis Details	Statistics Checklist Section				
Update Device Details?:	Background Information	Risk Assessment - Section A	Risk Assessment - Section B	Risk Assessment - Section C	Risk Assessment - Section D
	Date:	Severity:	Incidents in the last 12 months:	Manufacturer analysis:	

Copy Data From:							
Assessor:	Manufacturer documentation:	Incidents in last 24 months:	Manufacturer action:	ESTIMATED LEVEL OF INVESTIGATION:	FINAL LEVEL OF INVESTIGATION:		
Injured Party:	Device Recalls:	Incidents in last 36 months:	IVD status:	EXCEPTION TO INVESTIGATION LEVEL:			
Found Prior To Use:	Is AE covered by current recall:	Incidents Worldwide:	Number of potential contributing factors:				
Reusable:	Similar events (past 6 months):	Products supplied the last 12 months:	Specific factors identified:	ESTIMATED LEVEL OF PRIORITY:	FINAL LEVEL OF PRIORITY:		
	3 or more events - batch/model:	Products supplied last 24 months:	Number of potential sensitivities:	EXCEPTION TO PRIORITY LEVEL:			
	3 or more events - health district:	Products supplied last 36 months:	Specific sensitivities identified:				
	3 or more events - organisation:	Products supplied Worldwide:	Consultations during risk assessment:	Final Risk Assessment:			

Sponsor/Manufacturer Information Section

Search Sponsors:	Name:	Client #:
65620	Monash IVF Group	65620
Attention To:	Address 1:	Town/Suburb:
	Pelaco Building 1 Level 1 / 21-31 Goodwood Street	Richmond
State:	Postcode:	Fax:
VIC	3121	
Email:		

Investigation Information Section - Submitted by Sponsor/Manufacturer

Device Analysis Results:	Corrective/Preventative Actions:
	Ongoing
Details of Similar Events:	Additional Details (use for tables):
CAPA# Reference:	
Risk Assessment	
	
Frequency:	Severity:
Rating:	Type Cause and Outcome:
Expected Rate:	Number of Similar Events:
Actual Rate:	
Countries Similar Events Also Occurred:	
Completed Actions:	Planned Actions and Proposed Timelines:

Suspension of the test and working through communications with affected patients	Proposed re-validation of test under extended committee Appointment of external reviewer to completed RCA relating to validation
Additional Comments:	
<div></div>	

Reason for Level 1 Investigation

Details of Reasons	
Reason for Level 1 Investigation	
Sensitivity - requests for information	

Focus of Level 2 Investigation

Details of Focus	
Essential Principles	If 'Other' Selected

Sources of Evidence for Level 2

Details of Source			
Sources of Evidence	If 'Others' please specify here	Expected Sourcing Date	Date of Evidence Received

Evidence

Investigation Questions (Level 1 and Level 2):

Request for information made under Schedule 3 Part 6A.3 & 6A.4 of the Regulations - refer to TRIM D20-3634097

Potential Risks

Delays in response by product manufacturers:

☒

Delays in response by incident reporters:

☐

Delays in analysis within the TGA:

☒

Delays in reporting by other sources (e.g. clinical registries):

☐

Other Risks (which need to be specified):

Next Steps for Level 1 & Level 2 Investigations

Next Steps for Level 1 Investigation:

Next Steps for Level 2 Investigation:

Click [N] to begin a new Correspondence entry. Note that the Email address specified here will receive a notification if the Date Received is not filled in by the Date Expected.

Correspondence and Chronology Details										
Include?	Heading	Type L1	Type L2	Email	Sent	Expected	Received	Response	Notes	
<input type="checkbox"/>		Sponsor Routine Correspondence	Sponsor Request for Information Letter - s41JA		30/10/2020	26/11/2020	27/11/2020			

List of Problem Observed Codes - Click [N] to begin entering information.

Problem Observed Details				
Problem Observed (Level 1)	Problem Observed (Level 2)	Problem Observed (Level 3)	If 'Other' Selected	

Output Problem	Incorrect, Inadequate or Imprecise Result or Readings	False Positive Result	
----------------	---	-----------------------	--

Clinical signs symptoms and conditions

Details		
Level 1	Level 2	Level 3
Others	Insufficient Information	

Health Impact

Details		
Level 1	Level 2	Level 3
Delay to Treatment/ Therapy		

Investigation Findings

Finding Details			
Investigation Findings (Level 1)	Investigation Findings (Level 2)	Investigation Findings (Level 3)	If 'Other' Selected

Investigation Conclusion

Conclusion Details		
Investigation Conclusion (L1)	Investigation Conclusion (L2)	If Additional Conclusion Detail Requested

Investigation Outcomes

Outcome Details		
Outcome of Investigation (L1)	Outcome of Investigation (L2)	If Additional Conclusion Detail Requested

Investigation Summary

Latest Investigation (DII) where this DIR is the Primary DIR:	Latest Investigation (DII) where this DIR is a Related DIR:	Investigator:	Peer Review:	Peer Review Completed on Date:
			Yes	
Investigator's Notes:		Summary Findings:		Recall Number:
<p>NOTE - For Peer Review prior to closing After reviewing requested information - change to "Awaiting Final", send to OPR - Make a note to triage final report to [REDACTED]</p> <p>30 October 2020 - Request for information sent (Notice requiring information/documents to be provided under Schedule 3, Part 6A, Clause 6.3 and 6.4 of the Therapeutic Goods (Medical Devices) Regulations 2002) (D20-3641206) 27 November 2020 - Monash IVF responded (D20-3822434)</p> <p>Review of documents provided: 1. TGA Response to NI Questions v1.0 27Nov20 (Q1) - Provided a general description of the "cell-free" as compared "biopsy" pre-implantation genetic test (PGT) and the indication of results - (a) suitable for transfer (Euploid or "no abnormality detected"); (b) chromosomally abnormal & unsuitable for transfer (Aneuploid); (c) Inconclusive (unable to yield a result) - embryo treated in same way as an untested embryo. MIVF state cell-free PGT-A is a screening test only design to detect abnormal chromosome number (aneuploidy) for whole chromosome number which may lead to miscarriage or serious condition affecting pregnancy and live birth. WARNING to Patients:</p>				

Patients are advised: "This test is only a screening test and therefore cannot provide an absolute guarantee of the chromosome status of the embryo. In some embryos, the biopsied cell/s or culture media may not be representative of the whole embryo.

While every effort is made to ensure that the PGT-A test offered has the highest possible accuracy using the currently available technology, results are not 100% accurate. Therefore, prenatal diagnosis is highly recommended in an ensuing pregnancy."

(Q2) - Provided an updated procedure for reporting adverse event (their reporting procedure, at the time of the adverse event, was not provided) - Reviewed in section 2. Monash IVF Group Adverse Event Feedback Policy v5.0 9Nov20
- Their rationale for the 4 month delay in reporting their concerns regarding the significant increase in aneuploidy rates when they first became aware of the issue, in June 2020 as part of their routine surveillance program, and only notified TGA in October 2020 was:

(a) initial concern were raised in June 2020 which triggered a review of the validation data. This review revealed irregularities which led to the decision to suspend the test as of 25 September.

(b) The precautionary suspension of the test did not appear, to Monash IVF, to meet the Adverse Event Criteria listed on the TGA website.

(c) Confirmation that Monash IVF Group Adverse Events and Feedback Policy to the version provided (see (2) below) following discussion with [REDACTED] (TGA) to include reference and link to the NPAAC document: Requirements for the Development And Use Of In-House In Vitro Diagnostic Medical Devices (IVDs).

(Q3) Although Monash IVF identified increased rates of aneuploidy results in June 2020 between the surveillance population (n=805) compared to the validation population, they only suspended the test after their investigation identified irregularities in the validation data in September 2020 (when number of tests conducted was in the order of n=1300). Therefore between when they noticed concern and when they suspended use of the test, around an additional 500 tests had been conducted. Monash IVF claimed that the increased rates observed in June 2020 were NOT the grounds for suspending the test. The grounds for suspending the test was related to the IRREGULARITIES IN THE VALIDATION DATA which were identified in September 2020.

(Q4) Monash was asked to provide information on their documented procedures for reviewing experience gained in the post-product on phase and how corrective actions are applied to the design or production of such devices. Monash provided information for their genetics laboratory at Repromed advising they participate in ALL external QC programs available for their preimplantation genetic testing program but that no such external QC programs are available for the cell-free PGT.

An internal QC program was developed for the cell-free PGT.

Monash IVF claim that "results from our internal QC program have all been concordant since it was established", but did not provide any details or data.

Repromed also monitor false negative rates for all screening tests and that the false negative rate appeared to be in-line with the biopsy PGT test (false negative - calling an embryo healthy when it was not). The issue with the cell-free PGT-A test was false positive (calling an embryo unhealthy/aneuploid when it was healthy/euploid)

Monash IVF didn't actually answer the question asked. No details of the system for reviewing experience gained in the post-product on phase was provided. It was expected that copies of their QMS documents would be provided that demonstrated how they monitored and reviewed post-production data.

No information was provided on how corrective actions will be applied to the design or production of a device (IVD test), only that "learnings from our experience with this process will be incorporated into future design of any tests implemented across Monash IVF Group".

Could request additional information - please provide a copy of your documented procedures for reviewing experience gained in the post-product on phase and the process for how any necessary corrective actions will be applied to the design or production of a device.

Could request a copy of their internal QC program for the test.

(Q5) Referred to (3) below: IVD Validation Document - Next Generation Sequencing v1.0 for details of the design specification of the device, including relevant standards, risk analysis or other solutions adopted to ensure the device complies with the EPs.

(Q6) Referred to (4) below: Monash IVF Group validation of NGS summary for clinical evidence as required by the clinical evaluation procedures (Schedule 3, Part 8 of the Regs)

(Q7) Referred to (3 & 4) below: request for validation report.

2. Monash IVF Group Adverse Event Feedback Policy v5.0 9Nov20

Provides definitions of

Events - something that requires reporting in Riskman - incidents, near misses, feedback, suggestions for improvement, supplier issues, equipment failure, audit findings and hazards

3. IVD Validation Document - Next Generation Sequencing v1.0

Review for (Q5 above): details of the design specification of the device, including relevant standards, risk analysis or other solutions adopted to ensure the device complies with the EPs.

4. Monash IVF Group validation of NGS summary

Review for (Q6 above): clinical evidence as required by the clinical evaluation procedures (Schedule 3, Part 8 of the Regs)

Note: Letter generation buttons disabled if report not ready for website publication or risk analysis not completed.

Device Lookup

This section is used to match information provided via UDIR forms to ARTG informat on. You can select a Brand/Name from information provided in the 'Other Dev ces Involved' table below or enter informat on manually.

Other Dev ce (Entered):

Brand Name:

Manufacturer Name:

Dev ce ARTG #:

Other Devices

Dev ce ARTG No:

Manufacturer Name:

Sponsor/Supplier:

GMDN / UMDN Text:

Trade/Brand Name:

Serial #:

Model Number:

Batch #:

Lot #:

Expiry Date:

Related DIR Informat on - Click **New** to begin entering information.

Rec No

1

Samples Record - Cl ck **[N]** to begin entering informat on. **Note:** Sample # Generated on Save.

Rec No

1

Details

Sample Details

Additional Details

Date Entered:

LIMS #:

Sample Requested:

Sample Received:

Manufacturer:

GMDN:

Device Descript on:

Brand Name:

Serial Number:

Reason for Testing:

Samples from Reporter:

Samples from Sponsor:

Outcome of TGA's Testing:

Lot Number:

Batch Number:

Model Number:

Version Number:

Who sent the device to the TGA?:

Why does the TGA have the sample?:

Additional Patients

Click **[N]** to begin entering informat on.

Patient Details

Sex:

Weight:

Age:

Patient Focused Corrective Action Taken:

Patient History:

Injured - Extent of Injury:

Was dev ce directly linked to death?:

Was device directly linked to permanent disabilitiy?:

Consequence:

Other Consequence:

Describe any test (Lab, xray, etc.):

Additional Event Descript on:

Med cal Problem Dev ce Used For:

Additional Device Information

Where d d you get this device from?:

How reliant is the affected person on correct/safe operat on of this dev ce?:

Any other relevant informat on to a d assessing/investigating the inc dent?:

ileader.production.tga.gov.au/InformationLeaderAD/Forms/FormDetailPrint.aspx?sid=2008642495

9/10

Similar Events

Similar events - how many times?:

Date of Recent Report:

Event Reported To:

Reporter Reference Number:

Device Access - Alternate Device Contact Information Provided

Title:

First Name:

Last Name:

Phone:

Fax:

Email:

Incident Location Details

Occurred in Australia:

Organisation:

Address Line 1:

Address Line 2:

Town/Suburb:

State:

Postcode:

Flow Details DIR-REQ - Device Incident Request 283548

Request Details

ID	Type	Location	Status	Assigned By	Assigned To	Assigned On	Priority	Attachments
283548	DIR-REQ		Awaiting Final Report			23/09/2021	Normal	0

Signature Details

Role	IRIS Investigator	
User		
Signed At		
Comment		