



## Laboratories Branch

Operations Biotherapeutics Laboratories Operations Manual	
Procedure	Identification of the mRNA in modRNA BNT162b2 (1525) using RT-PCR Assay
Author	[REDACTED]
Owner	[REDACTED]
Authoriser	[REDACTED]
Revision #	2

TRIM link to SOP	D21-2063055	Date of assay	17/11/2021
TRIM link to PCR template	D21-2061428	TRIM link to assay data file	<a href="#">D21-3344781</a>
Operator	[REDACTED]	Checked by	[REDACTED]

**Reaction mixture preparation:** Prepare 15µL of reaction mixture per well plus 10% overage  
 Primers and probe are provided by Pfizer and are stored in aliquots at working concentration (x10)  
 Each assay requires a total of 12 wells for PCR and extraction controls plus 3 wells per test sample.

**Reporting:** Record the sample LIMS numbers in the table below prior to the assay. Once the assay and analysis are complete, record the assay validity criteria parameters on the worksheet below noting whether the assay is valid. For a valid assay, record the sample Ct values and whether identity has been confirmed in the table below.  
 Using the Quantstudio Design & Analysis Software, update the assay results file to include the LIMS numbers of each sample (replacing the placeholder letter designation), using the "Advanced Setup" pane of the "Plate" tab. Save the updated file in the assay specific data folder in TRIM (E21-219384).  
 Following completion of the assay, convert this worksheet into a .pdf and append to it the following: .pdf copies of amplification curve plots, a .pdf experiment report as generated by the Quantstudio Design & Analysis Software. Combine these into a single .pdf and file it in the assay specific data folder (E21-219384) along with the Quantstudio data file.

## Reagent information

Reagent	Manufacturer and Catalogue Number	Lot Number	Expiry	Notes
QIAamp Viral RNA Mini Kit	Qiagen, 52906	166026648	2022-03-15	
Buffer AVL	Qiagen, 52906	166025080	2022-03-15	
Buffers AW1 and AW2	Qiagen, 52906	166025919 & 166021395	N/A	
Ethanol 96-100%	Supelco	K50375083828	2023/07/31	
TE Buffer	Life Technologies, AM9849	01063124	N/A	
RT-PCR-grade water	Life Technologies, AM9935	2104121 & 2004106	N/A	
TaqPath 1-step RT-qPCR Master Mix, CG	Applied BioSystems, A15299	2293147	2022-01-30	
ModRNA1525 RT-PCR forward primer working stock		00711497-0162-M06	-	see last page for probe and primer details
ModRNA1525 RT-PCR reverse primer working stock		00711497-0162-M04	-	
ModRNA1525 RT-PCR probe working stock		00711497-0162-M02	-	

Record Details	Identification of the mRNA in modRNA BNT162b2 (1525) using RT-PCR Assay		
Last Editor	[REDACTED]	Edit Date	18/11/2021 11:02 AM
Print Date	18/11/2021 11:03 AM		Page 1 of 5

## Sample Preparation and Dilution

Positive PCR Control	Positive Extraction Control	Negative Extraction Control	Test Samples
<p>Dilution 1: 10µl DS RM (10µg/mL) + 90µL RT-PCR-grade water = 1µg/mL</p> <p>Dilution 2: 10µL Dilution 1 (1ug/mL) + 190µL RT-PCR-grade water = 50ng/mL = 50 pg/µL</p> <p>5µL Dilution 2 per RT-PCR well</p>	<p>Prepare with QIAamp Viral RNA Mini Kit. Lysis step composition: 140µL undiluted DP RM + 560 QIAamp Buffer AVL, without carrier RNA. Follow kit protocol Elute in 60µL Buffer</p> <p>Dilution 1: 10µL eluate + 990µL RT-PCR-grade water</p> <p>Dilution 2: 10µL Dilution 1+ 990µL RT-PCR-grade water</p> <p>5µL Dilution 2 per RT-PCR well</p>	<p>Prepare with QIAamp Viral RNA Mini Kit. Lysis step composition: 140µL TE buffer + 560 QIAamp Buffer AVL, without carrier RNA Follow kit protocol Elute in 60µL Buffer</p> <p>Dilution 1: 10µL eluate + 990µL RT-PCR-grade water</p> <p>Dilution 2: 10µL Dilution 1+ 990µL RT-PCR-grade water</p> <p>5µL Dilution 2 per RT-PCR well</p>	<p>Prepare with QIAamp Viral RNA Mini Kit. Lysis step composition: 140µL undiluted test sample + 560 QIAamp Buffer AVL, without carrier RNA Follow kit protocol Elute in 60µL Buffer</p> <p>Dilution 1: 10µL eluate + 990µL RT-PCR-grade water</p> <p>Dilution 2: 10µL Dilution 1+ 990µL RT-PCR-grade water</p> <p>5µL Dilution 2 per RT-PCR well</p>

## RT-PCR Reaction Mixture

Number of wells required	18	Volume of master mix per well	15µL	Total volume of master mix required (incl.10% overage)	283	Volume of sample per well	5µL
--------------------------	----	-------------------------------	------	--	-----	---------------------------	-----

Component	Volume per well	Total Volume (incl. 10% overage)
Mastermix	10 µL	<b>190</b>
Forward Primer	1 µL	<b>19</b>
Reverse Primer	1 µL	<b>19</b>
Probe	1 µL	<b>19</b>
RT-PCR grade water	2 µL	<b>36</b>
<b>Total volume of reaction mixture</b>		<b>283</b>

## Reaction mixture notes:

~~Mastermix used: TaqPath 1-step RT-qPCR Master Mix, CG~~  
4x TaqPath 1-step Master Mix, CG

Record Details

Identification of the mRNA in modRNA BNT162b2 (1525) using RT-PCR Assay

Last Editor

[REDACTED]

Edit Date

18/11/2021 11:02 AM

Print Date

18/11/2021 11:03 AM

Page 2 of 5

## PCR Plate layout:

	1	2	3	4	5	6	7	8	9	10	11	12
A	PCR Positive Control	PCR Positive Control	PCR Positive Control							Positive Extraction Control	Positive Extraction Control	Positive Extraction Control
B												
C	Sample A	Sample A	Sample A									
D	Sample B	Sample B	Sample B									
E	<del>Sample C</del>	<del>Sample C</del>	<del>Sample C</del>									
F	<del>Sample D</del>	<del>Sample D</del>	<del>Sample D</del>									
G												
H	No Template Control	No Template Control	No Template Control							Negative Extraction Control	Negative Extraction Control	Negative Extraction Control

Sample	Sample LIMS#:	Ct Value	Identity results ( <i>confirmed/not confirmed</i> ) (Ct must be < 32 for all sample replicates)
Sample A	2111004115	17.120, 16.998, 16.718	confirmed
Sample B	2111004218	16.447, 16.519, 16.427	confirmed
Sample C			
Sample D			

Record Details	Identification of the mRNA in modRNA BNT162b2 (1525) using RT-PCR Assay	Edit Date	18/11/2021 11:02 AM
Last Editor	██████████		
Print Date	18/11/2021 11:03 AM		Page 3 of 5

## Assay Validity Data

Criterion	Required value	Observed values	Validity
Ct of each Negative Extraction Control well	> 28 or undetermined	31.205, 32.264, 29.730	Valid
Ct of each No Template Control well	> 28 or undetermined	35.549, 35.872, 31.445	Valid
Ct of each PCR Positive Control well	< 28	13.277, 13.486, 13.483	Valid
Ct of each Positive Extraction Control well	< 28	18.392, 18.347, 18.568	Valid

## Sample interpretation

Sequence identity of the mRNA is considered confirmed if all test sample replicates show amplification curves with a Ct value of < 28.0000. If all test sample replicates show amplification curves with a Ct value of > 28.0000 the identity is considered not confirmed. If a mixture of results (with Ct values both greater and less than 28.0000) is found for a test sample it must be repeated.

Pipettes used: **Extraction step: 30006, 33190, 33166**  
**PCR step: 5646, 5653, 33087**

## Batch details for probes and primers:

For Primer: (Merck) Sigma, Ref# VC00021, SY21020241530-088, Lot# 3026595983-000020  
 Rev Primer: (Merck) Sigma, Ref# VC00021, SY21020241529-078, Lot# 3026595983-000030  
 100 µM stock prepared as per manufacturer's instructions, 20 µL of 100 µM added to 91 µL RT-PCR  
 Grade water to make 18 µM working stock, MJ 06Sep21  
 Probe: (Applied Biosystems) ThermoFisher, Lot# 7495294-1 C1, Ref# 4316034  
 100 µM stock prepared as per manufacturer's instructions; 5 µL 100 µM added to 95 µL RT-PCR Grade  
 Water to make 5 uM working stock, MJ 27Sep21

NOTE: Validity criteria for assay has been altered as per recommendation in D21-2294359.

Criteria used to assess validity of assay are as follows:

Ct of Negative Extraction Control (NEC) more than 8 Ct higher than Ct of Positive Extraction Control (PEC):

Ct of No Template Control (NTC) more than 8 Ct higher than Ct of PCR Positive Control:

Record Details

Identification of the mRNA in modRNA BNT162b2 (1525) using RT-PCR Assay

Last Editor

[REDACTED]

Edit Date


18/11/2021 11:02 AM

Print Date

18/11/2021 11:03 AM

Page 4 of 5

# Amplification Plot - BNT162b2 Identification assay - 17Nov21



■ 2111004115 ■ 2111004218 ■ NEC ■ NTC ■ PEC ■ Positive Control

# Experiment Results Report

## Experiment Summary

Experiment Name: 2021-01-11\_133749  
Experiment Type: Standard Curve  
Chemistry: TaqMan® Reagents  
BarCode:  
File Name: ID Pfizer\_21.eds  
Run Started: 11-17-2021 17:49:50 AEDT  
Run Finished: 11-17-2021 18:56:01 AEDT  
Run Duration: 66 minutes 10 seconds  
Date EDS Modified: 11-18-2021 10:43:02 AEDT  
Date EDS Created: 01-11-2021 15:21:26 AEDT  
User Name:  
Number of Wells Used: 18  
Number of Wells with Results: 18  
Instrument Name: QS3  
Instrument Type: QuantStudio™ 3 System  
Instrument Serial Number: 272322852  
Model/Block Type: QuantStudio™ 3 System / 96-Well 0.2-mL Block  
Quantification Cycle Setting: Ct  
Stage/Step for Analysis: Stage2, Step2  
Comments:



# Experiment Results Report

## Reagent Information

Type	Name	Part Number	Lot Number	Expiration Date
------	------	-------------	------------	-----------------



For Research Use Only, not for use in diagnostic procedures.

# Experiment Results Report

## Results Summary

Sample	Target	Qty Mean	Qty SD	Ct Mean	Ct SD
2111004115	Target 1			16.945	0.207
2111004218	Target 1			16.464	0.048
NEC	Target 1				
NTC	Target 1				
PEC	Target 1			18.436	0.117
Positive Control	Target 1			13.415	0.120





# Experiment Results Report

## Plate Layout

	1	2	3	4	5	6	7	8	9	10	11	12
A	Positive.. Target 1	Positive.. Target 1	Positive.. Target 1							PEC Target 1	PEC Target 1	PEC Target 1
B												
C	211100411 Target 1	211100411 Target 1	211100411 Target 1									
D	211100421 Target 1	211100421 Target 1	211100421 Target 1									
E												
F												
G												
H	NTC Target 1	NTC Target 1	<b>1</b> NTC Target 1					<b>2</b> NEC Target 1	<b>2</b> NEC Target 1		<b>2</b> NEC Target 1	



# Experiment Results Report

## Results Table

Well	Sample	Target	Task	Ct	Ct Mean	Ct SD	Qty	Qty Mean	Qty SD	Ct Threshold	Baseline Start	Baseline End	Cq Conf
A1	Positive Control	Target 1	U	13.277	13.415	0.120				0.076	3	9	0.938
A2	Positive Control	Target 1	U	13.486	13.415	0.120				0.076	3	9	0.943
A3	Positive Control	Target 1	U	13.483	13.415	0.120				0.076	3	9	0.957
H11	NEC	Target 1	N	29.730	31.066	1.273				0.076	3	25	0.974
A10	PEC	Target 1	U	18.392	18.436	0.117				0.076	3	14	0.952
A11	PEC	Target 1	U	18.347	18.436	0.117				0.076	3	14	0.949
A12	PEC	Target 1	U	18.568	18.436	0.117				0.076	3	14	0.947
C1	2111004115	Target 1	U	17.120	16.945	0.207				0.076	3	13	0.948
C2	2111004115	Target 1	U	16.998	16.945	0.207				0.076	3	13	0.949
C3	2111004115	Target 1	U	16.718	16.945	0.207				0.076	3	12	0.953
D1	2111004218	Target 1	U	16.447	16.464	0.048				0.076	3	12	0.947
D3	2111004218	Target 1	U	16.427	16.464	0.048				0.076	3	12	0.964
D2	2111004218	Target 1	U	16.519	16.464	0.048				0.076	3	12	0.946
H1	NTC	Target 1	N	35.549						0.076	3	31	0.924
H9	NEC	Target 1	N	32.264	31.066	1.273				0.076	3	28	0.982
H8	NEC	Target 1	N	31.205	31.066	1.273				0.076	3	26	0.967
H3	NTC	Target 1	N	31.445						0.076	3	27	0.965



# Experiment Results Report

Well	Sample	Target	Task	CT	CT Mean	CT SD	Qty	Qty Mean	Qty SD	CT Threshold	Baseline Start	Baseline End	Cq Conf
H2	NTC	Target 1	N	35.872						0.076	3	32	0.922

Legend: S = Standard, N = Non Template Control, U = Unknown, UND. = Undetermined



# Experiment Results Report

## QC Summary

Total Wells:96    Processed Wells:18    Targets Used:1    Well Setup:18    Flagged Wells:4    Samples Used:6

Flag	Description	Frequency	Locations
AMPNC	Amplification in negative control	4	H11, H9, H8, H3
BADROX	Bad passive reference signal	0	
BLFAIL	Baseline algorithm failed	0	
CQCONF	Low Cq confidence	0	
CTFAIL	Ct algorithm failed	0	
DRNMIN	Define acceptable delta Rn based on Ct range	0	
EXPFAIL	Exponential algorithm failed	0	
HIGHSD	High standard deviation in replicate group	3	H11, H9, H8
HIGHSD	High standard deviation in replicate group	0	
NOAMP	No amplification	0	
NOISE	Noise higher than others in plate	0	
NOSIGNAL	No signal in well	0	
OFFSCALE	Fluorescence is offscale	0	
OUTLIERRG	Outlier in replicate group	0	
PRFDROP	Passive reference signal changes near Ct	0	
PRFLOW	Low passive reference signal	0	
SPIKE	Noise spikes	0	
THOLDFAIL	Thresholding algorithm failed	0	

Customer is responsible for any validation of assays and compliance with regulatory requirements that pertain to their procedures and uses of reagents and/or instruments. ©2018 Thermo Fisher Scientific. All rights reserved. The trademarks mentioned herein are the property of Thermo Fisher Scientific or their respective owners.

