

### **3.2.P.3.5. VERIFICATION OF IN-PROCESS TEST METHODS [PUURS]**

Descriptions of the validation of the analytical procedures used for in-process control (IPT-C) are provided below.

#### **3.2.P.3.5.1. pH**

The pH is measured potentiometrically during preparation of the buffer using calibrated pH meters. This method meets compendial requirements.

#### **3.2.P.3.5.2. RNA Content by Fluorescence Assay**

A summary of the method validation for in-process testing is detailed in [Section 3.2.P.5.3 Fluorescence Assay](#).

#### **3.2.P.3.5.3. Bioburden**

The bioburden analytical procedure (performed following the principles described in USP <61>, Ph. Eur. 2.6.12, and JP 4.05, using membrane filtration) was verified for BNT162b2 drug product (DP) in-process samples. The method verification (challenge recovery test) challenges the test method to ensure that the test articles are non-inhibitory to the recovery of inoculated microorganisms. Verification of the pre-sterile filtration sample is detailed in this section.

The test article tested for the presence of bioburden is generally expected to be non-bacteriostatic/non-fungistatic or to be readily neutralized by membrane filtration. To demonstrate this, the bioburden assay was performed with the addition of a microbial inoculum (target titer of not more than 100 colony forming units (CFUs)) composed of compendial -recommended microorganisms and in-house microorganisms. Verification of the bioburden method consisted of challenge-recovery testing on 3 pre-sterile filtration test articles. The method verification consists of:

- Neutralizer efficacy; which is assessed by adding a low level of inoculum of a challenge organism to the test article and demonstrates that recovery is not inhibited by product residues on the filter.
- Neutralizer toxicity; which is assessed using a neutralizer control (without product) and demonstrates that recovery is not inhibited by the use of neutralization agent(s) and/or filter type.
- Negative controls (blanks); demonstrate that the materials used (e.g. media, filter, diluent) are free from contamination.
- Inoculum controls; demonstrate the micro-organisms being used are suitable for use.

Verification of the bioburden method consisted of challenge-recovery testing on 3 different test articles (material from bulk drug product lots EG5447, EH9979 and EK2538).

For the challenge recovery testing to be verified, the following acceptance criteria must be met:

- The inoculum control count of all test organisms on Tryptic Soy Agar (TSA) and Sabouraud Dextrose Agar + Chloramphenicol (SDAC) has to be between [REDACTED] CFU.
- The blank of the neutralizer control may not exhibit growth on TSA and SDAC.
- The blank of the product test on TSA and SDAC has to be [REDACTED] of the lowest inoculum control count of all test organisms on TSA and SDAC.
- Neutralizer efficacy microbial recovery from the inoculated product count has to be at least [REDACTED] compared to the inoculum control.
- Neutralizer toxicity microbial recovery from the neutralizer control count has to be at least [REDACTED] compared to the inoculum control.

#### 3.2.P.3.5.3.1. Bioburden Verification Data

For each verification run, 100 mL of the pre-sterile filtration sample was inoculated with not more than 100 CFUs of the challenge microorganism and filtered. This procedure was repeated for each challenge organism listed below:

- Bacterial test organisms: *Staphylococcus aureus* ATCC 6538, *Pseudomonas aeruginosa* ATCC 9027, *Bacillus subtilis* ATCC 6633, *Ralstonia picketti* (UP 71) in-house isolate common in water samples, *Stenotrophomonas maltophilia* (UP 07) in house isolate common in bioburden samples.
- Fungal test organisms: *Candida albicans* ATCC 10231, *Aspergillus brasiliensis* ATCC 16404.

The summary data that supports the verification of the pre-sterile filtration material is presented in the tables below.

**Table 3.2.P.3.5-1. Summary Data for Pre-sterile Filtration Bioburden Challenge Recovery: Trial 1 Using Lot EG5447**

Microorganism	Inoculum Control (CFU) (A)	Neutralizer Control Count (CFU) (B)	Inoculated Product Control Count (CFU) (C)	Neutralizer Toxicity (%) (B/A) × 100%	Neutralizer Efficacy (%) (C/A) × 100%	Acceptance Criteria
<b>Total Aerobic Microbial Count: Tryptic Soy Agar, 30-35 °C, ≤3 days</b>						
Blank	NA	0	0	NA	NA	B = 0 CFU C [redacted] of lowest inoculated product control count
<i>S. aureus</i>	[redacted]					A between [redacted] CFU
<i>R. pickettii</i>						(B/A) × 100% = [redacted] %
<i>S. maltophilia</i>						(C/A) × 100% = [redacted] %
<i>B. subtilis</i>						
<i>P. aeruginosa</i>						
<i>C. albicans</i>						
<i>A. brasiliensis</i>						
<b>Total Yeast and Mold Count: Sabouraud Dextrose Agar + Chloramphenicol, 20 -25 °C, ≤5 days</b>						
Blank	NA	0	0	NA	NA	B = 0 CFU C [redacted] of lowest inoculated product control count
<i>C. albicans</i>	[redacted]					A between [redacted] CFU
<i>A. brasiliensis</i>						(B/A) × 100% = [redacted] %
						(C/A) × 100% = [redacted] %

NA = not applicable

090177e196993238\Approved\Approved On: 26-Mar-2021 16:13 (GMT)

**Table 3.2.P.3.5-2. Summary Data for Pre-sterile Filtration Bioburden Challenge Recovery: Trial 2 Using Lot EH9979**

Microorganism	Inoculum Control (CFU) (A)	Neutralizer Control Count (CFU) (B)	Inoculated Product Control Count (CFU) (C)	Neutralizer Toxicity (%) (B/A) × 100%	Neutralizer Efficacy (%) (C/A) × 100%	Acceptance Criteria
<b>Total Aerobic Microbial Count: Tryptic Soy Agar, 30-35 °C, ≤3 days</b>						
Blank	NA	0	0	NA	NA	B = 0 CFU C [redacted] of lowest inoculated product control count
<i>S. aureus</i>	[redacted]					A between [redacted] CFU
<i>R. pickettii</i>						
<i>S. maltophilia</i>						(B/A) × 100% = [redacted] %
<i>B. subtilis</i>						
<i>P. aeruginosa</i>						(C/A) × 100% = [redacted] %
<i>C. albicans</i>						
<i>A. brasiliensis</i>						
<b>Total Yeast and Mold Count: Sabouraud Dextrose Agar + Chloramphenicol, 20 -25 °C, ≤5 days</b>						
Blank	NA	0	0	NA	NA	B = 0 CFU C [redacted] of lowest inoculated product control count
<i>C. albicans</i>	[redacted]					A between [redacted] CFU
<i>A. brasiliensis</i>						(B/A) × 100% = [redacted] %
						(C/A) × 100% = [redacted] %

NA = not applicable

090177e196993238\Approved\Approved On: 26-Mar-2021 16:13 (GMT)

**Table 3.2.P.3.5-3. Summary Data for Pre-sterile Filtration Bioburden Challenge Recovery: Trial 3 Using Lot EK2538**

Microorganism	Inoculum Control (CFU) (A)	Neutralizer Control Count (CFU) (B)	Inoculated Product Control Count (CFU) (C)	Neutralizer Toxicity (%) (B/A) × 100%	Neutralizer Efficacy (%) (C/A) × 100%	Acceptance Criteria
<b>Total Aerobic Microbial Count: Tryptic Soy Agar, 30-35 °C, ≤3 days</b>						
Blank	NA	0	0	NA	NA	B = 0 CFU C [redacted] of lowest inoculated product control count
<i>S. aureus</i>	[redacted]					A between [redacted] CFU
<i>R. pickettii</i>						
<i>S. maltophilia</i>						(B/A) × 100% = [redacted] %
<i>B. subtilis</i>						
<i>P. aeruginosa</i>						(C/A) × 100% = [redacted] %
<i>C. albicans</i>						
<i>A. brasiliensis</i>						
<b>Total Yeast and Mold Count: Sabouraud Dextrose Agar + Chloramphenicol, 20 -25 °C, ≤5 days</b>						
Blank	NA	0	0	NA	NA	B = [redacted] CFU C [redacted] of lowest inoculated product control count
<i>C. albicans</i>	[redacted]					A between [redacted] CFU
<i>A. brasiliensis</i>						(B/A) × 100% = [redacted] %
						(C/A) × 100% = [redacted] %

NA = not applicable

090177e196993238\Approved\Approved On: 26-Mar-2021 16:13 (GMT)

The results of the verification study (challenge recovery) demonstrate the bioburden method is effective in recovering microorganisms from the test article. The analytical procedure is verified for its intended use.