Device Incident Report : Medical Devices Branch - Device Vigilance and Monitoring

DIR : 44 - ID : 517161

Released by s22 on 18/03/2021 17:26:54

Report #:	Records Management #:	Reporter's Reference #:	Report Type:
68516	E21-295890	0172-21	Duplicate

ARTG: 320484

[Document Container URL](#)

Report Information Section

Enter the ARTG of the device or Duplicate DIR in the search box then press Tab to search for the DIR that this report duplicates. If an ARTG has already been entered below it will be entered in the search field for you.

Duplicate DIR No Search:	ARTG Search:	DIR #:
	320484	68783 - J-Lock PEEK Suture - Anchor, soft-tissue implanted, non-biodegradable

Report Status:	Sponsor's Reported Category:	Date of Adverse Event:	Date of Initial Report:
Closed	Death / Serious Injury	s22	05/03/2021
Date of Final Report:	Date of Initial TGA Action:	Reviewed by Team:	Date Response Received:
04/05/2021	05/03/2021	11/05/2021	
Date Completed:	Operator at Time of Event:	If 'Other' Operator Selected:	Reporter consents to contact by sponsor:
13/10/2021	Healthcare Professional		N/A
Source of Report:	If 'Other' Source Selected:	Type of Initial Action:	
Sponsor		For Team Meeting	

Event Description for Website Publication:

A part of the inserter handle broke and has remained in the bone with the anchor

Clinical Event Information:

The patient underwent a labral repair s22 Hip Arthroscopy. Two MectaLock (J-Lock) PEEK anchors were used to refixing the labrum to the acetabulum with no visual issues. Then the surgeon sent to the agent a picture of an X-ray asking if the image in the X-ray was part of the anchor inserter handle. It appears that a part of the inserter handle has remained in the bone with the anchor. The patient underwent revision surgery to remove the broken part.

The event was reported to the manufacturer on 29.01.2021, but only after the revision surgery, which took place on s22, we were sure that the piece that remained in the patient was part of the Medacta device.

Number of Incidents in Report:	Contact:	Alternative Person Title:	Alternative Person First Name:
1			
Alternative Person Surname:	Alternative Person Phone:	Alternative Person Fax:	Alternative Person Email:

Recorded Problems Observed

Recorded Problems Observed:

Material Integrity Problem -> Break -> Fracture

Clinical Signs, Symptoms and Conditions

Recorded Clinical Signs, Symptoms and Conditions:

Injury -> Foreign Body In Patient -> Device Embedded In Tissue or Plaque

Health Impact

Recorded Health Impacts:

Surgical Intervention -> Device Revision or Replacement ->

Patient Information

Sex:	Weight:	Age:
s22		s22
Patient Focused Corrective Action Taken:	Patient History:	

Patient Outcome/Consequences:

Describe any test (Lab, xray, etc.):

Injured - Extent of Injury:

Additional Patients Added:

0

Submitting Reporter Section

Search Reporter By Surname:

Reporter Title:

Position:

Address 1:

Country:

Mobile:

Reporter #:

First Name:

Address 2:

Postcode:

Email:

Additional Event Description:

Other medical devices currently using/implanted:

Medical Problem Device Used For:

Surname:

Company/Institution:

Town/Suburb:

Phone:

Last External Submission By:

Preferred Contact Method:

State:

Fax:

Initial Reporter Section

As Above?:

Search Reporter By Surname:

Title:

Position:

Address 1:

Postcode:

Mobile:

If No, fill out the following:

Initial Reporter #:

First Name:

Address 2:

Country:

Email:

Surname:

Company/Institution:

Town/Suburb:

Phone:

Allow the device company to contact you about the incident:

☐

Initial Reporter Confidential:

Preferred Contact Method:

State:

Fax:

Device Information Section

Product Exempt (Note: If not exempt, enter ARTG No):

Product Licence Category:

Brand Name:

Model #:

Purchase Date:

Search Device ARTG:

Device Class:

Initial Device Description:

Serial #:

Expiry Date:

Device ARTG #:

GMDN / UMDN Code:

Usage of Device:

Batch #:

Date of Implant:

Therapeutic Licence Type:

GMDN / UMDN Text:

Software Version:

Lot #:

Date of Explant:

		§22	§22 Document 1
Date of Initial Procedure:	Place of Implantation:	Reported Device Location:	Access Contact Title:
	§22	With Reporter	
Access Contact First Name:	Access Contact Surname:	Access Contact Phone:	Access Contact Fax:
Access Contact Email:	Licence Status:	Status Effective Date:	Additional Devices Added:
	A	19/07/2019	0

Manufacturer Information Section

Manufacturer Name:	Manufacturer Client Id:		Address 1:
Medacta International SA	45296		
Address 2:	Town/Suburb:	State/Province:	Country:
Postcode:	Phone:	Fax:	Email:
Manufacturer Informed:	Date Aware of Adverse Event:	Contact Title:	Contact First Name:
Yes	29/01/2021	Mr	§22
Contact Surname:	§22		

Supplier Information Section

Supplier Name:	Address 1:		Address 2:
Town/Suburb:	State:	Country:	Postcode:
Phone:	Fax:	Email:	Website:
Supplier Informed:	Date of Supplier Contact:	Contact Title:	Contact First Name:
Contact Surname:	Contact Phone:	Contact Fax:	Contact Email:

Report Status

For website publication:	Ready for Publication:	Investigated:	Investigation Reason:	Team Assignment:	Team Priority:
No	Yes	Yes	More information is required	Team B (IIB, Is & Im)	Routine

Team Review

Reviewed by Team:	Reason Sent To Meeting:	Outcome from team meeting:	
11/05/2021	Not enough information has been provided	Investigation (within DIR) is recommended	
Notes for Team meeting:			
1 other DIR, no recall/pmr, same device same issue currently under investigation. problem is with the disposable introducer in the kit, which has welding failure. it looks like the same issue with the other DIR 68783. But device in this DIR is from different lot, different implantation date. sending to meeting to investigate together.			
Outcomes from Team Meeting:			
lv 3 investigation - investigate together with DIR 68783.			

Initial Risk Analysis

Background Information	Risk Assessment - Section A	Risk Assessment - Section B	Risk Assessment - Section C	Risk Assessment - Section D
Date:	Severity:	Incidents in the last 12 months:	Manufacturer analysis:	Assessor:
			Manufacturer documentation:	

04/03/2021	5 - An illness/injury was resolved or prevented with treatment by a health professional			s22	Unknown - updated information from the manufacturer is required
Incidents in last 24 months:	Manufacturer action:	ESTIMATED LEVEL OF INVESTIGATION:	FINAL LEVEL OF INVESTIGATION:	Injured Party:	Device Recalls:
		Level 1 Investigation (to complete screening)		Patient	0. No recalls for similar incidents in Australia
Incidents in last 36 months:	IVD status:	EXCEPTION TO INVESTIGATION LEVEL:		Found Prior To Use:	Is AE covered by current recall:
				No	No
Incidents Worldwide:	Number of potential contributing factors:			Reusable:	Similar events (past 6 months):
	No			No	0 incidents
Products supplied the last 12 months:	Specific factors identified:	ESTIMATED LEVEL OF PRIORITY:	FINAL LEVEL OF PRIORITY:		3 or more events - batch/model:
		Routine			
Products supplied last 24 months:	Number of potential sensitivities:	EXCEPTION TO PRIORITY LEVEL:			3 or more events - health district:
	No				
Products supplied last 36 months:	Specific sensitivities identified:				3 or more events - organisation:
Products supplied Worldwide:	Consultations during risk assessment:	Final Risk Assessment:			
	I did none of the above incidents	No			

Additional Risk Analysis

Click 'N' to start a new risk analysis

Analysis Details	Statistics Checklist Section				
Update Device Details?:	Background Information	Risk Assessment - Section A	Risk Assessment - Section B	Risk Assessment - Section C	Risk Assessment - Section D
Yes	Date:	Severity:	Incidents in the last 12 months:	Manufacturer analysis:	
Copy Data From:	07/05/2021	5 - An illness/injury was resolved or prevented with treatment by a health professional		Yes	
Initial					
	Assessor:	Manufacturer documentation:	Incidents in last 24 months:	Manufacturer action:	ESTIMATED LEVEL OF INVESTIGATION:
	s22	Unknown - updated information from the manufacturer is required		No	Level 1 Investigation (to complete screening)
	Injured Party:	Device Recalls:	Incidents in last 36 months:	IVD status:	FINAL LEVEL OF INVESTIGATION:
	Patient	0. No recalls for similar incidents in Australia			Level 3 Investigation (for multiple DIRs/a DII)
	Found Prior To Use:	Is AE covered by current recall:	Incidents Worldwide:	Number of potential contributing factors:	EXCEPTION TO INVESTIGATION LEVEL:
	No	No		Yes - some potential factors (up to 3)	1 other DIR. no recall/pmr. same device same issue currently under investigation. problem is with the disposable introducer in the kit, which has welding failure. it looks like the same issue with the other DIR 68783. But device in this DIR is from different lot, different implantation date. sending to meeting to investigate together.
	Reusable:	Similar events (past 6 months):	Products supplied the last 12 months:	Specific factors identified:	ESTIMATED LEVEL OF PRIORITY:
	No	1 incidents		Production of device - manufacture/batch	Routine
		3 or more events - batch/model:	Products supplied last 24 months:	Number of potential sensitivities:	FINAL LEVEL OF PRIORITY:
				No	Routine
		3 or more events - health district:	Products supplied last 36 months:	Specific sensitivities identified:	EXCEPTION TO PRIORITY LEVEL:
		3 or more events - organisation:	Products supplied Worldwide:	Consultations during risk assessment:	Final Risk Assessment:
				I discussed issues with one of my peers	Yes

Sponsor/Manufacturer Information Section

Search Sponsors:	Name:	Client #:
45972	Medacta Australia Pty Ltd	45972

Attention To: s22	Address 1: Unit A1 / 16 Mars Road	Address 2: 	Town/Suburb: LANE COVE	Document 1
State: NSW	Postcode: 2066	Phone: 	Fax: 	
Email: s22@medacta.com.au				
Investigation Information Section - Submitted by Sponsor/Manufacturer				
Device Analysis Results: Batch review performed on 23-Feb-2021: Lot 1924147: 125 items manufactured and released on 24-Mar-2020. Expiration date: 2025-Feb-07. No anomalies found related to the problem. To date, 19 items of the same lot have been already sold without any other similar reported even. Clinical Evaluation performed by Medical Affairs Director: During an arthroscopic surgery of labrum refixation with a PEEK anchor, the tip of the positioning instrument broke, though the surgeon did not perceive the fracture. The tip of the instrument is visible in the postoperative X-rays, it remained in the cannulation of the anchor. We cannot determine the reasons for the fracture with the elements at hand. No clinical cause could be ascertained. We were informed that a different surgeon intervened to remove the broken instrument part; the future clinical pathway of the patient should be totally flawless and no further consequences of the event should be expected. Visual Inspection performed by R&D Project Manager: The component of the MectaLock PEEK inserter has been received and analyzed. It's evident that the laser welding failed and the tip component has been detached from the main shaft of the inserter. The proximal portion of the analyzed component resulted totally ruined and scratched, most likely due to the extraction, therefore is not possible to perform a detailed microscopic analysis of the failure mode of the welding. However, it's also evident that the component has been subjected to excessive lateral bending as the axis of the component results deviated from the original. This suggests that misuse could have occurred during the anchor's placement or that a not foreseeable condition of use required an excessive force, as the MectaLock PEEK anchors is supposed to be impacted axially only within the predrilled cylindrical pilot hole. Conclusion: The device history record review does not indicate any potentially related manufacturing issue. The component has been subjected to lateral bending as the axis of the component results deviated from the original. This suggests that suboptimal usage occurred during the anchor's placement or that a not foreseeable condition of use required an excessive force, as the MectaLock PEEK anchors are supposed to be impacted axially only within the predrilled cylindrical pilot hole.		Corrective/Preventative Actions: The event is an isolated case and no systematic action is currently indicated as applicable, however, it has been included in our active monitoring system for potential CAPA needs and opportunities.		
Details of Similar Events: 1. MectaLock instruments 05.10.001 MectaLock PEEK Ø2.4 – Short driver Lot. 1924147: this is the 1st event for the family, the first for the reference and lot. This is the first case reported. From the beginning till today, 03-Mar-2021 we have sold 201 pieces considering short and long drivers. The occurrence rate is 1/201=0.5%		Additional Details (use for tables): 		
CAPA# Reference: 				
Risk Assessment				
Frequency: 	Severity: 			
Rating: 		Type Cause and Outcome: 	Number of Similar Events: 	
Expected Rate: 	Actual Rate: 			
Countries Similar Events Also Occurred: 				
Completed Actions: 		Planned Actions and Proposed Timelines: 		
Additional Comments: DIR was closed on 13/10/2021.				

Reason for Level 1 Investigation

Details of Reasons

Reason for Level 1 Investigation	Document 1

Focus of Level 2 Investigation

Details of Focus	
Essential Principles	If 'Other' Selected
EP1-Use may compromise patient health and safety	
EP2-Design/construction may not conform with principles	

Sources of Evidence for Level 2

Details of Source			
Sources of Evidence	If 'Others' please specify here	Expected Sourcing Date	Date of Evidence Received
Information from Sponsors (e.g. 41JA, Questionnaires, Emails)		30/08/2021	28/08/2021

Evidence

Investigation Questions (Level 1 and Level 2):

1.Please provide a detailed manufacturer’s investigation analysis to confirm that it adequately demonstrates the conformity with EP2 design and construction principles.

2.Please send a copy of Instructions for Use (IFU).

Potential Risks

Delays in response by product manufacturers:

☒

Delays in response by incident reporters:

☐

Delays in analysis within the TGA:

☐

Delays in reporting by other sources (e.g. clinical registries):

☐

Other Risks (which need to be specified):

Next Steps for Level 1 & Level 2 Investigations

Next Steps for Level 1 Investigation:

Next Steps for Level 2 Investigation:

Case Closure

Click [N] to begin a new Correspondence entry. Note that the Email address specified here will receive a notification if the Date Received is not filled in by the Date Expected.

Correspondence and Chronology Details										
Include?	Heading	Type L1	Type L2	Email	Sent	Expected	Received	Response	Notes	
☐		Sponsor Routine Correspondence	Sponsor Notification Letter		14/05/2021					
☐		TGA Regulatory Action	Regulatory Information Request (s41)		30/07/2021		28/08/2021		s41JA request to the sponsor, Trim number: D21-2903364 s41JA response from the sponsor, Trim number: D21-3154881	
☐		Email (non-regulatory)	Information Request		09/09/2021		09/09/2021		Email sent for obtaining IFU, Trim number: D21-3071339 IFU sent by the sponsor, Trim number: D21-3075072	
☐		Sponsor Routine Correspondence	Sponsor DIR Closure Letter		13/10/2021				Trim number: D21-3211403	

List of Problem Observed Codes - Click [N] to begin entering information.

Problem Observed Details

Problem Observed (Level 1)	Problem Observed (Level 2)	Problem Observed (Level 3)	If 'Other' Selected
Material Integrity Problem	Break	Fracture	

Clinical signs symptoms and conditions

Details

Level 1	Level 2	Level 3
Injury	Foreign Body In Patient	Device Embedded In Tissue or Plaque

Health Impact

Details

Level 1	Level 2	Level 3
Surgical Intervention	Device Revision or Replacement	

Investigation Findings

Finding Details

Investigation Findings (Level 1)	Investigation Findings (Level 2)	Investigation Findings (Level 3)	If 'Other' Selected
No Findings Available			

Investigation Conclusion

Conclusion Details

Investigation Conclusion (L1)	Investigation Conclusion (L2)	If Additional Conclusion Detail Requested
Cause Not Established		

Investigation Outcomes

Outcome Details

Outcome of Investigation (L1)	Outcome of Investigation (L2)	If Additional Conclusion Detail Requested
Reviewed, No Further Action Required		

Investigation Summary

Latest Investigation (DII) where this DIR is the Primary DIR: 	Latest Investigation (DII) where this DIR is a Related DIR: 	Investigator: S22	Peer Review: No	
Investigator's Notes: Made a duplicate report after review of annual report and questionnaire response.		Summary Findings: A review of the manufacturer's analysis shows that the welding failed due to the component being subjected to stress of excessive lateral bending. The manufacturer notes that it is not possible to identify a clear or specific root cause for the mode of failure. A review of the warnings and precautions section in the Instructions for Use (IFU) has shown that the breakage of the MectaLock Suture Anchor can occur if any lateral loading is applied on the anchor/insertor while inserting the MectaLock anchor. The review of the risk management		Recall Number:

analysis document shows that all necessary actions have been implemented to minimize any potential risk associated with the laser welding integrity. In addition, this is supported by design verification and validation activities which confirmed the effectiveness of the device during its intended use. The sponsor notes that this is an isolated event and a review of the IRIS database has shown that there is only one similar adverse event which was investigated alongside with this DIR. No further investigation will occur at this time; however the TGA will continue to monitor the rate and pattern of occurrence and may re-open the file as appropriate.

Note: Letter generation buttons disabled if report not ready for website publication or risk analysis not completed.

Device Lookup

This section is used to match information provided via UDIR forms to ARTG information. You can select a Brand/Name from information provided in the 'Other Devices Involved' table below or enter information manually.

Other Device (Entered):	Brand Name:	Manufacturer Name:	Device ARTG #:

Other Devices					
Device ARTG No:	Manufacturer Name:	Sponsor/Supplier:	GMDN / UMDN Text:	Trade/Brand Name:	Serial #:
Model Number:	Batch #:	Lot #:	Expiry Date:		

Related DIR Information - Click **New** to begin entering information.

Rec No	
1	

Samples Record - Click **[N]** to begin entering information. **Note:** Sample # Generated on Save.

Rec No	Details	Sample Details			Additional Details				
1	Date Entered:	LIMS #:	Sample Requested:	Sample Received:	Manufacturer:	GMDN:	Device Description:	Brand Name:	Serial Number:
	Reason for Testing:	# Samples from Reporter:	# Samples from Sponsor:	Outcome of TGA's Testing:	Lot Number:	Batch Number:	Model Number:	Version Number:	
						Who sent the device to the TGA?:			Why does the TGA have the sample?:

Additional Patients

Click **[N]** to begin entering information.

Patient Details			
Sex:	Weight:	Age:	
Patient Focused Corrective Action Taken:		Patient History:	
Injured - Extent of Injury:	Was device directly linked to death?:	Was device directly linked to permanent disability?:	Consequence:

Other Consequence:	Describe any test (lab, xray, etc.):	Additional Event Description:	Medical Problem Device Used For: Document 1

Additional Device Information

Where did you get this device from?:	How reliant is the affected person on correct/safe operation of this device?:
Any other relevant information to aid assessing/investigating the incident?:	

Similar Events

Similar events - how many times?:	Date of Recent Report:

Device Access - Alternate Device Contact Information Provided

Title:	First Name:	Last Name:	Phone:
Fax:	Email:		

Incident Location Details

Occurred in Australia:	Organisation:	Address Line 1:	Address Line 2:
Town/Suburb:	State:	Postcode:	

Attachment(s) Details

Type	Open	Name	Size	Attached Within	Attached To
FILE		DIR 68516_Sponsor_Complete_Letter		169 Form	

Flow Details : DIR-REQ - Device Incident Request : 296794

Request Details

ID	Type	Location	Status	Assigned By	Assigned To	Assigned On	Priority	Attach
296794	DIR-REQ		Closed	S22	OPR Administration User	13/10/2021	Normal	0

Signature Details

Role	IRIS Investigator	
User	S22	
Signed At	24/02/2022 09:40:47	
Comment		

s22

From: IRIS <IRIS@health.gov.au>
Sent: Friday, 14 May 2021 3:02 PM
To: s22 @medacta.com.au
Subject: DIR 68416 - Sponsor Notification Letter [SEC=OFFICIAL]
Attachments: DIR 68516 - Sponsor Notification Letter.pdf

Dear s22

Please see the attached letter in relation to DIR 68516 (0172-21) - MectaLock instruments MectaLock PEEK Ø2.4 – Short driver- Ligament anchor

Kind regards,

s22

Device Support Team
Medical Device Incident Report Investigation Scheme (IRIS)
 Devices Post Market Monitoring | Medical Devices Surveillance Branch
 Therapeutic Goods Administration (TGA)
 Australian Government Department of Health

Email: IRIS@health.gov.au
 Phone: DST s22
 Phone: IRIS 1800 809 361
 Post: PO Box 100, Woden, ACT 2606
 Courier: 136 Narrabundah Lane, Symonston, ACT 2609

Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission. This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met

Please Note: Medical device sponsors and manufacturers are encouraged to submit adverse events via the online Medical Device Incident Reporting (MDIR) system. The MDIR allows sponsors and manufacturers to submit initial, follow-up and final reports and review reports already submitted to the TGA. The MDIR can be accessed via <https://apps.tga.gov.au/prod/mdir/MDIRSummary.aspx>, using a TGA eBusiness Service (eBS) user name and password. MDIR guidance documents and FAQ are available from the 'Training' section of the TGA eBS portal, or from the [TGA website](#).

MDIR Technical assistance is available via s22 or email (IRIS@health.gov.au). If you are a sponsor or manufacturer requiring a new account to gain access to the MDIR system or need to update your details with TGA eBS, please contact TBS Helpdesk via 1800 010 624 or email (eBS@health.gov.au).

Online reporting forms for Medical Device Healthcare Professionals/Users can be accessed via:

<https://apps.tga.gov.au/prod/mdir/udir03.aspx>

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

Australian Government
Department of Health
Therapeutic Goods Administration

**Australian Medical Device
Incident Report Investigation Scheme**

s22
Medacta Australia Pty Ltd
Unit A1 / 16 Mars Road
LANE COVE NSW 2066
Email: s22@medacta.com.au

File Reference: E21-295890
Sent by email

Dear s22

DEVICE INCIDENT REPORT DIR 68516 - ARTG # 320484 - Ligament anchor

An investigation into the incident reported to the Therapeutic Goods Administration concerning the above device is being carried out alongside DIR 68783.

A copy of the Medical Device Incident Report Investigation Scheme (IRIS) database entry, including is attached for your information.

Thank you for your support of the Medical Device Incident Report Investigation Scheme. Should you have any further queries concerning this report please contact our team on s22 or send an email to: IRIS@health.gov.au.

Yours sincerely

Signed electronically by

Administration Officer

Incident Report Investigation Scheme
Devices Post Market Monitoring Section
Therapeutic Goods Administration

14/05/2021

DIR 68516 - ARTG # 320484 - Ligament anchor

Reporter Reference #: 0172-21

Date of Adverse Event:

s22

Date of Initial Report:

05/03/2021

ARTG #:

320484

Brand Name:

MectaLock instruments MectaLock PEEK Ø2.4 – Short driver- Ligament anchor

Device Class:

Class IIb

Model #:

05.10.001

Serial #:

Software Version:

Batch #:

Lot #:

1924147

Manufacturer:

Medacta International SA [45296]

Sponsor:

Medacta Australia Pty Ltd [45972]

Unit A1 / 16 Mars Road

LANE COVE NSW 2066

Contact Name: s22

Phone:

Fax:

Email: s22@medacta.com.au

Reporter:

Mr s22

Complaint Coordinator

Medacta International SA

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Castel San Pietro TI 6874

Email: s22@medacta.ch

Confidential: N/A

Phone: s22

Fax:

Date of Implant:

s22

Date of Explant:

s22

Clinical Event Information:

The patient underwent a labral repair s22 Hip Arthroscopy. Two MectaLock (J-Lock) PEEK anchors were used to refixing the labrum to the acetabulum with no visual issues. Then the surgeon sent to the agent a picture of an X-ray asking if the image in the X-ray was part of the anchor inserter handle. It appears that a part of the inserter handle has remained in the bone with the anchor. The patient underwent revision surgery to remove the broken part.

The event was reported to the manufacturer on 29.01.2021, but only after the revision surgery, which took place on s22, we were sure that the piece that remained in the patient was part of the Medacta device.

Patient Outcome/Consequences:

Device Analysis Results:

Batch review performed on 23-Feb-2021:

Lot 1924147: 125 items manufactured and released on 24-Mar-2020. Expiration date: 2025-Feb-07. No anomalies found related to the problem. To date, 19 items of the same lot have been already sold without any other similar reported even.

Clinical Evaluation performed by Medical Affairs Director:

During an arthroscopic surgery of labrum refixation with a PEEK anchor, the tip of the positioning instrument broke, though the surgeon did not perceive the fracture. The tip of the instrument is visible in the postoperative X-rays, it remained in the cannulation of the anchor. We cannot determine the reasons for the fracture with the elements at hand.

No clinical cause could be ascertained. We were informed that a different surgeon intervened to remove the broken instrument part; the future clinical pathway of the patient should be totally flawless and no further consequences of the event should be expected.

Visual Inspection performed by R&D Project Manager:

The component of the MectaLock PEEK inserter has been received and analyzed. It's evident that the laser welding failed and the tip component has been detached from the main shaft of the inserter. The proximal portion of the analyzed component resulted totally ruined and scratched, most likely due to the extraction, therefore is not possible to perform a detailed microscopic analysis of the failure mode of the welding. However, it's also evident that the component has been subjected to excessive lateral bending as the axis of the component results deviated from the original. This suggests that misuse could have occurred during the anchor's placement or that a not foreseeable condition of use required an excessive force, as the MectaLock PEEK anchors is supposed to be impacted axially only within the predrilled cylindrical pilot hole.

Conclusion:

The device history record review does not indicate any potentially related manufacturing issue. The component has been subjected to lateral bending as the axis of the component results deviated from the original. This suggests that suboptimal usage occurred during the anchor's placement or that a not foreseeable condition of use required an excessive force, as the MectaLock PEEK anchors are supposed to be impacted axially only within the predrilled cylindrical pilot hole.

Corrective/Preventative Actions:

The event is an isolated case and no systematic action is currently indicated as applicable, however, it has been included in our active monitoring system for potential CAPA needs and opportunities.

Details of Similar Events:

1. MectaLock instruments 05.10.001 MectaLock PEEK Ø2.4 – Short driver Lot. 1924147: this is the 1st event for the family, the first for the reference and lot.

This is the first case reported. From the beginning till today, 03-Mar-2021 we have sold 201 pieces considering short and long drivers. The occurrence rate is $1/201=0.5\%$

Date Completed:

End of DIR 68516

s22

From: s22@Health.gov.au>
Sent: Wednesday, 13 October 2021 5:15 PM
To: s22@medacta.com.au
Subject: Sponsor Completion Letter for DIR 68516 [SEC=OFFICIAL]
Attachments: [D21-3211403] DIR 68516_Sponsor_Complete_Letter.PDF

Dear s22,

Good afternoon. Please find attached sponsor completion letter for DIR 68516. Thanking you.

Kind regards,

Dr s22
Departmental Officer
Devices Post Market Monitoring
Medical Devices Surveillance Branch
Email: s22@health.gov.au
Telephone number : s22

Therapeutic Goods Administration

Department of Health
PO Box 100
Woden ACT 2606 Australia

www.tga.gov.au

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met.

Important: *This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission.*

The Department of Health acknowledges the Traditional Custodians of Australia and their continued connection to land, sea and community. We pay our respects to them and their culture, and to all Elders both past and present.

P *Please consider the environment before printing this email.*

Australian Government
Department of Health
Therapeutic Goods Administration

**Australian Medical Device
Incident Report Investigation Scheme**

s22
Medacta Australia Pty Ltd
Unit A1 / 16 Mars Road
LANE COVE NSW 2066
Email: **s22**@medacta.com.au

File Reference: E21-295890
Sent by email

Dear **s22**

DEVICE INCIDENT REPORT DIR 68516 - ARTG # 320484 - Ligament anchor

An investigation into the incident reported to the Therapeutic Goods Administration concerning the above device is now complete.

A copy of the Medical Device Incident Report Investigation Scheme (IRIS) database entry, including the investigation summary is attached for your information.

Thank you for your support of the Medical Device Incident Report Investigation Scheme. Should you have any further queries concerning this report please contact our team on **s22** or

send an email to: IRIS@health.gov.au

Yours sincerely

Signed electronically by

Administration Officer

Incident Report Investigation Scheme
Devices Post Market Monitoring Section
Therapeutic Goods Administration

13/10/2021

DIR 68516 - ARTG # 320484 - Ligament anchor

Reporter Reference #: 0172-21

Date of Adverse Event:

s22

Date of Initial Report:

05/03/2021

ARTG #:

320484

Brand Name:

MectaLock instruments MectaLock PEEK Ø2.4 – Short driver- Ligament anchor

Device Class:

Class IIb

Model #:

05.10.001

Serial #:

Software Version:

Batch #:

Lot #:

1924147

Manufacturer:

Medacta International SA [45296]

Sponsor:

Medacta Australia Pty Ltd [45972]

Unit A1 / 16 Mars Road

LANE COVE NSW 2066

Contact Name: s22

Phone:

Fax:

Email: s22@medacta.com.au

Reporter:

Mr s22

Complaint Coordinator

Medacta International SA

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Castel San Pietro TI 6874

Email: s22@medacta.ch

Confidential: N/A

Phone: s22

Fax:

Date of Implant:

s22

Date of Explant:

s22

Clinical Event Information:

The patient underwent a labral repair s22 Hip Arthroscopy. Two MectaLock (J-Lock) PEEK anchors were used to refixing the labrum to the acetabulum with no visual issues. Then the surgeon sent to the agent a picture of an X-ray asking if the image in the X-ray was part of the anchor inserter handle. It appears that a part of the inserter handle has remained in the bone with the anchor. The patient underwent revision surgery to remove the broken part.

The event was reported to the manufacturer on 29.01.2021, but only after the revision surgery, which took place on s22, we were sure that the piece that remained in the patient was part of the Medacta device.

Patient Outcome/Consequences:

Device Analysis Results:

Batch review performed on 23-Feb-2021:

Lot 1924147: 125 items manufactured and released on 24-Mar-2020. Expiration date: 2025-Feb-07. No anomalies found related to the problem. To date, 19 items of the same lot have been already sold without any other similar reported even.

Clinical Evaluation performed by Medical Affairs Director:

During an arthroscopic surgery of labrum refixation with a PEEK anchor, the tip of the positioning instrument broke, though the surgeon did not perceive the fracture. The tip of the instrument is visible in the postoperative X-rays, it remained in the cannulation of the anchor. We cannot determine the reasons for the fracture with the elements at hand.

No clinical cause could be ascertained. We were informed that a different surgeon intervened to remove the broken instrument part; the future clinical pathway of the patient should be totally flawless and no further consequences of the event should be expected.

Visual Inspection performed by R&D Project Manager:

The component of the MectaLock PEEK inserter has been received and analyzed. It's evident that the laser welding failed and the tip component has been detached from the main shaft of the inserter. The proximal portion of the analyzed component resulted totally ruined and scratched, most likely due to the extraction, therefore is not possible to perform a detailed microscopic analysis of the failure mode of the welding. However, it's also evident that the component has been subjected to excessive lateral bending as the axis of the component results deviated from the original. This suggests that misuse could have occurred during the anchor's placement or that a not foreseeable condition of use required an excessive force, as the MectaLock PEEK anchors is supposed to be impacted axially only within the predrilled cylindrical pilot hole.

Conclusion:

The device history record review does not indicate any potentially related manufacturing issue. The component has been subjected to lateral bending as the axis of the component results deviated from the original. This suggests that suboptimal usage occurred during the anchor's placement or that a not foreseeable condition of use required an excessive force, as the MectaLock PEEK anchors are supposed to be impacted axially only within the predrilled cylindrical pilot hole.

Corrective/Preventative Actions:

The event is an isolated case and no systematic action is currently indicated as applicable, however, it has been included in our active monitoring system for potential CAPA needs and opportunities.

Details of Similar Events:

1. MectaLock instruments 05.10.001 MectaLock PEEK Ø2.4 – Short driver Lot. 1924147: this is the 1st event for the family, the first for the reference and lot.

This is the first case reported. From the beginning till today, 03-Mar-2021 we have sold 201 pieces considering short and long drivers. The occurrence rate is $1/201=0.5\%$

Number of Similar Events:

Rate of Similar Events:

Countries Similar Events Also Occurred:

Clinical Signs (Level 1)

Injury

Health Impacts (Level 1)

Surgical Intervention

Clinical Signs (Level 2)

Foreign Body In Patient

Health Impacts (Level 2)

Device Revision or Replacement

Type of Problem (Level 1)
Material Integrity Problem

Type of Problem (Level 2)
Break

Cause of Problem (Level 1)

Cause of Problem (Level 2)

No Findings Available

Outcome of Investigation

Reviewed, No Further Action Required

Summary of Investigation:

A review of the manufacturer's analysis shows that the welding failed due to the component being subjected to stress of excessive lateral bending. The manufacturer notes that it is not possible to identify a clear or specific root cause for the mode of failure. A review of the warnings and precautions section in the Instructions for Use (IFU) has shown that the breakage of the MectaLock Suture Anchor can occur if any lateral loading is applied on the anchor/insertor while inserting the MectaLock anchor. The review of the risk management analysis document shows that all necessary actions have been implemented to minimize any potential risk associated with the laser welding integrity. In addition, this is supported by design verification and validation activities which confirmed the effectiveness of the device during its intended use. The sponsor notes that this is an isolated event and a review of the IRIS database has shown that there is only one similar adverse event which was also investigated alongside with this DIR. No further investigation will occur at this time; however the TGA will continue to monitor the rate and pattern of occurrence and may re-open the file as appropriate.

Date Completed:

End of DIR 68516

Device Incident Report : Medical Devices Branch - Device Vigilance and Monitoring

DIR : 44 - ID : 517464

Released by §22 on 18/03/2021 17:26:54

Report #:	Records Management #:	Reporter's Reference #:	Report Type:
68783	E21-259665	§22 ID 69823	Final
ARTG: 320484	Document Container URL		

Report Information Section

Report Status:	Sponsor's Reported Category:	Date of Adverse Event:	Date of Initial Report:
Closed		§22	18/03/2021
Date of Final Report:	Date of Initial TGA Action:	Reviewed by Team:	Date Response Received:
18/03/2021	18/03/2021	23/03/2021	
Date Completed:	Operator at Time of Event:	If 'Other' Operator Selected:	Reporter consents to contact by sponsor:
06/10/2021	Doctor		No
Source of Report:	If 'Other' Source Selected:	Type of Initial Action:	
Hospital Administrator		Trend data only	

Event Description for Website Publication:

Anchor introducer snapped off during surgery and not found until follow-up appointment with surgeon.

Clinical Event Information:

Patient was admitted into surgery for right hip arthroscopy and femoral osteotomy and labral repair on the [redacted], Anchor introducer found retained in patient on surgeon appointment. Anchor introducer snapped off during surgery and not found until patient went for appointment with surgeon. Additional surgery required for removal of retained anchor introducer on §22

Number of Incidents in Report:	Contact:	Alternative Person Title:	Alternative Person First Name:
1	Reporter		
Alternative Person Surname:	Alternative Person Phone:	Alternative Person Fax:	Alternative Person Email:

Recorded Problems Observed

Recorded Problems Observed:

Material Integrity Problem -> Break -> Fracture

Clinical Signs, Symptoms and Conditions

Recorded Clinical Signs, Symptoms and Conditions:

Others -> Insufficient Information ->

Health Impact

Recorded Health Impacts:

Surgical Intervention -> Additional Surgery ->

Patient Information

Sex:	Weight:	Age:
§22		§22 old
Patient Focused Corrective Action Taken:		Patient History:
		NA
Patient Outcome/Consequences:		Additional Event Description:
Additional surgery required for removal of retained anchor introducer on §22		Anchor introducer snapped off during surgery and not found until patient went for appointment with surgeon.
Describe any test (Lab, xray, etc.):	Injured - Extent of Injury:	Consequence:
NA	Temporary Injury	Other
Other medical devices currently using/implanted:	Medical Problem Device Used For:	Additional Patients Added:
		Additional surgery to retrieve item

0

Submitting Reporter Section

Search Reporter By Surname:

§22

Reporter Title:

MS

Position:

§22

Address 1:

§22

Country:

Australia

Mobile:

Reporter #:

First Name:

§22

Address 2:

Postcode:

§22

Email:

§22

Surname:

§22

Company/Institution:

§22

Town/Suburb:

§22

Phone:

§22

Last External Submission By:

Preferred Contact Method:

Email

State:

§22

Fax:

Initial Reporter Section

As Above?:

Yes

Search Reporter By Surname:

Title:

Position:

Address 1:

Postcode:

Mobile:

If No, fill out the following:

Initial Reporter #:

First Name:

Address 2:

Country:

Email:

Surname:

Company/Institution:

Town/Suburb:

Phone:

Allow the device company
to contact you about the incident:☐

Initial Reporter Confidential:

Preferred Contact Method:

Device Information Section

Product Exempt (Note: If not exempt, enter ARTG No):

No

Product Licence Category:

Included

Brand Name:

Medacta J-Lock PEEK Suture - Anchor, soft-tissue implanted,
non-biodegradable

Model #:

Purchase Date:

Date of Initial Procedure:

§22

Access Contact First Name:

Search Device ARTG:

320484

Device Class:

Class IIb

Initial Device Description:

J-Lock PEEK Suture - Anchor, soft-tissue implanted, non-
biodegradable

Serial #:

Expiry Date:

17/07/2024

Place of Implantation:

Hip

Access Contact Surname:

Device ARTG #:

320484

GMDN / UMDN Code:

36174

Usage of Device:

Batch #:

Ref 05.11.001

Date of Implant:

§22

Reported Device Location:

Place of use

Access Contact Phone:

Therapeutic Licence Type:

Medical Device

GMDN / UMDN Text:

Ligament anchor

Software Version:

Lot #:

1954843

Date of Explant:

Access Contact Title:

Access Contact Fax:

Access Contact Email:

Licence Status:

A

Status Effective Date:

19/07/2019

Additional Devices Added:

0

Document 4

Manufacturer Information Section

Manufacturer Name:

Medacta International SA

Manufacturer Client Id:

45296

Address 1:

Address 2:

Town/Suburb:

State/Province:

Country:

Postcode:

Phone:

Fax:

Email:

Manufacturer Informed:

Date Aware of Adverse Event:

Contact Title:

Contact First Name:

Contact Surname:

Supplier Information Section

Supplier Name:

Medacta

Address 1:

Address 2:

Town/Suburb:

State:

Country:

Australia

Postcode:

Phone:

Fax:

Email:

Website:

Supplier Informed:

Yes

Date of Supplier Contact:

01/03/2021

Contact Title:

Contact First Name:

Contact Surname:

Contact Phone:

Contact Fax:

Contact Email:

Report Status

For website publication:

Yes

Ready for Publication:

Yes

Investigated:

Yes

Investigation Reason:

More information is required

Team Assignment:

Team B (IIB, Is & Im)

Team Priority:

Routine

Team Review

Reviewed by Team:

23/03/2021

Reason Sent To Meeting:

User report

Outcome from team meeting:

Investigation (within DIR) is recommended

Notes for Team meeting:

No DIRs, PMRs or Recalls. User report. Sent to meeting as unsure if this is the correct ARTG. Reviewed on-line IFU and found no mention of an introducer, could this possibly be referring to a component of a surgical kit under a different ARTG entry? (<https://www.medacta.com/EN/resources-med-reviewed-ifu>)

Outcomes from Team Meeting:

Level 2 investigation, send questionnaire

Initial Risk Analysis

Background Information	Risk Assessment - Section A	Risk Assessment - Section B	Risk Assessment - Section C	Risk Assessment - Section D	
Date:	Severity:	Incidents in the last 12 months:	Manufacturer analysis:	Assessor:	Manufacturer documentation:
18/03/2021	5 - An illness/injury was resolved or prevented with treatment by a health professional			\$22	Unknown - updated information from the manufacturer is required
Incidents in last 24 months:	Manufacturer action:	ESTIMATED LEVEL OF INVESTIGATION:	FINAL LEVEL OF INVESTIGATION:	Injured Party:	Device Recalls:
		Level 1 Investigation (to complete screening)	Level 1 Investigation (to complete screening)	Patient	0. No recalls for similar incidents in Australia

Incidents in last 36 months: 	IVD status: 	EXCEPTION TO INVESTIGATION LEVEL: Device not returned. Isolated incident, no other reports. May want to check if this is the correct ARTG entry - introducer may not be included in the ARTG entry may be a surgical kit component. https://www.medacta.com/EN/resources-med - reviewed IFU, no mention of "introducer"		Found Prior To Use: No	Is AE covered by Document 4: No
Incidents Worldwide: 	Number of potential contributing factors: Yes - some potential factors (up to 3)			Reusable: No	Similar events (past 6 months): 1 incidents
Products supplied the last 12 months: 	Specific factors identified: Use of device - according to IFUs, Use of device - interaction with other equipment	ESTIMATED LEVEL OF PRIORITY: Routine	FINAL LEVEL OF PRIORITY: Routine	3 or more events - batch/model: 	
Products supplied last 24 months: 	Number of potential sensitivities: No	EXCEPTION TO PRIORITY LEVEL: 		3 or more events - health district: 	
Products supplied last 36 months: 	Specific sensitivities identified: 			3 or more events - organisation: 	
Products supplied Worldwide: 	Consultations during risk assessment: I did none of the above incidents	Final Risk Assessment: Yes			

Sponsor/Manufacturer Information Section

Search Sponsors: 45972	Name: Medacta Australia Pty Ltd	Client #: 45972
Attention To: §22	Address 1: Unit A1 / 16 Mars Road	Town/Suburb: LANE COVE
State: NSW	Postcode: 2066	Fax:
Email: §22@medacta.com.au		

Investigation Information Section - Submitted by Sponsor/Manufacturer

Device Analysis Results: **Additional information received by Sponsor 03/06/2021** The component of the MectaLock PEEK inserter has been received and analysed. It's evident that the laser welding failed and the tip component has been detached from the main shaft of the inserter. The proximal portion of the analysed component was received scratched, most likely due to the extraction, therefore is not possible to perform a detailed microscopic analysis of the failure mode of the welding. A deeper analysis on the safety of the device has been carried out by means of further mechanical investigation on equivalent devices sourced from the same production batch, as representative of the component object of this complaint. In particular, the welding between the shaft and the inserter tip has been subject to stress limit tests which confirmed the safety and reliability of the device under load conditions that are commonly expected during a surgery. Therefore, it's not possible to identify a clear or specific root cause for the mode of failure, unless an outstanding condition required a not foreseeable usage during the surgical procedure.	Details of Similar Events:
Additional Details (use for tables): 	CAPA# Reference:
Risk Assessment	
Frequency: 	Severity:
Rating: 	
Type Cause and Outcome: 	Number of Similar Events:
Expected Rate: 	Actual Rate:
Countries Similar Events Also Occurred: 	
Completed Actions: 	Planned Actions and Proposed Timelines:

Additional Comments:

Amended final received 03/06/2021, DIR updated

DIR was closed on 06/10/2021.

Reason for Level 1 Investigation

Details of Reasons

Reason for Level 1 Investigation

Focus of Level 2 Investigation

Details of Focus

Essential Principles

If 'Other' Selected

EP13-Information may not have been provided

EP2-Design/construction may not conform with principles

Sources of Evidence for Level 2

Details of Source

Sources of Evidence

If 'Others' please specify here

Expected Sourcing Date

Date of Evidence Received

Information from Sponsors (e.g. 41JA, Questionnaires, Emails)

27/04/2021

26/04/2021

Evidence

Investigation Questions (Level 1 and Level 2):

1.Please provide the response of the standard questionnaire which contains a detailed manufacturer's investigation analysis to confirm that it adequately demonstrates the conformity with EP2 design and construction principles.

2.Please send a copy of Instructions for Use (IFU).

Potential Risks

Delays in response by product manufacturers:

Delays in response by incident reporters:

Delays in analysis within the TGA:

Delays in reporting by other sources (e.g. clinical registries):

Other Risks (which need to be specified):

Next Steps for Level 1 & Level 2 Investigations

Next Steps for Level 1 Investigation:

Next Steps for Level 2 Investigation:

Case Closure

Click [N] to begin a new Correspondence entry. Note that the Email address specified here will receive a notification if the Date Received is not filled in by the Date Expected.

Correspondence and Chronology Details

Include?	Heading	Type L1	Type L2	Email	Sent	Expected	Received	Response	Notes	Document 4
☐		Sponsor Routine Correspondence	Questionnaire - non-regulatory		26/03/2021	27/04/2021	26/04/2021			
☐		Reporter Routine Correspondence	Reporter Notification Letter		26/03/2021					
☐		TGA Regulatory Action	Regulatory Information Request (s41)		30/07/2021	30/08/2021	28/08/2021		s41JA request to the sponsor, Trim number: D21-2903364, s41JA response from the sponsor, Trim number: D21-3154881	
☐		Email (non-regulatory)	Information Request		09/09/2021	16/09/2021	09/09/2021		Email sent for obtaining IFU, TRIM number: D21-3071339 IFU sent by the sponsor, Trim number: D21-3075072	
☐		Reporter Routine Correspondence	Reporter DIR Closure Letter		06/10/2021				Trim number: D21-3178893	
☐		Sponsor Routine Correspondence	Sponsor DIR Closure Letter		06/10/2021				Trim number: D21-3178905	

List of Problem Observed Codes - Click **[N]** to begin entering information.

Problem Observed Details

Problem Observed (Level 1)	Problem Observed (Level 2)	Problem Observed (Level 3)	If 'Other' Selected
Material Integrity Problem	Break	Fracture	

Clinical signs symptoms and conditions

Details	
Level 1	Level 2
Others	Insufficient Information

Health Impact

Details	
Level 1	Level 2
Surgical Intervention	Additional Surgery

Investigation Findings

Finding Details	
Investigation Findings (Level 1)	Investigation Findings (Level 2)
No Findings Available	

Investigation Conclusion

Conclusion Details	
Investigation Conclusion (L1)	Investigation Conclusion (L2)
Cause Not Established	

Investigation Outcomes

Outcome Details		
Outcome of Investigation (L1)	Outcome of Investigation (L2)	If Additional Conclusion Detail Requested
Reviewed, No Further Action Required		

Investigation Summary

Latest Investigation (DII) where this DIR is the Primary DIR:	Latest Investigation (DII) where this DIR is a Related DIR:	Investigator:	Peer Review:	Peer Review Completion Date:
		s22	Yes	30/09/2021
Investigator's Notes:		Summary Findings:		Recall Number:
**To be Peer Reviewed by s22 **		A review of the manufacturer's analysis shows that the welding failed due to the component being subjected to stress of excessive lateral bending. The manufacturer notes that it is not possible to identify a clear or specific root cause for the mode of failure. A review of the warnings and precautions section in the Instructions for Use (IFU) has shown that the breakage of the MectaLock Suture Anchor can occur if any lateral loading is applied on the anchor/insertor while inserting the MectaLock anchor. The review of the risk management analysis document shows that all necessary actions have been implemented to minimize any potential risk associated with the laser welding integrity. In addition, this is supported by design verification and validation activities which confirmed the effectiveness of the device during its intended use. The sponsor notes that this is an isolated event and a review of the IRIS database has shown that there is only one similar adverse event. No further investigation will occur at this time; however the TGA will continue to monitor the rate and pattern of occurrence and may re-open the file as appropriate.		

Note: Letter generation buttons disabled if report not ready for website publication or risk analysis not completed.

Device Lookup

This section is used to match information provided via UDIR forms to ARTG information. You can select a Brand/Name from information provided in the 'Other Devices Involved' table below or enter information manually.

Other Device (Entered):	Brand Name:	Manufacturer Name:	Device ARTG #:

Other Devices

Device ARTG No:	Manufacturer Name:	Sponsor/Supplier:	GMDN / UMDN Text:	Trade/Brand Name:	Serial #:
Model Number:	Batch #:	Lot #:	Expiry Date:		

Related DIR Information - Click **New** to begin entering information.

Rec No	
1	

Samples Record - Click **[N]** to begin entering information. **Note:** Sample # Generated on Save.

Rec No	Details	Sample Details			Additional Details				
1	Date Entered:	LIMS #:	Sample Requested:	Sample Received:	Manufacturer:	GMDN:	Device Description:	Brand Name:	Serial Number:
	Reason for Testing:	# Samples from Reporter:	# Samples from Sponsor:	Outcome of TGA's Testing:	Lot Number:	Batch Number:	Model Number:	Version Number:	
		Who sent the device to the TGA?:					Why does the TGA have the sample?:		

Additional Patients

Click [N] to begin entering information.

Patient Details	
Sex:	Weight:
Age:	
Patient History:	
Patient Focused Corrective Action Taken:	
Injured - Extent of Injury:	Was device directly linked to death?:
Was device directly linked to permanent disability?:	
Consequence:	
Other Consequence:	Describe any test (Lab, xray, etc.):
Additional Event Description:	
Medical Problem Device Used For:	

Additional Device Information

Where did you get this device from?:	How reliant is the affected person on correct/safe operation of this device?:
Supplier	
Any other relevant information to aid assessing/investigating the incident?:	

Similar Events

Similar events - how many times?:	Date of Recent Report:	Event Reported To:	Reporter Reference Number:

Device Access - Alternate Device Contact Information Provided

Title:	First Name:	Last Name:	Phone:
Fax:	Email:		

Incident Location Details

Occurred in Australia:	Organisation:	Address Line 1:	Address Line 2:
Yes		159 Grey Street	
Town/Suburb:	State:	Postcode:	
East Melbourne	VIC	3002	

Attachment(s) Details

Type	Open	Name	Size	Attached Within	Attached To
FILE		DIR 68783 - Original User Report	431 Form		
FILE		R Amended DIR 68783 - Questionnaire Letter - R...	268 Form		
FILE		R Amended DIR 68783 - Questionnaire Letter - R...	1267 Form		

FILE		Peer Review Checklist	94	Form
FILE		DIR 68783_Reporter_Complete_Letter	158	Form
FILE		DIR 68783_Sponsor_Complete_Letter	162	Form

Document 4

Flow Details : DIR-REQ - Device Incident Request : 297807

Request Details

ID	Type	Location	Status	Assigned By	Assigned To	Assigned On	Priority	Attach
297807	DIR-REQ		Closed	S22	OPR Administration User	06/10/2021 ✓	Normal	0

Signature Details

Role	IRIS Investigator	
User	S22	
Signed At	06/10/2021 11:14:31	
Comment		

s22

From: IRIS <IRIS@health.gov.au>
Sent: Friday, 26 March 2021 12:26 PM
To: s22
Subject: DIR 68783 - Reporter Notification Letter [SEC=OFFICIAL, ACCESS=Personal-Privacy]
Attachments: DIR 68783 - Reporter Notification Letter.pdf

Dear s22

Please see the attached letter in relation to the Device Incident Report (DIR) that you submitted to the Therapeutic Goods Administration (TGA).

Kind regards,

s22

Device Support Team
Medical Device Incident Report Investigation Scheme (IRIS)
 Devices Post Market Monitoring | Medical Devices Surveillance Branch
 Therapeutic Goods Administration (TGA)
 Department of Health

Email: IRIS@health.gov.au
 Phone: DST s22
 Phone: IRIS 1800 809 361
 Post: PO Box 100, Woden, ACT 2606
 Courier: 136 Narrabundah Lane, Symonston, ACT 2609

***Important:** This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission. This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met*

Please Note: Medical device sponsors and manufacturers are encouraged to submit adverse events via the online Medical Device Incident Reporting (MDIR) system. The MDIR allows sponsors and manufacturers to submit initial, follow-up and final reports and review reports already submitted to the TGA. The MDIR can be accessed via <https://apps.tga.gov.au/prod/mdir/MDIRSummary.aspx>, using a TGA eBusiness Service (eBS) user name and password. MDIR guidance documents and FAQ are available from the 'Training' section of the TGA eBS portal, or from the [TGA website](#).

MDIR Technical assistance is available via s22 or email (IRIS@health.gov.au). If you are a sponsor or manufacturer requiring a new account to gain access to the MDIR system or need to update your details with TGA eBS, please contact TBS Helpdesk via 1800 010 624 or email (eBS@health.gov.au).

Online reporting forms for Medical Device Healthcare Professionals/Users can be accessed via:

<https://apps.tga.gov.au/prod/mdir/udir03.aspx>

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

Australian Government
Department of Health
 Therapeutic Goods Administration

**Australian Medical Device
 Incident Report Investigation Scheme**

MS s22
 Quality and Risk Consultant

File Reference: E21-259665
 Sent by email

s22
 s22
 s22

Email: s22

Dear MS s22,

**DEVICE INCIDENT REPORT DIR 68783 - Medacta J-Lock PEEK Suture - Anchor, soft-tissue
 implanted, non-biodegradable**

Thank you for your recent correspondence concerning a problem experienced with the above device.

The information you provided has been entered into the Medical Device Incident Report Investigation (IRIS) Database, where it will be evaluated against any previous incidents with the same or similar devices. This report is currently under investigation and a copy of the Device Incident Report (DIR) is attached for your reference.

Should you have any questions, please contact our team on s22 or email: IRIS@health.gov.au quoting the above DIR number. The TGA may contact you regarding this report.

Assistance and information regarding product refunds, returns or product compensation should be addressed with the Australian Competition and Consumer Commission (ACCC) via

<https://www.accc.gov.au/consumers/complaints-problems>

If you would like to lodge a complaint regarding the clinical practice of an Australian healthcare professional you should contact the Australian Health Practitioner Regulation Agency (AHPRA) via <https://www.ahpra.gov.au/notifications/make-a-complaint.aspx>, or the Health Care Complaints Commission within your State or Territory. More information on the respective State and Territory health complaint organisations can be accessed via

<https://www.ahpra.gov.au/notifications/further-information/health-complaints-organisations.asp>

Thank you for your support of the Medical Device Incident Report Investigation Scheme.

Yours sincerely
 Signed electronically by

Administrative Officer

Incident Report Investigation Scheme
 Devices Post Market Monitoring Section
 Therapeutic Goods Administration
 26/03/2021

DIR 68783 - Medacta J-Lock PEEK Suture - Anchor, soft-tissue implanted, non-biodegradable

Reporter Reference #: s22 ID 69823

Date of Adverse Event:

s22

Date of Final Report:

18/03/2021

ARTG #:
330998Brand Name:
Medacta J-Lock PEEK Suture - Anchor, soft-tissue implanted,
non-biodegradableDevice Class:
Class IIb

Model #:

Serial #:

Software Version:

Batch #:
Ref 05.11.001Lot #:
1954843Manufacturer:
Medacta International SASponsor:
Medacta Australia Pty Ltd
Unit A1 / 16 Mars Road
LANE COVE NSW 2066

Contact Name:

s22

Phone:

s22

Reporter:

Confidential: Yes

Date of Implant:

s22

Date of Explant:

Clinical Event Information:

Patient was admitted into surgery for right hip arthroscopy and femoral osteotomy and labral repair on the [redacted], Anchor introducer found retained in patient on surgeon appointment. Anchor introducer snapped off during surgery and not found until patient went for appointment with surgeon. Additional surgery required for removal of retained anchor introducer on s22

Patient Outcome/Consequences:

Additional surgery required for removal of retained anchor introducer on s22

Additional Event Description:

Anchor introducer snapped off during surgery and not found until patient went for appointment with surgeon.

End of DIR 68783

s22

From: IRIS <IRIS@health.gov.au>
Sent: Friday, 26 March 2021 12:30 PM
To: s22 @medacta.com.au
Subject: DIR 68783 - Questionnaire Letter - Response due COB 27/03/2021 [SEC=OFFICIAL]
Attachments: DIR 68783 - Questionnaire Letter.rtf

Dear s22

Please see the attached **Questionnaire Letter** for action.
***Please note the due date for a response is COB 27 March 2021.**

Kind regards,

s22

Device Support Team
Medical Device Incident Report Investigation Scheme (IRIS)
 Devices Post Market Monitoring | Medical Devices Surveillance Branch
 Therapeutic Goods Administration (TGA)
 Department of Health

Email: IRIS@health.gov.au
 Phone: DST s22
 Phone: IRIS 1800 809 361
 Post: PO Box 100, Woden, ACT 2606
 Courier: 136 Narrabundah Lane, Symonston, ACT 2609

***Important:** This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission. This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met*

Please Note: Medical device sponsors and manufacturers are encouraged to submit adverse events via the online Medical Device Incident Reporting (MDIR) system. The MDIR allows sponsors and manufacturers to submit initial, follow-up and final reports and review reports already submitted to the TGA. The MDIR can be accessed via <https://apps.tga.gov.au/prod/mdir/MDIRSummary.aspx>, using a TGA eBusiness Service (eBS) user name and password. MDIR guidance documents and FAQ are available from the 'Training' section of the TGA eBS portal, or from the [TGA website](#).

MDIR Technical assistance is available via s22 or email (IRIS@health.gov.au). If you are a sponsor or manufacturer requiring a new account to gain access to the MDIR system or need to update your details with TGA eBS, please contact TBS Helpdesk via 1800 010 624 or email (eBS@health.gov.au).

Online reporting forms for Medical Device Healthcare Professionals/Users can be accessed via:

<https://apps.tga.gov.au/prod/mdir/udir03.aspx>

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

Australian Government

Department of Health
Therapeutic Goods Administration

**Australian Medical Device
Incident Report Investigation Scheme**

s22

Medacta Australia Pty Ltd
Unit A1 / 16 Mars Road
LANE COVE NSW 2066
Email: s22@medacta.com.au

File Reference: E21-259665

Sent by email

Dear s22

**DEVICE INCIDENT REPORT DIR 68783 - ARTG # 330998 - Anchor, soft-tissue implanted,
non-biodegradable**

The Therapeutic Goods Administration has been advised of an incident involving the above product. A copy of the Device Incident Report (DIR) is attached.

To assist in the evaluation and resolution of this report, please provide the information requested in the attached questionnaire and return it to this office **within 20 working days of the date of this letter**, and no later than COB 27/04/2021.

Please send responses via email to IRIS@health.gov.au, referencing the DIR number.

If you are unable to respond with all the information requested by the due date please advise, **within the 20 days**, when a full response will be provided. Extensions of a reasonable time frame will be accepted depending on the seriousness of the complaint and the time requested.

Thank you for your cooperation. If you require further information please contact me on

s22 or email: IRIS@health.gov.au.

Yours sincerely
Signed electronically by

Administration Officer

Incident Report Investigation Scheme
Devices Post Market Monitoring Section
Therapeutic Goods Administration

26/03/2021

**- MEDICAL DEVICE INCIDENT REPORT INVESTIGATION SCHEME -
INITIAL REQUEST OF INFORMATION FROM LISTED SPONSOR**

Document 6

INITIAL REQUEST FOR INFORMATION FROM LISTED SPONSOR

Date: 26/03/2021

DIR: 68783

Manufacturer Name: Medacta International SA [45296]

Question/Requirement

1) Please confirm the device's Australian Register of Therapeutic Goods (ARTG) number

ARTG:

2) Do you currently supply or have you previously supplied this product, with the indicated Model/Serial/Batch/Lot numbers:

- a) To the Australian Market
- b) For Export

YES

NO

☐
☐☐
☐

3) How many of this model have been supplied

In Australia:

Worldwide:

4) How many of this batch (if applicable) have been supplied:

In Australia:

Worldwide:

5) Are you aware of this problem, as reported?

☐☐

6) If deemed necessary, is a sample of the mentioned device available for review and/or testing?

☐☐

7) Have you had any other reports of similar* problems with this product in Australia?

☐☐

If **YES**, how many:

If **YES**, please give details:

* Similar events are based on the clinical event description and not the cause of an event. Both the confirmed and unconfirmed rates for similar events are often required and beneficial to show the full outlook of the device and its use.

8) Is the manufacturer aware of reports of similar problems with this product?

☐

Document 6

☐

If **YES**, how many:

If **YES**, please give details:

9) Please provide details of any action you have taken, or intend to take, regarding this problem

10) Please provide details of the manufacturer's investigation to date, including expected **Manufacturer's investigation completion date**

Date:

11) When returning this response to the office of the Therapeutic Goods Administration, you are requested to attach the following (if checked):

☐

Sample of the product/device

☐

Operator's manual

☐

Product Specifications

☐

Technical Service Manual

☐

Descriptive product promotional documentation

☐

Clinical training manual in printed or video form

☐

Instructions for use, as supplied with the device

☐

In-house training documentation

☐

Device Packaging with printed instructions

☐

Evidence of compliance with the Essential Principles

☐

A summary of risk assessment activities performed by the manufacturer for the device, eg Risk Management Report required by Clause 8 of ISO 14971:2007.

☐

Service History, and Safety and Performance Test Results

☐

Configuration Information and Documentation

☐

HHE (Health Hazard Evaluation)

12) Additional Information required:

Document 6

13) If your device is an implantable pacemaker/defibrillator and/or associated lead you are asked to provide the following additional information:

1. Both published and unpublished clinical trial data where events of this type are analysed.
2. The number of reported events of ALL types (including unconfirmed events), the number of devices sold and the cumulative implant months for each device in this product family.
3. What material has been used to insulate, both internally and externally, the lead? (applicable to leads only)

Information Supplied By:

14)

Name		Phone	
Signature		Fax	
Position		Email	

This questionnaire and any appended documents should be returned to the TGA within **20 working days or as specified on page 1 of this letter.**

Electronic submission of all information is preferred.

Please send the requested information to email address: IRIS@health.gov.au

Document 6

For large size documents, please post a universal serial bus (USB), compact disc (CD), or digital versatile disc (DVD) via postal address:

Incident Report Investigation Scheme
Devices Post Market Monitoring Section
Therapeutic Goods Administration
PO Box 100
Woden ACT 2606

If you are sending a device/s to the TGA please follow the instructions via link:

<https://www.tga.gov.au/form/report-medical-device-adverse-event-medical-device-use>

Sponsors of products listed or registered on the Australian Register of Therapeutic Goods (ARTG) are reminded of their responsibilities under Section 31 and/or 41JA (as appropriate) of the Therapeutic Goods Act of 1989, to provide information relating to their product's formulation, composition, design specification, quality, method and place of manufacture, presentation, safety and efficacy, conformity to advertising regulations under the Act,

Reporter Reference #: s22 ID 69823

Date of Adverse Event:
s22Date of Final Report:
18/03/2021

ARTG #:	Brand Name:	
330998	Medacta J-Lock PEEK Suture - Anchor, soft-tissue implanted, non-biodegradable	
Device Class:	Model #:	Serial #:
Class IIb		
Software Version:	Batch #:	Lot #:
	Ref 05.11.001	1954843

Manufacturer:
Medacta International SA [45296]Sponsor:
Medacta Australia Pty Ltd [45972]
Unit A1 / 16 Mars Road
LANE COVE NSW 2066

Contact Name: s22

Phone: s22

Email: s22@medacta.com.au

Fax:

Reporter:

Confidential: Yes

Date of Implant:
s22

Date of Explant:

Clinical Event Information:

Patient was admitted into surgery for right hip arthroscopy and femoral osteotomy and labral repair on the [redacted], Anchor introducer found retained in patient on surgeon appointment.

Anchor introducer snapped off during surgery and not found until patient went for appointment with surgeon. Additional surgery required for removal of retained anchor introducer on s22

Patient Outcome/Consequences:

Additional surgery required for removal of retained anchor introducer on s22

Additional Event Description:

Anchor introducer snapped off during surgery and not found until patient went for appointment with surgeon.

Device Analysis Results:

Corrective/Preventative Actions:

Details of Similar Events:

Number of Similar Events:

Rate of Similar Events:

Countries Similar Events Also Occurred:

s22

From: IRIS <IRIS@health.gov.au>
Sent: Friday, 26 March 2021 12:33 PM
To: s22 @medacta.com.au
Subject: *Amended* DIR 68783 - Questionnaire Letter - Response due COB 27/04/2021
 [SEC=OFFICIAL]
Attachments: DIR 68783 - Questionnaire Letter.rtf

Dear s22

Please see the amended due date for a response is **COB 27 April 2021**.

Sorry for the confusion.

Kind regards,

s22

Device Support Team

Medical Device Incident Report Investigation Scheme (IRIS)

Devices Post Market Monitoring | Medical Devices Surveillance Branch
 Therapeutic Goods Administration (TGA)
 Department of Health

Email: IRIS@health.gov.au

Phone: DST s22

Phone: IRIS 1800 809 361

Post: PO Box 100, Woden, ACT 2606

Courier: 136 Narrabundah Lane, Symonston, ACT 2609

Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission. This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met

Please Note: Medical device sponsors and manufacturers are encouraged to submit adverse events via the online Medical Device Incident Reporting (MDIR) system. The MDIR allows sponsors and manufacturers to submit initial, follow-up and final reports and review reports already submitted to the TGA. The MDIR can be accessed via <https://apps.tga.gov.au/prod/mdir/MDIRSummary.aspx>, using a TGA eBusiness Service (eBS) user name and password. MDIR guidance documents and FAQ are available from the 'Training' section of the TGA eBS portal, or from the [TGA website](#).

MDIR Technical assistance is available via s22 or email (IRIS@health.gov.au). If you are a sponsor or manufacturer requiring a new account to gain access to the MDIR system or need to update your details with TGA eBS, please contact TBS Helpdesk via 1800 010 624 or email (eBS@health.gov.au).

Online reporting forms for Medical Device Healthcare Professionals/Users can be accessed via:

<https://apps.tga.gov.au/prod/mdir/udir03.aspx>

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

From: IRIS

Sent: Friday, 26 March 2021 12:30 PM

To: s22 [REDACTED]@medacta.com.au' s22 [REDACTED]@medacta.com.au>

Subject: DIR 68783 - Questionnaire Letter - Response due COB 27/03/2021 [SEC=OFFICIAL]

Dear s22 [REDACTED]

Please see the attached **Questionnaire Letter** for action.

***Please note the due date for a response is COB 27 March 2021.**

Kind regards,

s22 [REDACTED]

Device Support Team

Medical Device Incident Report Investigation Scheme (IRIS)

Devices Post Market Monitoring | Medical Devices Surveillance Branch

Therapeutic Goods Administration (TGA)

Department of Health

Email: IRIS@health.gov.au

Phone: DST s22 [REDACTED]

Phone: IRIS 1800 809 361

Post: PO Box 100, Woden, ACT 2606

Courier: 136 Narrabundah Lane, Symonston, ACT 2609

***Important:** This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission. This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met*

Please Note: Medical device sponsors and manufacturers are encouraged to submit adverse events via the online Medical Device Incident Reporting (MDIR) system. The MDIR allows sponsors and manufacturers to submit initial, follow-up and final reports and review reports already submitted to the TGA. The MDIR can be accessed via <https://apps.tga.gov.au/prod/mdir/MDIRSummary.aspx>, using a TGA eBusiness Service (eBS) user name and password. MDIR guidance documents and FAQ are available from the 'Training' section of the TGA eBS portal, or from the [TGA website](#).

MDIR Technical assistance is available via s22 [REDACTED] or email (IRIS@health.gov.au). If you are a sponsor or manufacturer requiring a new account to gain access to the MDIR system or need to update your details with TGA eBS, please contact TBS Helpdesk via 1800 010 624 or email (eBS@health.gov.au).

Online reporting forms for Medical Device Healthcare Professionals/Users can be accessed via:

<https://apps.tga.gov.au/prod/mdir/udir03.aspx>

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

Australian Government

Department of Health
Therapeutic Goods Administration

**Australian Medical Device
Incident Report Investigation Scheme**

s22

Medacta Australia Pty Ltd
Unit A1 / 16 Mars Road
LANE COVE NSW 2066
Email: s22@medacta.com.au

File Reference: E21-259665

Sent by email

Dear s22

**DEVICE INCIDENT REPORT DIR 68783 - ARTG # 330998 - Anchor, soft-tissue implanted,
non-biodegradable**

The Therapeutic Goods Administration has been advised of an incident involving the above product. A copy of the Device Incident Report (DIR) is attached.

To assist in the evaluation and resolution of this report, please provide the information requested in the attached questionnaire and return it to this office **within 20 working days of the date of this letter**, and no later than COB 27/04/2021.

Please send responses via email to IRIS@health.gov.au, referencing the DIR number.

If you are unable to respond with all the information requested by the due date please advise, **within the 20 days**, when a full response will be provided. Extensions of a reasonable time frame will be accepted depending on the seriousness of the complaint and the time requested.

Thank you for your cooperation. If you require further information please contact me on

s22 or email: IRIS@health.gov.au.

Yours sincerely
Signed electronically by

Administration Officer

Incident Report Investigation Scheme
Devices Post Market Monitoring Section
Therapeutic Goods Administration

26/03/2021

**- MEDICAL DEVICE INCIDENT REPORT INVESTIGATION SCHEME -
INITIAL REQUEST OF INFORMATION FROM LISTED SPONSOR**

Document 7

INITIAL REQUEST FOR INFORMATION FROM LISTED SPONSOR

Date: 26/03/2021

DIR: 68783

Manufacturer Name: Medacta International SA [45296]

Question/Requirement

1) Please confirm the device's Australian Register of Therapeutic Goods (ARTG) number

ARTG:

2) Do you currently supply or have you previously supplied this product, with the indicated Model/Serial/Batch/Lot numbers:

- a) To the Australian Market
- b) For Export

YES

NO

☐
☐☐
☐

3) How many of this model have been supplied

In Australia:

Worldwide:

4) How many of this batch (if applicable) have been supplied:

In Australia:

Worldwide:

5) Are you aware of this problem, as reported?

☐☐

6) If deemed necessary, is a sample of the mentioned device available for review and/or testing?

☐☐

7) Have you had any other reports of similar* problems with this product in Australia?

☐☐

If **YES**, how many:

If **YES**, please give details:

* Similar events are based on the clinical event description and not the cause of an event. Both the confirmed and unconfirmed rates for similar events are often required and beneficial to show the full outlook of the device and its use.

8) Is the manufacturer aware of reports of similar problems with this product?

☐

Document 7

☐

If **YES**, how many:

If **YES**, please give details:

9) Please provide details of any action you have taken, or intend to take, regarding this problem

10) Please provide details of the manufacturer's investigation to date, including expected **Manufacturer's investigation completion date**

Date:

11) When returning this response to the office of the Therapeutic Goods Administration, you are requested to attach the following (if checked):

☐

Sample of the product/device

☐

Operator's manual

☐

Product Specifications

☐

Technical Service Manual

☐

Descriptive product promotional documentation

☐

Clinical training manual in printed or video form

☐

Instructions for use, as supplied with the device

☐

In-house training documentation

☐

Device Packaging with printed instructions

☐

Evidence of compliance with the Essential Principles

☐

A summary of risk assessment activities performed by the manufacturer for the device, eg Risk Management Report required by Clause 8 of ISO 14971:2007.

☐

Service History, and Safety and Performance Test Results

☐

Configuration Information and Documentation

☐

HHE (Health Hazard Evaluation)

12) Additional Information required:

Document 7

13) If your device is an implantable pacemaker/defibrillator and/or associated lead you are asked to provide the following additional information:

1. Both published and unpublished clinical trial data where events of this type are analysed.
2. The number of reported events of ALL types (including unconfirmed events), the number of devices sold and the cumulative implant months for each device in this product family.
3. What material has been used to insulate, both internally and externally, the lead? (applicable to leads only)

Information Supplied By:

14)

Name		Phone	
Signature		Fax	
Position		Email	

This questionnaire and any appended documents should be returned to the TGA within **20 working days or as specified on page 1 of this letter.**

Electronic submission of all information is preferred.

Please send the requested information to email address: IRIS@health.gov.au

Document 7

For large size documents, please post a universal serial bus (USB), compact disc (CD), or digital versatile disc (DVD) via postal address:

Incident Report Investigation Scheme
Devices Post Market Monitoring Section
Therapeutic Goods Administration
PO Box 100
Woden ACT 2606

If you are sending a device/s to the TGA please follow the instructions via link:

<https://www.tga.gov.au/form/report-medical-device-adverse-event-medical-device-use>

Sponsors of products listed or registered on the Australian Register of Therapeutic Goods (ARTG) are reminded of their responsibilities under Section 31 and/or 41JA (as appropriate) of the Therapeutic Goods Act of 1989, to provide information relating to their product's formulation, composition, design specification, quality, method and place of manufacture, presentation, safety and efficacy, conformity to advertising regulations under the Act,

Reporter Reference #: s22 ID 69823

Date of Adverse Event:
s22Date of Final Report:
18/03/2021

ARTG #:	Brand Name:	
330998	Medacta J-Lock PEEK Suture - Anchor, soft-tissue implanted, non-biodegradable	
Device Class:	Model #:	Serial #:
Class IIb		
Software Version:	Batch #:	Lot #:
	Ref 05.11.001	1954843

Manufacturer:
Medacta International SA [45296]Sponsor:
Medacta Australia Pty Ltd [45972]
Unit A1 / 16 Mars Road
LANE COVE NSW 2066

Contact Name: s22

Phone: s22

Email: s22@medacta.com.au

Fax:

Reporter:

Confidential: Yes

Date of Implant:
s22

Date of Explant:

Clinical Event Information:

Patient was admitted into surgery for right hip arthroscopy and femoral osteotomy and labral repair on the [redacted], Anchor introducer found retained in patient on surgeon appointment.

Anchor introducer snapped off during surgery and not found until patient went for appointment with surgeon. Additional surgery required for removal of retained anchor introducer on s22

Patient Outcome/Consequences:

Additional surgery required for removal of retained anchor introducer on s22

Additional Event Description:

Anchor introducer snapped off during surgery and not found until patient went for appointment with surgeon.

Device Analysis Results:

Corrective/Preventative Actions:

Details of Similar Events:

Number of Similar Events:

Rate of Similar Events:

Countries Similar Events Also Occurred:

s22

From: s22@medacta.com.au>
Sent: Friday, 26 March 2021 12:52 PM
To: IRIS
Subject: Re: *Amended* DIR 68783 - Questionnaire Letter - Response due COB 27/04/2021 [SEC=OFFICIAL]

REMINDER: Think before you click! This email originated from outside our organisation. Only click links or open attachments if you recognise the sender and know the content is safe.

Thank you for your correspondence. We will attend to this as soon as possible

s22
s22
Medacta Australia
s22

On 26 Mar 2021, at 12:33, IRIS <IRIS@health.gov.au> wrote:

Dear s22

Please see the amended due date for a response is **COB 27 April 2021**.

Sorry for the confusion.

Kind regards,

s22

Device Support Team

Medical Device Incident Report Investigation Scheme (IRIS)

Devices Post Market Monitoring | Medical Devices Surveillance Branch
Therapeutic Goods Administration (TGA)
Department of Health

Email: IRIS@health.gov.au

Phone: DST s22

Phone: IRIS 1800 809 361

Post: PO Box 100, Woden, ACT 2606

Courier: 136 Narrabundah Lane, Symonston, ACT 2609

<image001.gif>

Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission. This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met

Please Note: Medical device sponsors and manufacturers are encouraged to submit adverse events via the online Medical Device Incident Reporting (MDIR) system. The MDIR allows sponsors and manufacturers to submit initial, follow-up and final reports and review reports already submitted to the TGA. The MDIR can be accessed via

<https://apps.tga.gov.au/prod/mdir/MDIRSummary.aspx>, using a TGA eBusiness Service (eBS) user name and password. MDIR guidance documents and FAQ are available from the 'Training' section of the TGA eBS portal, or from the [TGA website](#).

MDIR Technical assistance is available via **s22** or email (IRIS@health.gov.au). If you are a sponsor or manufacturer requiring a new account to gain access to the MDIR system or need to update your details with TGA eBS, please contact TBS Helpdesk via 1800 010 624 or email (eBS@health.gov.au).

Online reporting forms for Medical Device Healthcare Professionals/Users can be accessed via:

<https://apps.tga.gov.au/prod/mdir/udir03.aspx>

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

From: IRIS

Sent: Friday, 26 March 2021 12:30 PM

To: **s22** @medacta.com.au' **s22** @medacta.com.au>

Subject: DIR 68783 - Questionnaire Letter - Response due COB 27/03/2021 [SEC=OFFICIAL]

Dear **s22**

Please see the attached **Questionnaire Letter** for action.

***Please note the due date for a response is COB 27 March 2021.**

Kind regards,

s22

Device Support Team

Medical Device Incident Report Investigation Scheme (IRIS)

Devices Post Market Monitoring | Medical Devices Surveillance Branch

Therapeutic Goods Administration (TGA)

Department of Health

Email: IRIS@health.gov.au

Phone: **DST** **s22**

Phone: **IRIS** 1800 809 361

Post: PO Box 100, Woden, ACT 2606

Courier: 136 Narrabundah Lane, Symonston, ACT 2609

<image001.gif>

Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission. This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met

Please Note: Medical device sponsors and manufacturers are encouraged to submit adverse events via the online Medical Device Incident Reporting (MDIR) system. The MDIR allows sponsors and manufacturers to submit initial, follow-up and final reports and review reports already submitted to the TGA. The MDIR can be accessed via

<https://apps.tga.gov.au/prod/mdir/MDIRSummary.aspx>, using a TGA eBusiness Service (eBS) user name and password. MDIR guidance documents and FAQ are available from the 'Training' section of the TGA eBS portal, or from the [TGA website](#).

MDIR Technical assistance is available via **s22** or email (IRIS@health.gov.au). If you are a sponsor or manufacturer requiring a new account to gain access to the MDIR system or need to update your details with TGA eBS, please contact TBS Helpdesk via 1800 010 624 or email (eBS@health.gov.au).

Online reporting forms for Medical Device Healthcare Professionals/Users can be accessed via:

<https://apps.tga.gov.au/prod/mdir/udir03.aspx>

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

"Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission."

Message protected by MailGuard: e-mail anti-virus, anti-spam and content filtering.

<http://www.mailguard.com.au>

[Report this message as spam](#)

<DIR 68783 - Questionnaire Letter.rtf>

Message protected by MailGuard: e-mail anti-virus, anti-spam and content filtering.

<http://www.mailguard.com.au>

s22

From: s22@medacta.ch>
Sent: Monday, 26 April 2021 6:40 PM
To: IRIS
Cc: s22; s22; s22
Subject: R: *Amended* DIR 68783 - Questionnaire Letter - Response due COB 27/04/2021 [SEC=OFFICIAL]
Attachments: DIR 68783 - Questionnaire Letter signed.pdf; ANNEX 1 (Analysis).pdf

REMINDER: Think before you click! This email originated from outside our organisation. Only click links or open attachments if you recognise the sender and know the content is safe.

Dear all,
 as per your request you will find attached the Questionnaire Letter related to the DIR 68783.
 The complaint was registered as Medacta number 0172-21 and already communicated to the TGA with the DIR 68516.
 The item involved in the event is available and it will be sent to your attention as soon as possible.
 Do not hesitate to contact me if you have any issues,
 Kind regards

s22 s22 <i>Regulatory Affairs Department</i> s22@medacta.ch www.medacta.ch	Medacta International SA Strada Regina CH-6874 Castel San Pietro Switzerland Pho: s22 Mob: s22 Fax: +41 91 696 60 66
---	--

From: IRIS <IRIS@health.gov.au>
Date: 26 March 2021 at 12:33:06 AEDT
To: s22@medacta.com.au>
Subject: *Amended* DIR 68783 - Questionnaire Letter - Response due COB 27/04/2021 [SEC=OFFICIAL]

Dear s22
 Please see the amended due date for a response is **COB 27 April 2021**.
 Sorry for the confusion.
 Kind regards,
 s22
Device Support Team
Medical Device Incident Report Investigation Scheme (IRIS)
 Devices Post Market Monitoring | Medical Devices Surveillance Branch
 Therapeutic Goods Administration (TGA)
 Department of Health
 Email: IRIS@health.gov.au
 Phone: DST s22
 Phone: IRIS 1800 809 361
 Post: PO Box 100, Woden, ACT 2606
 Courier: 136 Narrabundah Lane, Symonston, ACT 2609

Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission. This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met

Please Note: Medical device sponsors and manufacturers are encouraged to submit adverse events via the online Medical Device Incident Reporting (MDIR) system. The MDIR allows sponsors and manufacturers to submit initial, follow-up and final reports and review reports already submitted to the TGA. The MDIR can be accessed via <https://apps.tga.gov.au/prod/mdir/MDIRSummary.aspx>, using a TGA eBusiness Service (eBS) user name and password. MDIR guidance documents and FAQ are available from the Training section of the TGA eBS portal, or from the [TGA website](#).

MDIR Technical assistance is available via [s22](#) or email (IRIS@health.gov.au). If you are a sponsor or manufacturer requiring a new account to gain access to the MDIR system or need to update your details with TGA eBS, please contact TBS Helpdesk via 1800 010 624 or email (eBS@health.gov.au).

Online reporting forms for Medical Device Healthcare Professionals/Users can be accessed via:

<https://apps.tga.gov.au/prod/mdir/udir03.aspx>

For ongoing information and updates please subscribe to the TGAs [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

From: IRIS

Sent: Friday, 26 March 2021 12:30 PM

To: [s22](#) @medacta.com.au' [s22](#) @medacta.com.au>

Subject: DIR 68783 - Questionnaire Letter - Response due COB 27/03/2021 [SEC=OFFICIAL]

Dear [s22](#)

Please see the attached **Questionnaire Letter** for action.

***Please note the due date for a response is COB 27 March 2021.**

Kind regards,

[s22](#)

Device Support Team

Medical Device Incident Report Investigation Scheme (IRIS)

Devices Post Market Monitoring | Medical Devices Surveillance Branch

Therapeutic Goods Administration (TGA)

Department of Health

Email: IRIS@health.gov.au

Phone: **DST** [s22](#)

Phone: **IRIS** 1800 809 361

Post: PO Box 100, Woden, ACT 2606

Courier: 136 Narrabundah Lane, Symonston, ACT 2609

Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission. This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met

Please Note: Medical device sponsors and manufacturers are encouraged to submit adverse events via the online Medical Device Incident Reporting (MDIR) system. The MDIR allows sponsors and manufacturers to submit initial, follow-up and final reports and review reports already submitted to the TGA. The MDIR can be accessed via <https://apps.tga.gov.au/prod/mdir/MDIRSummary.aspx>, using a TGA eBusiness Service (eBS) user name and password. MDIR guidance documents and FAQ are available from the Training section of the TGA eBS portal, or from the [TGA website](#).

MDIR Technical assistance is available via [s22](#) or email (IRIS@health.gov.au). If you are a sponsor or manufacturer requiring a new account to gain access to the MDIR system or need to update your details with TGA eBS, please contact TBS Helpdesk via 1800 010 624 or email (eBS@health.gov.au).

Online reporting forms for Medical Device Healthcare Professionals/Users can be accessed via:

<https://apps.tga.gov.au/prod/mdir/udir03.aspx>

For ongoing information and updates please subscribe to the TGAs [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

"Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission."

Message protected by MailGuard: e-mail anti-virus, anti-spam and content filtering.
<http://www.mailguard.com.au>

[Report this message as spam](#)

Message protected by MailGuard: e-mail anti-virus, anti-spam and content filtering.
<http://www.mailguard.com.au>

**- MEDICAL DEVICE INCIDENT REPORT INVESTIGATION SCHEME -
INITIAL REQUEST OF INFORMATION FROM LISTED SPONSOR**

Document 9

INITIAL REQUEST FOR INFORMATION FROM LISTED SPONSOR

Date: 26/04/2021

DIR: 68783

Manufacturer Name: Medacta International SA [45296]

Question/Requirement

1) Please confirm the device's Australian Register of Therapeutic Goods (ARTG) number

ARTG: 320484

2) Do you currently supply or have you previously supplied this product, with the indicated Model/Serial/Batch/Lot numbers:

- a) To the Australian Market
- b) For Export

YES

NO

☒

☐

☒

☐

3) How many of this model have been supplied

In Australia:

s47

Worldwide:

s47

4) How many of this batch (if applicable) have been supplied:

In Australia:

s47

Worldwide:

s47

5) Are you aware of this problem, as reported?

☒

☐

6) If deemed necessary, is a sample of the mentioned device available for review and/or testing?

☒

☐

7) Have you had any other reports of similar* problems with this product in Australia?

☐

☒

If **YES**, how many:

If **YES**, please give details:

* Similar events are based on the clinical event description and not the cause of an event. Both the confirmed and unconfirmed rates for similar events are often required and beneficial to show the full outlook of the device and its use.

8) Is the manufacturer aware of reports of similar problems with this product?

☐

Document 9

☒

If **YES**, how many:

If **YES**, please give details:

9) Please provide details of any action you have taken, or intend to take, regarding this problem

The event is an isolated case and no systematic action is currently indicated as applicable, however it has been included in our active monitoring system for potential CAPA needs and opportunities.

10) Please provide details of the manufacturer's investigation to date, including expected **Manufacturer's investigation completion date**

Date: 19-May-2021

See the annex 1 attached (Analysis)

11) When returning this response to the office of the Therapeutic Goods Administration, you are requested to attach the following (if checked):

☒

Sample of the product/device

☐

Operator's manual

☐

Product Specifications

☐

Technical Service Manual

☐

Descriptive product promotional documentation

☐

Clinical training manual in printed or video form

☐

Instructions for use, as supplied with the device

☐

In-house training documentation

☐

Device Packaging with printed instructions

☐

Evidence of compliance with the Essential Principles

☐

A summary of risk assessment activities performed by the manufacturer for the device, eg Risk Management Report required by Clause 8 of ISO 14971:2007.

☐

Service History, and Safety and Performance Test Results

☐

Configuration Information and Documentation

☐

HHE (Health Hazard Evaluation)

12) Additional Information required:

Document 9

The complaint was registered as Medacta number 0172-21 and communicated to the TGA with the DIR 68516.

Comparing the data received from you with the DIR 68783, you will find here below the incongruence we have found:

Correct data:

ARTG number: 320484

Primary surgery date s22

Lot: 1924147

Correct data:

ARTG number: 330998

Primary surgery date s22

Lot: 1954843

The ARTG number specified in the document (330998) is not correct because related to all-suture anchors while the event is related to PEEK anchors (ARTG n° 320484).

13) If your device is an implantable pacemaker / defibrillator and/or associated lead you are asked to provide the following additional information:

1. Both published and unpublished clinical trial data where events of this type are analysed.
2. The number of reported events of ALL types (including unconfirmed events), the number of devices sold and the cumulative implant months for each device in this product family.
3. What material has been used to insulate, both internally and externally, the lead? (applicable to leads only)

Information Supplied By:

14)

Name s22

Phone s22

Signature

Fax s22

Position s22

Email s22@medacta.ch

This questionnaire and any appended documents should be returned to the TGA within **20 working days or as specified on page 1 of this letter.**

Electronic submission of all information is preferred.

Please send the requested information to email address: IRIS@health.gov.au

Document 9

For large size documents, please post a universal serial bus (USB), compact disc (CD), or digital versatile disc (DVD) via postal address:

Incident Report Investigation Scheme
Devices Post Market Monitoring Section
Therapeutic Goods Administration
PO Box 100
Woden ACT 2606

If you are sending a device/s to the TGA please follow the instructions via link:

<https://www.tga.gov.au/form/report-medical-device-adverse-event-medical-device-use>

Sponsors of products listed or registered on the Australian Register of Therapeutic Goods (ARTG) are reminded of their responsibilities under Section 31 and/or 41JA (as appropriate) of the Therapeutic Goods Act of 1989, to provide information relating to their product's formulation, composition, design specification, quality, method and place of manufacture, presentation, safety and efficacy, conformity to advertising regulations under the Act,

Reporter Reference #: s22 ID 69823

Date of Adverse Event:

s22

Date of Final Report:

18/03/2021

ARTG #:	Brand Name:	
330998	Medacta J-Lock PEEK Suture - Anchor, soft-tissue implanted, non-biodegradable	
Device Class:	Model #:	Serial #:
Class IIb		
Software Version:	Batch #:	Lot #:
	Ref 05.11.001	1954843

Manufacturer:

Medacta International SA [45296]

Sponsor:

Medacta Australia Pty Ltd [45972]

Unit A1 / 16 Mars Road

LANE COVE NSW 2066

Contact Name: s22

Phone: s22

Email: s22@medacta.com.au

Fax:

Reporter:

Confidential: Yes

Date of Implant:

s22

Date of Explant:

Clinical Event Information:

Patient was admitted into surgery for right hip arthroscopy and femoral osteotomy and labral repair on the [redacted], Anchor introducer found retained in patient on surgeon appointment.

Anchor introducer snapped off during surgery and not found until patient went for appointment with surgeon. Additional surgery required for removal of retained anchor introducer on s22

Patient Outcome/Consequences:

Additional surgery required for removal of retained anchor introducer on s22

Additional Event Description:

Anchor introducer snapped off during surgery and not found until patient went for appointment with surgeon.

Device Analysis Results:

Corrective/Preventative Actions:

Details of Similar Events:

Number of Similar Events:

Rate of Similar Events:

Countries Similar Events Also Occurred:

ANNEX 1 (Analysis)**Clinical Evaluation performed by S22 :**

During an arthroscopic surgery of labrum refixation with a PEEK anchor, the tip of the positioning instrument broke, though the surgeon did not perceive the fracture. The tip of the instrument is visible in the postoperative Xrays, it remained in the cannulation of the anchor. We cannot determine the reasons for the fracture with the elements at hand. No clinical cause could be ascertained. We were informed that a different surgeon intervened to remove the broken instrument part; the future clinical pathway of the patient should be totally flawless and no further consequences of the event should be expected.

Visual Inspection performed by S22 :

The component of the MectaLock PEEK inserter has been received and analysed. It's evident that the laser welding failed and the tip component has been detached from the main shaft of the inserter. The proximal portion of the analysed component resulted totally ruined and scratched, most likely due to the extraction, therefore is not possible to perform a detailed microscopic analysis of the failure mode of the welding. However, it's also evident that the component has been subjected to excessive lateral bending as the axis of the component results deviated from the original. This suggests that a misuse occurred during the anchor's placement as the MectaLock PEEK anchors is supposed to be impacted axially only within the predrilled cylindrical pilot hole.

Batch review performed on 23-Feb-2021

Lot 1924147: 125 items manufactured and released on 24-Mar-2020. Expiration date: 2025-Feb-07. No anomalies found related to the problem. To date, 19 items of the same lot have been already sold without any other similar reported event

s22

From: s22@Health.gov.au>
Sent: Wednesday, 5 May 2021 6:00 PM
To: s22@svha.org.au
Subject: Sample received DIR 68783 [SEC=OFFICIAL]

Dear MS s22,

Thank you for the sample you have sent to the TGA in regards to your Device Incident report (DIR) 68783 for the (J-Lock PEEK Suture - Anchor, soft-tissue implanted, non-biodegradable) . The TGA has limited testing capabilities upon receipt of a device and the best way to have a device tested is to send it back to the manufacturer for analysis. The TGA is seeking your permission to send the device onto the manufacturer for extensive testing to assist us with the investigation of the DIR. Without this analysis a root cause would not be able to be determined.

If you do not wish for us to send it to the sponsor for this analysis we are able to send it back to you should you wish to seek independent testing. Again, thanking you.

Kind regards,

Dr s22
 Departmental Officer
 Devices Post Market Monitoring
 Medical Devices Surveillance Branch
 Email: s22@health.gov.au
 Telephone number : s22

Therapeutic Goods Administration

Department of Health
 PO Box 100
 Woden ACT 2606 Australia

www.tga.gov.au

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met.

Important: *This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission.*

The Department of Health acknowledges the Traditional Custodians of Australia and their continued connection to land, sea and community. We pay our respects to them and their culture, and to all Elders both past and present.

P *Please consider the environment before printing this email.*

s22

From: s22@Health.gov.au>
Sent: Wednesday, 12 May 2021 5:06 PM
To: s22
Subject: RE: DIR 68783 [SEC=OFFICIAL]

Hi s22

Thanks a lot for getting back to me. I think we don't need to take any photos of the device and thanks for suggesting though. Please, let me know when you have sent samples and I will send an email to the sponsor letting know about that. Again, thanking you.

Kind regards,

Dr s22
 Departmental Officer
 Devices Post Market Monitoring
 Medical Devices Surveillance Branch
 Email: s22@health.gov.au
 Telephone number : s22

Therapeutic Goods Administration

Department of Health
 PO Box 100
 Woden ACT 2606 Australia

www.tga.gov.au

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met.

Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission.

The Department of Health acknowledges the Traditional Custodians of Australia and their continued connection to land, sea and community. We pay our respects to them and their culture, and to all Elders both past and present.

P Please consider the environment before printing this email.

From: s22@health.gov.au>
Sent: Wednesday, 12 May 2021 1:21 PM
To: s22@Health.gov.au>
Subject: RE: DIR 68783 [SEC=OFFICIAL]

Hi s22

No Problem we can send it on.
 Would you like us to take any photos of the device.
 If not I will organise for it to out tomorrow.

Happy to do whatever is required.

Kind regards

s22

From: s22@Health.gov.au>
Sent: Wednesday, 12 May 2021 12:38 PM
To: s22@health.gov.au>
Subject: DIR 68783 [SEC=OFFICIAL]

Dear s22

Good afternoon. We recently received the sample from the reporter (s22) regarding DIR 68783 and the sample has been passed to the lab team as it is biohazard. After getting the permission from the reporter, we would like to send the sample to the sponsor as a part of our investigation process.

The address details of the sponsor is given below:

Name: s22
Sponsor's name: Medacta Australia Pty Ltd
Unit A1/ 16 Mars Road LANE COVE
NSW 2066
Contact number : s22

I would really appreciate if the sample is sent to the sponsor by the lab. Please, let me know if you need any further information regarding this. Again, thanking you.

Kind regards,

Dr s22
Departmental Officer
Devices Post Market Monitoring
Medical Devices Surveillance Branch
Email: s22@health.gov.au
Telephone number : s22

Therapeutic Goods Administration

Department of Health
PO Box 100
Woden ACT 2606 Australia

www.tga.gov.au

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met.

Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission.

The Department of Health acknowledges the Traditional Custodians of Australia and their continued connection to land, sea and community. We pay our respects to them and their culture, and to all Elders both past and present.

P Please consider the environment before printing this email.

s22

From: s22@medacta.com.au>
Sent: Tuesday, 22 June 2021 11:37 AM
To: s22
Cc: s22; s22
Subject: RE: DIR 68783 [SEC=OFFICIAL]
Attachments: final_report_0172-21.pdf

REMINDER: Think before you click! This email originated from outside our organisation. Only click links or open attachments if you recognise the sender and know the content is safe.

Dear Dr s22

Please find attached the internal analysis your require in relation to DIR 68783.

Please do not hesitate to contact me if I can be of further assistance

Yours faithfully

s22

s22

s22

s22@medacta.com.au
www.medacta.com

Medacta Australia
 Unit A1, 16 Mars Road
 Lane Cove, NSW 2066
 Australia
 Pho: s22
 Mob: s22
 Fax: +612 9420 2578

IMPORTANT NOTICE

The content of this email is confidential and is for exclusive use of the intended recipient(s). If you are not the intended recipient(s) please note that any form of distribution, copying or use of this communication or the information in it is strictly forbidden and may be unlawful. Please note that, although this email has been virus-checked, neither Medacta International SA nor the sender accepts any responsibility for viruses, and it is your responsibility to scan the email and attachments (if any). If you believe you have received this communication in error, please return it (with the Object "received in error") to the sender, and then delete the email and destroy any copies of it.

We would be grateful if you would also copy the communication to customerservice@medacta.com.au

 Please consider the environment before printing this email

From: s22@Health.gov.au>
Sent: Friday, 11 June 2021 3:43 PM
To: s22@medacta.com.au>
Subject: RE: DIR 68783 [SEC=OFFICIAL]

Hi s22

Good afternoon. We contacted you previously to provide us the detailed device analysis of the specimen for DIR 68783 (J-Lock PEEK Suture - Anchor, soft-tissue implanted, non-biodegradable). Please provide us the detailed device analysis of the specimen received from s22 by COB 2 July, 2021. Again, thanking you.

Kind regards,

Dr s22

Departmental Officer

Devices Post Market Monitoring

Medical Devices Surveillance Branch

Email: s22@health.gov.au

Telephone number : s22

Therapeutic Goods Administration

Department of Health

PO Box 100

Woden ACT 2606 Australia

www.tga.gov.au

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met.

Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission.

The Department of Health acknowledges the Traditional Custodians of Australia and their continued connection to land, sea and community. We pay our respects to them and their culture, and to all Elders both past and present.

P Please consider the environment before printing this email.

From: s22

Sent: Tuesday, 25 May 2021 9:46 AM

To: s22@medacta.com.au' s22@medacta.com.au>

Subject: DIR 68783 [SEC=OFFICIAL]

Dear s22,

Good morning. We recently received the sample regarding DIR 68783 (J-Lock PEEK Suture - Anchor, soft-tissue implanted, non-biodegradable) from the reporter (s22). The sample was sent by TGA Lab after seeking permission from the reporter on last week by courier. Could you please provide a date analysis of this specimen for our review which is a part of our investigation ? Again, thanking you.

Kind regards,

Dr s22

Departmental Officer

Devices Post Market Monitoring

Medical Devices Surveillance Branch

Email: s22@health.gov.au

Telephone number : s22

Therapeutic Goods Administration

Department of Health

PO Box 100

Woden ACT 2606 Australia

www.tga.gov.au

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met.

Important: *This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission.*

The Department of Health acknowledges the Traditional Custodians of Australia and their continued connection to land, sea and community. We pay our respects to them and their culture, and to all Elders both past and present.

P *Please consider the environment before printing this email.*

"Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission."

Message protected by MailGuard: e-mail anti-virus, anti-spam and content filtering.

<http://www.mailguard.com.au>

[Report this message as spam](#)

Message protected by MailGuard: e-mail anti-virus, anti-spam and content filtering.

<http://www.mailguard.com.au>

Confidential – Do not forward to third parties without Medacta International SA's authorization

Complaint nr:	0	1	7	2	-	2	1
	COMPLAINT NUMBER				YEAR		

1. EVENT NOTIFICATION

Reported by	s22	Company name	MEDACTA AUSTRALIA Pty. Ltd
Country	Australia	Date of complaint	January 29, 2021
Surgeon	Mr s22	Hospital	s22
Surgeon of the primary surgery	Mr s22	Date of primary surgery	s22
Surgeon of the revision surgery	Mr s22	Date of revision surgery	s22

2. DEVICES INVOLVED

Product Family	Reference n.	Product Description	Lot n.	SN	Qty
MectaLock instruments	05.10.001	MectaLock PEEK Ø2.4 – Short driver	1924147	NA	1

3. EVENT DESCRIPTION

The patient underwent a labral repair (s22), Hip Arthroscopy. Two MectaLock (J-Lock) PEEK anchors were used for refixing the labrum to the acetabulum with no visual issues. Then the surgeon sent to the agent a picture of an X-ray asking if the image in the X-ray was part of the anchor inserter handle. It appears that a part of the inserter handle has remained in the bone with the anchor. The patient underwent revision surgery to remove the broken part.

The event was reported to the manufacturer on 29.01.2021, but only after the revision surgery, which took place on s22, we were sure that the piece that remained in the patient was part of the Medacta device.

4. INVESTIGATION RESULTS

Clinical Evaluation performed by Medical Affairs Director:

During an arthroscopic surgery of labrum refixation with a PEEK anchor, the tip of the positioning instrument broke, though the surgeon did not perceive the fracture. The tip of the instrument is visible in the postoperative X-rays, it remained in the cannulation of the anchor. We cannot determine the reasons for the fracture with the elements at hand.

No clinical cause could be ascertained. We were informed that a different surgeon intervened to remove the broken instrument part; the future clinical pathway of the patient should be totally flawless and no further consequences of the event should be expected.

Confidential – Do not forward to third parties without Medacta International SA's authorization

Complaint nr:	0	1	7	2	-	2	1
	COMPLAINT NUMBER				YEAR		

Visual Inspection performed by R&D Project Manager:

The component of the MectaLock PEEK inserter has been received and analyzed. It's evident that the laser welding failed and the tip component has been detached from the main shaft of the inserter. The proximal portion of the analyzed component was received scratched, most likely due to the extraction, therefore is not possible to perform a detailed microscopic analysis of the failure mode of the welding. A deep analysis of the device has been carried out by means of mechanical testing of equivalent devices sourced from the same production batch as representative of the component involved in this complaint. In particular, the welding between the shaft and the inserter tip has been subjected to stress limit tests which confirmed the safety and reliability of the device under load conditions that are commonly expected during surgery. Therefore, it's not possible to identify a clear or specific root cause for the mode of failure, unless an outstanding condition required an unforeseeable usage during the surgical procedure.

Batch review performed on 23-Feb-2021

Lot 1924147: 125 items manufactured and released on 24-Mar-2020. Expiration date: 2025-Feb-07. No anomalies found related to the problem. To date, 19 items of the same lot have been sold without any similar reported event.

5. ROOT CAUSE

The device history record review does not indicate any potentially related manufacturing issue. Mechanical testing has confirmed reliability of the device under foreseeable loading conditions. It could be that anomalous conditions of application may have resulted in the breakage but no confirmation can be made in this case.

6. CORRECTIVE/PREVENTIVE ACTIONS

The event is an isolated case and no systematic action is currently indicated as applicable, however it has been included in our active monitoring system for potential CAPA needs and opportunities.

Confidential – Do not forward to third parties without Medacta International SA's authorization

Complaint nr:	0	1	7	2	-	2	1
	COMPLAINT NUMBER				YEAR		

	Function	Name	Signature	Date
Prepared by:	s22	s22	Electronic legal signature	June 01, 2021
Completed by:	s22	s22	Electronic legal signature	June 02, 2021
Verified by: (only if needed)	s22	s22	Electronic legal signature	June 02, 2021

s22

From: s22@medacta.ch>
Sent: Thursday, 3 June 2021 1:16 AM
To: IRIS
Cc: s22; s22
Subject: R: *Amended* DIR 68783 - Questionnaire Letter - Response due COB 27/04/2021 [SEC=OFFICIAL]
Attachments: Internal analysis - Complaint 0172-21 DIR 68516.pdf

REMINDER: Think before you click! This email originated from outside our organisation. Only click links or open attachments if you recognise the sender and know the content is safe.

Dear all,

considering the Medacta complaint number **0172-21**, already reported to the TGA with the **DIR 68516**, we have performed some additional tests that you will find summarized in the attachment. Given that the final report (**DIR 68516**) was already sent to the TGA I kindly ask you to update the Visual Inspection previously reported as follows:

The component of the MectaLock PEEK inserter has been received and analysed. It's evident that the laser welding failed and the tip component has been detached from the main shaft of the inserter. The proximal portion of the analysed component was received scratched, most likely due to the extraction, therefore is not possible to perform a detailed microscopic analysis of the failure mode of the welding.

A deeper analysis on the safety of the device has been carried out by means of further mechanical investigation on equivalent devices sourced from the same production batch, as representative of the component object of this complaint. In particular, the welding between the shaft and the inserter tip has been subject to stress limit tests which confirmed the safety and reliability of the device under load conditions that are commonly expected during a surgery. Therefore, it's not possible to identify a clear or specific root cause for the mode of failure, unless an outstanding condition required a not foreseeable usage during the surgical procedure.

In addition, as per your below request named **DIR 68783**, you should have received the item involved in the Medacta complaint number **0172-21** that was sent to your attention. Please, let us know if any additional information is requested.

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Dear all,

as per your request you will find attached the Questionnaire Letter related to the DIR 68783.

The complaint was registered as Medacta number 0172-21 and already communicated to the TGA with the DIR 68516.

The item involved in the event is available and it will be sent to your attention as soon as possible.

Do not hesitate to contact me if you have any issues,
Kind regards

s22 [REDACTED] Complaint Coordinator <i>Regulatory Affairs Department</i> s22 [REDACTED] @medacta.ch www.medacta.ch	Medacta International SA Strada Regina CH-6874 Castel San Pietro Switzerland Pho: s22 [REDACTED] Mob: s22 [REDACTED] Fax: +41 91 696 60 66
---	--

From: IRIS <IRIS@health.gov.au>

Date: 26 March 2021 at 12:33:06 AEDT

To: s22 [REDACTED]@medacta.com.au>

Subject: *Amended* DIR 68783 - Questionnaire Letter - Response due COB 27/04/2021
[SEC=OFFICIAL]

Dear s22 [REDACTED]

Please see the amended due date for a response is **COB 27 April 2021**.

Sorry for the confusion.

Kind regards,

s22 [REDACTED]

Device Support Team

Medical Device Incident Report Investigation Scheme (IRIS)

Devices Post Market Monitoring | Medical Devices Surveillance Branch
Therapeutic Goods Administration (TGA)
Department of Health

Email: IRIS@health.gov.au

Phone: **DST** s22 [REDACTED]

Phone: **IRIS** 1800 809 361

Post: PO Box 100, Woden, ACT 2606

Courier: 136 Narrabundah Lane, Symonston, ACT 2609

Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission. This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met

Please Note: Medical device sponsors and manufacturers are encouraged to submit adverse events via the online Medical Device Incident Reporting (MDIR) system. The MDIR allows sponsors and manufacturers to submit initial, follow-up and final reports and review reports already submitted to the TGA. The MDIR can be accessed via <https://apps.tga.gov.au/prod/mdir/MDIRSummary.aspx>, using a TGA eBusiness Service (eBS) user name and password. MDIR guidance documents and FAQ are available from the 'Training' section of the TGA eBS portal, or from the [TGA website](#).

MDIR Technical assistance is available via [s22](#) or email (IRIS@health.gov.au). If you are a sponsor or manufacturer requiring a new account to gain access to the MDIR system or need to update your details with TGA eBS, please contact TBS Helpdesk via 1800 010 624 or email (eBS@health.gov.au).

Online reporting forms for Medical Device Healthcare Professionals/Users can be accessed via:

<https://apps.tga.gov.au/prod/mdir/udir03.aspx>

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

From: IRIS

Sent: Friday, 26 March 2021 12:30 PM

To: [s22](#) @medacta.com.au' [s22](#) @medacta.com.au>

Subject: DIR 68783 - Questionnaire Letter - Response due COB 27/03/2021 [SEC=OFFICIAL]

Dear [s22](#)

Please see the attached **Questionnaire Letter** for action.

***Please note the due date for a response is COB 27 March 2021.**

Kind regards,

[s22](#)

Device Support Team

Medical Device Incident Report Investigation Scheme (IRIS)

Devices Post Market Monitoring | Medical Devices Surveillance Branch

Therapeutic Goods Administration (TGA)

Department of Health

Email: IRIS@health.gov.au

Phone: **DST** [s22](#)

Phone: **IRIS** 1800 809 361

Post: PO Box 100, Woden, ACT 2606

Courier: 136 Narrabundah Lane, Symonston, ACT 2609

Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission. This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met

Please Note: Medical device sponsors and manufacturers are encouraged to submit adverse events via the online Medical Device Incident Reporting (MDIR) system. The MDIR allows sponsors and manufacturers to submit initial, follow-up and final reports and review reports already submitted to the TGA. The MDIR can be accessed via <https://apps.tga.gov.au/prod/mdir/MDIRSummary.aspx>, using a TGA eBusiness Service (eBS) user name and password. MDIR guidance documents and FAQ are available from the 'Training' section of the TGA eBS portal, or from the [TGA website](#).

MDIR Technical assistance is available via [s22](#) or email (IRIS@health.gov.au). If you are a sponsor or manufacturer requiring a new account to gain access to the MDIR system or need to update your details with TGA eBS, please contact TBS Helpdesk via 1800 010 624 or email (eBS@health.gov.au).

Online reporting forms for Medical Device Healthcare Professionals/Users can be accessed via:

<https://apps.tga.gov.au/prod/mdir/udir03.aspx>

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

"Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission."

Message protected by MailGuard: e-mail anti-virus, anti-spam and content filtering.
<http://www.mailguard.com.au>

[Report this message as spam](#)

Message protected by MailGuard: e-mail anti-virus, anti-spam and content filtering.
<http://www.mailguard.com.au>

Analysis definition

The receipt of the broken component (inserter tip only) initiated an investigation procedure in place within the internal Medacta Quality System in order to assess the reliability of the parts currently on stock and to get the confirmation of the reliability of the current design of the device (fully assembled component: J-Lock inserter - Ø2.4, ref. 66.03.020).

Figure 1 - Received component


Figure 2 - Tip assembly


Figure 3 - Complete device J-Lock inserter

The analysis has been carried out within two main activities:

1. Documental inspection – Manufacturing process

The goal of this analysis is to assess the conformity to specifications of all the parts manufactured within the same batch, according to the inspection documents and the sampling plan in force.

2. Design confirmation – welding mechanical investigation

The goal of this analysis is to assess the welding resistance by means of two mechanical test setups, in order to mimic the most common loads that could occur during the device use. The two mechanical tests are:

- 3-points flexion test
- Axial instability test


Results

1. Documental inspection – Manufacturing process



The documental analysis assessed the conformity of the entire batch according to the planned control inspection in place within the existing manufacturing process (welding inspected on 100% of the parts, material specifications respected)

2. Design confirmation – welding mechanical investigation

3-points flexion test


- Pure flexion test, focused on the welding to simulate lateral load
- Mechanical worst case vs. simple bending against surrounding structures



The welding supported max loads of 54Kg, not showing physical separation between the components. This load is considerably out of the normal scale of loads values occurring during a surgery. In a simple leverage with the inserter, the surrounding structures would fail prior to reach that limit.

Axial instability test

- Axial load to simulate stress on the welding during hammering
- Mechanical worst case: constrained hard stop vs. pilot hole

The shaft supported max loads of 85Kg for short and 60Kg for long version. These values exceed a typical insertion in pure bone (no pilot hole). The curvature occurs along the central portion of the shaft, keeping the welding in a safe portion not affected by the resulting flexion.

No torsion test has been considered on the subject inserter since it's evident that in case of excessive torsion stress, within the threaded connection between the PEEK anchor and the inserter's tip, the PEEK thread would fail much before any other metal structure, leaving the inserter unaffected.

Conclusions

A deeper analysis on the safety of the device has been carried out by means of further mechanical investigation on equivalent devices sourced from the same production batch, as representative of the component object of this complaint.

The manufacturing process analysis and the mechanical tests performed provided a confirmation that the actual design can be considered suitable and safe for the expected application.

In particular, the welding between the shaft and the inserter tip has been subject to stress limit tests which confirmed the safety and reliability of the device under load conditions that are commonly expected during a surgery. Therefore, it's not possible to identify a clear or specific root cause for the mode of failure, unless an outstanding condition required a not foreseeable usage during the surgical procedure.

s22

s22

01 Jul 2021

s22

From: s22
Sent: Friday, 10 September 2021 11:44 AM
To: s22
Cc: s22 @medacta.com.au; s22
Subject: RE: s41JA request for DIR 68783 & 68516 [SEC=OFFICIAL]

Dear s22

Good morning. Thank you so much for providing a copy of IFU.

Kind regards,

Dr s22
 Departmental Officer
 Devices Post Market Monitoring
 Medical Devices Surveillance Branch
 Email: s22 @health.gov.au
 Telephone number : s22

Therapeutic Goods Administration

Department of Health
 PO Box 100
 Woden ACT 2606 Australia

www.tga.gov.au

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met.

Important: *This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission.*

The Department of Health acknowledges the Traditional Custodians of Australia and their continued connection to land, sea and community. We pay our respects to them and their culture, and to all Elders both past and present.

P *Please consider the environment before printing this email.*

From: s22 @medacta.ch>
Sent: Thursday, 9 September 2021 5:55 PM
To: s22 @Health.gov.au>
Cc: s22 @medacta.com.au; s22 @medacta.ch>
Subject: R: s41JA request for DIR 68783 & 68516 [SEC=OFFICIAL]

REMINDER: Think before you click! This email originated from outside our organisation. Only click links or open attachments if you recognise the sender and know the content is safe.

Dear s22

You will find attached the MectaLock suture anchors IFU (Ref: 75.09.045 rev. 08) in force.
 At the page number 4, "**Breakage of the MectaLock Suture Anchor can occur if:**", is reported as follows:

- Any lateral loading is applied on the anchor/insertion while inserting the MectaLock anchor

Do not hesitate to contact me if you have any issues.

Kind regards,

s22 s22 Regulatory Affairs Department s22 @medacta.ch www.medacta.ch	Medacta International SA Strada Regina CH-6874 Castel San Pietro Switzerland Pho: s22 Mob: s22 Fax: +41 91 696 60 66
--	--

Da: **s22** [@Health.gov.au](mailto:s22@Health.gov.au)>
 Inviato: giovedì 9 settembre 2021 06:15
 A: **s22** [@medacta.ch](mailto:s22@medacta.ch)>
 Cc: **s22** [@medacta.com.au](mailto:s22@medacta.com.au); **s22** [@medacta.ch](mailto:s22@medacta.ch)>
 Oggetto: RE: s41JA request for DIR 68783 & 68516 [SEC=OFFICIAL]

Dear **s22**

Good afternoon. A review of the manufacturer's investigation shows that the component has been subjected to excessive lateral bending as the axis of the component results deviated from the original. Please provide a copy of Instructions for Use (IFU) for this device and show where it is mentioned in IFU that not doing this 'excessive lateral bending' in a warning or use space. Thanking you.

Kind regards,
Dr s22
 Departmental Officer
 Devices Post Market Monitoring
 Medical Devices Surveillance Branch
 Email: **s22** [@health.gov.au](mailto:s22@health.gov.au)
 Telephone number : **s22**

Therapeutic Goods Administration

Department of Health
 PO Box 100
 Woden ACT 2606 Australia
www.tga.gov.au

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met.

Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission.

The Department of Health acknowledges the Traditional Custodians of Australia and their continued connection to land, sea and community. We pay our respects to them and their culture, and to all Elders both past and present.

P Please consider the environment before printing this email.

From: s22
 Sent: Wednesday, 1 September 2021 5:02 PM
 To: s22 @medacta.ch>
 Cc: s22 @medacta.com.au; s22 @medacta.ch>
 Subject: RE: s41JA request for DIR 68783 & 68516 [SEC=OFFICIAL]

Dear s22

Good afternoon. Thank you so much for providing the requested information. We are currently reviewing it. Thanking you.

Kind regards,

Dr s22
 Departmental Officer
 Devices Post Market Monitoring
 Medical Devices Surveillance Branch
 Email: s22 @health.gov.au
 Telephone number : s22

Therapeutic Goods Administration

Department of Health
 PO Box 100
 Woden ACT 2606 Australia

www.tga.gov.au

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met.

Important: *This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission.*

The Department of Health acknowledges the Traditional Custodians of Australia and their continued connection to land, sea and community. We pay our respects to them and their culture, and to all Elders both past and present.

P *Please consider the environment before printing this email.*

From: s22 @medacta.ch>
 Sent: Saturday, 28 August 2021 1:03 AM
 To: s22 @Health.gov.au>
 Cc: IRIS <IRIS@health.gov.au>; s22 @medacta.com.au; s22 @medacta.ch>
 Subject: R: s41JA request for DIR 68783 & 68516 [SEC=OFFICIAL]

REMINDER: Think before you click! This email originated from outside our organisation. Only click links or open attachments if you recognise the sender and know the content is safe.

Dear Dr s22

you will find here below the answers to your requests related to DIR 68783 and 68516:

1. Please provide a description of the adverse events and the rate for those events. Separate the events for Australia and the rest of the world and provide each in a percentage format. Provide a list of complaints and numbers for those similar complaints.

This is the first event reported, no additional instrument breakages have been communicated considering both MectaLock PEEK Ø2.4 short and long drivers. In total, we have sold s47 short and long drivers (MectaLock PEEK Ø2.4).

2. In your detailed device analysis response dated 22 June 2021, it states “the laser welding failed, and the tip component has been detached from the main shaft of the inserter”. Please provide details of your risk mitigation strategies in order to avoid the laser welding failure to ensure the device safety. You will find attached a dedicated rationale.

Do not hesitate to contact me if you have any issues.

Kind regards,

s22 s22 Regulatory Affairs Department s22 @medacta.ch www.medacta.ch	Medacta International SA Strada Regina CH-6874 Castel San Pietro Switzerland Pho: s22 Mob: s22 Fax: +41 91 696 60 66
---	--

From: s22 @Health.gov.au>
Sent: Friday, 30 July 2021 2:34 PM
To: s22 @medacta.com.au>
Subject: s41JA request for DIR 68783 & 68516 [SEC=OFFICIAL]

Dear s22,

Good afternoon. Please find attached section 41JA request for DIR 68783 & 68516. Please note that the requested information must be provided within 20 working days from the date of this letter and no later than **30 August 2021**.
 Thanking you.

Kind regards,
Dr s22
 Departmental Officer
 Devices Post Market Monitoring
 Medical Devices Surveillance Branch
 Email: s22 @health.gov.au
 Telephone number : s22

Therapeutic Goods Administration
 Department of Health
 PO Box 100
 Woden ACT 2606 Australia
www.tga.gov.au

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met.

Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission.

The Department of Health acknowledges the Traditional Custodians of Australia and their continued connection to land, sea and community. We pay our respects to them and their culture, and to all Elders both past and present.

P *Please consider the environment before printing this email.*

"Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission."

Message protected by MailGuard: e-mail anti-virus, anti-spam and content filtering.

<http://www.mailguard.com.au>

[Report this message as spam](#)

"Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission."

Australian Government
Department of Health
Therapeutic Goods Administration

File Reference: E21-259665
Sent by email

s22
Medacta Australia Pty Ltd
Unit A1 / 16 Mars Road
LANE COVE NSW 2066

Email: **s22** [@medacta.com](mailto:s22@medacta.com)

Dear **s22**,

DEVICE INCIDENT REPORT DIR 68783 & 68516
Notice requiring information/documents to be provided
under section 41JA of the Therapeutic Goods Act

I have made a decision under subsection 41JA(1) of the *Therapeutic Goods Act 1989* (the Act) to require you to provide information/documents about Medacta J-Lock PEEK Suture - Anchor, soft-tissue implanted, non-biodegradable. I am the delegate of the Secretary of the Department of Health for the purposes of section 41JA of the Act. The information is required for the purposes of the investigation of the above incident report.

You are requested to provide the information/documents listed below including details of any action you have taken in relation to this matter:

1. Please provide a description of the adverse events and the rate for those events. Separate the events for Australia and the rest of the world and provide each in a percentage format. Provide a list of complaints and numbers for those similar complaints.
2. In your detailed device analysis response dated 22 June 2021, it states "the laser welding failed, and the tip component has been detached from the main shaft of the inserter". Please provide details of your risk mitigation strategies in order to avoid the laser welding failure to ensure the device safety.

The requested information must be provided within **20** working days from the date of this letter and no later than **30 August 2021**.

Address for sending the requested information

Please submit all documents in electronic format only. The requested information should be sent to Dr **s22** via:

- Email to: IRIS@health.gov.au

For large size documents, please post a universal serial bus (USB), compact disc (CD), or digital versatile disc (DVD) via:

- post to:

Devices Post Market Monitoring Section
Medical Surveillance Branch
Therapeutic Goods Administration
PO Box 100
Woden ACT 2606

- courier to:

Devices Post Market Monitoring Section
Medical Devices Surveillance Branch
Therapeutic Goods Administration
136 Narrabundah Lane
Symonston ACT 2609

Please quote file reference E21-259665 in all correspondence relating to this request. This information is essential for the TGA to track and identify all pieces of data in relation to this device.

Please note that failure to comply with this request for information may attract a range of sanctions. These are outlined in Attachment A.

Should you wish to seek a review of my decision to require you to provide the information/documents about Osseo-Frame - Dental implant system your rights of review are outlined in Attachment B.

Please note that in the event that medical devices are cancelled from the ARTG following a notice under subsection 41JA(1) of the Act requiring information/documents to be provided and a proposal to cancel, the Secretary is required under 41GP of the Act to publish in the Commonwealth Gazette or the Department's website particulars of the cancellation of the entry from the ARTG as soon as practicable after the cancellation. The TGA's practice is to publish information on its website about that cancellation including the name of the sponsor and product and the grounds on which the goods were cancelled. This information is updated in the event that the cancellation is subject to internal review or review by the Administrative Appeals Tribunal, or other form of legal review. [For more information see <http://www.tga.gov.au/about/compliance.htm>]

Please contact the undersigned or Dr **s22** on **s22** or email: IRIS@health.gov.au if any aspect of this request for information requires clarification.

Yours sincerely

Signed electronically by

s22

Delegate of the Secretary
Devices Post Market Monitoring Section
Medical Devices Surveillance Branch
Therapeutic Goods Administration

30 July 2021

ATTACHMENT A

41JA Secretary may require information or documents

(1) The Secretary may, by written notice given to a person:

(a) who is an applicant for a conformity assessment certificate that would relate to a kind of medical device; or

(b) who holds a conformity assessment certificate, or an Australian conformity assessment body certificate, that relates to a kind of medical device; or

(ba) who held, at any time during the notice period under subsection (2), a conformity assessment certificate, or an Australian conformity assessment body certificate, that related to a kind of medical device; or

(c) who is an applicant for the inclusion of a kind of medical device in the Register; or

(d) in relation to whom a kind of medical device is included in the Register; or

(da) in relation to whom a kind of medical device was, at any time during the notice period under subsection (2), included in the Register;

require the person to give to the Secretary information or documents, relating to devices of that kind, that are relevant to one or more of the following:

(e) whether the devices comply with the essential principles;

(f) whether the conformity assessment procedures have been applied to the devices or whether requirements, comparable to those procedures, have been applied to the devices;

(g) whether the devices comply with conditions (if any) imposed on a conformity assessment certificate issued in respect of the device or the inclusion of the device in the Register;

(h) whether either of the following has not been complied with in relation to the devices:

(i) an applicable provision of the Therapeutic Goods Advertising Code;

(ii) any other requirement relating to advertising applicable under Part 5-1 or under the regulations;

(i) if the kind of medical device is included in the Register in relation to the person—whether medical devices of that kind are being:

(i) supplied in Australia; or

(ii) imported into Australia; or

(iii) exported from Australia;

(ia) the safety and efficacy of the devices for the purposes for which they are to be used;

(ib) the regulatory history of the devices in another country;

(j) any other matter prescribed by the regulations for the purposes of this paragraph.

(1AA) If a notice is given under subsection (1) to a person covered by paragraph (1)(ba), then paragraphs (1)(e) to (j) (to the extent to which they are relevant) apply in relation to the period the person held the certificate.

(1AB) If a notice is given under subsection (1) to a person covered by paragraph (1)(da), then paragraphs (1)(e) to (j) (to the extent to which they are relevant) apply in relation to the period the kind of medical device was included in the Register.

(1A) The Secretary may, by written notice given to a person who is an applicant for a conformity assessment certificate, require the person to give to the Secretary such further information concerning the application as is specified in the notice.

(1B) Requirements under subsections (1) and (1A) may be included in the same notice.

(1C) The Secretary may, by written notice given to a person who holds a conformity assessment certificate, require the person to give to the Secretary specified information to be used by the Secretary in deciding whether to suspend the certificate under section 41EM, or to revoke the certificate under section 41ET, in relation to the circumstances referred to in paragraph 41ET(1)(e).

(1D) Requirements under subsections (1) and (1C) may be included in the same notice.

(1E) The Secretary may, by written notice given to an Australian corporation that has been an Australian conformity assessment body require the corporation to give to the Secretary specified information, or specified documents, relating to:

- (a) the certification-related activities carried on by the corporation while the corporation was an Australian conformity assessment body; or
- (b) the conditions referred to in subsection 41EWA(5) that applied while the corporation was an Australian conformity assessment body.

(2) For the purposes of paragraphs (1)(ba) and (da), the notice period is the period:

- (a) of the length specified in the regulations; and
- (b) ending on the day before the Secretary gives the notice under subsection (1).

(3) Nothing in this section affects the operation of Part VIIC of the *Crimes Act 1914* (which includes provisions that, in certain circumstances, relieve persons from the requirement to disclose spent convictions and require persons aware of such convictions to disregard them).

41JB Complying with the Secretary's requirements

(1) The person must give the information or documents to the Secretary:

- (a) within such reasonable time, being not less than 10 working days from the day on which the notice is given, as is specified in the notice; and
- (b) in such form as is specified in the notice.

(2) The form may require or permit information to be given in accordance with specified software requirements:

- (a) on a specified kind of data processing device; or
- (b) by way of a specified kind of electronic transmission.

Offence for failing to comply with a notice

(3) A person commits an offence if:

- (a) the person is given a notice under section 41JA; and
- (aa) the person is covered by paragraph 41JA(1)(b), (ba), (d) or (da) or subsection 41JA(1E); and
- (b) the person fails to comply with the notice.

Penalty: 500 penalty units.

Note: Failure to comply with the notice might also lead to suspension or revocation of a conformity assessment certificate (see Divisions 3 and 4 of Part 4-4) or suspension or cancellation of the entry of a kind of medical device in the Register (see Part 4-6).

(3A) Subsection (3) does not apply if the person has a reasonable excuse.

Note: A defendant bears an evidential burden in relation to the matter in subsection (3A): see subsection 13.3(3) of the *Criminal Code*.

(3B) A person commits an offence if:

(a) the person is given a notice under section 41JA; and

(b) the person is covered by paragraph 41JA(1)(b), (ba), (d) or (da) or subsection 41JA(1E); and

(d) the person fails to comply with the notice.

Penalty: 100 penalty units.

(3C) An offence against subsection (3B) is an offence of strict liability.

(3D) Subsection (3B) does not apply if the person has a reasonable excuse.

Note: A defendant bears an evidential burden in relation to the matter in subsection (3D): see subsection 13.3(3) of the *Criminal Code*.

Offences for giving false or misleading information in purported compliance with a notice

(4) A person commits an offence if:

(a) the person is given a notice under section 41JA in relation to a kind of medical device; and

(b) the person gives information in purported compliance with the notice; and

(c) the information is false or misleading in a material particular; and

(d) either:

(i) the use of the kind of medical device has resulted in, will result in, or is likely to result in, harm or injury to any person; or

(ii) the use of the kind of medical device, if the kind of medical device were used, would result in, or would be likely to result in, harm or injury to any person.

Penalty: Imprisonment for 5 years or 4,000 penalty units, or both.

Note 1: A jury may acquit a person of an offence against this subsection and may convict the person of an offence against subsection (7) instead: see section 53A.

Note 2: For the liability of an executive officer of a body corporate, see sections 54B and 54BA.

ATTACHMENT B

Request for reconsideration of an initial decision

This decision is a reviewable initial decision under section 60 of the Act. Under section 60, a person whose interests are affected by a 'reviewable' initial decision, can seek reconsideration of the initial decision.

As this document constitutes written notice of the making of an initial decision being given by the Secretary, a request for reconsideration of this initial decision must be given to the Minister within 90 days and be accompanied by any information that you wish to have considered. A request for reconsideration given to the Minister outside the statutory 90 day reconsideration period cannot be accepted.

The Minister may either personally undertake a request for reconsideration of an initial decision or delegate to an officer of the Department with the appropriate delegation.

Under section 60(3A) of the Act, the Minister (or the Minister's delegate) is not able to consider any information provided after the notification is made of a request for reconsideration of an initial decision unless the information is provided in response to a request from the Minister (or the Minister's delegate), or it is information that indicates that the quality, safety or efficacy of the relevant therapeutic goods is unacceptable.

Guidelines for requesting reconsideration of an initial decision

A request for reconsideration should be made in writing, signed and dated by the person requesting reconsideration, should be titled "<insert person/company name> - Request for Reconsideration Under Section 60 of the *Therapeutic Goods Act 1989*" and should include the following:

- a copy of the initial decision notification letter (or other evidence of notification);
- identify, and describe with as much specificity as possible, which component(s) of the initial decision should be reconsidered and set out the reasons why reconsideration is requested;
- any information/documentation in support of the request, clearly labelled to correspond with (any or each of) the reasons why reconsideration is requested; and
- an email address nominated for the purposes of receiving correspondence in relation to the request for reconsideration.

All requests for reconsideration should be given to the Minister by email:

Email: '**minister.hunt.DLO@health.gov.au**' and copied to '**decision.review@health.gov.au**'

Requests for reconsideration that include dossiers (or similar bulk material) that cannot easily be attached to the request given first by email, may then be submitted on a USB drive or CD sent by express post or registered mail to:

Mail: **Minister for Health**
Suite M1 40
c/- Parliament House
CANBERRA ACT 2600

Subject to the *Administrative Appeals Tribunal Act 1975* (AAT Act), if you are dissatisfied with the decision upon reconsideration by the Minister (or the Minister's delegate), you can apply to the Administrative Appeals Tribunal (AAT) for a review of that decision upon reconsideration.

NOTE: This initial decision remains in effect unless and until it is revoked or revoked and substituted by the Minister (or the Minister's delegate) as a result of a request for reconsideration under section 60 of the Act OR is set aside, varied or remitted by the AAT or is otherwise overturned or stayed.

Further information on Section 60 and the Reviewable Initial Decisions can be found at the following TGA webpages:

Section 60 - <https://www.tga.gov.au/publication/guidance-requesting-reconsideration-initial-decision>

Reviewable Initial Decisions - <https://www.tga.gov.au/publication/reviewable-initial-decisions-under-therapeutic-goods-act-regulations-0>

s22

From: s22
Sent: Thursday, 9 September 2021 5:55 PM
To: s22
Cc: s22 @medacta.com.au; s22
Subject: R: s41JA request for DIR 68783 & 68516 [SEC=OFFICIAL]
Attachments: IFU MectaLock PEEK and MectaLock TI EN 75.09.045 rev.08.pdf

REMINDER: Think before you click! This email originated from outside our organisation. Only click links or open attachments if you recognise the sender and know the content is safe.

Dear s22

You will find attached the MectaLock suture anchors IFU (Ref: 75.09.045 rev. 08) in force.
 At the page number 4, ***“Breakage of the MectaLock Suture Anchor can occur if:”***, is reported as follows:

- Any lateral loading is applied on the anchor/insertor while inserting the MectaLock anchor

Do not hesitate to contact me if you have any issues.

Kind regards,

s22 s22 Regulatory Affairs Department s22 @medacta.ch www.medacta.ch	Medacta International SA Strada Regina CH-6874 Castel San Pietro Switzerland Pho: s22 Mob: s22 Fax: +41 91 696 60 66
--	--

Da: s22 @Health.gov.au>
Inviato: giovedì 9 settembre 2021 06:15
A: s22 @medacta.ch>
Cc: s22 @medacta.com.au; s22 @medacta.ch>
Oggetto: RE: s41JA request for DIR 68783 & 68516 [SEC=OFFICIAL]

Dear s22

Good afternoon. A review of the manufacturer's investigation shows that the component has been subjected to excessive lateral bending as the axis of the component results deviated from the original. Please provide a copy of Instructions for Use (IFU) for this device and show where it is mentioned in IFU that not doing this 'excessive lateral bending' in a warning or use space. Thanking you.

Kind regards,

Dr s22
 Departmental Officer
 Devices Post Market Monitoring
 Medical Devices Surveillance Branch
 Email: s22 @health.gov.au
 Telephone number : s22

Therapeutic Goods Administration

Department of Health
PO Box 100
Woden ACT 2606 Australia
www.tga.gov.au

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met.

Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission.

The Department of Health acknowledges the Traditional Custodians of Australia and their continued connection to land, sea and community. We pay our respects to them and their culture, and to all Elders both past and present.

P Please consider the environment before printing this email.

From: s22
Sent: Wednesday, 1 September 2021 5:02 PM
To: s22 @medacta.ch>
Cc: s22 @medacta.com.au; s22 @medacta.ch>
Subject: RE: s41JA request for DIR 68783 & 68516 [SEC=OFFICIAL]

Dear s22

Good afternoon. Thank you so much for providing the requested information. We are currently reviewing it.
Thanking you.

Kind regards,

Dr s22
Departmental Officer
Devices Post Market Monitoring
Medical Devices Surveillance Branch
Email: s22 @health.gov.au
Telephone number s22

Therapeutic Goods Administration

Department of Health
PO Box 100
Woden ACT 2606 Australia
www.tga.gov.au

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met.

Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission.

The Department of Health acknowledges the Traditional Custodians of Australia and their continued connection to land, sea and community. We pay our respects to them and their culture, and to all Elders both past and present.

P Please consider the environment before printing this email.

From: s22 [REDACTED]@medacta.ch>
 Sent: Saturday, 28 August 2021 1:03 AM
 To: s22 [REDACTED]@Health.gov.au>
 Cc: IRIS <IRIS@health.gov.au>; s22 [REDACTED]@medacta.com.au; s22 [REDACTED]@medacta.ch>
 Subject: R: s41JA request for DIR 68783 & 68516 [SEC=OFFICIAL]

REMINDER: Think before you click! This email originated from outside our organisation. Only click links or open attachments if you recognise the sender and know the content is safe.

Dear Dr s22 [REDACTED]

you will find here below the answers to your requests related to DIR 68783 and 68516:

1. Please provide a description of the adverse events and the rate for those events. Separate the events for Australia and the rest of the world and provide each in a percentage format. Provide a list of complaints and numbers for those similar complaints.
 This is the first event reported, no additional instrument breakages have been communicated considering both MectaLock PEEK Ø2.4 short and long drivers. In total, we have sold s47 [REDACTED] short and long drivers (MectaLock PEEK Ø2.4).
2. In your detailed device analysis response dated 22 June 2021, it states "the laser welding failed, and the tip component has been detached from the main shaft of the inserter". Please provide details of your risk mitigation strategies in order to avoid the laser welding failure to ensure the device safety. You will find attached a dedicated rationale.

Do not hesitate to contact me if you have any issues.

Kind regards,

s22 [REDACTED] s22 [REDACTED] <i>Regulatory Affairs Department</i> s22 [REDACTED]@medacta.ch www.medacta.ch	Medacta International SA Strada Regina CH-6874 Castel San Pietro Switzerland Pho: s22 [REDACTED] Mob: s22 [REDACTED] Fax: +41 91 696 60 66
--	--

From: s22 [REDACTED]@Health.gov.au>
 Sent: Friday, 30 July 2021 2:34 PM
 To: s22 [REDACTED]@medacta.com.au>
 Subject: s41JA request for DIR 68783 & 68516 [SEC=OFFICIAL]

Dear s22 [REDACTED],

Good afternoon. Please find attached section 41JA request for DIR 68783 & 68516. Please note that the requested information must be provided within 20 working days from the date of this letter and no later than **30 August 2021**. Thanking you.

Kind regards,
 Dr s22 [REDACTED]
 Departmental Officer
 Devices Post Market Monitoring

Medical Devices Surveillance Branch
Email: s22 [REDACTED]@health.gov.au
Telephone number : s22 [REDACTED]

Therapeutic Goods Administration

Department of Health
PO Box 100
Woden ACT 2606 Australia

www.tga.gov.au

For ongoing information and updates please subscribe to the TGA's [Medical Devices Information](#) and [IVDs Information](#) email subscription services.

This response is general information given to you without prejudice; it is not binding on the TGA and you should get your own independent legal advice to ensure that all of the legislative requirements are met.

Important: *This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission.*

The Department of Health acknowledges the Traditional Custodians of Australia and their continued connection to land, sea and community. We pay our respects to them and their culture, and to all Elders both past and present.

P *Please consider the environment before printing this email.*

"Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission."

Message protected by MailGuard: e-mail anti-virus, anti-spam and content filtering.

<http://www.mailguard.com.au>

[Report this message as spam](#)

"Important: This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission."

Ref: 75.09.045 rev. 08

CE 2797
2020

Strada Regina
CH - 6874 Castel San Pietro - Switzerland
Phone +41 91 696 60 60 - Fax +41 91 696 60 66
info@medacta.ch - www.medacta.com

Manufactured by:
Medacta International SA
CH-6874 Castel San Pietro - Switzerland

MECTALOCK SUTURE ANCHORS

CAUTION

For U.S. specific instructions for use, please go to www.medacta.com or call toll free 800-901-7836.

SYMBOLS / SIMBOLI / PICTOGRAMMES / PIKTOGRAMME / PICTOGRAMAS / PICTOGRAMAS / PICTOGRAMMES / ΕΙΚΟΝΟΣΥΜΒΟΛΑ

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ENGLISH: MECTALOCK SUTURE ANCHORS - INSTRUCTIONS FOR USE

Important notice: the device(s) can be prescribed and implanted only by a doctor legally authorized to perform this type of surgery.

1. GENERAL

Before any surgery, the surgeon must be familiar with the product literature and operative technique and must read carefully these instructions for use. Patient selection is as important as implant placement or positioning. Unsuitable functional requirements may contribute to reduce the implant life. The warnings must be heeded, and the instructions for use must be strictly followed.

2. PRODUCT DESCRIPTION

These instructions for use are intended for all products described here below.

CATALOGUE NUMBER	ANCHOR DESCRIPTION	COMPATIBLE REUSABLE DRILL N.
05.11.001	MectaLock PEEK Ø2.4 Long driver	05.11.10.0035
05.11.002	MectaLock PEEK Ø2.9 Long driver	05.11.10.0036
05.11.003	MectaLock PEEK Ø3.4 Long driver	05.11.10.0037
05.10.001	MectaLock PEEK Ø2.4 Short driver	05.10.10.0015
05.10.002	MectaLock PEEK Ø2.9 Short driver	05.10.10.0016
05.10.003	MectaLock PEEK Ø3.4 Short driver	05.10.10.0017
05.10.011	MectaLock TI Ø5.0	Not needed
05.10.012	MectaLock TI Ø6.5	Not needed

Medacta MectaLock PEEK (knotless) and MectaLock TI (preloaded) are manufactured respectively in PolyEtherEtherKeton (PEEK) or Titanium alloy. The MectaLock Suture Anchors are provided with sutures directly preloaded into the implant (MectaLock TI) or placed directly into the final packaging (MectaLock PEEK).

CATALOGUE NUMBER	ANCHOR FAMILY	INCLUDED SUTURES
05.11.001	MectaLock PEEK	1x USP#2/EP#5 HS Fiber 39" UHMWPE suture
05.11.002		
05.11.003		
05.10.001		
05.10.002		
05.10.003		
05.10.011	MectaLock TI	2x USP#2/EP#5 HS Fiber 39" preloaded sutures
05.10.012		

Material Specifications

Anchors: PolyEtherEtherKeton (PEEK) [MectaLock PEEK]; Ti6Al4V ELI ISO 5832-3 [MectaLock TI]

Driver shaft/body: Stainless steel

Driver Handle: Polycarbonate medical grade

Suture: non-absorbable USP#2/EP#5 HS Fiber 39" Ultra-High molecular weight polyethylene (average diameter: 0.500-0.599mm)

3. INTENDED USE / INDICATIONS

The non-absorbable suture anchors are intended to:

- MectaLock TI Suture anchors: soft tissue refixation within the shoulder joint (i.e.: cuff rotator repair, biceps tenodesis)
- MectaLock PEEK Suture anchors: treatment of the hip and shoulder instability by means of the refixation of the acetabular or glenoidal labrum to the bone

The device is indicated for use in both arthroscopic and open surgical procedures. After the suture is anchored to the bone, it is used to reattach soft tissues to the bone, including the following system specific indications:

- MectaLock TI Suture anchors: soft tissue refixation within the shoulder joint
- MectaLock PEEK Suture anchors: treatment of the hip and shoulder instability by means of the refixation of the acetabular or glenoidal labrum to the bone

The suture anchor system in that way stabilizes the damaged soft tissue, in combination with appropriate postoperative immobilization, throughout the healing period.

4. CONTRAINDICATIONS

- Pathological conditions in the soft tissue to be repaired or reconstructed which would adversely affect suture fixation
- Pathological conditions of bone which adversely affect the MectaLock Suture Anchors
- Physical condition that would eliminate, or tend to eliminate, adequate implant support or retard healing
- Conditions which tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period
- Attachment of artificial ligament or other implants
- Foreign body sensitivity, known or suspected allergies to implant and/or instrument materials

5. WARNINGS AND PRECAUTIONS

The success of the operation depends on compliance with the operative technique supplied, and the proper use of the dedicated instruments specially designed for that range of implants. Surgeon must choose proper implant size based on specific procedure and patient history. Malpositioning may reduce implant longevity and lead to early implant failure. The disposable driver is single use and must be disposed of according to hospital policy and procedure.

One Medacta label for each device must be properly applied to the Medacta Implant Card provided, and delivered to the patient.

Ensure the MectaLock Suture Anchors are used with the corresponding compatible reusable drill bit. Use the corresponding compatible reusable drill bit to create the pilot hole to create the pilot socket. Maintain proper alignment during insertion of the anchor and prevent disengagement from the driver. The risk of MectaLock Suture Anchor breakage during implantation is reduced by following the specified instructions for use listed below. Proper selection and placement of the implant are important considerations in the successful utilization of this device. Proper orientation and alignment of instruments is important during implantation of the MectaLock Suture Anchor to minimize possible breakage of the anchor. Breakage of the MectaLock Suture Anchor can occur if:

- The pilot hole, where expected, is not drilled to the proper depth
- The MectaLock Suture Anchor is not properly aligned with the pilot hole, where expected
- Any lateral loading is applied on the anchor/insertor while inserting the MectaLock anchor
- The MectaLock Suture Anchor driver is used for prying

Do not disengage the sutures from the disposable driver while inserting the MectaLock TI Suture Anchor before the final positioning depth is reached. Do not use any other suture other than non-absorbable USP#2/EP#5 UHMWPE suture.

MRI Compatibility

MectaLock TI:

Non-clinical testing has demonstrated titanium implants are MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions. Static magnetic field of 1.5 Tesla and 3 Tesla, with:

- Maximum spatial field gradient of 12,800 G/cm (128 T/m)
- Maximum force product of 231,000,000 G²/cm (231 T²/m)

Theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode). Under the scan conditions defined above, titanium implants expected to produce a maximum temperature rise of less than:

- 1.5°C (2 W/kg, 1.5 Tesla) RF-related temperature increase with a background temperature increase of ≈1.3°C (2 W/kg, 1.5 Tesla)
- 3.0°C (2 W/kg, 3 Tesla) RF-related temperature increase with a background temperature increase of ≈2.3°C (2 W/kg, 3 Tesla)

after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 20.5 mm from the titanium implant when imaged with a gradient echo pulse sequence and a 3 Tesla MR system.

The MectaLock PEEK Suture Anchors are MR Safe.

6. MEDACTA INTERNATIONAL IMPLANTS

Medacta International is not responsible for the use of its implant components in combination with a component from another manufacturer (unless otherwise specified by Medacta International in the surgical steps), therefore we advise against such a use. The MectaLock Suture Anchors should never be reimplanted. While an implant may appear undamaged, microscopic imperfections may occur and cause implant failure. The operating surgeon must be aware that even a very small superficial damage, caused for instance by a sharp tool or electrocautery, can have an influence on the endurance of the device and can lead to failure. Medacta MectaLock Suture Anchors are non-pyrogenic.

7. RISK FACTORS

The following risk factors, individually or together, may result in poor clinical outcomes:

- Inadequate bone quality, (e.g. osteoporosis, osteopenia)
- Systemic diseases or metabolic disorders
- History of infections or recurrent falls
- Drug dependence and abuse of alcohol and medicaments
- Mental incapacity of patient to understand the instructions of the physician and to comply with them
- Local bone tumors

8. INSTRUCTIONS FOR USE

PREOPERATIVE PHASE

The surgeon should verify possible patient physical limitations and mental deficiencies and he should also discuss with the patient all the details of the procedure and implant. The discussion should consider the limitations of the procedure and the constraints imposed by the selected implant. The factors which could limit the performance and

stability of the implant, e.g. level of activity, should be set out to improve the patient's chances to avoid complications. The necessity to follow the postoperative instructions given by the surgeon should be fully understood by the patient. A stock of sterile implants of suitable sizes must be available and checked by the operator before surgery.

9. HANDLING

To avoid scratching or damaging the implants, these should be handled with the utmost care by qualified personnel and in an environment where conditions of hygiene are controlled. The implants should be kept in their undamaged packages.


10. SURGICAL TECHNIQUE

The surgeon should be fully familiar with the surgical technique. Supplementary information about the surgical techniques (brochure and video) and products are available on request. Careful preoperative planning is essential.

MECTALOCK PEEK SUTURE ANCHOR

- Pass the USP#2/EP#5 HS Fiber 39" UHMWPE suture through soft tissue with surgeon's preferred suture-passing technique; let the two strands of the suture come out from the patient through the arthroscopic portal
- Find the desired location for the placement of the MectaLock PEEK anchor; insert a percutaneous guidewire (max Ø1.5mm) through the proper arthroscopic portal to reach the desired location
- A dedicated aimer can be used as a drill guide and to aid in pointing the location of the pilot hole in the bone. Insert the Medacta

Arthroscopic Aimer coupled with a cannulated obturator over the guidewire; the surgeon can choose between different aimer's tip according to his preference and the anatomical structures to be treated. Aimers and obturators are provided in two different lengths:

Aimers/obturators and drills can be coupled according to the following table:

AIMER		OBTURATOR		REUSABLE DRILL BIT	
CATALOGUE N.	DESCRIPTION	CATALOGUE N.	DESCRIPTION	CATALOGUE N.	DESCRIPTION
05.11.10.0030	Crown tip Long	05.11.10.0033	Obturator Long	05.11.10.0035	Drill Ø2.4 Long
05.11.10.0031	Fishmouth tip Long			05.11.10.0036	Drill Ø2.9 Long
05.11.10.0032	Fork tip Long			05.11.10.0037	Drill Ø3.4 Long
05.10.10.0010	Crown tip Short	05.10.10.0013	Obturator Short	05.10.10.0015	Drill Ø2.4 Short
05.10.10.0011	Fishmouth tip Short			05.10.10.0016	Drill Ø2.9 Short
05.10.10.0012	Fork tip Short			05.10.10.0017	Drill Ø3.4 Short

- Using the corresponding compatible reusable drill bit, create a hole in the bone in the desired location for the placement of the MectaLock PEEK anchor
- Make sure to drill until the distal depth mark of the drill is aligned with the distal depth mark of the aimer or until the drill is mechanically blocked against the aimer

Load the MectaLock PEEK anchor with the sutures coming out of the portal: while firmly holding the two strands of the suture with one hand, align them with the slot of the anchor's tip. Gently push the anchor against the suture to let it clip into the tip's eyelet.

Rotate the anchor clockwise by 90°


Align the main axis of the driver with the suture (1) and then, while still holding the suture gently tensioned, slide the anchor over the suture (2) towards the patient through the arthroscopic portal, until the pilot hole is reached.

Once reaching the pilot hole, regulate the suture tension according to the surgeon's preference and tap the plastic handle of the driver until the MectaLock PEEK anchor is correctly placed. The correct depth is achieved when the distal laser marking is flush with the bone surface.

- Unscrew the disposable driver
- Palpate the repaired soft tissue using an arthroscopic probe to ensure that the repair is secure
- Cut and remove sutures in excess


MECTALOCK TI SUTURE ANCHOR

- Find the desired location for the placement of the MectaLock TI anchor; insert the preloaded driver with the MectaLock TI anchor through the proper arthroscopic portal to reach the desired location.
- Gently tap on the back handle of the disposable driver to let the MectaLock TI anchor's tip pierce the bone in the chosen spot; while keeping a gentle axial force on the disposable driver, screw the MectaLock TI anchor into the bone until the distal laser marking on the driver is flush with the bone surface.

Unlock the sutures from the plastic handle and gently disengage the driver's tip from the MectaLock TI anchor, pulling the driver out through the arthroscopic portal.

- Manage the 2x USP#2/EP#5 HS Fiber 39" UHMWPE sutures through soft tissue with surgeon's preferred suture-passing and knotting technique
- Palpate the repaired soft tissue using an arthroscopic probe to ensure that the repair is secure
- Cut and remove suture in excess

11. POSTOPERATIVE CARE AND FOLLOW-UP

The surgeon should caution the patients to control their level of activity and avoid excessive loads on the operated limb. Moreover, the surgeon should make the patients aware of the precautions to be taken in terms of exercises, treatments and limitations on activities, any limitations reported on the label, as well as exposure to magnetic fields. The patient must be told that implants can affect the results of computer tomography (CT) or magnetic resonance imaging (MRI) scans. Periodic follow-up are recommended to make comparisons with the immediate postoperative condition and anticipate implant related complications. Excessive physical activity, and operated limb traumas may cause early failure of the implant. If the case occurs, it is necessary to place the patient under supervision, evaluate the possible progression of the deterioration, and weigh the benefit of early revision.

12. ADVERSE EFFECTS AND COMPLICATIONS

Adverse effects that can occur in the reattachment of soft tissue to bone in orthopedic surgical procedures:

- Infection, both deep and superficial
- Allergies, mild inflammatory and foreign body reactions to implant material

Some adverse effects can ultimately lead to death.

General complications include:

- Venous thrombosis with/without pulmonary embolism
- Cardiovascular or pulmonary disturbances
- Haematomas
- Systemic allergic reactions
- Systemic pain

13. PACKAGING

The MectaLock Suture Anchors are supplied in single-use packages. For components delivered sterile, the sterilization method is indicated on the label. The expiration date must be checked on the label as well as the package integrity to ensure that sterility of the contents has not been compromised. If the package is damaged or has been previously opened, do not use the component. Do not re-sterilize.

14. STORAGE

The packages must be stored in a cool, dry place, away from light.

15. TRADEMARKS

Medacta is registered trademark of Medacta International SA, Castel San Pietro, Switzerland.

The reference text is the English text.

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s22

From: s22
Sent: Saturday, 28 August 2021 1:03 AM
To: s22
Cc: IRIS; s22@medacta.com.au; s22
Subject: R: s41JA request for DIR 68783 & 68516 [SEC=OFFICIAL]
Attachments: TGA response (s41JA request for DIR 68783 & 68516).pdf

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1. Please provide a description of the adverse events and the rate for those events. Separate the events for Australia and the rest of the world and provide each in a percentage format. Provide a list of complaints and numbers for those similar complaints.
 This is the first event reported, no additional instrument breakages have been communicated considering both MectaLock PEEK Ø2.4 short and long drivers. In total, we have sold s47 short and long drivers (MectaLock PEEK Ø2.4).
2. In your detailed device analysis response dated 22 June 2021, it states "the laser welding failed, and the tip component has been detached from the main shaft of the inserter". Please provide details of your risk mitigation strategies in order to avoid the laser welding failure to ensure the device safety. You will find attached a dedicated rationale.

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From: s22@Health.gov.au>
Sent: Friday, 30 July 2021 2:34 PM
To: s22@medacta.com.au>
Subject: s41JA request for DIR 68783 & 68516 [SEC=OFFICIAL]

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[Report this message as spam](#)

Event description

The patient underwent a labral repair s22 hip arthroscopy. Two MectaLock (J-Lock) PEEK anchors were used for refixing the labrum to the acetabulum with no visual issues. Then the surgeon sent to the agent a picture of an X-ray asking if the image in the X-ray was part of the anchor inserter handle. It appears that a part of the inserter handle has remained in the bone with the anchor. The patient underwent revision surgery to remove the broken part.

The event was reported to the manufacturer on 29 Jan 2021, but only after the revision surgery, which took place on s22, we were sure that the piece that remained in the patient was part of the Medacta device.

TGA question

In your detailed device analysis response dated 22 June 2021, it states “the laser welding failed, and the tip component has been detached from the main shaft of the inserter”. Please provide details of your risk mitigation strategies in order to avoid the laser welding failure to ensure the device safety.

Reply

The laser welding on the MectaLock inserter component is a technical specification associated to a critical process which is performed within the manufacturing process of the device. In particular, the residual failure risk associated to the welding is evaluated within the assembly process of the component as reported in the Risk analysis (Annex 2 of TF LI-SR-01, ITEM 1).

Item (per M07.26.xxxx)	Hazard (potential source of harm)	Harm (physical injury or damage to health of people, or damage to property or the environment)	Rationale for Probability of Harm	Severity of Harm	Probability of Harm	Residual Risk
B.02.3.10	Assembly process is not adequate to provide a safe and effective final product	Implant failure before the expected lifetime with consequent revision surgery. Another operation required with anaesthesia or septic risks	An adequate and controlled manufacturing process flow is implemented for the assembling of the “MectaLock PEEK Suture anchor” final product. (See Annex 20 of TF LI-SR-01). Validation activities and mechanical resistance tests are performed on products that are conform to all the specifications.	Very Serious	Very Low	Acceptable
B.02.3.10	Assembly process is not adequate to provide a safe and effective final product	Implant failure before the expected lifetime with consequent revision surgery. Another operation required with anaesthesia or septic risks	An adequate and controlled manufacturing process flow is implemented for the assembling of the “MectaLock PEEK Suture anchor” final product. (see Annex 20 of TF LI-SR-01). Validation activities and mechanical resistance tests are performed on products that are conform to all the specifications.	Very Serious	Very Low	Acceptable

As indicated within the manufacturing process flow of the device (FSR 0049, see extract below), the laser welding is one of the steps subject to a 100% Full Verification performed by the supplier and by Medacta, in order to verify that the specifications reported on the technical drawings are fully respected.

Analysis of Complaint 0172-21

4.3.6 SM Mectalock – Assembly (Cycle C0086)

Step #	Step Description	Planning Documents	Medacta / Supplier	Process FMEA	Qualification and Validation Documents	Deviation	Results	Other documents [e.g. rationales, additional information]
Start	Machining, laser marking, laser welding, cleaning of anchor inserter	SPC-01.009.000 rev.4	Maire-Didier	RA-01.009.000 rev.1	Full verification FVR-01.009.003	NONE	PASS	N/A

The supplier performs laser welding following a process that has been validated internally, defining specific technical parameters of the machines reported within the following Validation documents:

	IQ	OQ
Laser welding (internal MAIRE-DIDIER)	QI0016	QO/RQ00004

Moreover, Medacta performs inspections on the subject device according to a defined quality plan, specified within the inspection document CA51059. In particular, the laser welding is checked on 100% of the components in accordance with the internal procedure IC 08.55, which specifies the related acceptance criteria.

Design validation activities have been also carried out with final products to assess the required mechanical performance during all conditions of use and to confirm the safety and integrity of the device (reported in the Technical File Annex 8 of TF LI-SR-01, ITEM 1).

Conclusions

All necessary actions have been implemented to minimize any potential risk associated with the laser welding integrity. This is supported by design verification and validation activities which confirmed the effectiveness of the device during its intended use. The total residual risk of the whole system is therefore acceptable, with consideration of the individual risk and all risks from a cumulative perspective. It is therefore confirmed that the benefits, as correlated to the user needs, outweigh the collective and individual risks.

Description	Prepared by [CP]	Verified by [AQ]	Approved by [R]	Date
Reply to TGA question #2 – Complaint 0172-21	s22	s22	s22 	27 Aug 21

IRIS DIR Peer Review Procedure Checklist

DIR No.	68783	TRIM file No:	E21-259665
Lead investigator:	s22		
Review start date:	23/09/2021		

Instructions:

Save this document into the correct TRIM DIR number file

Complete **each** checklist question by ticking the appropriate checkbox.

Record any comments or requests to the investigator at the end of the document.

Return the review checklist and all files / documents / samples to the lead investigator

Discuss any comments with the lead investigator

Review of INITIAL requests to the sponsor / reporter

	Yes	No	NA
1.1 Are the questions asked or the information requested of the sponsor and/or reporter in line with the Team Meeting recommendations for the investigation?	☒	☐	☐
1.2 Are the questions asked or the information requested of the sponsor and/or reporter appropriate to the issue to be investigated?	☒	☐	☐
1.3 Can you suggest any additional information to request from the sponsor/reporter? If YES, list these in Review Comments section.	☐	☒	☐

Review of the Investigation

	Yes	No	NA
2.1 Has the investigator followed the appropriate WI for the level of investigation (see DPMMS WI 1.3 or 1.10)	☒	☐	☐
2.2 Are the administration details recorded on the DIR accurate for the device (eg ARTG number)?	☒	☐	☐
2.3 Have all relevant test results/reports, manufacturer/sponsor reports been included in the investigator's final outcomes/recommendation?	☒	☐	☐
2.4 Do the investigator's outcomes/recommendations address the issue as reported?	☒	☐	☐
2.5 Are the problems observed, Clinical signs and symptoms and health Impact consistent with the event clinical description?	☒	☐	☐
2.6 Are the Investigation findings, conclusions and OUTCOME decisions consistent with the investigation conducted?	☒	☐	☐

Review Comments

Comments/Recommendations:

- Correspondence and Chronology page has not been completed correctly. It should include a catalogue of all correspondence and evidence awkward (e.g. samples collected, and s41JAs sent/received).
- The evidence page has not been filled out correctly. This page should only be filled out between a team meeting and sending your questionnaire/RFI via DST. It gives us a snapshot of what you aim to complete in your investigation, including what potential EPs are in breach, and what questions you need answered. It is not an ongoing page that you continuously update as your investigation progresses, that is the purpose of the correspondence and chronology page.
- The investigation summary is great. However, I could not find any of the risk management analysis documents referenced within TRIM. Did you obtain this document? Or did you just fail to file it into TRIM?

Investigator's response:

Thanks for your valuable recommendations.

1. I have completed correspondence and chronology page putting the date and trim number for s41JA sent/received, IFU request/ received.
2. The evidence page has been filled based on your recommendation.
3. The risk management analysis document was not trimmed earlier and now I added into the trim container.

Reviewer Name: **S22**

Completion Date: 30 September 2021

Signature Signed Electronically

(Do not sign until the Review comment/recommendations have been agreed and finalised with the Investigator.)

s22

From: s22@Health.gov.au>
Sent: Wednesday, 6 October 2021 11:00 AM
To: s22
Subject: Reporter Completion Letter for DIR 68783 [SEC=OFFICIAL]
Attachments: [D21-3178893] DIR 68783_ Reporter_Complete_Letter.PDF

Dear MS s22

Good morning. Please find attached reporter completion letter for DIR 68783. Thanking you.

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MS **s22**
Quality and Risk Consultant
s22
s22
s22
Email: **s22**

File Reference: E21-259665
Sent by email

Dear MS **s22**

**DEVICE INCIDENT REPORT DIR 68783 - Medacta J-Lock PEEK Suture - Anchor, soft-tissue
implanted, non-biodegradable**

An investigation into the incident you reported to the Therapeutic Goods Administration concerning the above device is now complete.

A copy of the Incident Report Investigation Scheme (IRIS) database entry, including the investigation summary is attached for your information.

Thank you for your support of the Medical Device Incident Report Investigation Scheme. Should you have any questions regarding this report please contact our team on **s22** or send an email to: IRIS@health.gov.au

Yours sincerely,
Signed electronically by

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Investigation Summary:

A review of the manufacturer's analysis shows that the welding failed due to the component being subjected to stress of excessive lateral bending. The manufacturer notes that it is not possible to identify a clear or specific root cause for the mode of failure. A review of the warnings and precautions section in the Instructions for Use (IFU) has shown that the breakage of the MectaLock Suture Anchor can occur if any lateral loading is applied on the anchor/insertor while inserting the MectaLock anchor. The review of the risk management analysis document shows that all necessary actions have been implemented to minimize any potential risk associated with the laser welding integrity. In addition, this is supported by design verification and validation activities which confirmed the effectiveness of the device during its intended use. The sponsor notes that this is an isolated event and a review of the IRIS database has shown that there is only one similar adverse event. No further investigation will occur at this time; however the TGA will continue to monitor the rate and pattern of occurrence and may re-open the file as appropriate.

Date Completed:

End of DIR 68783

From: s22
To: s22 @medacta.com.au
Subject: Sponsor Completion Letter for DIR 68783 [SEC=OFFICIAL]
Date: Wednesday, 6 October 2021 11:12:33 AM
Attachments: [D21-3178905] DIR 68783_Sponsor_Complete_Letter.PDF

Dear s22,

Good morning. Please find attached sponsor completion letter for DIR 68783. Thanking you.

Kind regards,

Dr s22
Departmental Officer
Devices Post Market Monitoring
Medical Devices Surveillance Branch
Email: s22 @health.gov.au
Telephone number : s22

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***Important:** This transmission is intended only for the use of the addressee and may contain confidential or legally privileged information. If you are not the intended recipient, you are notified that any use or dissemination of this communication is strictly prohibited. If you receive this transmission in error please notify the author immediately and delete all copies of this transmission.*

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File Reference: E21-259665
Sent by email

Dear **s22**

DEVICE INCIDENT REPORT DIR 68783 - ARTG # 320484 - Ligament anchor

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Device Analysis Results:

Additional information received by Sponsor 03/06/2021

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Cause of Problem (Level 2)

Summary of Investigation:

A review of the manufacturer's analysis shows that the welding failed due to the component being subjected to stress of excessive lateral bending. The manufacturer notes that it is not possible to identify a clear or specific root cause for the mode of failure. A review of the warnings and precautions section in the Instructions for Use (IFU) has shown that the breakage of the MectaLock Suture Anchor can occur if any lateral loading is applied on the anchor/insertor while inserting the MectaLock anchor. The review of the risk management analysis document shows that all necessary actions have been implemented to minimize any potential risk associated with the laser welding integrity. In addition, this is supported by design verification and validation activities which confirmed the effectiveness of the device during its intended use. The sponsor notes that this is an isolated event and a review of the IRIS database has shown that there is only one similar adverse event. No further investigation will occur at this time; however the TGA will continue to monitor the rate and pattern of occurrence and may re-open the file as appropriate.

Date Completed:

End of DIR 68783

**Device Incident Report : Medical Devices Branch - Device Vigilance and Monitoring**

DIR : 46 - ID : 571181

Released by s22 on 28/09/2022 14:44:45

Report #:	Records Management #:	Reporter's Reference #:	Report Type:
89036		1273-23	Final

ARTG: 320484	Document Container URL
--------------	--

Report Information Section

Report Status:	Sponsor's Reported Category:	Date of Adverse Event:	Date of Initial Report:
Closed	Other	s22	12/07/2023
Date of Final Report:	Date of Initial TGA Action:	Reviewed by Team:	Date Response Received:
05/10/2023	12/07/2023		
Date Completed:	Operator at Time of Event:	If 'Other' Operator Selected:	Reporter consents to contact by sponsor:
16/10/2023	Healthcare Professional		N/A
Source of Report:	If 'Other' Source Selected:	Type of Initial Action:	
Sponsor		Trend data only	

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The tip of the inserter was broken.

Clinical Event Information:

During anchor insertion, the surgeon felt like a lot of force was required to mallet into bone. Once the anchor was fully implanted, the introducer came away from the anchor immediately (without need for counterclockwise turns). The broken tip of the inserter was visible in the intraoperative xray.

It is feasible that, during the malleting phase, the anchor and the inserter were not fully axially aligned with the pilot hole and that the malleting was not fully axially aligned with the handle of the inserter. Slotted Aimer was not used. The tip was left inside the patient.

Suture fixation was adequate and surgeon used 2 more JLOCK anchors successfully following and 1 prior to the breakage. No revision surgery is planned.

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Recorded Problems Observed

Recorded Problems Observed:

Appropriate Term/Code Not Available -> -> ;
Material Integrity Problem -> Break -> Material Fragmentation

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Health Impact

Recorded Health Impacts:

Change in Therapeutic Response -> Unexpected Therapeutic Effects ->

Patient Information

Sex:	Weight:	Age:
s22		s22 old
Patient Focused Corrective Action Taken:	Patient History:	
Patient Outcome/Consequences:	Additional Event Description:	

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Medical Problem Device Used For: Document 21

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Date of Explant:

Access Contact Title:

	§22	With Manufacturer	Document 21
Access Contact First Name:	Access Contact Surname:	Access Contact Phone:	Access Contact Fax:
Access Contact Email:	Licence Status:	Status Effective Date:	Additional Devices Added:
	A	19/07/2019	0

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Manufacturer Name:		Manufacturer Client Id:		Address 1:	
Medacta International SA		45296			
Address 2:		State/Province:		Country:	
Postcode:		Fax:		Email:	
Manufacturer Informed:		Contact Title:		Contact First Name:	
Yes		Mr.		§22	
Contact Surname:					
§22					

Supplier Information Section

Supplier Name:		Address 1:		Address 2:	
Town/Suburb:		Country:		Postcode:	
Phone:		Email:		Website:	
Supplier Informed:		Contact Title:		Contact First Name:	
Contact Surname:		Contact Fax:		Contact Email:	

Report Status

For website publication:	Ready for Publication:	Investigated:	Investigation Reason:	Team Assignment:	Team Priority:
Yes	Yes	No	Event determined to be an isolated one	Team B (Iib, Is & Im)	Not Investigated

Team Review

Reviewed by Team:	Reason Sent To Meeting:	Outcome from team meeting:	
Notes for Team meeting:			
Outcomes from Team Meeting:			

Initial Risk Analysis

Background Information	Risk Assessment - Section A	Risk Assessment - Section B	Risk Assessment - Section C	Risk Assessment - Section D	
Date:	Severity:	Incidents in the last 12 months:	Manufacturer analysis:	Assessor:	Manufacturer documentation:
16/10/2023	0 - Insufficient information to determine level of severity		Yes	§22	Unknown - updated information from the manufacturer is required
Incidents in last 24 months:	Manufacturer action:	ESTIMATED LEVEL OF INVESTIGATION:	FINAL LEVEL OF INVESTIGATION:	Injured Party:	Device Recalls:

	No	Level 1 Investigation (to complete screening)	Screening only	Patient	0. No recalls for similar incidents in Australia
Incidents in last 36 months:	IVD status:	EXCEPTION TO INVESTIGATION LEVEL:		Found Prior To Use:	Is AE covered by current recall:
		Isolated event		No	
Incidents Worldwide:	Number of potential contributing factors:			Reusable:	Similar events (past 6 months):
	No			No	0 incidents
Products supplied the last 12 months:	Specific factors identified:	ESTIMATED LEVEL OF PRIORITY:	FINAL LEVEL OF PRIORITY:		3 or more events - batch/model:
		Routine	Routine		
Products supplied last 24 months:	Number of potential sensitivities:	EXCEPTION TO PRIORITY LEVEL:			3 or more events - health district:
	No	Isolated event			
Products supplied last 36 months:	Specific sensitivities identified:				3 or more events - organisation:
Products supplied Worldwide:	Consultations during risk assessment:	Final Risk Assessment:			
	I did none of the above incidents	Yes			

Additional Risk Analysis

Click 'N' to start a new risk analysis

Analysis Details	Statistics Checklist Section					
Update Device Details?:	Background Information	Risk Assessment - Section A	Risk Assessment - Section B	Risk Assessment - Section C	Risk Assessment - Section D	
	Date:	Severity:	Incidents in the last 12 months:	Manufacturer analysis:		
Copy Data From:						
	Assessor:	Manufacturer documentation:	Incidents in last 24 months:	Manufacturer action:	ESTIMATED LEVEL OF INVESTIGATION:	FINAL LEVEL OF INVESTIGATION:
	Injured Party:	Device Recalls:	Incidents in last 36 months:	IVD status:	EXCEPTION TO INVESTIGATION LEVEL:	
	Found Prior To Use:	Is AE covered by current recall:	Incidents Worldwide:	Number of potential contributing factors:		
	Reusable:	Similar events (past 6 months):	Products supplied the last 12 months:	Specific factors identified:	ESTIMATED LEVEL OF PRIORITY:	FINAL LEVEL OF PRIORITY:
		3 or more events - batch/model:	Products supplied last 24 months:	Number of potential sensitivities:	EXCEPTION TO PRIORITY LEVEL:	
		3 or more events - health district:	Products supplied last 36 months:	Specific sensitivities identified:		
		3 or more events - organisation:	Products supplied Worldwide:	Consultations during risk assessment:	Final Risk Assessment:	

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Search Sponsors:	Name:		Client #:
45972	Medacta Australia Pty Ltd		45972
Attention To:	Address 1:	Address 2:	Town/Suburb:
	Unit A1 / 16 Mars Road		LANE COVE
State:	Postcode:	Phone:	Fax:
NSW	2066		

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Device Analysis Results:

Batch review performed on 5 July 2023
Lot 2227487: 75 items manufactured and released on 06-Feb-2023. Expiration date: 2027-Dec-12. No anomalies found related to the problem. To date, 16 items of the same lot have been sold with no similar reported event during the period of review.

Visual Inspection performed by R&D Project Manager

The breakage of the inserter occurred in correspondence to the change of section between the distal cylindrical portion and the thread, leaving the most distal pin (thread + cone) within the PEEK anchor. The breakage in this section suggests that the assembled inserter and anchor have been subjected to considerable lateral loads. As confirmed by the event report, the direction of insertion was not correspondant with the pilot hole axis, suggesting a misalignment between the inserter and the bone. Hammering the inserter in this configuration may result in an unforeseeable lateral hitting load condition, bringing to the weakest section breakage. A preliminary microscopic analysis of the broken section suggests an instantaneous failure, compatible with a hammering action, therefore excluding a cyclic load condition.
The constrained condition of the pilot hole prevented the anchor from failing (internally also supported by the metal pin) concentrating the loads on the interface section between the anchor and inserter.

Clinical Evaluation performed by Medical Affairs Department

While inserting an anchor in hip arthroscopy surgery, the surgeon reported the breakage of the introducer tip, after feeling hard insertion. As shown by the intraoperative X-ray at our disposal, the metallic tip of the instrument remained seated in the bone. Even though the material of the instrument is medical-grade stainless steel it does not belong to the series approved for long-term implantation and therefore we cannot guarantee that leaving the fragment in situ may not cause secondary problems. However, given that in this case the fragment is completely surrounded by bone, it is extremely unlikely that a secondary adverse effect would ever be triggered. Retrieving the fragment would undoubtedly be a very invasive procedure and therefore the surgeon has understandably decided to leave it in place.
It is difficult to determine the cause of this breakage, but the report indicates that the patient may have had sclerotic bone which made insertion difficult. With the information at hand, we cannot make a definitive conclusion.

Additional Details (use for tables):

Details of Similar Events:

MectaLock instruments 05.11.002 MectaLock PEEK Ø2.9 – Long driver Lot. 2227487: 1st case on the ref and lot.

Considering these and similar references, this is the 2nd case reported due to inserter breakage. However, the breakage of the previous complaint is different, it was a breakage at the welding point and lead to an additional surgery to remove the part. The occurrence rate is 0.26%.

Lot 2227487: 75 items manufactured and released on 06-Feb-2023. Expiration date: 2027-Dec-12. No anomalies found related to the problem. To date, 16 items of the same lot have been sold with no similar reported event during the period of review.

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Completed Actions:

Planned Actions and Proposed Timelines:

Additional Comments:

DIR closed 16/10/2023

Reason for Level 1 Investigation

Details of Reasons

Reason for Level 1 Investigation

Focus of Level 2 Investigation

Details of Focus

Essential Principles					If 'Other' Selected					Document 21	

Sources of Evidence for Level 2

Details of Source											
Sources of Evidence				If 'Others' please specify here			Expected Sourcing Date			Date of Evidence Received	

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Click **[N]** to begin a new Correspondence entry. Note that the Email address specified here will receive a notification if the Date Received is not filled in by the Date Expected.

Correspondence and Chronology Details										
Include?	Heading	Type L1	Type L2	Email	Sent	Expected	Received	Response	Notes	

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Problem Observed Details					
Problem Observed (Level 1)		Problem Observed (Level 2)		Problem Observed (Level 3)	If 'Other' Selected
Material Integrity Problem		Break		Material Fragmentation	
Appropriate Term/Code Not Available				Foreign body left in patient	

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Change in Therapeutic Response	Unexpected Therapeutic Effects		Document 21
--------------------------------	--------------------------------	--	-------------

Investigation Findings

Finding Details

Investigation Findings (Level 1)	Investigation Findings (Level 2)	Investigation Findings (Level 3)	If 'Other' Selected	
Patient Sample Problem				

Investigation Conclusion

Conclusion Details

Investigation Conclusion (L1)	Investigation Conclusion (L2)	If Additional Conclusion Detail Requested	
Cause Traced to User	Unintended Use Error Caused or Contributed to Event		

Investigation Outcomes

Outcome Details

Outcome of Investigation (L1)	Outcome of Investigation (L2)	If Additional Conclusion Detail Requested	
Reviewed, for Trending Purposes Only			

Investigation Summary

Latest Investigation (DII) where this DIR is the Primary DIR: 	Latest Investigation (DII) where this DIR is a Related DIR: 	Investigator: 	Peer Review: No	
Investigator's Notes: 	Summary Findings: No further investigation will occur at this time, however the TGA will continue to monitor the rate and pattern of occurrence and may re-open the file as appropriate.			Recall Number:

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Other Device (Entered): 	Brand Name: 	Manufacturer Name: 	Device ARTG #: 	
--	----------------------------	-----------------------------------	-------------------------------	--

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Samples Record - Click [N] to begin entering information. **Note:** Sample # Generated on Save.

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Patient Details							
Sex:		Weight:		Age:			
Patient Focused Corrective Action Taken:				Patient History:			
Injured - Extent of Injury:		Was device directly linked to death?:		Was device directly linked to permanent disability?:		Consequence:	
Other Consequence:		Describe any test (Lab, xray, etc.):		Additional Event Description:		Medical Problem Device Used For:	

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Flow Details : DIR-REQ - Device Incident Request : 392939

Request Details

ID	Type	Location	Status	Assigned By	Assigned To	Assigned On /	Priority	Attach
392939	DIR-REQ		Closed	\$22	OPR Administration User	16/10/2023	Normal	0

Signature Details

Role	IRIS Investigator	
User	\$22	
Signed At	16/10/2023 16:02:33	
Comment		

From: s22@medacta.ch>
Sent: Thursday, 10 February 2022 10:18 PM
To: s22
Cc: IRIS; s22; s22; s22@medacta.com.au; s22; s22
Subject: R: RESPONSE OVERDUE RE: Review of Medacta Australia Annual Report - response due COB 21 December 2021 [SEC=OFFICIAL]
Attachments: R: *Amended* DIR 68783 - Questionnaire Letter - Response due COB 27/04/2021 [SEC=OFFICIAL]

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Dear s22,
 I'm s22, s22 of Medacta International.
 I've personally followed the event involved in your communication and I would like to answer to the questions raised to Mr s22.
 As per the e-mail attached, the event and the incongruency between the DIR 68516 and DIR 68783 was reported to the Competent Authority in April 2021.
 As reported in your table, the event was communicated to the TGA by the sponsor as DIR 68516 and by the user as DIR 68783. They are related to the same event.
 The correct item involved in the event is the MectaLock PEEK Ø2.4 – Short driver reference 05.10.001 Lot. 1924147 that was wrongly reported by the user, as lot 1954843, in the DIR 68783.
 As per the above clarification, I would suggest to consider DIR 68783 as duplicate of the DIR 68516 that reports all the correct information related to this event.
 Do not hesitate to contact me if you have any issues.
 Kind regards,

s22 s22 <i>Regulatory Affairs Department</i> s22@medacta.ch www.medacta.ch	Medacta International SA Strada Regina CH-6874 Castel San Pietro Switzerland Pho: s22 Mob: s22 Fax: +41 91 696 60 66
---	--

Da: s22@health.gov.au>
Inviato: giovedì, 10 febbraio 2022 07:17
A: s22@medacta.com.au>
Cc: s22@medacta.com.au; s22@medacta.ch; s22@medacta.ch>
Oggetto: RE: RESPONSE OVERDUE RE: Review of Medacta Australia Annual Report - response due COB 21 December 2021 [SEC=OFFICIAL]
 Dear s22,
 Thank you for your email along with the attachment. The lot number and date of implant is different in the 2 DIRs as shown in the table below:

DIR Number	Report type	Lot number	Model number	Date of Implant
68783	User	1954843	Unknown	s22
68516	Sponsor	1924147	05.10.001	s22

Please confirm if the DIR 68516 is a duplicate of DIR 68783 so we can make it a duplicate in our system? If these are 2 different incidents with 2 different devices the DIR 68516 will not be considered a duplicate.

Please provide the response by **COB 17 February 2022**.

Thanks & Regards,

s22
s22

Devices Post Market Monitoring Section
Medical Devices Surveillance Branch

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Cc: s22@medacta.com.au>; s22@medacta.ch>; s22

s22@medacta.ch>

Subject: RE: RESPONSE OVERDUE RE: Review of Medacta Australia Annual Report - response due COB 21 December 2021 [SEC=OFFICIAL]

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Dear s22

A response to this enquiry was provided on 13/12/21. Please see attachment.

I will call shortly to clarify

Best regards

s22

From: s22

Sent: Thursday, 13 January 2022 3:29 PM

To: s22; s22

Cc: s22; s22; s22

Subject: RESPONSE OVERDUE RE: Review of Medacta Australia Annual Report - response due COB 21 December 2021 [SEC=OFFICIAL]

Dear s22

A response to below information request is overdue.

Please provide a response to this request soon as possible, or further actions may be taken in the next step.

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 **Please consider the environment before printing this email.**

From: s22 @health.gov.au>

Sent: Tuesday, 7 December 2021 4:15 PM

To: s22 @medacta.com.au>

Cc: s22 @medacta.com.au>; s22 @medacta.ch>; s22 @medacta.ch>

Subject: Review of Medacta Australia Annual Report - response due COB 21 December 2021 [SEC=OFFICIAL]

Dear s22

Thank you for your 2020-21 Annual Report submission.

I am reviewing your report.

A search in the IRIS database has revealed 2 DIRs (DIR 68783 and 68516) received for ARTG 320484 for the reporting period 1 July 2020 – 30 June 2021. I have noted only DIR 68516 is included within your annual report.

Action required by close of business 21 December 2021:

1. Please provide a justification on the discrepancy observed in DIR numbers for ARTG 320484.
2. Based on your answer to question 1, please advise if an amendment to reflect the correct adverse events (and their corresponding DIRs) number for ARTG 320484 is needed. If so, please provide the amended annual report.

Kind regards,

s22
 Departmental Officer
 Devices Post Market Monitoring
 Medical Devices Surveillance Branch

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 Please consider the environment before printing this email.

From: s22 [redacted] <[redacted]@medacta.com.au>

Sent: Tuesday, 28 September 2021 9:54 AM

To: IRIS <IRIS@health.gov.au>; IRIS <IRIS@health.gov.au>

Cc: s22 [redacted] <[redacted]@medacta.com.au>; s22 [redacted] <[redacted]@medacta.ch>; s22 [redacted]

s22 [redacted] <[redacted]@medacta.ch>

Subject: Medacta Australia Annual Report

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Dear Sir/madam,

Please find attached the 2021 Annual report for Medacta Australia.

Should you require further information please do not hesitate to contact me .

Yours faithfully

s22 [redacted]

For and on behalf of Medacta Australia Pty Ltd

IMPORTANT NOTICE

The content of this email is confidential and is for exclusive use of the intended recipient(s). If you are not the intended recipient(s) please note that any form of distribution, copying or use of this communication or the information in it is strictly forbidden and may be unlawful.

Please note that, although this email has been virus-checked, neither Medacta International SA nor the sender accepts any responsibility for

viruses, and it is your responsibility to scan the email and attachments (if any).

If you believe you have received this communication in error, please return it (with the Object "received in error") to the sender, and then delete the email and destroy any copies of it.

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<https://www.mailguard.com.au>

s22

From: s22@medacta.com.au>
Sent: Thursday, 22 September 2022 8:33 PM
To: IRIS
Cc: s22; s22; s22; s22
Subject: Medacta Australia Annual Report
Attachments: 2022 Annual Report - Medacta Australia.xlsx

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Dear Sir/Madam,

Please find attached the annual report of Medacta Australia Pty Ltd.

Please do not hesitate to contact me should you require further information

Kind Regards

s22

s22

s22

s22@medacta.com.au
www.medacta.com

Medacta Australia

Unit A1, 16 Mars Road
Lane Cove, NSW 2066
Australia

Pho: s22
Mob: s22
Fax: +612 9420 2578

IMPORTANT NOTICE

The content of this email is confidential and is for exclusive use of the intended recipient(s). If you are not the intended recipient(s) please note that any form of distribution, copying or use of this communication or the information in it is strictly forbidden and may be unlawful.

Please note that, although this email has been virus-checked, neither Medacta International SA nor the sender accepts any responsibility for viruses, and it is your responsibility to scan the email and attachments (if any).

If you believe you have received this communication in error, please return it (with the Object "received in error") to the sender, and then delete the email and destroy any copies of it.

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<https://www.mailguard.com.au>


Class	Start Date	ARTG ID	Sponsor	Device Category	GMDN	Term	Intended Purpose	Sponsor ID	Address (line 1)	Address (line 2)	Town	State	Postcode	Manufacturer ID	Manufacturer Name
s22															
s22	19/07/2019	320484	Medacta Australia Pty Ltd	Non-active implantable devices	36174	Ligament anchor	The Meralock PEEK Suture anchors are intended for use in arthroscopic or open surgical approaches for fixation of su	45972	Unit 61 / 15 Mars Road		LANE COVE	NSW	2066	45296	Medacta International SA

ARTG Number	Product name (description)	Model Number (or Name)	Quantity supplied in Australia	Quantity supplied Worldwide (incl Australia)	Number of complaints (Aust / WW)	Number of Adverse Events (Aust / WW)	Comment relevant qualifications relating to data
Data from the 1st July 2020 to the 30th June 2021 - WW includes Australia							
320484	Ligament anchor	MectaLock PEEK Suture Anchors	S4 /		2/3	0/0	

ARTGNumber	Type of adverse events	Number Total (AU included)	Number AU	Percentage world Wide	Percentage in Australia	TGA DIR Number	Regulatory action	Comments
	(description and code)					if relevant	regulatory/corrective action/notification by the manufacturer	relevant qualifications relating to type of issue and actions

s22

ARTG Number	Type of complaint	Number Total (AU included)	Number AU	Percentage world Wide	Percentage in Australia	TGA DIR Number	Regulatory action	Comments
	(description and code)					if relevant	regulatory/corrective action/notification by the manufacturer	relevant qualifications relating to type of issue and actions
320484	Failure of press-fit	1	0	0.75%	0.00%	/	Nil	no related defect found
	Implant breakage	3	2	2.24%	3.77%	/	Nil	no related defect found

s22

From: s22 @health.gov.au>
Sent: Thursday, 10 February 2022 5:17 PM
To: s22
Cc: s22; s22; s22
Subject: RE: RESPONSE OVERDUE RE: Review of Medacta Australia Annual Report - response due COB 21 December 2021 [SEC=OFFICIAL]
Attachments: RE: Review of Medacta Australia Annual Report - response due COB 21 Decem... (264 KB)

Dear s22

Thank you for your email along with the attachment. The lot number and date of implant is different in the 2 DIRs as shown in the table below:

DIR Number	Report type	Lot number	Model number	Date of Implant
68783	User	1954843	Unknown	s22
68516	Sponsor	1924147	05.10.001	s22

Please confirm if the DIR 68516 is a duplicate of DIR 68783 so we can make it a duplicate in our system? If these are 2 different incidents with 2 different devices the DIR 68516 will not be considered a duplicate.

Please provide the response by **COB 17 February 2022**.

Thanks & Regards,

s22
s22
Devices Post Market Monitoring Section
Medical Devices Surveillance Branch

Phone: s22
Email: s22 @health.gov.au

Therapeutic Goods Administration
Department of Health
PO Box 100
Woden ACT 2606
www.tga.gov.au

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From: s22 @medacta.com.au>
Sent: Thursday, 13 January 2022 3:46 PM
To: s22 @health.gov.au>
Cc: s22 @medacta.com.au>; s22 @medacta.ch>; s22 @medacta.ch>
Subject: RE: RESPONSE OVERDUE RE: Review of Medacta Australia Annual Report - response due COB 21 December 2021 [SEC=OFFICIAL]

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Dear s22
 A response to this enquiry was provided on 13/12/21. Please see attachment.
 I will call shortly to clarify
 Best regards
 s22

From: s22
Sent: Thursday, 13 January 2022 3:29 PM
To: s22 ; s22
Cc: s22 ; s22 ; s22
Subject: RESPONSE OVERDUE RE: Review of Medacta Australia Annual Report - response due COB 21 December 2021 [SEC=OFFICIAL]

Dear s22
 A response to below information request is overdue.
 Please provide a response to this request soon as possible, or further actions may be taken in the next step.
 Kind regards,
 s22
 Departmental Officer
 Devices Post Market Monitoring
 Medical Devices Surveillance Branch

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 Please consider the environment before printing this email.

From: s22 [REDACTED] <[REDACTED]@health.gov.au>

Sent: Tuesday, 7 December 2021 4:15 PM

To: s22 [REDACTED] <[REDACTED]@medacta.com.au>

Cc: s22 [REDACTED] <[REDACTED]@medacta.com.au>; s22 [REDACTED] <[REDACTED]@medacta.ch>; s22 [REDACTED] <[REDACTED]@medacta.ch>

Subject: Review of Medacta Australia Annual Report - response due COB 21 December 2021 [SEC=OFFICIAL]

Dear s22 [REDACTED]

Thank you for your 2020-21 Annual Report submission.

I am reviewing your report.

A search in the IRIS database has revealed 2 DIRs (DIR 68783 and 68516) received for ARTG 320484 for the reporting period 1 July 2020 – 30 June 2021. I have noted only DIR 68516 is included within your annual report.

Action required **by close of business 21 December 2021**:

1. Please provide a justification on the discrepancy observed in DIR numbers for ARTG 320484.
2. Based on your answer to question 1, please advise if an amendment to reflect the correct adverse events (and their corresponding DIRs) number for ARTG 320484 is needed. If so, please provide the amended annual report.

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 Please consider the environment before printing this email.

From: s22 [REDACTED] <[REDACTED]@medacta.com.au>

Sent: Tuesday, 28 September 2021 9:54 AM

To: IRIS <IRIS@health.gov.au>; IRIS <IRIS@health.gov.au>

Cc: s22 [REDACTED] <[REDACTED]@medacta.com.au>; s22 [REDACTED] <[REDACTED]@medacta.ch>; s22 [REDACTED] <[REDACTED]@medacta.ch>

Subject: Medacta Australia Annual Report

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Dear Sir/madam,

Please find attached the 2021 Annual report for Medacta Australia.

Should you require further information please do not hesitate to contact me .

Yours faithfully

s22

For and on behalf of Medacta Australia Pty Ltd

s22 le

s22

s22 @medacta.com.au

www.medacta.com**Medacta Australia**

Unit A1, 16 Mars Road

Lane Cove, NSW 2066

Australia

Pho: s22

Mob: s22

Fax: +612 9420 2578

IMPORTANT NOTICE

The content of this email is confidential and is for exclusive use of the intended recipient(s). If you are not the intended recipient(s) please note that any form of distribution, copying or use of this communication or the information in it is strictly forbidden and may be unlawful.

Please note that, although this email has been virus-checked, neither Medacta International SA nor the sender accepts any responsibility for

viruses, and it is your responsibility to scan the email and attachments (if any).

If you believe you have received this communication in error, please return it (with the Object "received in error") to the sender, and then delete the email and destroy any copies of it.

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<https://www.mailguard.com.au>

s22

From: s22@medacta.com.au>
Sent: Monday, 13 December 2021 1:57 PM
To: s22
Cc: s22; s22; s22
Subject: RE: Review of Medacta Australia Annual Report - response due COB 21 December 2021 [SEC=OFFICIAL]
Attachments: R: *Amended* DIR 68783 - Questionnaire Letter - Response due COB 27/04/20... (215 KB)

Dear Miss s22

1. DIR 68783 and 68516 refer to the same adverse event # 0172-21. In summary, we had originally reported the event as DIR 68516, and in parallel the TGA also received a notification from the hospital and registered the event as DIR 68783. Attempt to clarify the discrepancy had subsequently been made to the TGA by Medacta. (please see attached correspondence dated 26/4/21) On the TGA website DIR 68783 is the only available reference via the Medacta account.
2. Based on the response provided in point 1 above, Medacta do not propose any amendments to the annual report

Should further clarification be required, please do not hesitate to contact me

Thank you and best regards

s22

s22	Medacta Australia
s22	Unit A1, 16 Mars Road
	Lane Cove, NSW 2066
	Australia
s22@medacta.com.au	Pho: s22
www.medacta.com	Mob: s22
	Fax: +612 9420 2578

IMPORTANT NOTICE

The content of this email is confidential and is for exclusive use of the intended recipient(s). If you are not the intended recipient(s) please note that any form of distribution, copying or use of this communication or the information in it is strictly forbidden and may be unlawful.

Please note that, although this email has been virus-checked, neither Medacta International SA nor the sender accepts any responsibility for viruses, and it is your responsibility to scan the email and attachments (if any).

If you believe you have received this communication in error, please return it (with the Object "received in error") to the sender, and then delete the email and destroy any copies of it.

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From: s22
Sent: Tuesday, 7 December 2021 4:15 PM
To: s22
Cc: s22 ; s22 ; s22
Subject: Review of Medacta Australia Annual Report - response due COB 21 December 2021 [SEC=OFFICIAL]

Dear s22

Thank you for your 2020-21 Annual Report submission.
I am reviewing your report.

A search in the IRIS database has revealed 2 DIRs (DIR 68783 and 68516) received for ARTG 320484 for the reporting period 1 July 2020 – 30 June 2021. I have noted only DIR 68516 is included within your annual report.

Action required by close of business 21 December 2021:

1. Please provide a justification on the discrepancy observed in DIR numbers for ARTG 320484.
2. Based on your answer to question 1, please advise if an amendment to reflect the correct adverse events (and their corresponding DIRs) number for ARTG 320484 is needed. If so, please provide the amended annual report.

Kind regards,

s22
Departmental Officer
Devices Post Market Monitoring
Medical Devices Surveillance Branch

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 Please consider the environment before printing this email.

From: s22 [REDACTED] <[REDACTED]@medacta.com.au>
Sent: Tuesday, 28 September 2021 9:54 AM
To: IRIS <IRIS@health.gov.au>; IRIS <IRIS@health.gov.au>
Cc: s22 [REDACTED] <[REDACTED]@medacta.com.au>; s22 [REDACTED] <[REDACTED]@medacta.ch>; s22 [REDACTED] <[REDACTED]@medacta.ch>
Subject: Medacta Australia Annual Report

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Dear Sir/madam,

Please find attached the 2021 Annual report for Medacta Australia.

Should you require further information please do not hesitate to contact me .

Yours faithfully

s22 [REDACTED]
For and on behalf of Medacta Australia Pty Ltd

s22 [REDACTED]	Medacta Australia
s22 [REDACTED]	Unit A1, 16 Mars Road
	Lane Cove, NSW 2066
	Australia
s22 [REDACTED] <[REDACTED]@medacta.com.au>	Pho: s22 [REDACTED]
www.medacta.com	Mob: s22 [REDACTED]
	Fax: +612 9420 2578

IMPORTANT NOTICE

The content of this email is confidential and is for exclusive use of the intended recipient(s). If you are not the intended recipient(s) please note that any form of distribution, copying or use of this communication or the information in it is strictly forbidden and may be unlawful.

Please note that, although this email has been virus-checked, neither Medacta International SA nor the sender accepts any responsibility for viruses, and it is your responsibility to scan the email and attachments (if any).

If you believe you have received this communication in error, please return it (with the Object "received in error") to the sender, and then delete the email and destroy any copies of it.

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[Report this message as spam](#)

s22

From: s22@medacta.com.au>
Sent: Tuesday, 28 September 2021 9:54 AM
To: IRIS; IRIS
Cc: s22; s22; s22
Subject: Medacta Australia Annual Report
Attachments: 2021 Annual Report - Medacta Australia.xlsx

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Dear Sir/madam,

Please find attached the 2021 Annual report for Medacta Australia.

Should you require further information please do not hesitate to contact me .

Yours faithfully

s22

For and on behalf of Medacta Australia Pty Ltd

s22

Medacta Australia

Unit A1, 16 Mars Road
Lane Cove, NSW 2066
Australia

s22

Pho: s22

s22@medacta.com.au

Mob: s22

www.medacta.com

Fax: +612 9420 2578

IMPORTANT NOTICE

The content of this email is confidential and is for exclusive use of the intended recipient(s). If you are not the intended recipient(s) please note that any form of distribution, copying or use of this communication or the information in it is strictly forbidden and may be unlawful.

Please note that, although this email has been virus-checked, neither Medacta International SA nor the sender accepts any responsibility for viruses, and it is your responsibility to scan the email and attachments (if any).

If you believe you have received this communication in error, please return it (with the Object "received in error") to the sender, and then delete the email and destroy any copies of it.

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<http://www.mailguard.com.au>

[illegible]

ARTGNumber	Product name (description)	Model Number (or Name)	Quantitysuppl iedinAustralia	Quantitysuppl iedWorldwide	Number of complaints	Number of Adverse Events	Comment
				(inclAustralia)	(Aust / WW)	(Aust / WW)	relevant qualifications relating to data
Datafromthe1stjuly2020tothe30thJune2021- WWincludesAustralia							
320484	Ligament anchor	MectaLock PEEK Suture Anchors	S4 /		0/6	1/1	


ARTGNumber	Type of adverse events	Number Total (AU included)	Number AU	Percentage world Wide	Percentage in Australia	TGA DIR Number	Regulatory action	Comments
	(description and code)					if relevant	regulatory/corrective action/notification by the manufacturer	relevant qualifications relating to type of issue and actions

s22

320484	<i>Instrument breakage</i>	<i>1</i>	<i>1</i>	<i>0.82%</i>	<i>1.43%</i>	<i>68516</i>	<i>Nil</i>	<i>no action - no trend observed</i>
---------------	----------------------------	----------	----------	--------------	--------------	--------------	------------	--------------------------------------

s22


s22

ARTG Number	Type of complaint	Number Total (AU included)	Number AU	Percentage world Wide	Percentage in Australia	TGA DIR Number	Regulatory action	Comments
	(description and code)					if relevant	regulatory/corrective action/notification by the manufacturer	relevant qualifications relating to type of issue and actions
S22								
320484	Implant breakage	4	0	3.28%	0.00%	/	Nil	no related defect found
	Implant deformation	2	0	1.64%	0.00%	/	Nil	no related defect found

s22

s22

s22

From: s22@medacta.com.au>
Sent: Monday, 28 September 2020 11:17 AM
To: IRIS
Cc: s22; s22; s22
Subject: Medacta Australia 2020 Annual Report
Attachments: 2020 Annual Report - Medacta Australia.xlsx

REMINDER : Think before you click! This email originated from outside our organisation. Only click links or open attachments if you recognise the sender and know the content is safe.

Dear Sir/madam,

Please find attached the 2020 Annual report for Medacta Australia Pty Ltd.

Should you require further information please do not hesitate to contact me .

Yours faithfully

s22

s22

s22

s22@medacta.com.au
www.medacta.com

Medacta Australia
 Unit A1, 16 Mars Road
 Lane Cove, NSW 2066
 Australia
 Pho: s22
 Mob: s22
 Fax: +612 9420 2578

IMPORTANT NOTICE

The content of this email is confidential and is for exclusive use of the intended recipient(s). If you are not the intended recipient(s) please note that any form of distribution, copying or use of this communication or the information in it is strictly forbidden and may be unlawful.

Please note that, although this email has been virus-checked, neither Medacta International SA nor the sender accepts any responsibility for viruses, and it is your responsibility to scan the email and attachments (if any).

If you believe you have received this communication in error, please return it (with the Object "received in error") to the sender, and then delete the email and destroy any copies of it.

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
ARTG Number	Product name (description)	Model Number (or Name)	Quantity supplied in Australia	Quantity supplied Worldwide (incl Australia)	Number of complaints (Aust / WW)	Number of Adverse Events (Aust / WW)	Comment relevant qualifications relating to data
Data from the 1st July 2019 to the 30th June 2020 - WW includes Australia							
s22							
320484	Ligament anchor	MectaLock PEEK Suture Anchors	s4 /		0/4	0/0	
s22							

ARTG Number	Type of adverse events (description and code)	Number Total (AU included)	Number AU	Percentage world Wide	Percentage in Australia	TGA DIR Number if relevant	Regulatory action regulatory/correc tive action/notification by the manufacturer	Comments relevant qualifications relating to type of issue and actions
----------------	--	-------------------------------	-----------	--------------------------	----------------------------	-----------------------------------	--	--


s22

--	--	--	--	--	--	--	--	--

s22

s22

320484	Deformation	1	0	0.694%	0.000%	/	Nil	no related defect found
	Implant breakage	1	0	0.694%	0.000%	/	Nil	no related defect found
	Implant deformation	1	0	0.694%	0.000%	/	Nil	no related defect found
	Implants assembling	1	0	0.694%	0.000%	/	Nil	no related defect found