

Serial Number Location Chart

| Example | Serial Number Type | Interpretation |
|------------|--------------------------|--|
| 11DE003250 | Invacare - Serial Number | 1st 2 digits represent year, first letter represents month |

| Product | Serial #/Date Code | Location |
|---------------------------|---------------------|---|
| Nutron Series | Invacare - Serial # | On the front of the frame - right side. |
| At'm (Discontinued) | Invacare - Serial # | On the front of the frame - left side. |
| Power 9000 (Discontinued) | Invacare - Serial # | On the front of the frame - right side. |
| Power Tiger | Invacare - Serial # | On the bottom right crossbar. |
| P9000XDT | Invacare - Serial # | On the front of the frame - right side. |
| Pronto M41 | Invacare - Serial # | Under arm pad right side. |
| Pronto M71 (Discontinued) | Invacare - Serial # | underneath seat in rear, should not have to remove seat |
| Pronto M50/M51 | Invacare - Serial # | On the right side of the battery box |
| Pronto M61 | Invacare - Serial # | On rear of frame by the headtube-right side. |
| Pronto M91 | Invacare - Serial # | On the inside of the right rear frame |
| Pronto M94 | Invacare - Serial # | On the inside of the right rear frame. |
| Storm RWD | Invacare - Serial # | On the rear of the frame - right side. |
| TDX 3/4/5 (Discontinued) | Invacare - Serial # | Between rear casters/remove shroud/look for tag near floor |
| TDXSI-HD | Invacare - Serial # | Inside of right rear swing arm, near the caster. |
| TDX-SI | Invacare - Serial # | Inside of right rear swing arm. |
| TDX-SC | Invacare - Serial # | Inside of right rear swing arm. |
| TDX-SR | Invacare - Serial # | Inside of right rear swing arm. |
| TDX-SP | Invacare - Serial # | Inside of right rear swing arm. |
| TDX-Spree | Invacare - Serial # | Inside of right rear swing arm. |
| Xterra (Discontinued) | Invacare - Serial # | Under the seat on the left side as you sit in the chair. |
| FDX (Front Wheel Drive) | Invacare - Serial # | On the right side of the battery box (behind the drive wheel) |

