



**CERTIFICATE OF ANALYSIS**  
**FITkit® TESTING SERVICES**

No QA090128

Customer:

**Therapeutic Goods Administration**  
**Office of Laboratories and Scientific Services**  
136 Narrabundah Lane, Symonston ACT 2609, Australia

Description of sample: syringe (inner surface)  
Arrival date: 19 October 2009  
Testing method: Capture enzyme immunoassay method (EIA; IEMA)  
Compliant with: EN 455-3:2006; ASTM D7427 - 08

**Note:** On request of customer the extraction step of testing process was altered. Extraction solution (0.5 mL of PBS solution, +4°C) was withdrawn into syringe and ejected immediately. Extraction solutions of three syringes were mixed into one sample.

**Test Results**

Icosagen code	Sample code	Hev b 1 µg/L	Hev b 3 µg/L	Hev b 5 µg/L	Hev b 6.02 µg/L	Total µg/L
LAT09106	sample 1	UD	UD	UD	UD	UD

UD - Undetectable

Detection limits:

Hev b 1 < 0.050 µg/g; Hev b 3 < 0.050 µg/g; Hev b 5 < 0.025 µg/g; Hev b 6.02 < 0.025 µg/g

30 October 2009

  
-----  
Ranno Rätsep  
Proxy Quality Manager

As of 02 March 2009, Icosagen AS is the new business name of the former Quattromed AS.

**Note:**

The results related only to the tested item.

The test certificate shall not be reproduced except in full, without written approval of Icosagen AS laboratory.

Icosagen AS

Nooruse 9, 50411 Tartu, Estonia

Phone +372 737 7070, Fax +372 737 7077, info@icosagen.com

Reg No: 10565884, VAT No: EE100573055

www.icosagen.com

