

Status

Row Labels

14/16 Metasul Heads
14/16 Taper Head - 14/16 Taper Head
1st Q Add On Sulcus Intraocular Lenses - 1st Q Add On Sulcus IOL
2-butyl cyanoacrylate (Variclose)
3-0 Supramid Extra Suture
4350-45 Intramuscular Electrode (Interstim)
ABI electrode and placing system
Active Knee System
Ad Tech Spencer Probe Depth Electrode - Ad Tech Spencer Probe Depth Electrode
Advantim Knee System - Advantim Knee System
Advantim tibial insert unconstrained
Aequalis Perform+ Glenoid
AGC Knee Tradition Bearings
Agluna Femoral Shaft
allergEAZE patch test chambers
AMK Knee
Amplia CRT-D - Amplia CRT-D
Aniridia Intraocular Lens model 68
APR Liner - APR Liner
Aptus Heli FX system - Aptus
ArthroSurface HemiCAP
ArthroSurface Patello Femoral
Assianchor - Hanita Lenses
Atriclep Pro
AtriClip®; Left Atrial Appendage Exclusion System
Attain Stability LV Pacing Leads
Axioma Modular Metal Back
Bardex IC Anti-Infective Latex Foley Catheter - Bardex IC Anti-Infective Latex Foley Catheter
Beflex MRI safe leads - Beflex
Biodegradable co Polymer - InSpace System
Biodegradable Temporising Matrix
Biodesign Enterocutaneous Fistula Plug
Biomet Comprehensive Shoulder Titanium Versa-Dial Humeral Head
Biotechni Filler Stem
BlastGen - Origio
Bollard Skeletal Anchor
Bone Plate - Jiel Medical Corp
Boston Scientific IntellaNAV - Boston Scientific IntellaNAV
Braun Alcohol Concentrate 95% (Ethanol 95%)
Calcium Chloride Dihydrate
Calcium Hydroxylapatite - Radiesse voice
Calcium Ionophore - GM508 CultActive
Cardiomed Cooling Jacket
CARILLION Mitral Contour System
CeraRoot Dental Implants
COAPT Complete Control - COAPT
CollaPlug

CollaPlug - CollaPlug
Collatape
Collatape - CollaPlug - CollaCote
Comprehensive Reverse Shoulder Titanium Glenosphere
CoreValve Evolut R system
CX Infinion Lead
Cystistat
Deep Dish Classic Bearings - LCS
Deflux
Dehydrated Alcohol Injection - Ethanol 96%
Dehydrated Alcohol Injection (aka Ethanol 98%, aka Alcohol 98%, aka Ethyl Alcohol 98%, Alcohol 100%)
DePuy Acetabular Re-inforcement Rings
Dixi Medical Cortical Electrodes
Dixi Medical Depth Electrodes and Instruments
Double Array Cochlear Implant
Edora 8 HF-T QP - Edora 8 HF-T QP
Edwards sapien 3 THV
eFlow Electronic Nebuliser
EmbryoGen
EmbryoGen Ver 2
Emperion Revision Stem
Endotine Forehead/Brow
Endotine Midface
Endotine Mid-face
Endovascular Stent Graft System - Aorfix
Endrys Transseptal Needle
Enterra Gastric Unipolar Leads #4351
Epi Clear Corneal Epithelium Removal Device
ESOP Hip Implant
Ethanol 100%
Ethicon PDS Tape
Evolution Adaptive Inserts
Fermatron (Hyaluronan)
Fin II Insert XLPE 36mm 5 56/58 (Part # 110361356)
Flexitouch PD32-U - Tactile Medical
Freestyle Libre Metre
FSP Studio software
Genesis 1 Articular Insert
Global Steptech Glenoid & Instruments
Gold Eyelid Weight (Meddev)
Graft Jacket
Guidezilla Guide Extension Catheter
Harris femoral head
Heli Fx Endoanchor
HemiCAP PF Wave
Hexloc acetabular liner
Hintegra Total Ankle Replacement Polyethylene liner insert
Human Dermal Allograft - Arthroflex Dermal Allograft
Hyaluronic Acid (Cystistat)
Hypertonic saline contact lens

ICX- Implant
Ilex Skin Protectant Paste
Implantable PHAKIC Contact Lens
In-Pact Falcon Paclitaxel Eluting Balloon
InSpace System
IntellaMap Orion
Intellijoint Hip
Intellijoint Hip System
Interax ISA Tibial Mobile Bearing Polyethylene Insert
Intergard Synergy Vascular Graft
Intra-lock dental implant system
Intraocular Contact Lens (Implantable Contact Lens ICL)
IntraSpin System Single 220volts / 50 HZ - IntraSpin System Single 220volts / 50 HZ
Iodine-125 Seed (Echoseed)
Jupiter Sternal Protection Device
Juvenile cartilaginous allograft tissue - DeNovo NT
Kinemax Plus
Kitazato Cryotop Safety Kit
KX+ DURAT LS insert L 8mm, 10mm
KX+ DURAT LS insert M 10mm, 12mm &
Lactosorb Midface Distractor - Paediatric
Legion revision tibial baseplate
M8 Lead Apdaptor
Macrolane
Magnum pressfit acetabular cup
MatriDerm
Maxim and AGC tibial bearing
Maxim tibial bearing
MBK poly liner
Melody Transcatheter Pulmonary Valve & Delivery System
MG Unicompartmental Knee Articular Surface
Mi1200 SYNCHRONY PIN (Flex 28)
MIAR Cannula 11014L
Micra Transcatheter Pacing System
MiRCSP Cannula
Modular Acetabular System - Ringloc 40mm 10deg liner
Monoblock stem
Morcher Aniridia Implant 67B
Morcher Aniridia IOL
Morcher aniridia IOL Model 68
Morcher Type 7D PMMA Sulcus Lens
Morcher XtraFocus Intraocular Implant
Motiva Implant Matrix
MP3 Osteobiol
Nagor Mammary Implants Saline filled textured
Nagor Type SFX-HP Saline filled textured - high profile
Nellix
Neuroform Atlas Stent System
Nexgen CRA System
Nickel sensitive hinge

Nucleus 24 Double Array Cochlear Implant
Ologen
Opti-fix acetabular cup liner - Opti-fix acetabular cup liner
Optifix Stem
OSS IM Total Femoral Rod
Osteogen
Osteograf N
Osteolock
Osteonics Series II
Ova Cryo Kit Type M
Oxinium Femoral Head 14/16 Taper
paclitaxel
Paclitaxel - Eluting Balloon (In-Pact Pacific)
Pantera LUX
Passeo-18 Lux Paclitaxel - Releasing PTA Balloon
PCA tibial insert implant - Duracon
Persona TKR
PFC Curved Tibial Insert
PFC Total Knee Replacement
Platinum VR 1210 (Single Chamber Implantable Cardioverter Defibrillator (ICD))
PMMA Implantable Intraocular Lens KL-100 (Anterior) - Mediphacos PMMA Intraocular Lens
Porex (TiNbN) Coated Femoral Component (SLED)
Portico Transcatheter Aortic Valve Implant and SJM TAVI Transfemoral Delivery System
Precision M8 Adapter
Precision Novi
Prime Implant Analog
ProOsteon
Prostalac (Prosthesis Antibiotic Loaded Acrylic Cement) Hip
Proximal Loading Revision Hip System - PLR
PyC Aequalis Humeral Head
QUINTIP Individual Skin Test Device
Radiesse Injectable Implant
Randelli short body and liner
Re Motion Total Wrist Implant System Xray templates
RENEXUS
RESOLUTE ONYX STENT
Revision Articular Insert BCS
ROCC tibial implants - ROCC tibial implants
RT-Plus Tibial Insert
SafeSept Transseptal Guidewire - Model No. SS-135 - SafeSept
SensoriMotor Active Rehabilitation Training Arm - SMART
Sheathless Eucath
Sic Invent Dental Implant
Silikon 1000 (purified polydimethylsiloxane)
SMR - L2 Liners
SMR Metal Back Liner 6mm
SMR Stemless
SMR TT Metal Back Uncemented Glenoid - Lima Orthopaedics
St Jude Medical Portico Valve
Staar Implantable Contact Lens (ICL) ICM

Staar Toric (ICL)
Stainless Steel Omega Shaped Spring
Stryker PCA Knee Insert - Stryker PCA Knee Insert
Stryker Vitalock Liner
SuperCross Microcatheter - Vascular Solutions
Surgiform
Synovasure
Tapered Precision Impression Transfer
Thoraflex Hybrid
Tigris Stent
Titanium Breathe Implant - a Wengen
Total Ankle Replacement Agility Tibial Insert
Trabecular Metal Glenoid Augment Device
TYRX Absorbable Antibacterial Envelope
Uni elbow and rHead implant system xray templates
VAIOS Shoulder Replacement - VAIOS Shoulder Replacement
Vercise M8 adapter
Vibrant Soundbridge System
Vicryl Tape
Visian Implantable Contact Lens
Visionist X4
VNS Therapy Aspire SR Pulse Generator
Xen Gel Stent
Xenform Soft Tissue Repair Matrix - Mesh, Surgical
Zenker's diverticulum overtube - Zenker's diverticulum overtube (ZDO-22-30)
Zimmer MGII - Knee - Tibial Bearing - Zimmer MGII
Zimmer trilogy locking ring
ZSI 100 FTM Penile Implant

Grand Total

Approved