



**Australian Government**

**Department of Health and Ageing**  
Therapeutic Goods Administration

# ECRI codes

ECRI Universal Medical Device Nomenclature  
System cross referenced to Australian Device  
Groups (ADG)

March 1998

**TGA** Health Safety  
Regulation



## About the Therapeutic Goods Administration (TGA)

- The TGA is a division of the Australian Government Department of Health and Ageing, and is responsible for regulating medicines and medical devices.
- TGA administers the *Therapeutic Goods Act 1989* (the Act), applying a risk management approach designed to ensure therapeutic goods supplied in Australia meet acceptable standards of quality, safety and efficacy (performance), when necessary.
- The work of the TGA is based on applying scientific and clinical expertise to decision-making, to ensure that the benefits to consumers outweigh any risks associated with the use of medicines and medical devices.
- The TGA relies on the public, healthcare professionals and industry to report problems with medicines or medical devices. TGA investigates reports received by it to determine any necessary regulatory action.
- To report a problem with a medicine or medical device, please see the information on the TGA website.

### **Copyright**

© Commonwealth of Australia 1998

This work is copyright. Apart from any use as permitted under the Copyright Act 1968, no part may be reproduced by any process without prior written permission from the Commonwealth. Requests and inquiries concerning reproduction and rights should be addressed to the Commonwealth Copyright Administration, Attorney General's Department, National Circuit, Barton ACT 2600 or posted at <http://www.ag.gov.au/cca>

---

**ECRI UMDNS DEVICE TERM****ECRI IMDC**    **ADG CODE****A**

Accessories, Endoscopes	*****	ENDOAC
Acupuncture Kits	10014	ACPTEQ
Adapter Cables, Pressure Transducer	16969	TRDSNS
Adapter/Connectors, Angioplasty Balloon Dilation Catheters	17787	ANGPDV
Adaptometers	10024	OPHDEQ
Adenotomes	10025	SGINAC
Adhesive Strips, Butterfly	10278	BDDRST
Adhesive Strips, General Purpose	16769	BDDRCL
		BDDRST
Adhesive Strips, Hypoallergenic	10028	BDDRST
Adhesive Strips, Waterproof	10029	BDDRCL
		BDDRST
Adhesives, Aerosol	10035	ADHMED
Adhesives Enhancers	*****	ADHMED
Adhesives, Liquid	10036	ADHMED
Aerosol Generators	10046	ATINNB
Air Detectors, Heart-Lung Bypass	15010	CPBYEQ
Airway Obstruction Removal Devices, Emergency	10058	SGINAC
Airways, Esophageal Obturating	10056	AWPSDA
Airways, Esophageal Obturating/Tracheal	17986	AWPSDA
Airways, Nasopharyngeal	10057	AWPSDA
Airways, Oropharyngeal	10059	AWPSDA
Alarms	15550	IVASAC
Alarms, Enuresis	11588	INCOEQ
Alarms, Oxygen Depletion	12879	MGASEQ
Alloys, Dental	10077	MTDTRS
Ambulation Support Units	10086	AMBEQP
Amniocentesis Kit	15274	AMNOKT
Amnioscopes	10088	ENDOSC
Amniotomes	10089	SGINAC
Analgesia Units	15014	AGUNAC
Analgesia Units, Audio	16954	AGUNAC
Analgesia Units, Cryogenic	17140	AGUNAC
Analgesia Units, Inhalation	16953	AGUNAC
Anchor, Fixation	*****	GRFIMP
Anesthesia Kits, Brachial Plexus	10125	ANESKT
Anesthesia Kits, Caudal	10126	ANESKT
Anesthesia Kits, Epidural	10127	ANESKT
Anesthesia Kits, Glossopharyngeal	10128	ANESKT
Anesthesia Kits, Laryngeal-Tracheal	16827	ANESKT
Anesthesia Kits, Other	15202	ANESKT
Anesthesia Kits, Paracervical	10129	ANESKT

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Anesthesia Kits, Pudendal	10130	ANESKT
Anesthesia Kits, Saddle Block	16776	ANESKT
Anesthesia Kits, Spinal	10131	ANESKT
Anesthesia Unit Absorbers, Carbon Dioxide	10140	ANEQAC
Anesthesia Unit Absorbers, Halogenated Vapour	17615	ANEQAC
Anesthesia Unit Carbon Dioxide Absorbents	17509	ANEQAC
Anesthesia Unit Gas Scavengers	10142	ANEQAC
Anesthesia Unit Vaporizers	10144	ANEQAC
Anesthesia Unit Ventilators	10145	ANEQAC
Anesthesia Units	10134	ANEQAC
Anesthesimeters	10148	ANEQAC
Angiographic Kits, Digital	16546	RADACC
Angiographic Kits, Special Procedure	16547	RADACC
Angioplasty System	17558	ANGPDV
Angioplasty Systems, Atherectomy	17559	ANGPDV
Angioplasty Systems, Atherotomy	17560	ANGPDV
Angioplasty Systems, Laser	17561	ANGPDV
Angioplasty Systems, Radio-Frequency Thermal	17562	ANGPDV
Angioscopes	16996	ENDOSC
Angiotomes	10152	SGINAC
Annuloplasty Rings	16039	GRFIMP
Anomaloscopes	16342	ENDOSC
Anoscopes	10156	ENDOSC
Anticoagulants	*****	BLOCOL
Antrotomes	17426	SGINAC
Apex Cardiographs	10164	CARDIO
Apex Locators, Endodontic	16355	DENTEQ
Apheresis Units	16405	BLOCOL
Apnea Monitors	12575	RESPMA
Apnea Monitors, Neonatal	17172	FENEMA
Apnea Monitors, Recording	17885	RESPMA
Applanometers	10168	OPHDEQ
Applicators, Antiseptic	10173	APPMED
Applicators, Antiseptic, Towelette	17988	SWABCL
Applicators, Antiseptic, Swab	13913	SWABCL
Applicators, Clip	10174	SGINAC
Applicators, Eye	10175	APPMED
Applicators, Tracheal	10177	APPMED
Applicators, Vaginal	16397	APPMED
Applicators, Other	15689	APPMED
Appliers, Aneurysm Clip	16443	SGINAC
Appliers, Cerclage, Orthopedic	16444	ORTHEQ
Appliers, Hemostatic Clip, Disposable	17684	SGINAC
Appliers, Hemostatic Clip, Reusable	16445	SGINAC
Aprons, Surgical	16398	APPRO

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Aqueous/Vitreous Humor Replacement Medium	16844	IOFLUD
Aqueous/Vitreous Humor Replacement Medium Kits	17642	IOFLUD
Argon Beam Coagulation Systems	17739	COAGTS
Arterioscopes	10197	ENDOSC
Arthrogram Kits	15316	RADACC
Arthroscopes	10198	ENDOSC
Arthroscopic Shaver System	17918	ESUNAC
Aspirators	10208	ASDVAC
Aspirators, Dental	10212	DENTEQ
Aspirators, Emergency	15016	ASDVAC
Aspirators, Endocervical	16227	GYNOBE
Aspirators, Infant	10214	ASDVAC
Aspirators, Liposuction	17103	ASDVAC
Aspirators, Low Volume	10215	ASDVAC
Aspirators, Nasal	10216	ASDVAC
Aspirators, Surgical	10217	ASDVAC
Aspirators, Thoracic	10218	ASDVAC
Aspirators, Tracheal	10219	ASDVAC
Aspirators, Ultrasonic	15756	ASDVAC
Aspirators, Uterine	10222	GYNOBE
Aspirators, Wound	10223	ASDVAC
Atmocauteries	10225	COAGTS
Atomizers	10226	ATINNB
Audiometers	10228	AUDIEQ
Audiometric Booths	10229	AUDIEQ
Audiometric Evoked-Potential Units	16303	EEGEQA
Auditory Function Screening Devices	17601	AUDIEQ
Autotransfusion Units, Blood Processing	17537	BLOCOL
Autotransfusion Units, Whole Blood Recovery	17538	BLOCOL
Awls	15275	SGINAC

**B**

Ballistocardiographs	10266	CARDIO
Balloons, Gastric, Appetite Suppressing	17202	GRFIMP
Balloons, Nasal	12699	CATVAR
Bandages/Dressings, Adhesive	10276	BDDRCL
		BDDRST
Bandages, Elastic	10279	BDDRCL
Bandages, Gauze	10281	BDDRST
Bandages, Pressure	10284	BDDRCL
Bandages, Scultetus	10286	BDDRCL
Bandages, Self-Adherent	10288	BDDRCL
Bandages, Traction	10289	BDDRCL
Bandages, Tubular	10291	BDDRCL

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Bandages, Other	15203	BDDRCL
Bands, Spinal	*****	GRFIMP
Barrier Protectants	*****	EMOLPR
Baths, Paraffin, Physical Therapy	12956	BATHTH
Baths, Water, Shaker	16772	BATHTH
Baths, Whirlpool	14450	BATHTH
Baths, Dialysate Warming	15825	HDHFEQ
Beds, Air-Fluidized	16889	PREREA
Beds, Birthing	15732	FURNIP
Beds, Circle Electric	10345	FURNIP
Beds, Electric	10347	FURNIP
Beds, Electric, Obese	15760	FURNIP
Beds, Flotation Therapy	10348	PREREA
Beds, Flotation Therapy, Neonatal	16168	PREREA
Beds, Fluoroscopic	10351	FURNIP
Beds, Hydraulic	10353	FURNIP
Beds, Labor	10356	FURNIP
Beds, Orthopedic	10360	FURNIP
Beds, Rocking	10363	FURNIP
Behavioral Therapy Units	10369	BIOSYS
Belts, Electrode	16226	ELECTR
Belts, Rib	13397	SGINAC
Bile Collection Bags	10386	DRKTAC
Biofeedback Systems	10396	BIOSYS
Biopsy Kits	10401	CELTCE
Biopsy Kits, Bone Marrow	16833	CELTCE
Biopsy Kits, Renal	16834	CELTCE
Biopsy Kits, Soft Tissue	16835	CELTCE
Biopsy Kits, Ultrasonic Aspiration	10403	CELTCE
Bistouries, Tracheal	16409	TTCEQU
Bladder Support, Incontinence	*****	INCOEQ
Blades	15558	SGDRIL
Blades, Bone Cutting	10456	SGINAC
Blades, Cast Cutter	15774	SGINAC
Blades, Knife	12234	SGINAC
Blades, Saw	13447	SGINAC
Blankets	10415	BLANMD
Blood Administration Sets	10421	IVASAC
Blood Cholesterol Monitors	*****	DIAHMU
Blood Collection Tube Adaptors	17814	FITADP
Blood Donor Sets	10426	BLOCOL
Blood-Flow Detectors, Ultrasonic	10429	FLOWDV
Blood Gas Monitors, Carbon Dioxide	12588	BLGSMA
Blood Gas Monitors, Oxygen	12590	BLGSMA
Blood Gas/pH Analyzers	15709	BLGSMA

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Blood Gas/pH Monitors, Extracorporeal	17680	BLGSMA
Blood Cell Processors	16857	BLOCOL
Blood Cell, Washer Bags	16752	BLOCOL
Blood Glucose Monitors, Portable	16488	DIAHMU
Blood Lactate Monitor	*****	DIAHMU
Body Exhaust Systems	16176	PHYSEV
Body Fat/Lean Ratio Analyzers	17417	PHYSEV
Bolts, Bone	16077	GRFIMP
Bone Matrix Implants	17600	MATERL
Bone Matrix Implants, Artificial	17751	MATERL
Bone Matrix Implants, Biological	17756	AUGPSM
Bone Mill	*****	SGINAC
Bone Wax	10459	GRFIMP
Bottles, Drainage	11303	DRKTAC
Bottles, Evacuated	15809	DRKTAC
Bougies	15560	AWPSDA
		SGINAC
Bougies, Ear	11369	SGINAC
Bougies, Esophageal	11597	SGINAC
Bougies, Urethral	14286	SGINAC
Brachytherapy Applicators Afterloading	17732	RADTUN
Brachytherapy Sources	17518	RADTUN
Brachytherapy Systems, Remote Afterloading	17517	RADTUN
Breathing Bags	10486	VBDVAC
Breathing Circuit Adapters	10123	VBDVAC
Breathing Circuits, Anesthesia	10139	ANEQAC
Breathing Circuits, Oxygen Administration	16987	MGASEQ
Breathing Circuits, Ventilator	15003	VBDVAC
Breathing Devices, Self-Contained	17089	VBDVAC
Bronchoscopes, Flexible	15073	ENDOSC
Bronchoscopes, Flexible, Video	17662	ENDOSC
Bronchoscopes, Rigid	15074	ENDOSC
Brushes, Cytology	15018	CELTCE
Brushes, Eye	10505	SGINAC
Brushes, Intramedullary	15867	SGINAC
Brushes, Surgical Scrub	10507	SWABCL
Bubble Detectors	10510	CPBYEQ
Burn Kits	10516	BDDRST
Burs, Corneal	16413	SGDRIL
Burs, Cranial	10520	SGDRIL
Burs, Dental	10521	DENTEQ
Burs, Microsurgical	15883	SGINAC
Burs, Oral Surgery	11341	SGDRIL
Burs, Other	15204	SGDRIL
Burs, Podiatric	10522	SGDRIL

**ECRI UMDNS DEVICE TERM****ECRI IMDC**    **ADG CODE****C**

Cables/Leads	16312	PRMSEQ
Cables/Leads, Catheter	*****	CATKEQ
Cables/Leads, Electrosurgical Unit	11496	ESUNAC
Cables/Leads, ECG	15754	ELECTR
Cables/Leads, ECG, Radiolucent	17469	ELECTR
Cables/Leads, EEG	16313	ELECTR
Cables/Leads, EMG	16314	EDGEQA
Cables/Leads, Nuclear	*****	NUCMEQ
Cables/Leads, Other	16316	STIMAC
Cables/Leads, TENS	16315	STIMAC
Calipers, Ophthalmic	16415	SGINAC
Calipers, Orthopedic	16416	SGINAC
Calipers, Other	16417	SGINAC
Calipers, Skinfold, Electronic	16267	PHYSEV
Caloric Irrigation Units	10548	ENTEQU
Cameras	10549	DENTEQ
Cameras, Fundus	10551	OPHDEQ
Cameras, Gamma, Mobile	16891	NUCMEQ
Cameras, Gamma, Stationary	16892	NUCMEQ
Canisters, Aspirator Collection	10211	ASDVAC
Cannulae, Aortic	15768	CANNUL
Cannulae, Arterial	10564	CANNUL
Cannulae, Arthroscopy Drainage	15982	CANNUL
Cannulae, Aspirating	10566	CANNUL
Cannulae, Brain	10567	CANNUL
Cannulae, Bronchial	16106	CANNUL
Cannulae, Coronary Artery	10569	CANNUL
Cannulae, Coronary Sinus	17644	CANNUL
Cannulae, Ear	10571	CANNUL
Cannulae, Epidural	10572	CANNUL
Cannulae, Eye	10573	CANNUL
Cannulae, Eye, Cyclodialysis	10574	CANNUL
Cannulae, Eye, Irrigating/Aspirating	17899	CANNUL
Cannulae, Eye, Lacrimal	10575	CANNUL
Cannulae, Femoral	10576	CANNUL
Cannulae, Hemodialysis	15021	CANNUL
Cannulae, Liposuction	17640	CANNUL
Cannulae, Nasal	16422	CANNUL
Cannulae, Nasal Oxygen	12700	CANNUL
Cannulae, Nasal Oxygen, Continuous Positive Airway Pressure	12701	CANNUL
Cannulae, Sinus	10578	CANNUL
Cannulae, Suprapubic	16420	CANNUL



**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Cannulae, Thoracic	17991	CANNUL
Cannulae, Tracheostomy	14089	CANNUL
Cannulae, Uterine	10580	CANNUL
Cannulae, Vena Caval	10581	CANNUL
Cannulae, Venous	10582	CANNUL
Cannulae, Ventricular	16783	CANNUL
Cannulae, Other	15206	CANNUL
Caps, Bone	16082	GRFIMP
Caps, Cervical	16423	CONBAR
Caps, Hypothermia	16219	THERCO
Caps, Nerve	16083	GRFIMP
Caps, Ophthalmic, Quarter Globe	16064	GRFIMP
Carbon Dioxide Analyzers	10588	BLGSMA
Carbon Dioxide Monitors, Exhaled Gas	16938	BLGSMA
Carbon Dioxide Monitors, Transcutaneous	15970	BLGSMA
Carbon Monoxide Analyzers	15093	BLGSMA
Cardiac Output Units, Direct Fick	10614	CAOUEQ
Cardiac Output Units, Dye Dilution	16177	CAOUEQ
Cardiac Output Units, Impedance	17496	CAOUEQ
Cardiac Output Units, Radioisotope	15795	CAOUEQ
Cardiac Output Units, Thermal Dilution	10615	CAOUEQ
Cardiac Output Units, Ultrasonic	17190	CAOUEQ
Cardiac Output Units, Other	15207	CAOUEQ
Cardiodynameters	10617	CARDMA
Cardiointegrators	12799	CARDIO
Cardiokymographs	16522	CARDIO
Cardioplegia Solution Administration Kits	16163	CPBYEQ
Cardioverters Automatic, Implantable	*****	DEFIMP
Cartridges, Freon, Cryosurgical	16636	MGASEQ
Cartridges, Freon, Pneumatic Splint	16637	MGASEQ
Cartridges, Freon, Pneumatic Tourniquet	16638	MGASEQ
Carts, Resuscitation	10647	RESEMA
Cast Cutter Vacuums	10656	SGINAC
Cast Cutters, Electric	16340	SGINAC
Cast Cutters, Manual	16341	SGINAC
Cast Material	10665	MTC SPL
Cataract Extraction Units	10673	OPHINS
Cataract Extraction Units, Phacoemulsification	17596	OPHINS
Catheter Care Kits	10676	CATKEQ
Catheter Flush Valves, Continuous	15320	CATKEQ
Catheter Injection Ports	16858	IJIFPT
Catheter Introducers	10678	CATKEQ
Catheter Introducers-Hemostasis Valve	17578	CATKEQ
Catheter Plugs	13066	CATKEQ
Catheterization Kits	15564	CATKEQ

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Catheterization Kits, Cardiac	10598	CATKEQ
Catheterization Kits, Central Venous	16615	CATKEQ
Catheterization Kits, Intermittent	16321	CATKEQ
Catheterization Kits, Intravenous	12161	CATKEQ
Catheterization Kits, Peritoneal Dialysis	16992	CATKEQ
Catheters, Anesthetic Conduction	16430	CATVAR
Catheters, Angiography	10688	CATVAS
Catheters, Angioplasty	17183	CATVAS
Catheters, Angioplasty, Atherectomy	17519	CATVAS
Catheters, Angioplasty, Atherotomy	17565	CATVAS
Catheters, Angioplasty, Balloon Dilatation	17184	CATVAS
Catheters, Angioplasty, Balloon Dilatation, Coronary Perfusion	17521	CATVAS
Catheters, Angioplasty, Laser	17185	CATVAS
Catheters, Angioplasty, Radio-Frequency Thermal	17566	CATVAS
Catheters, Arterial	10689	CATVAS
Catheters, Arterial Micro-Flow	10691	CATVAS
Catheters, Arthroscopy	15983	CATVAR
Catheters, Barium Enema	10694	ENMAKT
Catheters, Bartholin's Gland	10695	CATVAR
Catheters, Biliary	10696	CATVAR
Catheters, Biliary Manometric	17799	CATVAR
Catheters, Bronchial Balloon	17827	CATVAR
Catheters, Cardiac	10698	CATVAS
Catheters, Cardiac, Ventricular	17613	CATVAS
Catheters, Cardiac Ablation	16655	CATVAS
Catheters, Cardiac Acoustic	10699	CATVAS
Catheters, Cardiac Balloon, Oxygenation Monitoring	16656	CATVAS
Catheters, Cardiac Balloon, Pacing Electrode	16654	CATVAS
Catheters, Cardiac Balloon, Pulmonary Artery	16777	CATVAS
Catheters, Cardiac Balloon, Pulmonary Artery, Oximetric	17902	CATVAS
Catheters, Cardiac Balloon, Thermal Dilution	10754	CATVAS
Catheters, Cardiac Balloon, Other	16778	CATVAS
Catheters, Cardiac Bypass	10701	CATVAS
Catheters, Cardiac Percutaneous	10702	CATVAS
Catheters, Cardiac Transseptal	17478	CATVAS
Catheters, Central Venous, General Purpose	16746	CATVAS
Catheters, Central Venous, Long Term	16586	CATVAS
Catheters, Central Venous, Peripherally Inserted	18017	CATVAS
Catheters, Cerebrospinal	10704	CATVAR
Catheters, Cholangiography	16429	CATVAR
Catheters, Coronary Artery Infusion	17871	CATVAS
Catheters, Coronary Perfusion	10706	CATVAS
Catheters, DePezzer	10709	CATURO
Catheters, Dye Dilution	15885	CATVAS
Catheters, Embolectomy	10714	CATVAS

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Catheters, Epidural	10717	CATVAR
Catheters, Esophageal Balloon	10716	CATVAR
Catheters, Eustachian	10718	CATVAR
Catheters, Femoral	15023	CATVAS
Catheters, Hemodialysis	15022	CATVAR
Catheters, Hemostatic	10723	CATVAR
Catheters, Hydrocephalic	16133	CATVAR
Catheters, Intracranial, Ventricular	10769	CATVAR
Catheters, Intra-Aortic Balloon	10725	CATVAS
Catheters, Intrauterine	16431	CATVAR
Catheters, Intravascular, Guiding	17846	CATVAS
Catheters, Intravenous, Peripheral	10727	CATVAS
Catheters, Irrigation	10730	CATVAR
Catheters, Jejunostomy	10731	CATVAR
Catheters, Mushroom	10732	CATURO
Catheters, Nasal Oxygen	12702	CATVAR
Catheters, Nasopharyngeal	16432	CATVAR
Catheters, Nelaton	10734	CATURO
Catheters, Nephrostomy	10735	CATURO
Catheters, Occlusion	10736	CATVAS
Catheters, Olive Tip	10737	CATURO
Catheters, Perfusion	10739	CATVAR
Catheters, Pericardium Drainage	10741	CATVAR
Catheters, Peritoneal	16433	CATVAR
Catheters, Peritoneal Dialysis	10742	CATVAR
Catheters, Peritoneal Dialysis, Long Term, Indwelling	16109	CATVAR
Catheters, Rectal	10746	CATVAR
Catheters, Septostomy	10747	CATVAR
Catheters, Single-Needle Hemodialysis	10748	CATVAR
Catheters, Snare	17927	CATVAR
Catheters, Suction/Insufflation	17795	CATVAR
Catheters, Suprapubic	10751	CATURO
Catheters, Thrombectomy	10756	CATVAS
Catheters, Torque	*****	CATVAS
Catheters, Tracheal Oxygen	17450	CATVAR
Catheters, Transtracheal Oxygen	17940	CATVAR
Catheters, Umbilical	10759	CATVAS
Catheters, Umbrella	10760	CATVAS
Catheters, Ureteral	10761	CATURO
Catheters, Urethral	10762	CATURO
Catheters, Urethral, Diagnostic	10763	CATURO
Catheters, Urinary, Condom	10765	CATURO
Catheters, Urinary, Foley	10720	CATURO
Catheters, Urinary, Intermittent	17683	CATURO
Catheters, Urinary, Irrigation	10767	CATURO

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Catheters, Valvuloplasty, Balloon Dilatation	17453	CATVAS
Catheters, Vessel Dilatation Balloon	15307	CATVAS
Catheters, Other	15209	CATVAR
Cavity Liners	16182	MTDTRS
Cavity Varnish	16183	MTDTRS
Cell Separation System	*****	EXCOTH
Cement Dispensers, Orthopedic	16472	ORTHEQ
Cement, Dental, Ethoxy Benzoic Acid	16703	MTDTRS
Cement, Dental, Glass Ionomer	16704	MTDTRS
Cement, Dental, Polycarboxylate	16705	MTDTRS
Cement, Dental, Reinforced Zinc-Oxide Eugenol	16706	MTDTRS
Cement, Dental, Resin Composite	16707	MTDTRS
Cement, Dental, Silicosphosphate	16708	MTDTRS
Cement, Dental, Zinc Phosphate	16710	MTDTRS
Cement, Dental, Zinc-Oxide Eugenol	16709	MTDTRS
Cement, Dental, Other	16711	MTDTRS
Cement, Orthodontic	17943	MTDTNR
Cement, Orthopedic	12830	MATERL
Cement, Ostomy Appliance	16434	OSTSUB
Cerclage, Orthopedic Fixation	16435	ORTHEQ
Cerebral Blood-Flow Measurement Units, Xenon Clearance	16237	FLOWDV
Cerebral Function Monitors	12597	NEURDE
Chairs, Birthing	15024	FURNIP
Chairs, Cardiac	10790	FURNIP
Chairs, Dental	10792	DENTEQ
Chairs, Examination, Treatment	16437	FURNIP
Chairs, Examination/Treatment, ENT General	10794	FURNIP
Chairs, Flotation Therapy	10795	PREREA
Chairs, Rotating Vestibular Function	17217	NEURDE
Chambers, Hyperbaric	12061	CHAMPR
Chambers, Hypobaric	12076	CHAMPR
Chisels, Bone	10825	SGINAC
Chisels, Mastoid	10827	SGINAC
Chisels, Middle Ear	16400	SGINAC
Chisels, Nasal	10828	SGINAC
Chisels, Orthopedic	10829	SGINAC
Chisels, Other	15211	SGINAC
Choledochoscopes	10831	ENDOSC
Chondrotomes	16805	SGINAC
Chronaximeters	10833	NEURDE
Circulatory Assist Units, Cardiac	10840	CRASUN
Circulatory Assist Units, Intra-Aortic Balloon	10846	CRASUN
Circulatory Assist Units, Ventricular	10847	CRASUN
Circulatory Assist Units, Venous Return	10849	CRASUN
Circumcision Trays, Disposable	10859	SGINAC

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Clamps, Aorta	10864	SGINAC
Clamps, Artery	10865	SGINAC
Clamps, Bone	10866	SGINAC
Clamps, Bronchus	10867	SGINAC
Clamps, Bulldog	10868	SGINAC
Clamps, Cannula	16449	SGINAC
Clamps, Carotid	15878	SGINAC
Clamps, Circumcision	10869	SGINAC
Clamps, Coarctation	15879	SGINAC
Clamps, Hemorrhoidal	10870	SGINAC
Clamps, Intestinal	10871	SGINAC
Clamps, Ophthalmic	16450	SGINAC
Clamps, Otologic	16451	SGINAC
Clamps, Patent Ductus	10872	SGINAC
Clamps, Pediatric	15881	SGINAC
Clamps, Penile	16452	SGINAC
Clamps, Peripheral Vascular	10873	SGINAC
Clamps, Rectal	15671	SGINAC
Clamps, Skull	15025	SGINAC
Clamps, Tubing	10875	SGINAC
Clamps, Umbilical	10876	SGINAC
Clamps, Uterine	16453	SGINAC
Clamps, Vascular	15882	SGINAC
Clamps, Other	15212	SGINAC
Cleaners, Contact Lens	10995	COLECP
Cleaners, Room Air	16629	FILTER
Cleaning Solutions, Contact Lens	17430	COLECP
Clips	10898	IMPABS
Clip Appliers	10894	SGINAC
Clip Removers	10895	SGINAC
Clips, Aneurysm	10905	SGINAC
		GRFIMP
Clips, Drape	10902	SGINAC
Clips, Fallopian Tube	14174	CONTDV
Clips, Hemostatic	10903	SGINAC
Clips, Iris	10906	SGINAC
Clips, Lens Implantation	16457	SGINAC
Clips, Ligating	16052	GRFIMP
Clips, Penis	10908	SGINAC
Clips, Scalp	10909	SGINAC
Clips, Skin	10910	SGINAC
Clips, Suture	10911	SGINAC
Clips, Tube	15736	SGINAC
Clips, Vascular	16028	SGINAC
		GRFIMP

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Clips, Vena Caval	10915	SGINAC GRFIMP
Clips, Wound	10916	SGINAC
Clips, Other	15213	SGINAC
Coagulation Analyzers, Whole Blood	16749	DIAHMU
Coils/Pouches, Blood Warming	16845	BLOODW
Collectors Sputum	10943	ASDVAC
Collagen, Injectable	17853	AUGPSM
Collimators, Radiographic	16389	RADGEQ
Colonoscopes	10950	ENDOSC
Coloscopes	10951	ENDOSC
Colposcopes	10960	ENDOSC
Composite Restorative Material Kits, Dental, Light-Cured	16736	MTDTRS
Composite Restorative Material Kits, Dental, Other	16740	MTDTRS
Composite Restorative Materials, Dental, Self-Cured	16725	MTDTRS
Composite Restorative Materials, Dental, Other	16731	MTDTRS
Compresses, Gauze	10966	BDDRCL BDDRST
Compression Devices, Radiographic	16559	RADGEQ
Compression Units, Intermittent	10969	AEMBDV
Compression Units, Intermittent, Cryotherapy	17226	CRYOEQ
Compression Units, Sequential	16837	AEMBDV
Compressors, Air	10972	ATINNB
Compressors, Orbital	10975	SGINAC
Computer Information Systems, Bedside	17223	FURNIP
Computers, Cardiac Catheterization Laboratory	10980	CATKEQ
Computers, ECG Interpretation	10981	CARDIO
Computers, Nuclear Medicine	15678	NUCMEQ
Computers, Pulmonary Function Laboratory	10984	RESPMA
Computers, Radiotherapy Planning System	13281	RADTUN
Computers, Radiotherapy Record and Verify	17532	RADTUN
Computers, Stress Exercise	15027	PHYSEV
Condoms	10985	CONDOM
Condoms, Female	18079	CONBAR
Conductivity Meters	16280	ACPTEQ
Conformers, Ophthalmic	16065	OPHINS
Contact Lens Care Kits	*****	COLECP
Contact Lenses, Therapeutic	17524	CONTLN
Contact Lenses, Vision Corrective	17525	CONTLN
Continuous Positive Airway Pressure Units	11001	AWPSDA
Contraceptive, Barrier	30094	CONBAR
Contraceptive Drug Delivery Units, Implantable	16964	CONTDV
Controllers	11009	CPBYEQ
Controllers, Temperature, Humidifier	15664	THERCO
Controllers, Temperature, Programmable	16539	THERCO

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Controls, Liquid Level, Cryogenic Converter	16542 *****	AGUNAC SGINAC
Corneal Storage Media	17606	CELTCE
Corneas, Eye	17108	CIPHAO
Cornish Hemostatic Wool	11017	HEMOMT
Correlators, Retinal	11023	OPHDEQ
Cotton Balls	11028	BDDRCL
Cotton Rolls	13414	BDDRCL
Counters, Gamma	16381	NUCMEQ
Counters, Intravenous Drop	12168	IVASAC
Counters, Scintillation	15184	NUCMEQ
Cover, Probe	*****	ULTRDV
Covers	15571	APPRO
Covers, Bur Hole	15749	GRFIMP
Covers, Shoe, Conductive	13575	APPRO
Covers, Thermometer, Electronic Probe	13116	BODYTP
Cricothyrotomes	11052	TTCEQU
Cricothyrotomy Kits	15028	TTCEQU
Crimpers	16463	SGINAC
Crushers, Fascia	11060	SGINAC
Cryometers	11066	CRYOEQ
Cryosurgical Units	18051	CRYOEQ
Cryosurgical Units, General-Purpose	11067	CRYOEQ
Cryosurgical Units, Ophthalmic	11068	CRYOEQ
Cuff Inflators	11071	IVASAC
Cuffs, Blood Pressure	11073	BLPRDV
Cuffs, Tourniquet, Disposable	16907	TORQAC
Cuffs, Tourniquet, Reusable	17230	TORQAC
Cuffs, Tracheal Tube	14082	AWPSDA
Cuffs, Tracheostomy Tube	14094	TTCEQU
Cuffs, Vascular Catheter	17470	CATVAS
Culdocentesis Kits	15313	SGINAC
Culdoscopes	11079	ENDOSC
Curets	11084	SGINAC
Current Limiters, Patient Leads	11089	DEFEXT
Cutters	15573	SGINAC
Cutters, Bone	10455	SGINAC
Cutters, Meniscus	17117	SGINAC
Cutters, Pin and Wire	13041	STAPAC
Cutters, Suture	16224	STAPAC
Cuvette	17749	BLGSMA
Cystitomes	11113	SGINAC
Cystotomes	16002	SGINAC
Cystometers	11111	URODEQ
Cystoscopes, Flexible	17144	ENDOSC

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Cystoscopes, Rigid	17145	ENDOSC
Cystourethroscopes	11114	ENDOSC

**D**

Data Analysis Systems, Ultrasound	17203	ULTRDV
Data Analysis Systems, Ultrasound, Cardiac	17204	ULTRDV
Data Management Systems	17762	PDMSYS
Data Management Systems, Anesthesia	17504	PDMSYS
Data Management Systems, EGC	17205	PDMSYS
Data Management Systems, Hemodynamic	17766	PDMSYS
Data Management Systems, Obstetrical	17939	PDMSYS
Data Management Systems, Patient	16361	PDMSYS
Decapitation Instruments	11122	SGINAC
Decompression Units, Abdominal	11126	DECPEQ
Defibrillator Accessories	*****	DEFACC
Defibrillator Accessories, Implantable	*****	DEFIMP
Defibrillator Analyzers	11127	DEFEXT
Defibrillator Paddle Pads	11130	DEFEXT
Defibrillator/Monitor/Pacemakers	17882	DEFEXT
Defibrillator/Monitors	11129	DEFEXT
Defibrillator/Monitors, Line-Powered	15029	DEFEXT
Defibrillator/Monitors, Transtelephonic	17579	DEFEXT
Defibrillators Automatic, Implantable	16652	DEFIMP
Defibrillators, Automatic, External	17116	DEFEXT
Defibrillators, Battery-Powered	11134	DEFEXT
Defibrillators, Internal	11133	DEFEXT
Defibrillators, Line-Powered	11137	DEFEXT
Delivery Kits	11141	GYNOBE
Delivery Kits, for Drugs	*****	KTPWOD
Densitometers, Radiographic	16548	RADGEQ
Dental Amalgam	16367	MTDTRS
Dental Amalgam, Capsule	16712	MTDTRS
Dental Amalgam, Powder	16714	MTDTRS
Dental Amalgam, Tablet	16713	MTDTRS
Dental Delivery Units	11165	DENTEQ
Dental Engines	11156	DENTEQ
Dental Etching Liquids	17737	MTDTRS
Dental Foil, Gold	11159	MTDTRS
Dental Hand Instruments, Endodontic	16662	DENTEQ
Dental Hand Instruments, Operative	16663	DENTEQ
Dental Hand Instruments, Orthodontic	16664	DENTEQ
Dental Hand Instruments, Periodontal	16665	DENTEQ
Dental Hand Instruments, Surgical	16666	DENTEQ
Dental Hand Instruments, Other	16667	DENTEQ



**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Dental Protectants	*****	MTDTNR
Dental Syringes, Anesthetic	17111	DENTEQ
Dental Syringes, Irrigating	17112	DENTEQ
Depressors, Orbital	16465	SGINAC
Depressors, Vessel	14367	SGINAC
Dermabrasion Brushes	16344	SGINAC
Dermabrasion Units	11177	SGINAC
Dermatoscopes	18021	ENDOSC
Dermoscopes	18193	ENDOSC
Dermoscopes, Facial	18194	ENDOSC
Dermatomes	11179	SGINAC
Desiccators	15130	ESUNAC
Device Kits with Drugs	*****	KTPCOD
Devices, Albumin Coated	*****	MATHAO
Dialysate Conductivity Meters	11206	HDHFEQ
Dialysate Delivery Systems, Multipatient	11211	HDHFEQ
Dialysate Delivery Systems, Single-Patient	11213	HDHFEQ
Dialyzer Reprocessing Units	16018	HDHFEQ
Dialyzers, Hemodialysis, Coil	11233	HDHFEQ
Dialyzers, Hemodialysis, Hollow-Fiber	11234	HDHFEQ
Dialyzers, Hemodialysis, Parallel-Plate	11235	HDHFEQ
Dialyzers, High-Permeability	16485	HDHFEQ
Dialyzers, Serum/Urine	15131	HDHFEQ
Diathermy Unit Analyzers, Shortwave	16362	DIATEQ
Diathermy Units	11244	OPHINS
Diathermy Units, Microwave	11245	DIATEQ
Diathermy Units, Ophthalmic	15874	DIATEQ
Diathermy Units, Shortwave	11246	DIATEQ
Digital Angiography Systems	16560	RADGEQ
Dilatometers	11253	URODEQ
Dilators	11254	GYNOBE
Dilators, Cardiac	11255	SGINAC
Dilators, Cervical	16469	SGINAC
Dilators, Common Duct	11256	SGINAC
Dilators, Esophageal	11257	SGINAC
Dilators, Lacrimal	11259	SGINAC
Dilators, Nasal	11260	SGINAC
Dilators, Pulmonary	11261	SGINAC
Dilators, Rectal	11262	SGINAC
Dilators, Salivary Duct	16550	SGINAC
Dilators, Tracheal	11263	TTCEQU
Dilators, Tubal	11264	SGINAC
Dilators, Ureteral	16468	SGINAC
Dilators, Urethral	11265	SGINAC
Dilators, Uterine	11266	SGINAC

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Dilators, Vaginal	11267	SGINAC
Dilators, Vascular	11268	SGINAC
Dilators, Other	15215	SGINAC
Directoscopes	11274	ENTEQU
Disinfecting units, Liquid, Flexible Endoscope	11279	STERIL
Dissectors	11290	SGINAC
Dissector, Electrosurgical	*****	ESUNAC
Dolorimeters	11293	NEURDE
Domes, Pressure Transducer	15293	TRDSNS
Dosimeters, Radiation	11295	RADACC
Douche Kits	11297	IRILAV
Douche Nozzles	11298	IRILAV
Drainage Bags	11301	DRKTAC
Drainage Systems, Pleural	10817	DRKTAC
Drains, Sump	15270	DRKTAC
Drains, T	11307	DRKTAC
Drains, Thoracic	11308	DRKTAC
Drains, Vent	11310	DRKTAC
Drains, Wound	11305	DRKTAC
Drapes, Angiographic	16006	APPPRO
Drapes, Instrument	15643	APPPRO
Drapes, Laser-Resistant	17244	SGINAC
Drapes, Laser Unit	17243	APPPRO
Drapes, Microscopes	16474	APPPRO
Drapes, Surgical, Disposable	15646	APPPRO
Dressing Kits	11314	BDDRST
Dressing Pads	11326	BDDRST
Dressing, Fluff	11751	BDDRST
Dressings, Aerosol	11321	BDDRCL
Dressings, Burn	11322	BDDRST
Dressings, Cavity Wound	17926	BDDRST
Dressings, Foam	11323	BDDRST
Dressings, Germicidal	11324	BDDRST
Dressings, Nonadherent	11325	BDDRST
Dressings, Occlusive	15704	BDDRST
Dressings, Other	15216	BDDRCL
		BDDRST
Dressings, Periodontal	15806	MTDTNR
Dressings, Resorbable Antiadhesion	17670	BDDRST
Dressings, Roller Gauze	11327	BDDRST
Dressings, Skin Graft, Donor Site	15319	BDDRST
Dressings, Tracheostomy Tube	15624	BDDRST
Dressings, Transparent Adhesive	17428	BDDRST
Dressings, Universal	11328	BDDRST
Drills, Biopsy	11330	SGDRIL

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Drills, Bone	11331	SGDRIL
Drills, Cannulated	11332	SGDRIL
Drills, Cervical	11333	SGDRIL
Drills, Fingernail	11337	SGDRIL
Drills, Intramedullary	11338	SGDRIL
Drills, Middle Ear Surgery	11339	SGDRIL
Drills, Perforator	11342	SGDRIL
Driver/Extractors	11345	ORTHEQ
Driver/Extractors, Bone Nail-Pin	11346	SGINAC
Driver/Extractors, Bone Plate	11347	SGINAC
Drivers, Bone Staple	10458	SGINAC
Drivers, Nail	12398	SGINAC
Duodenoscopes	11359	ENDOSC
Dynamometers, Squeeze	11367	PHYSEV

**E**

Echoencephalographs	11387	ULTRDV
ECG, Software	*****	CARDIO
EEG, Software	*****	EEGEQA
Eikonometers	11391	OPHDEQ
Ejectors, Saliva	11393	DENTEQ
Elastomer Driven Drug Infusor	*****	INFDEV
Elastomers, Silicone Rubber	11394	MATERL
Electrical Conversion Kits, Mechanical Bed	16917	FURNIP
Electrooculographs	11480	EOPHRC
Electrocardiocontourographs	11402	CARDIO
Electrocardiographs, Interpretive	16231	CARDIO
Electrocardiographs, Multichannel	11411	CARDIO
Electrocardiographs, Signal-Averaging	17687	CARDIO
Electrocardiographs, Single-Channel	11413	CARDIO
Electrocautery Units, Battery-Powered	11421	COAGTS
Electrocautery Units, Line-Powered	11422	COAGTS
Electrocochleographs	17495	EDGEQA
Electroconvulsive Therapy Units	11484	ECTHUN
Electrode Caps, Electroencephalographic	17554	ELECTR
Electrode Holders, Electrosurgical	11497	ESUNAC
Electrode Monitoring Systems, Electrosurgical, Return	16010	ESUNAC
Electrodes, pH	15137	ELECTR
Electrodes, Blood Gas, pH	17686	ELECTR
Electrodes, Blood Gas, Carbon Dioxide	11432	ELECTR
Electrodes, Blood Gas, Oxygen	11433	ELECTR
Electrodes, Catheter Tip	11434	ELECTR
Electrodes, Cystoscopic	11436	ELECTR
Electrodes, Defibrillation	15033	ELECTR

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Electrodes, Dorsal Column	17214	ELECTR
Electrodes, Electrooculogram	11443	ELECTR
Electrodes, Electroanalgesic TENS	17191	ELECTR
Electrodes, Electrocardiographic	11439	ELECTR
Electrodes, Electrocardiographic, Neonatal	17460	ELECTR
Electrodes, Electroencephalographic	11440	ELECTR
Electrodes, Electromyographic	11441	ELECTR
Electrodes, Electronystagmographic	11442	ELECTR
Electrodes, Electroschock	11444	ELECTR
Electrodes, Electrosurgical	15579	ESUNAC
Electrodes, Electrosurgical, Active	16860	ESUNAC
Electrodes, Electrosurgical, Active, Foot-Controlled	16206	ESUNAC
Electrodes, Electrosurgical, Active, Hand-Controlled	11499	ESUNAC
Electrodes, Electrosurgical, Return	11500	ESUNAC
Electrodes, Esophageal	11446	ELECTR
Electrodes, Fetal Scalp	11447	ELECTR
Electrodes, Impedance Analyzers	16999	ELECTR
Electrodes, Intracardiac	16995	ELECTR
Electrodes, Intracardiac, Conduction Study	17612	ELECTR
Electrodes, Nasopharyngeal	11452	ELECTR
Electrodes, Neurological	11453	ELECTR
Electrodes, Neuromuscular Stimulator	11454	ELECTR
Electrodes, Pacemaker, Noninvasive Transcutaneous	17154	ELECTR
Electrodes, Pacemaker, Transthoracic	16127	ELECTR
Electrodes, Electroretinographic	11460	ELECTR
Electrodes, Sweat Test	11461	ELECTR
Electrodes, Other	15217	ELECTR
Electroencephalographs	11467	EDGEQA
Electroglottographs	16034	EDGEQA
Electromanometers	11473	PRMSEQ
Electromyographs	11474	EDGEQA
Electronystagmographs	11479	EOPHRC
Electroretinographs	11482	EOPHRC
Electrorheographs	11483	EDGEQA
Electrosleep Units	11485	STIMAC
Electrosurgical Unit Adapters, Cable	11493	ESUNAC
Electrosurgical Unit Adapters, Hand Control	11494	ESUNAC
Electrosurgical Unit Analyzers	11489	ESUNAC
Electrosurgical Units, General-Purpose	16137	ESUNAC
Electrosurgical Units, Specialty	16138	ESUNAC
Electrotherapeutic Units	11503	ELTHEQ
Elevators, Corneal	16479	SGINAC
Elevators, ENT	11506	SGINAC
Elevators, Neurosurgical	11507	SGINAC
Elevators, Orthopedic	11508	SGINAC

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Elevators, Uterine	15677	SGINAC
Elevators, Other	15218	SGINAC
Embolization Prostheses, Intravascular	15034	GRFIMP
Emergency Kits, Anaphylactic	16160	FAEMKT
Emergency Kits, CPR	16161	FAEMKT
Encephaloscopes	11526	ENDOSC
Endodontic Filling Materials	17611	MTDTRS
Endodontic Medicaments	*****	MTDTNR
Endoradiosondes	11530	TELEME
Endoscopic Kit	*****	KTPWOD
Endoscopic Power Supplies, Battery Operated	11577	ENDOSC
Enema Bags	11581	ENMAKT
Enema Kits	11582	ENMAKT
Enema Kits, Barium	16514	ENMAKT
Enema Tips	11583	ENMAKT
Enteral Feeding Bags	11675	EFEDKT
Enteral Feeding Kits	11677	EFEDKT
Entoptoscopes	15288	OPHDEQ
Enucleators	11587	SGINAC
Epilators	11589	OPHINS
Ergometers, Bicycle	10383	PHYSEV
Ergometers, Other	15219	PHYSEV
Esophageal Coolers	11599	THERCO
Esophageal Motility Analyzers	11600	GIDIEQ
Esophagoscopes	11603	ENDOSC
Esthesiometers	11604	NEURDE
Evoked-Potential Units	11614	EEGEQA
Excavators, Endaural	11617	SGINAC
Exchange Transfusion Kits	11618	BLOCOL
Exercisers, Arm	11625	EXERCS
Exercisers, Back	17201	EXERCS
Exercisers, Continuous Passive Motion, Hand	17137	EXERCS
Exercisers, Continuous Passive Motion, Lower Limb	17138	EXERCS
Exercisers, Continuous Passive Motion, Upper Limb	17139	EXERCS
Exercisers, Finger	11628	EXERCS
Exercisers, Hand	11629	EXERCS
Exercisers, Jaw	17802	EXERCS
Exercisers, Leg and Ankle	11630	EXERCS
Exercisers, Neck	18008	EXERCS
Exercisers, Trapeze	11638	EXERCS
Exercisers, Other	15220	EXERCS
Exophthalmometers	16345	OPHDEQ
Expanders, Skin	16045	SGINAC
Explorers	16482	SGINAC
Export only in vitro diagnostics	*****	DIAEXO

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

External Pacemaker Analysers	17721	PMKRAC
Extractors, Metal, Magnetic	12521	SGINAC
Extractors, Nail	12696	SGINAC
Eye Movement Monitors	12609	OPHDEQ
Eye Pads	11661	BDDRST
ECG Monitor Analyzers	11379	CARDIO
ECG Monitors	12599	CARDIO
ECG Monitors, Arrhythmia	12601	CARDIO
ECG Monitors, Telemetric	13988	TELEME
EEG Monitors	12602	EDGEQA
EEG Monitors, Spectral Display	17412	EDGEQA
EEC Monitors, Telemetric	15719	TELEME
ENT Treatment Units	11585	ENTEQU

**F**

Falloposcope	*****	ENDOSC
Fertility Analysers, Thermometric	17645	BODYTP
Fetal Heart Detectors, Phono	11695	FENEMA
Fetal Heart Detectors, Ultrasonic	11696	FENEMA
Fetal Monitors	12610	FENEMA
Fetal Monitors, Ultrasonic	12613	FENEMA
Fibrillators	11700	DEFEXT
Files, Bone	11702	SGINAC
Filiforms	11704	SGINAC
Filters, Air	11711	FILTER
Filters, Bacteria	11712	FILTER
Filters, Blood	11713	FILTER
Filters, Gas	15649	FILTER
Filters, Heart-Lung Bypass Priming	17580	CPBYEQ
Filters, Hemoconcentration	17598	FILTER
Filters, Intravenous Tubing	11716	FILTER
Filters, Kidney Stone	11717	FILTER
Filters, Pheresis	17850	FILTER
Filters, Shunt, Hydrocephalic	15688	SHUNTS
Filters, Syringe	15283	FILTER
Filters, Vena Cava	11718	FILTER
Filters, Ventilator	14352	FILTER
Filtration Systems	16904	FILTER
First Aid Kits	11723	FAEMKT
First Aid Kits, Dental	*****	FAEMKT
Fittings/Adapters	11726	FITADP
Fittings/Adapters, Luer	11729	FITADP
Fittings/Adapters, Pneumatic, DISS	11727	FITADP
Fittings/Adapters, Pneumatic, Quick-Connect	11731	FITADP

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Fittings/Adapters, Pneumatic, Other	16797	FITADP
Flow Mapping Dopplers	16595	FLOWDV
Flow Timers	11744	FLOWDV
Flow Totalizers	11745	FLOWDV
Flowmeters, Blood, Electromagnetic	10431	FLOWDV
Flowmeters, Blood, Laser	16903	FLOWDV
Flowmeters, Blood, Ultrasonic	10432	FLOWDV
Flowmeters, Gas	11748	MGASEQ
Flowmeters, Intravenous	12171	FLOWDV
Flowmeters, Pulmonary Function	15965	RESPMA
Fluids, Biological	*****	AUGPSM
Fluids, Mechanical, Therapeutic Use Only	*****	GASMEC
Fluoroscopic Units, Hand-Held	17533	RADGEQ
Forceps, Biopsy	11775	CELTCE
Forceps, Bone	15670	SGINAC
Forceps, Dissecting	15800	SGINAC
Forceps, Electrosurgical	11502	ESUNAC
Forceps, Epilation	11780	SGINAC
Forceps, Fixation	11781	SGINAC
Forceps, Gallbladder	11782	SGINAC
Forceps, Hemostatic	11784	SGINAC
Forceps, Intestinal	11785	SGINAC
Forceps, Kidney	16519	SGINAC
Forceps, Lung	11787	SGINAC
Forceps, Obstetrical	11788	SGINAC
Forceps, Ophthalmic	16209	OPHINS
Forceps, Specimen	11790	SGINAC
Forceps, Sponge	11791	SGINAC
Forceps, Stone Manipulation	11793	SGINAC
Forceps, Suction	11794	SGINAC
Forceps, Tissue	11797	SGINAC
Forceps, Tonsil	15672	SGINAC
Forceps, Other	15221	SGINAC
Frames, Turning	14256	SGINAC

**G**

Gait Analyzers	16334	PHYSEV
Gasometers	11850	BLGSMA
Gas Delivery Units, Argon Beam Coagulation	17738	COAGTS
Gas Sampling Units	11848	MGASEQ
Gases, Mechanical, Therapeutic, Use Only	*****	GASMEC
Gastroduodenoscopes	11853	ENDOSC
Gastrographs	11854	GIDIEQ
Gastrointestinal Motility Analyzers	15036	GIDIEQ

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Gastroplasty Bands	17649	GRFIMP
Gastrosopes, Flexible	11857	ENDOSC
Gastrosopes, Flexible, Video	17663	ENDOSC
Gastrosopes, Rigid	11858	ENDOSC
Gauze	11859	BDDRCL BDDRST
Gel, Electrode	11425	GELMED
Gel, Resorbable Antiadhesion	*****	MATERL
Gel, Ultrasonic Coupling	15321	GELMED
Gel, Wound	*****	BDDRST
Gels, Resorbable Antiadhesion	*****	AADIMP
Generators, Ultrasonic	*****	ULTRDV
Gingival Retraction Cords	16672	MTDTNR
Gingival Retraction Kits	16352	MTDTNR
Glaucoma Filtration Implants	17618	STENTS
Gloves, Examination/Treatment	11882	GLOVES
Gloves, Pressure	16500	BDDRCL
Gloves, Surgical	11883	GLOVES
Gloves, Surgical, Hypoallergenic	11884	GLOVES
Goniometers	11893	PHYSEV
Gouges	11895	SGINAC
Gowns, Operating Room, Disposable	11901	APPPRO
Gowns, Other	15223	APPPRO
Gowns, Patient, Disposable	11904	APPPRO
Grafts, Bone	11910	GRFIMP
Grafts, Bone, Synthetic	16966	GRFIMP
Grafts, Dura, Biological	17462	AUGPSM
Grafts, Skin	11912	GRFIMP
Guide Wires	11925	SGINAC
Guides	11927	OPHINS
Guides, Drill	11930	SGINAC
Guides, Gigli Saw	11931	SGINAC
Guides, Suturing	17688	STAPAC
Guides, Other	15224	SGINAC
Guillotines, Adenoid	11934	SGINAC
Guillotines, Rib	15038	SGINAC
Guillotines, Tonsil	11935	SGINAC
Gynaecological Units	*****	GYNOBE

**H**

Halogenated Anesthetics Analyzers	11941	BLGSMA
Handpieces	17953	OPHINS
Handpieces, Dental	11161	DENTEQ
Handpiece, Surgical Drill, Electric/Pneumatic	*****	SGDRIL



**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Handpieces, Surgical	17949	SGDRIL
Head Halters	11958	SGINAC
Hearing Aids	11967	HEARAD
Hearing Aids, Programmable	17666	HEARAD
Heart Rate Monitors, R-Wave	12618	CARDMA
Heart Rate Monitors, Other	15233	CARDMA
Heart-Lung Bypass Unit Level Detectors	15011	CPBYEQ
Heart-Lung Bypass Unit Tubing Sets	15012	CPBYEQ
Heart-Lung Bypass Units	11969	CPBYEQ
Heat Exchangers, Heart-Lung Bypass	11973	CPBYEQ
Heat Exchangers, Regional Perfusion	13318	RPUNAC
Heat Exchanging Jackets, Cardiac	17510	CPBYEQ
Heaters, Insufflators	*****	INSUFF
Heating Blocks	15105	THERCO
Heating Pads, Electric	11989	THERCO
Heating Units	17849	ATINNB
Heating Units, Hot Pack	16509	THERCO
Helium Analyzers	15892	BLGSMA
Helmets, Surgical	11995	SGINAC
Hemodialysis Tubing Sets	11225	HDHFEQ
Hemodialysis Unit Air Bubble-Foam Detectors	11219	HDHFEQ
Hemodialysis Unit Blood Line Clamps	16136	HDHFEQ
Hemodialysis Unit Blood Pumps	11220	HDHFEQ
Hemodialysis Unit Conductivity Monitors	11222	HDHFEQ
Hemodialysis Unit Single-Needle Controllers	11224	HDHFEQ
Hemodialysis Unit Single-Pass Converters	15031	HDHFEQ
Hemodialysis Units	11218	HDHFEQ
Hemofiltration Units	15039	HDHFEQ
Hemoperfusion Systems	12001	HDHFEQ
Hemostatic Materials, Animal Origin	*****	HEMOAO
Hemostatic Media	17944	HEMOMT
Hemostatic Media, Powder	17945	HEMOMT
High Level Disinfectant	*****	DISIWC
His Bundle Detectors	12009	CARDIO
Hoods, Oxygen	12027	MGASEQ
Hooks, Cordotomy	12029	SGINAC
Hooks, Incus	12031	SGINAC
Hooks, IUD Removal	12032	CONTDV
Hooks, Rectal	15673	SGINAC
Hooks, Rhinoplastic	12033	SGINAC
Hooks, Scleral Fixation	12034	SGINAC
Hooks, Skin	12035	SGINAC
Hooks, Strabismus	12036	SGINAC
Hooks, Sympathectomy	12037	SGINAC
Hooks, Tracheal	12038	SGINAC

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Hooks, Other	15225	SGINAC
Hospital Grade	*****	DISINF
Hospital Grade Disinfectant	*****	DISINF DISIWC
Hot Dry Air/Fluidized Medium Units	16982	HUMIDS
Hot Moist Air Therapy Units	17110	HUMIDS
Household/Commercial Grade Disinfectant	*****	DISIWC
Humidification Chambers	*****	HUMIDS
Humidifiers, Heat/Moisture Exchange	15645	HUMIDS
Humidifiers, Heated	12050	HUMIDS
Humidifiers, Nonheated	12051	HUMIDS
Hydroxyapatite	*****	MATERL
Hyperthermia Units, Circulating-Air	17700	THERCO
Hyperthermia Units, Circulating-Fluid	17648	THERCO
Hyperthermia Units, Forced-Air	17950	THERCO
Hyperthermia Units, Radio-Frequency	12066	THERCO
Hyperthermia Units, Ultrasonic	17571	THERCO
Hypo/Hyperthermia Blankets	12067	THERCO
Hypo/Hyperthermia Units, Central	12074	THERCO
Hypo/Hyperthermia Units, Mobile, General-Purpose	12075	THERCO
Hypo/Hyperthermia Units, Mobile, Heart-Lung Bypass	17206	THERCO
Hypothermia Units	12078	THERCO
Hysteroscopes	12081	ENDOSC

**I**

Image Intensifiers	15963	RADGEQ
Immobilizers, Cervical	15642	SGINAC
Immobilizers, Infant	12097	SGINAC
Impedance Cardiographs	121033	CARDIO
Impedance Pneumograph Monitors	12621	RESPMA
Implantable Pacemaker Analyzers	17637	PMKRAC
Implantable Pacemaker Analyzers, Transtelephonic	17638	PMKRAC
Incision and Drainage Kits	12108	DRKTAC
Incontinence Clamps	12109	INCOEQ
Incontinence Control Unit, Fecal	17505	INCOEQ
Incubators, Infant, Mobile	17432	INCINF
Incubators, Infant, Transport	12114	INCINF
Inflators, Angioplasty Balloon	17541	ANGPDV
Inflators, Catheters, Non-Angioplasty	*****	CATKEQ
Information Systems, Radiology	17175	NEURDE
Infusion Controller Administration Sets	17984	IVASAC
Infusion Controllers	11010	DRGINF IVASAC
Infusion Pump Administration Sets	16579	DRGINF

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Infusion Pumps, Ambulatory, Insulin	17159	DRGINF
Infusion Pumps, Ambulatory, Other	16721	DRGINF
Infusion Pumps, Ambulatory, Syringe	15782	DRGINF
Infusion Pumps, General Purpose	13215	DRGINF
Infusion Pumps, Implantable	16143	DDSIPA
		DRGINF
Infusion Pumps, Micro	16722	DRGINF
Infusion Pumps, Multichannel	17634	DRGINF
Infusion Pumps, Patient-Controlled Analgesic	16924	DDSIPA
		DRGINF
Infusion Pumps, Syringe	13217	DRGINF
Infusion Pumps, Syringe, Oxytocin	16167	DRGINF
Infusion Systems, Insulin, Closed Loop	15781	DRGINF
Infusion-Withdrawal Pumps	13220	DRGINF
Infusion Wires	*****	SGINAC
Inhalers, Aerosol	12128	ATINNB
Injection/Infusion Ports, Implantable	16854	IJIFPT
Injectors, Angiographic	15284	SYRINJ
Injectors, Contrast Media	17968	SYRINJ
Injectors, Contrast Media, Computed Tomography	17969	SYRINJ
Injectors, Contrast Media, Lymphangiography	15285	SYRINJ
Injectors, Hand-Held	16594	SYRINJ
Injectors, Insulin	*****	SYRINS
Injectors, Medication/Vaccine	12504	SYRINJ
Injectors, Pigmentation	17660	SYRINJ
Injectors, Medication/Vaccine, Syringe	12132	SYRINJ
Injectors, Thermal Dilution	12134	CAOUEQ
Injectors, Uterine	15622	SYRINJ
Insect Sting Emergency Kits	12135	FAEMKT
Inserters	17543	SGINAC
Inserters, Intraocular Lens	17544	SGINAC
Insufflators	12144	INSUFF
Insufflators, Arthroscopic	17216	ORTHEQ
Insufflators, Hysteroscopic	16848	INSUFF
Insufflators, Laparoscopic	16849	INSUFF
Insufflators, Ophthalmic	17924	INSUFF
Insulation Pads, Cardiac	17120	THERCO
Insulin Pump Administration Sets	16588	DRGINF
Intermediate Level Disinfectant	*****	DISIWC
Intracorneal Rings	18103	GRFIMP
Intrauterine Balloons	12155	CATVAR
Intrauterine Devices	12156	IUDEVC
Intravenous Administration Sets, Buret	12159	IVASAC
Intravenous Administration Sets, General-Purpose	16649	IVASAC
Intravenous Administration Sets, Nonabsorbing	16631	IVASAC

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Intravenous Catheter - Needle Immobilizers	12162	IVASAC
Intravenous Extension Tubing Kits	12170	IVASAC
Intravenous Flow Regulators, Dial-Calibrated	16789	IVASAC
Intravenous Fluid Containers	12172	IVASAC
Intubation Scopes, Flexible Fiberoptic	16331	ENDOSC
In Vitro Fertilisation Instrumentation	17248	ARTFER
Ionized Air Generators	12184	ATINNB
Iontophoresis Units	12185	IONPHO
Irrigation/Distension Systems	17676	IRILAV
Irrigation/Distension Systems, Arthroscopic	17233	IRILAV
Irrigation/Distension Systems, Hysteroscopic	17677	IRILAV
Irrigation Kits	15585	IRILAV
Irrigation Kits, Bladder	10406	IRILAV
Irrigation Kits, Ear	11371	IRILAV
Irrigation Kits, Eye	11655	IRILAV
Irrigation Kits, Perineal	12996	IRILAV
Irrigation Kits, Urinary	14301	IRILAV
Irrigation Kits, Wound	14462	IRILAV
IMV Supplemental Gas Circuits	16803	VBDVAC

**J**

Jaw Fracture Fixation Devices	15801	ORTHEQ
-------------------------------	-------	--------

**K**

Keratomes	12222	OPHINS SGINAC
Keratoscopes	16347	OPHDEQ
Kinesthesiometers	12229	NEURDE
Knife Handles	12235	SGINAC
Knife Needles	12237	NEEDAC
Knives, Adenotome	16583	SGINAC
Knives, Amputating	12240	SGINAC
Knives, Cataract	12241	SGINAC
Knives, Commissurotomy	12242	SGINAC
Knives, Dermatome	12243	SGINAC
Knives, Dural Hook	12244	SGINAC
Knives, Ear	12245	SGINAC
Knives, Keratome	12247	SGINAC
Knives, Meniscus	12248	SGINAC
Knives, Scalpel	12252	SGINAC
Knives, Scalpel Haemostatic	17705	ESUNAC
Knives, Skin Grafting	12253	SGINAC
Knives, Sternum	12254	SGINAC

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Knives, Tonsil	12255	SGINAC
Knives, Other	15226	SGINAC
Knot-Tying Instrumentation, Endoscopic	17971	ENDOSC

**L**

Laminaria Tents	12270	GYNOBE
Lancets, Blood	10440	BLOCOL
Lancing Devices, Blood	16380	BLOCOL
Laparoscope Holder	16240	SGINAC
Laparoscopes	12291	ENDOSC
Laparotomy Rings	15052	SGINAC
Laryngoscopes, Flexible	15075	ENDOSC
Laryngoscopes, Rigid	15076	ENDOSC
Laryngostroboscopes	12294	ENDOSC
Larynx, Artificial	10204	SPTHUN
Laser Contact Probes	17193	LASERS
Laser Beam Analysis Systems	17973	LASERS
Laser Delivery System, Fiberoptic	17807	LASERS
Laser Delivery Systems	17881	LASERS
Laser Delivery Systems, Waveguide	17866	LASERS
Lasers, Diode, Surgical	18183	LASERS
Lasers, Ophthalmic, Argon	16945	LASERS
Lasers, Ophthalmic, Argon/Krypton	17773	LASERS
Lasers, Ophthalmic, Diode	17808	LASERS
Lasers, Ophthalmic, Dye	17481	LASERS
Lasers, Ophthalmic, Excimer	17702	LASERS
Lasers, Ophthalmic, Krypton	16946	LASERS
Lasers, Ophthalmic, Nd:YAG	16947	LASERS
Lasers, Ophthalmic, Other	17446	LASERS
Lasers, Photocoagulator	12297	LASERS
Lasers, Surgical, Argon	16941	LASERS
Lasers, Surgical, Carbon Dioxide	16942	LASERS
Lasers, Surgical, Carbon Dioxide, General-Purpose	17771	LASERS
Lasers, Surgical, Carbon Dioxide, Specialty	17772	LASERS
Lasers, Surgical, Copper Vapor	17959	LASERS
Lasers, Surgical, Dye	17482	LASERS
Lasers, Surgical, Erbium:YAG	17768	LASERS
Lasers, Surgical, Excimer	17161	LASERS
Lasers, Surgical, Holmium:YAG	17769	LASERS
Lasers, Surgical, Nd:YAG	16943	LASERS
Lasers, Surgical, Nd:YAG/Carbon Dioxide	17774	LASERS
Lasers, Surgical, Ruby	17815	LASERS
Lasers, Surgical, Other	17447	LASERS
Lasers, Therapeutic	12299	LASERS

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Lasers, Other	16948	LASERS
Lavage Units, Colonic	12301	IRILAV
Lavage Units, Corneal	12302	IRILAV
Lavage Units, Dental-Oral	12304	IRILAV
Lavage Units, Otorhinolaryngological	12305	IRILAV
Lavage Units, Gastric	17189	IRILAV
Lavage Units, Surgical	12306	IRILAV
Lavage Units, Uterine	17520	IRILAV
Lead Introducer	*****	STIMAC
Leads	16882	AISTIM
Leads, Dorsal Column, Implantable	17215	AISTIM
Leads, Implantable Defibrillator	16653	LDDFIP
Leads, Implantable Defibrillator/Pacemaker	17964	LDDFIP
Leads, Pacemaker	16128	PMKRAC
Leads, Pacemaker, Implantable Endocardial	12900	LDPMIE
Leads, Pacemaker, Implantable Myocardial	11458	LDPMIO
Leads, Pacemaker, Temporary Endocardial	16129	PMKRAC
Leads, Pacemaker, Temporary Myocardial	16989	PMKRAC
Leads, Pacemaker, Transesophageal	17511	PMKRAC
Lens Expressors	12318	SGINAC
Lens Loops	12319	SGINAC
Lens Spoons	12321	SGINAC
Lenses, Intraocular, Anterior Chamber	16068	IOLATC
Lenses, Intraocular, Iridocapsular Fixation	16069	IOLICF
Lenses, Intraocular, Iris Fixation	16070	IOLIFX
Lenses, Intraocular, Posterior Chamber	16071	IOLPMM IOLPTC
Lesion Generators	12327	ESUNAC
Lifts, Patient Transfer	12330	AMBEQP
Ligators	12232	IMPABS
Ligators, Hemorrhoid	12335	SGINAC
Ligators, Umbilical	15041	SGINAC
Ligature Carriers	12336	SGINAC
Lights, Dental Activator	16386	DENTEQ
Lights, Dental, Intraoral	12352	DENTEQ
Lights, Infrared	12278	PHOTTH
Lights, Phototherapy	16003	PHOTTH
Lights, Surgical	12282	SGINAC
Lights, Ultraviolet, Germicidal	12284	PHOTTH
Lights, Ultraviolet, Other	15762	PHOTTH
Lights, Ultraviolet, UVA	15731	PHOTTH
Linear Accelerators	12364	RADTUN
Liquid Oxygen Units, Individual	16853	MGASEQ
Lithotripters	16228	LTPLTT
Lithotripters, Electrohydraulic	16229	LTPLTT

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Lithotripters, Extracorporeal	16758	LTPLTT
Lithotripters, Extracorporeal, Gallstone	17454	LTPLTT
Lithotripters, Extracorporeal, Renal	17455	LTPLTT
Lithotripters, Laser	17468	LASERS
Lithotripters, Ultrasonic	16230	LTPLTT
Lithotrites	12379	LTPLTT
Locking Mechanisms, IV Wye Line	17501	IVASAC
Loops, Endarterectomy	12394	SGINAC
Loops, Tracheal	12396	SGINAC
Loops, Vascular	12397	SGINAC
Low Level Disinfectant	*****	DISIWC
Lubricants, Contact Lens	*****	COLECP
Lubricants, Instrument	15772	LUBRIC
Lubricating Jellies	12401	LUBRIC
Lumbar Puncture Kits	12404	LBPTKT
Lumboscopes	15790	ENDOSC
Lymphangiographic Kits	12409	RADACC

**M**

Magnetic Resonance Imaging (MRI) Units	16260	MRIMUN
Magnetocardiographs	12416	CARDIO
Magnetometers	12417	NEURDE
Magnets, Eye	15269	OPHINS
Magnifiers, Operating	15653	MICROP
Mallets, Bone	12423	SGINAC
Mallets, Other	15229	SGINAC
Manifolds, Gas	15042	STOPCK
Manifolds, Liquid	15043	STOPCK
Manifolds, Stopcock	16387	STOPCK
Manometer Sets, Cerebrospinal Fluid	13107	LBPTKT
Manometer Sets, Venous, Central/Peripheral	10776	CARDMA
Manometers, Aneroid	12432	PRMSEQ
Markers, Implantable Radiographic	12441	RADACC
Markers, Ocular	12442	OPHINS
Masks	12447	RESPMA
Masks, Air-Oxygen, Aerosol Administration	12449	ATINNB
Masks, Air-Oxygen, Nonrebreathing	12450	MGASEQ
Masks, Air-Oxygen, Partial Rebreathing	12451	MGASEQ
Masks, Air-Oxygen, Venturi	12452	MGASEQ
Masks, Air-Oxygen, Other	15231	MGASEQ
Masks, Analgesia	16194	AGUNAC
Masks, Anesthesia	12453	ANEQAC
Masks, Eye, Phototherapy	11659	PHOTTH
Masks, Mouth-to-Mask, Resuscitation Training	17550	RESCAQ

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Masks, Other	15230	AWPSDA
Masks, Resuscitator	17170	RESCAQ
Masks, Surgical	12458	APPPRO
Masks, Tracheostomy	12459	TTCEQU
Masks, Tracheostomy, Aerosol Administration	12460	ATINNB
Mediastinoscopes	12499	ENDOSC
Medicated Pads	12934	BDDRST
Medical Gas Cylinders	16501	MGASEQ
Mesh Bags, Pacemaker	16038	PMKRAC
Mesh, Metallic	16048	SGINAC
Mesh, Polyglycolic Acid	17247	GRFIMP
Mesh, Polymeric	12510	GRFIMP
Metabolic Carts	10297	BSMTUN
Metabolic Carts, Stress Exercise	17895	BSMTUN
Metal Detectors, Magnetic	12516	LOSEFI
Metal Detectors, Photographic	12518	LOSEFI
Metal Detectors, Ultrasonic	12519	LOSEFI
Metreurynters	12524	GYNOBE
Microelectrode Amplifiers	12529	ELECTR
Microelectrodes	12531	ELECTR
Microscopes	12536	MICROP
Microscopes, Corneal	12537	MICROP
Microscopes, Dermatologic	17719	MICROP
Microscopes, Ear	12538	MICROP
Microscopes, Operating	12539	MICROP
Microsurgical Instruments	15621	SGINAC
Middle Ear Analyzers	15634	ENTEQU
Mid Stream Collection Kits	12542	UROCOL
Monitor, Urea	*****	HDHFEQ
Monitors, CPR	*****	RESCAQ
Motion Analysis Systems	17929	PHYSEV
Moulds, Vaginal	14322	SGINAC
Multi Dose Drug Vial Adapters	17896	FITADP
Multiple Medical Gas Monitors	17443	BLGSMA
Multiple Medical Gas Monitors, Respired	17444	BLGSMA
Multiple Medical Gas Monitors, Respired/Anesthetic	17445	BLGSMA
Musculoskeletal Analysis System	17763	PHYSEV
Myelogram Kits	12689	RADACC
Myeloscopes	12690	ENDOSC
Myographs	12692	EDGEQA

**N**

Nails, Bone	16078	GRFIMP
Nasal Septum Buttons	12704	SGINAC



**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Nasogastric Tube Holder	17514	TUBEST
Nasopharyngoscopes	12709	ENDOSC
Nasoscopes	12710	ENDOSC
Nebulizers, Heated	12716	ATINNB
Nebulizers, Nonheated	15045	ATINNB
Nebulizers, Ultrasonic	12719	ATINNB
Needle Adapters	16829	NEEDAC
Needle Guides	15868	SGINAC
Needle Holders	12726	SGINAC
Needle Introducers	12727	SGINAC NEEDAC
Needles, Acupuncture	12730	NEEDAC
Needles, Airway	12732	NEEDAC
Needles, Anesthesia	12733	NEEDAC
Needles, Aspirating	16774	NEEDAC
Needles, Biopsy	12734	NEEDAC
Needles, Blood-Collecting	12736	NEEDAC
Needles, Blunt	12737	NEEDAC
Needles, Cardiac	12738	NEEDAC
Needles, Cholangiographic	12739	NEEDAC
Needles, Dental	12740	DENTEQ
Needles, Dialysis	12741	NEEDAC
Needles, Endoscopic Sclerotherapy	17569	NEEDAC
Needles, Epidural	17834	NEEDAC
Needles, Eye	12744	NEEDAC
Needles, Filter	16266	NEEDAC
Needles, Hypodermic	12745	NEEDAC
Needles, Intra-Arterial	12747	NEEDAC
Needles, Intravenous	12748	NEEDAC
Needles, Pneumoperitoneal	12750	NEEDAC
Needles, Radiographic	12751	NEEDAC
Needles, Scalp-Vein Type	12752	NEEDAC
Needles, Spinal	12753	NEEDAC
Needles, Subcutaneous Injection/Infusion Port	17180	NEEDAC
Needles, Suture	12754	NEEDAC
Needles, Tine Test	15679	NEEDAC
Needles, Transfer	16627	NEEDAC
Needles, Vented	16628	NEEDAC
Needles, Other	15234	NEEDAC
Nephroscopes, Flexible	15289	ENDOSC
Nephroscopes, Rigid	15290	ENDOSC
Nerve Function Monitors	17582	NEURDE
Nerve Locators	16253	LOSEFI
Nerve Locators, Facial	16254	LOSEFI
Nippers, Malleus	12766	SGINAC

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Nitrous Oxide Analyzers	12768	BLGSMA
Non-absorbant menstrual collector	*****	MECLDV

**O**

Obturators	12789	GYNOBE
Occluders	17730	AADIMP
Occluders, Vascular	17731	GRFIMP
Oculoplethysmographs	12791	OPHDEQ
Ophthalmic Accessories	*****	OPHINS
Ophthalmic Applicators	*****	OPHINS
Ophthalmic Coagulations Brushes	*****	OPHINS
Ophthalmic Fixation Units	16466	OPHDEQ
Ophthalmic Foil	16066	SGINAC
Ophthalmodiaphanosopes	12804	OPHDEQ
Ophthalmodiastimeters	12805	OPHDEQ
Ophthalmographs	12809	OPHDEQ
Ophthalmoleukoscopes	12810	OPHDEQ
Ophthalmometers	12811	OPHDEQ
Ophthalmometrosopes	12813	OPHDEQ
Ophthalmoscopes	12815	OPHDEQ
Ophthalmoscopes, Direct	12817	OPHDEQ
Ophthalmoscopes, Indirect	12818	OPHDEQ
Ophthalmotropometers	12820	OPHDEQ
Organ Preservation Solutions	18004	EXCOTH
Orthodontic Instruments	12829	DENTEQ
Orthopedic External Fixation Systems, Fracture	15767	ORTHEQ
Orthopedic External Fixation Systems, Spinal	16046	ORTHEQ
Orthopedic External Fixation Systems, Other	16059	ORTHEQ
Orthopedic Internal Fixation Systems, Fracture	12833	GRFIMP
Orthopedic Internal Fixation Systems, Spinal	15766	GRFIMP
Orthopedic Internal Fixation Systems, Other	16031	GRFIMP
Orthoses, Elbow	17910	ORTHOS
Orthoses, Wrist	17911	ORTHOS
Oscillotonometers	12841	BLPRDV
Osteotomes, Orthopedic	12845	SGINAC
Osteotomes, Plastic Surgery	12846	SGINAC
Osteotomes, Other	15235	SGINAC
Ostomy Appliances	12847	OSTAPL
Ostomy Collection Bags	16459	OSTAPL
Otosopes	12849	ENDOSC
Oximeters, Cerebral	17942	BLGSMA
Oximeters, In Vitro	15201	BLGSMA
Oximeters, In Vitro, Dual-Wavelength	17655	BLGSMA
Oximeters, In Vitro, Multiwavelength	16618	BLGSMA

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Oximeters, Intracardiac	15200	BLGSMA
Oximeters Pulse	17148	BLGSMA
Oxygen Administration Kits	12855	MGASEQ
Oxygen Analyzers, Continuous Controlling	12861	BLGSMA
Oxygen Analyzers, Sampling	15015	BLGSMA
Oxygen Concentrators	12873	MGASEQ
Oxygen Generators	12883	MGASEQ
Oxygen Meters	15046	BLGSMA
Oxygen Monitors	12863	BLGSMA
Oxygen Monitors, Transcutaneous	12592	BLGSMA
Oxygen-Air Proportioners	12876	MGASEQ
Oxygenators, Extracorporeal Bubble	11974	CPBYEQ
Oxygenators, Extracorporeal Membrane	17643	CPBYEQ
Oxygenators, Regional Perfusion	13319	RPUNAC

**P**

pH Monitors	12634	EDGEQA
Pacemaker Accessories	*****	PMKRAC
Pacemaker Lead Adapters, Implantable	17635	LDADPM
Pacemaker Programmers	15993	PMKRAC
Pacemakers, Cardiac, External Invasive	12912	PMKRAC
Pacemakers, Cardiac, External Noninvasive	16516	PMKRAC
Pacemakers, Cardiac, Implantable	12913	PMRIMP
Pacemakers, Cardiac, Transesophageal	17494	PMKRAC
Pacers, Counting, CPR	15771	RESCAQ
Pachymeters	16330	OPHDEQ
Pacing System Analyzers, Intraoperative	17639	PMKRAC
Packers, Gauze	15047	SGINAC
Pads, Circulating-Fluid, Localized Heating/Cooling	11988	THERCO
Pain Gauges, Cutaneous	12950	NEURDE
Papillotomes	15625	SGINAC
Parturiometers	12962	GYNOBE
Patches, Myocardial	13165	LDDFIP
Pelviscopes	12981	ENDOSC
Penile Vacuum Devices	17744	TORQAC
Perforators, Amniotic Membrane	12990	SGINAC
Perforators, Cranial, Automatic	11334	SGDRIL
Perforators, Other	15236	SGINAC
Perimeters, Automatic	16918	OPHDEQ
Peritoneal Dialysis Catheter Adapters	16963	PDEQAC
Peritoneal Dialysis Sterile Connection Units	17434	PDEQAC
Peritoneal Dialysis Unit Tubing Sets	11227	PDEQAC
Peritoneal Dialysis Units	11226	PDEQAC
Peritoneoscopes	13003	ENDOSC

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Personal Protection Packs	17232	APPPRO
Pessaries	13004	GYNOBE
Pessaries, Diaphragm	13005	DIAPHR
Phaneroscopes	13006	ENDOSC
Pharyngoscopes	13011	ENDOSC
Pheresis Sets	16811	BLOCOL
Phlebographs	13012	CARDMA
Phonocardiographs	13017	CARDIO
Phonocatheters	13018	CARDIO
Phonostethographs	13020	CARDIO
Photocatheters	13022	CATVAR
Photocoagulators	13024	COAGTS
Photofluoroscopes	13027	RADGEQ
Phototherapy Units, Hyperbilirubinemia	17515	PHOTTH
Phototherapy Units, Other	17516	PHOTTH
Physiologic Monitoring Systems, Acute Care	12647	CARDMA
Physiologic Monitoring Systems, Cardiac Catheterization	12648	CATKEQ
Physiologic Monitoring Systems, Impotence	17441	EDGEQA
Physiologic Monitoring Systems, Neonatal	15791	FENEMA
Physiologic Monitoring Systems, Stress Exercise	12649	PHYSEV
Physiologic Monitoring Systems, Stress Exercise, Cardiac	17723	PHYSEV
Physiologic Monitoring Systems, Stress Exercise, Pulmonary	17724	PHYSEV
Physiologic Monitoring Systems, Telemetric	13987	TELEME
Physiologic Recording Systems	16306	EDGEQA
Picks, Middle Ear	13039	SGINAC
Picks, Vein	*****	SGINAC
Pins, Bone	16085	GRFIMP
Pins, Bone, Absorbable	*****	IMPABS
Pins, Dental Retention	16700	DENTEQ
Plates, Bone	13050	GRFIMP
Plates, Cranioplasty	13052	GRFIMP
Pledgets	13055	SPONGE
Plethysmographs	13056	PLETHS
Plethysmographs, Body	13059	PLETHS
Plugs	15847	AUGPSM
		GRFIMP
Plugs, Enuresis	*****	INCOEQ
Pneumatic Counterpressure Suits	13069	PREREA
Pneumo Sleeve	*****	SGINAC
Pneumographs	13078	RESPMA
Pneumotachographs	13080	RESPMA
Pneumothorax Units	13081	PNUMTX
Point, Paper, Endodontic	*****	MTDTNR
Polysomnographs	17458	EDGEQA
Position Emission Tomography (PET) Systems	16375	NUCMEQ

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Positioning Aids	16223	OPHINS
Positioning Aids, Radiographic/Radiotherapeutic	13260	RADTUN
Power Meters, Ultrasonic, Diagnostic	16842	ULTRDV
Power Meters, Ultrasonic, Therapeutic	16843	ULTRDV
Prep Kits	13097	ANESKT
Pressure Alarms, Airway	14351	AWPSDA
Pressure Infusor, Powered	*****	INFDEV
Pressure Infusors	13100	IVASAC
Pressure Measuring Unit	13102	PRMSEQ
Pressure Measuring Units, Compartmental	17196	PRMSEQ
Pressure Measuring Units, Esophageal	13108	PRMSEQ
Pressure Measuring Units, Inspiratory	13109	PRMSEQ
Pressure Measuring Units, Intracranial	13110	PRMSEQ
Pressure Monitoring Extension Tubing Sets	16621	PRMSEQ
Pressure Monitoring Tubing Sets	15640	PRMSEQ
Pressure Monitors, Airway	17425	RESPMA
Pressure Monitors, Blood, General/Invasive	16764	CARDMA
Pressure Monitors, Intracranial	16763	PRMSEQ
Pressure Pads, Alternating, Disposable	16144	PREREA
Pressure Pads, Alternating, Reusable	16145	PREREA
Pressure Pads, Gel	12944	PREREA
Pressure Reducers, Intraocular	15048	OPHDEQ
Probes	13117	LTPLTT
Probes, Arthroscopy	16980	ORTHEQ
Probes, Common Duct	13118	SGINAC
Probes, Eye	13119	OPHINS
Probes, Eye, Lacrimal	13120	OPHINS
Probes, Eye, Other	15238	OPHINS
Probes, Fistula	13121	SGINAC
Probes, Gynaecological	*****	GYNOBE
Probes, Microlaryngeal	15877	SGINAC
Probes, Nuclear	*****	NUCMEQ
Probes, Other	15237	SGINAC
Probes, Oximeter	17594	BLGSMA
		TRDSNS
Probes, Sinus	13124	SGINAC
Probes, Thermometer	13125	BODYTP
Probes, Ultrasonic	16272	ULTRDV
Procedure Trays, Aspiration	17716	KTPWOD
Procedure Trays, Surgical	17168	KTPWOD
Proctoscopes	13126	ENDOSC
Proctosigmoidoscopes	13127	ENDOSC
Prophylactic Units, Dental, Ultrasonic	16693	DENTEQ
Prophylactic Units, Dental, Water	16694	DENTEQ
Prostatic Hypertrophy Treatment, Microwave	17905	THERCO

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Prostheses, Arm	13134	GRFIMP
Prostheses, Artificial Heart, Temporary	16977	GRFIMP
Prostheses, Bile Duct	13135	GRFIMP
Prostheses, Bladder Support	13136	GRFIMP
Prostheses, Blood Vessel, Artificial	13177	GRFIMP
Prostheses, Blood Vessel, Biological	11913	PROVAO
Prostheses, Bone, Freeze Dried	16862	AUGPSM
Prostheses, Bowel Sphincter	*****	GRFIMP
Prostheses, Breast, Internal	16111	PROBIP
Prostheses, Cardiac Valve, Artificial	15869	PHVART
Prostheses, Cardiac Valve, Biological	15870	PHVHAO
Prostheses, Carpal	16089	GRFIMP
Prostheses, Cerclage, Bone	16140	GRFIMP
Prostheses, Cerclage, Cervix	13139	GRFIMP
Prostheses, Chin	13140	GRFIMP
Prostheses, Cochlear	16035	GRFIMP
Prostheses, Dental, Implantable	16744	DENTEQ
Prostheses, Ear	13143	GRFIMP
Prostheses, Esophageal	13145	GRFIMP
Prostheses, Eustachian Tube	15324	GRFIMP
Prostheses, Evisceration	16067	GRFIMP
Prostheses, Eyelid	15049	GRFIMP
Prostheses, Fallopian Tube	13148	GRFIMP
Prostheses, Femoral Head	16125	GRFIMP
Prostheses, Femoral, Upper	13149	GRFIMP
Prostheses, Finger/Thumb	16094	GRFIMP
Prostheses, Foot	13150	GRFIMP
Prostheses, Foot Arch	16135	GRFIMP
Prostheses, Gastroesophageal Antireflux	16165	GRFIMP
Prostheses, Hair	16611	GRFIMP
Prostheses, Hand	13151	GRFIMP
Prostheses, Heart	13152	GRFIMP
Prostheses, Joint	13158	GRFIMP
Prostheses, Joint, Ankle	16132	GRFIMP
Prostheses, Joint, Ankle, Tibial Component	16088	GRFIMP
Prostheses, Joint, Elbow	16105	GRFIMP
Prostheses, Joint, Hip	16150	GRFIMP
Prostheses, Joint, Hip Acetabular Component	16084	GRFIMP
Prostheses, Joint, Hip Femoral Component	16095	GRFIMP
Prostheses, Joint, Knee	16072	GRFIMP
Prostheses, Joint, Knee, Femoral Component	16097	GRFIMP
Prostheses, Joint Knee, Total	16096	GRFIMP
Prostheses, Joint, Mandible	16962	GRFIMP
Prostheses, Joint, Shoulder	16099	GRFIMP
Prostheses, Joint, Toe	16134	GRFIMP

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Prostheses, Joint, Wrist	13168	GRFIMP
Prostheses, Joint, Temporomandibular, Total	17479	GRFIMP
Prostheses, Larynx	13159	GRFIMP
Prostheses, Leg	13160	GRFIMP
Prostheses, Ligament	16242	GRFIMP
Prostheses, Mammary, Internal, Gel-Filled	17855	PROBIP
Prostheses, Mammary, Internal, Inflatable	17854	PROBIP
Prostheses, Mandibular	16047	GRFIMP
Prostheses, Maxillofacial	13162	GRFIMP
Prostheses, Meningeal	15961	GRFIMP
Prostheses, Muscle	13164	GRFIMP
Prostheses, Nasal, Dorsal	16036	GRFIMP
Prostheses, Nerve Sheath	16057	GRFIMP
Prostheses, Nipple	16130	GRFIMP
Prostheses, Orbital Rim	16113	GRFIMP
Prostheses, Ossicular, Incus and Stapes	16107	GRFIMP
Prostheses, Ossicular, Total	16108	GRFIMP
Prostheses, Patellar	16121	GRFIMP
Prostheses, Pericardial	17800	GRFIMP
Prostheses, Penile	13167	GRFIMP
Prostheses, Penile Inflatable	17997	GRFIMP
Prostheses, Radial Head	17101	GRFIMP
Prostheses, Radial/Ulna Head	16259	GRFIMP
Prostheses, Rhinoplasty	16115	GRFIMP
Prostheses, Tendon	13171	GRFIMP
Prostheses, Testicle	13172	GRFIMP
Prostheses, Tibial	13173	GRFIMP
Prostheses, Tracheal	15050	GRFIMP
Prostheses, Tracheobronchial	13174	GRFIMP
Prostheses, Urethral	13175	GRFIMP
Prostheses, Urethral Sphincter	13176	GRFIMP
Prostheses, Zygomatic	16116	GRFIMP
Prosthesis Implantation Instruments, Orthopedic	13180	SGINAC
Protectors, Waterproofing, Cast-Bandage	17583	MTC SPL
Pulmonary Function Analyzers, Adult	13182	RESPMA
Pulmonary Function Analyzers, Pediatric	17699	RESPMA
Pulsatile Pressure Generators	13188	CPBYEQ
Pulsatile Pressure Generators, Heart-Lung Bypass Unit	11975	CPBYEQ
Pulse Monitors	12657	CARDMA
Pumps, Alternating Pressure Pad	16076	PREREA
Pumps, Breast	10485	ASDVAC
Pumps, Dye Dilution	13201	CAOUEQ
Pumps, Enteral Feeding	13209	EFEDKT
Pumps, Extracorporeal Perfusion	13203	CPBYEQ
Pumps, Insufflating	16850	INSUFF

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Pumps, Obstetrical Decompression	13224	DECPEQ
Pumps, Perfusion Manometer	*****	PRMSEQ
Pumps, Urinary Collection Bag	16790	UROCOL
Punches, Aortic	13229	SGINAC
Punches, Biopsy	13230	SGINAC
Punches, Bone	13231	SGINAC
Punches, Corneoscleral	13232	SGINAC
Punches, Dermal	13233	SGINAC
Punches, Spinal	16939	SGINAC
Punches, Other	15239	SGINAC
Punctum Plugs, Lacrimal	17958	GRFIMP
Pupillographs	13235	OPHDEQ
Pupillometers	13236	OPHDEQ
Pyeloscopes	13237	ENDOSC

**R**

Radioaerosol Administration Sets	16784	NUCMEQ
Radiographic Units, Chest	10822	RADGEQ
Radiographic Units, Cystoscopy	15946	RADGEQ
Radiographic Units, Dental	13269	DENTEQ
Radiographic Units, General-Purpose	13271	RADGEQ
Radiographic Units, Head	15947	RADGEQ
Radiographic Units, Mammographic	12425	RADGEQ
Radiographic Units, Mobile	13272	RADGEQ
Radiographic Units, Pneumoencephalographic	16555	RADGEQ
Radiographic Units, Podiatric	17669	RADGEQ
Radiographic Units, Polytomographic	16556	RADGEQ
Radiographic Units, Skeletal	16557	RADGEQ
Radiographic Units, Tomographic	14059	RADGEQ
Radiographic/Fluoroscopic Table Systems	16885	RADGEQ
Radiographic/Fluoroscopic Table Systems, Urological	16212	FURNIP
Radiographic/Fluoroscopic Units, Angiographic	16597	RADGEQ
Radiographic/Fluoroscopic Units, Cardiac Catheterization	17192	RADGEQ
Radiographic/Fluoroscopic Units, General-Purpose	11757	RADGEQ
Radiographic/Fluoroscopic Units, Mobile	11758	RADGEQ
Radiographic/Fluoroscopic Units, Special Procedure	15951	RADGEQ
Radiographic/Tomographic Table Systems	16884	RADGEQ
Radiometers, Laser	17147	LASERS
Radiotherapy Beam Block Shaping System	17502	RADTUN
Radiotherapy, Delivery System	*****	RADTUN
Radiotherapy Units, Cobalt	16972	RADTUN
Radiotherapy Units, Orthovoltage	16973	RADTUN
Raspatories	13287	SGINAC
Rasps, Bone	15741	SGINAC



**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Rasps, Middle Ear	13289	SGINAC
Rasps, Nasal	13290	SGINAC
Rasps, Plastic Surgery	13291	SGINAC
Rasps, Other	15240	SGINAC
Razors, Skin Prep	16580	SGINAC
Reagent, Calibrator	17030	DIAHMU
Reagent, Quality Control, Coagulation	17059	DIAHMU
Reagent, Quality Control Other	17069	DIAHMU
Reagent Strips, Ophthalmological	*****	OPHDEQ
Reagents, Hepatitis C Antibody	17914	DIAHCV
Reagents, Hepatitis C, Nucleic Acid	*****	DIAHCV
Reagents, Hepatitis Surface Antigen	17360	DIAHCV
Reagents, Human Immunodeficiency Virus Antibody	17367	DIAHIV
Reagents, Human Immunodeficiency Virus Antigens	*****	DIAHIV
Reagents, Human Immunodeficiency Virus, Nucleic Acid	*****	DIAHIV
Reagents, Occult Blood, Fecal	17410	DIAHMU
Reagents, Ovulation Prediction	17485	DIAHMU
Reagents, Pregnancy Testing	17484	DIAHMU
Reagents, Strip	17418	OPHDEQ
Reagents, Strip, Cholesterol	17656	DIAHMU
Reagents, Strip, Glucose	17419	DIAHMU
Reagents, Strip, Lactate	*****	DIAHMU
Reagents, Strip, Ophthalmological	*****	OPHDEQ
Reagents, Strip, Urinalysis	17480	DIAHMU
Reamers, Orthopedic	13297	SGINAC
Reconstructions Blocks, Plastic Surgery	16043	MATERL
Reconstruction Sheeting, Plastic Surgery	16118	MATERL
Recorders, Echocardiographic	11385	ULTRDV
Recorders, Long-Term, pH	16888	EDGEQA
Recorders, Long Term, Blood Pressure, Portable	12386	BLPRDV
Recorders, Long Term, ECG	12387	CARDIO
Recorders, Long Term, ECG, Portable	12388	CARDIO
Recorders, Long-Term, EEG	12389	EEGEQA
Recorders, Long-Term, Oxygen	12390	EDGEQA
Recorders, Long-Term, Respiration	16305	RESPMA
Recorders, Long-Term, Temperature	17572	EDGEQA
Recorders, Long-Term, Trend	12391	EDGEQA
Recorders, Sphygmomanometers, Automatic	16271	BLPRDV
Recorders, Tape, ECG	16009	CARDIO
Recorders, Tape, Physiologic	16008	EDGEQA
Recorders, Temperature	16543	BODYTP
Refractors, Ophthalmologic	13313	OPHDEQ
Regional Perfusion Units	13317	RPUNAC
Regulators, High-Pressure Gas	13323	MGASEQ
Regulators, Suction, Low Volume	13329	ASDVAC

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Regulators, Suction, Surgical	15051	ASDVAC
Regulators, Suction, Thoracic	13332	ASDVAC
Regulators, Suction, Tracheal	13333	ASDVAC
Resectoscopes	13335	ENDOSC
Reservoirs, Cerebral Spinal Fluid	15974	SHUNTS
Reservoirs, Cardiotomy	13338	CPBYEQ
Resinous Compounds, Cranial	16131	MATERL
Respiration Monitors	12662	VBDVAC
Restorative Material Kits, Dental	16732	MTDTRS
Restorative Material Kits, Dental, Crown and Bridge	16733	MTDTRS
Restorative Material Kits, Dental, Unfilled	16739	MTDTRS
Restorative Materials, Dental, Crown and Bridge	16723	MTDTRS
Restorative Materials, Dental, Other	17619	MTDTRS
Restrictors, Orthopedic Cement	16126	ORTHEQ
Resuscitators, Cardiac	13361	RESCAQ
Resuscitators, Cardiopulmonary	13363	RESCAQ
Resuscitators, CPR	*****	RESCAQ
Resuscitators, Pulmonary, Exhaled Air	17141	RESCAQ
Resuscitators, Pulmonary, Gas-Powered	13366	RESCAQ
Resuscitators, Pulmonary, Manual, Disposable	17592	RESCAQ
Resuscitators, Pulmonary, Manual, Reusable	17591	RESCAQ
Retainers, Visceral	13371	SGINAC
Retinoscopes	13372	OPHDEQ
Retractors, Abdominal	13375	SGINAC
Retractors, Bladder	13376	SGINAC
Retractors, Brain	13377	SGINAC
Retractors, Brain Decompression	13378	SGINAC
Retractors, Cardiac	13379	SGINAC
Retractors, Dental	13380	DENTEQ
Retractors, Eye	13381	SGINAC
Retractors, Fiberoptic-Illuminated	15635	SGINAC
Retractors, Hand Operating	13383	SGINAC
Retractors, Laminectomy	13384	SGINAC
Retractors, Mastoid	15675	SGINAC
Retractors, Orbital	13385	SGINAC
Retractors, Orthopedic	13386	SGINAC
Retractors, Plastic Surgery	13387	SGINAC
Retractors, Rectal	13388	SGINAC
Retractors, Rib	13389	SGINAC
Retractors, Self-Retaining Operating	13390	SGINAC
Retractors, Thoracic	13391	SGINAC
Retractors, Vaginal	13392	SGINAC
Retractors, Vessel	13393	SGINAC
Retractors, Other	15242	SGINAC
Retrieval Baskets, Biliary Stone	15629	SGINAC

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Retrieval Baskets, Urological Stone	17574	SGINAC
Rhinomanometers	15737	ENTEQU
Rhinoscopes	13396	ENDOSC
Rings, Fallopian Tube	*****	CONTDV
Rings, Tourniquets	*****	TORQAC
Rods	*****	IMPABS
Rongeurs, Dental	15713	DENTEQ
Rongeurs, Intervertebral Disc	15669	SGINAC
Rongeurs, Neurosurgical	13418	SGINAC
Rongeurs, Rib	13419	SGINAC
Rongeurs, Other	15243	SGINAC
Root Canal Posts	16202	DENTEQ
Routers, Bone Cutting	13425	SGINAC
Rubber Dams, Dental	11155	DENTEQ

**S**

Saliva, Artificial	*****	ARTBYF
Sampling Kits, Blood	10446	BLOCOL
Sampling Kits, Blood, Fetal Scalp	16873	BLOCOL
Sampling Kits, Blood Gas	10438	BLOCOL
Sampling Kits, Endometrial	11531	SGINAC
Saws, Bone, Electrical	13450	SGINAC
Saws, Bone, Manual	13451	SGINAC
Saws, Bone, Micro	16716	SGINAC
Saws, Bone, Pneumatic	13452	SGINAC
Saws, Bone, Sternal	16718	SGINAC
Saws, Other	15244	SGDRIL
Scalers	13454	DENTEQ
Scalers, Dental, Ultrasonic	17497	DENTEQ
Scanners	13468	LASERS
Scanners, Computed Tomography, Cine	16899	CTSACC
Scanners, Computed Tomography, Full Body	15956	CTSACC
Scanners, Computed Tomography, Head	15955	CTSACC
Scanners, Computed Tomography, Mobile	16766	CTSACC
Scanners, Long Term Recording, ECG	15295	CARDIO
Scanners, Long-Term Recording, EEG	16307	EEGEQA
Scanners, Long Term-Recording, Respiration	16304	RESPMA
Scanners, Ultrasonic, Abdominal	16241	ULTRDV
Scanners, Ultrasonic, Cardiac	17422	ULTRDV
Scanners, Ultrasonic, Compound B	16598	ULTRDV
Scanners, Ultrasonic, Dedicated Linear Array	16747	ULTRDV
Scanners, Ultrasonic, General-Purpose	15976	ULTRDV
Scanners, Ultrasonic, Mammographic	15656	ULTRDV
Scanners, Ultrasonic, Obstetric/Gynecologic	15657	ULTRDV

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Scanners, Ultrasonic, Ophthalmic	11389	ULTRDV
Scanners, Ultrasonic, Other	15659	ULTRDV
Scanners, Ultrasonic, Vascular	15957	ULTRDV
Scarifiers	13472	OPHINS
Scissors, Cardiovascular	13482	SGINAC
Scissors, Dissecting	13492	SGINAC
Scissors, Eye, Corneal	13486	SGINAC
Scissors, Eye, Enucleation	13487	SGINAC
Scissors, Eye, Iris	13488	SGINAC
Scissors, Eye, Suture	13489	SGINAC
Scissors, Eye, Tenotomy	13490	SGINAC
Scissors, Eye, Other	15246	SGINAC
Scissors, Fistula	15977	SGINAC
Scissors, Gynecological	13493	GYNOBE
Scissors, Nasal	13495	SGINAC
Scissors, Neurosurgical	13496	SGINAC
Scissors, Orthopedic	13498	SGINAC
Scissors, Pediatric	13499	SGINAC
Scissors, Plastic Surgery	13500	SGINAC
Scissors, Rectal	13501	SGINAC
Scissors, Thoracic	13504	SGINAC
Scissors, Tonsil	16520	SGINAC
Scissors, Other	15245	SGINAC
Scleral Buckling Devices	15794	SGINAC
Sclerotomes	13507	SGINAC
Scoops, Common Duct	15597	SGINAC
Scoops, Gallstone	13509	SGINAC
Scooters, Powered	17733	AMBEQP
Scrapers, Cytology	13512	CELTCE
Screwdrivers, Bone	13517	SGINAC
Screws	16055	IMPABS
Screws, Bone	16101	GRFIMP
Screws, Cranioplasty Plate	16056	GRFIMP
Searchers, Mastoid	13532	LOSEFI
Semisolid, Wound Dressings	*****	DRSWSS
Sensitometers, Radiographic	16563	RADGEQ
Sensor, Radiation	*****	RADACC
Sensors	13536	TRDSNS
Sensors, Oxygen	13538	TRDSNS
Separators	13540	DENTEQ
Separators, Dural	13542	SGINAC
Separators, Plasma	16901	BLOCOL
Separators, Pylorus	13543	SGINAC
Separators, Stapes	15054	SGINAC
Sheets, Cranioplasty	16058	GRFIMP

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Sheets, Drapes, Disposable	15705	APPRO
Sheets, Examination Table, Disposable	15738	APPRO
Sheets, Operating Room, Disposable	15707	APPRO
Shields, Eye	11663	SGINAC
Shields, Eye, Collagen	17652	CIPHAO
Shields, X-ray	14490	RDPRAP
Shields, X-ray, Apron	14491	RDPRAP
Shields, X-ray, Eye	15871	RDPRAP
Shields, X-ray, Glove	14493	RDPRAP
Shields, X-ray, Gonadal	14494	RDPRAP
Shields, X-ray, Mask	16496	RDPRAP
Shields, X-ray, Other	15266	RDPRAP
Shields, X-ray, Portable	15872	RDPRAP
Shields, X-ray, Throat	16564	RDPRAP
Shunt Adapters	13584	SHUNTS
Shunts, Arteriovenous	13586	SHUNTS
Shunts, Carotid Artery	17797	GRFIMP
Shunts, Endolymphatic	16033	SHUNTS
Shunts, Hydrocephalic	16244	SHUNTS
Shunts, Peritoneovenous	13589	SHUNTS
Shunts, Pleuro Peritoneal	17218	SHUNTS
Shunts, Pulmonary Artery	17811	GRFIMP
Sigmoidoscopes, Flexible	15057	ENDOSC
Sigmoidoscopes, Flexible, Video	17664	ENDOSC
Sigmoidoscopes, Rigid	15058	ENDOSC
Silicone Rubber Sheeting	13598	MATNIM
Simulators, Radiotherapy	13280	RADTUN
Single line listing only is required (non-itemised)	*****	ORTHCO
Skin Graft Cutters	13613	SGINAC
Skin Graft Cutters, Expanded Mesh	13614	SGINAC
Sleeves	16146	GRFIMP
Sleeves, Eye Muscle	16074	SGINAC
Sleeves, Topical Oxygen Therapy	13620	RPUNAC
Slings, Arm	13622	BDDRCL
Slings, Patient Lift	17661	FURNIP
Slit Lamps	12281	OPHDEQ
Smoke Evacuation Systems, Surgical	16262	LASERS
Snake Bite Kits	13629	FAEMKT
Snares, Nasal	15676	SGINAC
Snares, Polyp	13631	SGINAC
Snares, Rectal	15988	SGINAC
Snares, Tonsil	13632	SGINAC
Snares, Other	15247	SGINAC
Sounds, Urethral	13643	SGINAC
Sounds, Uterine	13644	SGINAC

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Sounds, Other	15248	SGINAC
Spacer	18084	ATINNB
Spacers, Tendon	16102	SGINAC
Spatulas, Brain	13646	SGINAC
Spatulas, Eye	13648	SGINAC
Spatulas, Lung	13649	SGINAC
Spatulas, Microsurgical	13650	SGINAC
Spatulas, Middle Ear	13651	SGINAC
Spatulas, Other	15249	SGINAC
Specula, Ear	13662	SGINAC
Specula, Eye	13663	SGINAC
Specula, Nasal	13664	SGINAC
Specula, Rectal	13665	SGINAC
Specula, Vaginal	13666	SGINAC
Specula, Other	15250	SGINAC
Speech Amplifiers	15747	SPTHUN
Speech Therapy Units	13668	SPTHUN
Spheres, Eye	13146	GRFIMP
Sphincteroscopes	15787	ENDOSC
Sphygmocardiographs	13670	CARDIO
Sphygmographs	13671	CARDIO
Sphygmomanometers, Aneroid	16156	BLPRDV
Sphygmomanometers, Counterpressure	16986	BLPRDV
Sphygmomanometers, Electronic, Automatic	16173	BLPRDV
Sphygmomanometers, Electronic, Manual	16174	BLPRDV
Sphygmomanometers, Mercury	16158	BLPRDV
Spirometers, Bronchial	10492	RESPMA
Spirometers, Diagnostic	13680	RESPMA
Spirometers, Monitoring	13679	RESPMA
Splints, Other	15251	ORTHEQ
Splints, Pneumatic	13687	ORTHEQ
Splints, Vacuum	13689	ORTHEQ
Sponges, Dissector	13698	SPONGE
Sponges, Filmated	13699	SPONGE
Sponges, Gauze	13700	SPONGE
Sponges, Germicidal	13701	CONTDV
Sponges, Hemostatic, Absorbable	16791	HEMOMT
Sponges, Laparotomy	15085	SPONGE
Sponges, Neuro	13702	SPONGE
Sponges, Ophthalmic	13704	SPONGE
Sponges, Other	15252	SPONGE
Sponges, Rayon Cellulose	13703	SPONGE
Sponges, Scrub	15281	SWABCL
Sponges, X-ray Detectable	13705	SPONGE
Sprays, Refrigerant	*****	SPRAYS

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Spreaders, Rib	13709	SGINAC
Spreaders, Vein	13710	SGINAC
Spreaders, Other	15253	SGINAC
Spring Infusor, Ambulatory	*****	INFDEV
Spuds, Eye	16025	SGINAC
Spuds, Eye, Electric	15692	SGINAC
Stabilizer	*****	SGINAC
Stabilizer, Vascular	*****	SGINAC
Standing Frames, Mobile	16912	AMBEQP
Staples	16771	IMPABS
Staple Removers	16787	STAPAC
Staples, Bone	16103	GRFIMP
		SGINAC
Staples, Mandibular	16203	SGINAC
Staples, Tissue	15964	STAPAC
Stents	15784	IMPABS
		STENTS
Stents, Biliary	17672	STENTS
Stents, Bronchial	17957	STENTS
Stents, Esophageal	17928	STENTS
Stents, Pancreatic	17717	STENTS
Stents, Tracheal	17448	STENTS
Stents, Ureteral	16040	STENTS
Stents, Urethral	17889	STENTS
Stents, Vaginal	16041	STENTS
Stents, Vascular	17461	STENTS
Stereotactic, Biopsy Systems, Mammographic	17833	RADGEQ
Stereotactic Headframes	13727	SGINAC
		STCRSS
Stereotactic Radiosurgical Systems	18053	STCRSS
Stereotactic Radiosurgical Systems - Electrosurgical	*****	STCRSS
Stereotactic Radiosurgical Systems, Gamma	17641	STCRSS
Stereotactic Radiosurgical Systems, Linear Accelerator	18054	STCRSS
Sterilant	17920	DISIWC
Sterilising Units	13737	STERIL
Sterilising Units, Dry Heat	13739	STERIL
Sterilising Units, Steam, Tabletop	16142	STERIL
Stethoscopes, Electronic	13754	STETHS
Stethoscopes, Esophageal	11601	STETHS
Stethoscopes, Mechanical	13755	STETHS
Stethoscopes, Ultrasonic	13756	STETHS
Stimulator Analyzers, Transcutaneous Nerve	15631	STIMAC
Stimulators	13762	STIMAC
Stimulators, Acupuncture	13763	STIMAC
Stimulators, Antihidrosis	17421	AISTIM

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Stimulators, Auditory	13764	STIMAC AISTIM
Stimulators, Antiseizure	17913	STIMAC
Stimulators, Aversion Conditioning	13765	BIOSYS
Stimulators, Bone Growth	15087	AISTIM STIMAC
Stimulators, Carotid Sinus Nerve	13768	AISTIM
Stimulators, Cerebellar	16122	AISTIM
Stimulators, Cranial Electrotherapy	17758	STIMAC
Stimulators, Diagnostic Cardiological	17121	STIMAC
Stimulators, Diagnostic Neurological	16263	STIMAC
Stimulators, Diaphragmatic/Phrenic Nerve	16060	AISTIM
Stimulators, Electroanalgesic, Dorsal Column	13769	STIMAC AISTIM
Stimulators, Electroanalgesic, Intracerebral/Subcortical	16061	AISTIM
Stimulators, Electroanalgesic, Spinal Cord	17241	STIMAC AISTIM
Stimulators, Electroanalgesic, Transcutaneous Electrical Nerve	13782	STIMAC
Stimulators, Incontinence Management	17784	STIMAC
Stimulators, Neuromuscular, Functional	16250	STIMAC AISTIM
Stimulators, Neuromuscular, Therapeutic	16255	STIMAC AISTIM
Stimulators, Peripheral-Nerve-Block Monitor	16252	STIMAC
Stimulators, Scoliosis Treatment	17240	AISTIM
Stimulators, Somatosensory	16256	STIMAC
Stimulators, Tactile	13780	STIMAC
Stimulators, Urinary	16029	STIMAC AISTIM
Stockinettes	13785	BDDRCL
Stockings, Compression	13789	AEMBDV
Stockings, Pneumatic	13793	AEMBDV
Stone Dislodgers	16471	SGINAC
Stone Dislodgers, Ureteral	16261	SGINAC
Stopcocks	13803	STOPCK
Strain Detectors	13805	TRDSNS
Strain Gauges	13806	TRDSNS
Stress Tables, Nuclear Cardiology	16024	FURNIP
String, Gastrointestinal	13822	SGINAC
String, Medical-Surgical	17586	STAPAC
Strippers, Artery	13824	SGINAC
Strippers, Tendon	13827	SGINAC
Strippers, Vein	13828	SGINAC
Strippers, Other	15254	SGINAC
Strips, Nasal	*****	RESPMA



**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Stylets	13833	SGINAC
Stylets, Bronchial	13834	SGINAC
Stylets, Intravenous	13835	SGINAC
Stylets, Needle	13836	SGINAC
Stylets, Peritoneal	13837	SGINAC
Stylets, Tracheal Tube	14084	SGINAC
Stylets, Urethral	13839	SGINAC
Suction Irrigators	13845	IRILAV
Suction Kits	13846	ASDVAC
Suction Tips, Electrosurgical	15314	ESUNAC
Suction Tips, Fiberoptic Illuminated	15636	ASDVAC
Suction Tips, General Purpose	16883	ASDVAC
Surface Coils, Magnetic Resonance Imaging (MRI)	17542	MRIMUN
Surface Spray Disinfectant	*****	DISIWC
Surgical Drill System	*****	SGDRIL
Surgical Head Coverings	13882	APPRO
Surgical Packings	13886	BDDRCL BDDRST
Surgical Packs	13887	KTPWOD
Suture Bolsters	13890	STAPAC
Suture Buttons	13891	STAPAC
Suture Kits	13892	STAPAC
Suture Removal Kits	13894	STAPAC
Suture Rings	13895	STAPAC
Suture Units, Automatic	15065	STAPAC
Suture Units, Other	16264	STAPAC
Sutures, Catgut	13899	SUTABS
Sutures, Collagen	13898	SUTABS
Sutures, Cotton	13900	SUTNAB
Sutures, Iodized	15873	SUTNAB
Sutures, Linen	13902	SUTNAB
Sutures, Nylon	13905	SUTNAB
Sutures, Polybutester	17245	SUTNAB
Sutures, Polydioxanone	16584	SUTABS
Sutures, Polyester	13906	SUTNAB
Sutures, Polyethylene	13907	SUTNAB
Sutures, Polyglactin	17471	SUTABS
Sutures, Polyglycolic Acid	13908	SUTABS
Sutures, Polyglyconate	17246	SUTABS
Sutures, Polypropylene Monofilament	13909	SUTNAB
Sutures, Polytetrafluoroethylene	17467	SUTNAB
Sutures, Silk	13910	SUTNAB
Sutures, Silver	15972	SUTNAB
Sutures, Stainless Steel Monofilament	15971	SUTSSS
Sutures, Stainless Steel Multifilament	13904	SUTSSS

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Sutures, Other	15255	SUTABS SUTNAB
Swabs, Cotton	15066	BDDRCL
Swabs, Oral Care	16619	APPMED
Swabs, Specimen Collecting	13914	CELTCE
Symblepharon Rings	16073	SGINAC
Synchronizers	15604	VBDVAC
Synchronizers, Computed Axial Tomography	17238	CTSACC
Synchronizers, Magnetic Resonance Imaging	17239	MRIMUN
Synovial Fluid, Artificial	*****	ARTBYF
Syringe/Needle Adapters	13922	SYRINJ
Syringes, Anesthesia	13931	SYRINJ
Syringes, Angiographic	15286	SYRINJ
Syringes, Aspirating	16650	SYRINJ
Syringes, Blood Gas	16785	SYRINJ
Syringes, Cartridge	16585	SYRINJ
Syringes, Catheter	13936	SYRINJ
Syringes, Ear	11376	IRILAV
Syringes, Hypodermic	13940	SYRINJ
Syringes, Insulin	13941	SYRINS
Syringes, Intraligamentary	16702	SYRINJ
Syringes, Irrigating	13942	IRILAV
Syringes, Metered Delivery	16824	SYRINJ
Syringes, Micro	16878	SYRINJ
Syringes, Tuberculin	13945	SYRINJ
Syringes, Other	15256	SYRINJ

**T**

Tables, Examination/Treatment	13958	FURNIP
Tables, Obstetrical	13960	FURNIP
Tables, Operating	13961	FURNIP
Tables, Orthopedic	13962	FURNIP
Tables, Special Procedure	16601	SGINAC
Tables, Urological	13969	FURNIP
Tacks, Sacculotomy	16080	SGINAC
Tampons	16164	TAMPON
Tapes, Adhesive, General Purpose	16866	ADHSVT
Tapes, Adhesive, Hypoallergenic	10031	ADHSVT
Tapes, Adhesive, Waterproof	10032	ADHSVT
Tapes, Gauze, Self-Adhering	16913	BDDRCL
Tapes, Umbilical	13983	BDDRST
Tears, Artificial	*****	ARTBYF
Telemetry Units, Multichannel,	15722	TELEME
Telemetry Units, Physiologic, Other	15294	TELEME

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Teleradiology Systems	16994	TELEME
Television Systems, Endoscope	15367	ENDOSC
Temperature Alarms, Airway	17427	AWPSDA
Temperature Monitors	12672	BODYTP
Temperature Monitors, Airway	16717	VBDVAC
Tenaculums, Thyroid	13997	SGINAC
Tenaculums, Uterine	13998	SGINAC
Tenaculums, Other	15258	SGINAC
Tenotomes	13999	SGINAC
Tents, Mist	12554	HUMIDS
Tents, Mist, Face	12557	HUMIDS
Tents, Oxygen	12893	MGASEQ
Testers, Pulp	13187	DENTEQ
Thermographic Devices, Infrared	14022	BODYTP
Thermographic Devices, Liquid-Crystal	14023	BODYTP
Thermometer Kits	14025	BODYTP
Thermometers, Colour-Indicating	16804	BODYTP
Thermometers, Electronic, Continuous	14034	BODYTP
Thermometers, Electronic, Intermittent	14035	BODYTP
Thermometers, Infrared	14036	BODYTP
Thermometers, Infrared, Ear	17887	BODYTP
Thermometers, Infrared, Skin	17888	BODYTP
Thermometers, Mercury	14038	BODYTP
Thermometers, Phosphorescent/Fiber Optic	16960	BODYTP
Thorascopes	14047	ENDOSC
Thyroid Function Units	15143	NUCMEQ
Thyroid Uptake Systems	16567	NUCMEQ
Tinnitus Masker Analyzers	15632	HEARAD
Tinnitus Therapy Units	14056	HEARAD
Tissue Morcellator	*****	SGINAC
Tissue Reconstructive Material	17875	AUGPSM
Tissue Reconstructive Materials, Solid	16153	GRFIMP
Tissue Storage Media	*****	CELTCE
Tongs, Bone	14063	SGINAC
Tongs, Skull	14064	SGINAC
Tonographs	14067	OPHDEQ
Tonometers, Ophthalmic	16809	OPHDEQ
Tonosscopes	14069	OPHDEQ
Tonsillectomes	14070	SGINAC
Toothbrushes, Vacuum Supplemented	17597	DENTEQ
Topographic Display Units, EEG/Evoked-Potential	17413	EEGEQA
Torque Device	*****	SGINAC
Tourniquets, Cardiovascular	14075	TORQAC
Tourniquets, Pneumatic	14074	TORQAC
Tourniquets, Pneumatic, Automatic Rotating	10234	TORQAC

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Towels, Surgical	14077	APPPRO
Tracheal Tube Adapters	14080	AWPSDA
Tracheal Tube Cuff Inflators	14081	TUBEST
Tracheal Tube Holders	16489	AWPSDA
Tracheal Tube Introducers	17694	TUBEST
Tracheostomy Buttons	14093	TTCEQU
Tracheostomy Care Kits	14090	TTCEQU
Tracheostomy Tube Adapters	14092	TTCEQU
Tracheotomy Kits	14099	TTCEQU
Traction Units, Intermittent	14106	ORTHEQ
Traction Units, Intermittent, Mobile	14108	ORTHEQ
Traction Units, Static, Bed	14110	ORTHEQ
Traction Units, Static, Chair	14111	ORTHEQ
Traction Units, Static, Other	14109	ORTHEQ
Transducers	14115	TRDSNS
Transducers, Flow	14117	TRDSNS
Transducers, Flow, Catheter Tip	17473	TRDSNS
Transducers, Force	14118	TRDSNS
Transducers, Penile Tumescence	17442	TRDSNS
Transducers, Pressure	14119	TRDSNS
Transducers, Pressure, Catheter Tip	15071	TRDSNS
Transducers, Pressure, Disposable	16909	TRDSNS
Transducers, Ultrasonic	14121	TRDSNS
Transfer Sets	17164	BLOCOL
Transfer Sets, Continuous Ambulatory Peritoneal Dialysis	17163	PDEQAC
Transfer Sets, Fluid	16610	IVASAC
Transfusion Kits	14126	BLOCOL
Transmitter/Receiver Systems, Telephone, ECG, Multichannel	15685	TELEME
Transmitter/Receiver Systems, Telephone, ECG, Single-Channel	15684	TELEME
Transmitter/Receiver Systems, Telephone, Fetal Monitor	16507	TELEME
Traps, Bubble	10511	CPBYEQ
Treadmills	14141	PHYSEV
Trephines, Bone	14147	SGINAC
Trephines, Corneal	14148	SGINAC
Trephines, Kidney	14149	SGINAC
Trephines, Liver	14150	SGINAC
Trephines, Lung	14151	SGINAC
Trephines, Skull	14152	SGINAC
Trephines, Other	15259	SGINAC
Trocars, Abdominal	14155	SGINAC
Trocars, Amniotic	14156	SGINAC
Trocars, Antrum	14157	SGINAC
Trocars, Cardiovascular	14158	SGINAC
Trocars, Gallbladder	14159	SGINAC
Trocars, Laryngeal	14160	SGINAC

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Trocars, Suprapubic	16421	SGINAC
Trocars, Thoracic	14161	SGINAC
Trocars, Tracheal	14162	CANNUL
Trocars, Other	15260	SGINAC
Tropometers	14163	OPHDEQ
Tubes	14175	ATINNB
Tube Protector	*****	TTCEQU
Tubes, Blood Collection	14183	BLOCOL
Tubes, Blood, Microcollection	15315	BLOCOL
Tubes, Bronchial	15322	TUBEST
Tubes, Colon	14187	TUBEST
Tubes, Connecting	14188	TUBEST
Tubes, Decompression	14190	TUBEST
Tubes, Drainage	14191	TUBEST
Tubes, Enema	14194	ENMAKT
Tubes, Esophageal, Blakemore	14196	TUBEST
Tubes, Esophageal, Replogle	14197	TUBEST
Tubes, Esophageal, Sengstaken	14198	TUBEST
Tubes, Esophageal, Other	15261	TUBEST
Tubes, Feeding, Enterostomy	16799	EFEDKT
Tubes, Feeding, Nasoenteral	16798	EFEDKT
Tubes, Gastrointestinal Decompression, Baker Jejunostomy	15998	TUBEST
Tubes, Gastrointestinal Decompression, Cantor	14211	TUBEST
Tubes, Gastrointestinal Decompression, Dennis	14212	TUBEST
Tubes, Gastrointestinal Decompression, Miller-Abbott	14213	TUBEST
Tubes, Gastrointestinal Decompression, Wangensteen	14209	TUBEST
Tubes, Gastrointestinal Decompression, Other	15263	TUBEST
Tubes, Gastrointestinal, Einhorn	14203	TUBEST
Tubes, Gastrointestinal, Harris	14204	TUBEST
Tubes, Gastrointestinal, Homer Smathers	14205	TUBEST
Tubes, Gastrointestinal, Kaslow	14206	TUBEST
Tubes, Gastrointestinal, Rehfus	14207	TUBEST
Tubes, Gastrointestinal, Ryles	14208	TUBEST
Tubes, Gastrointestinal, Other	15262	TUBEST
Tubes, Gastrostomy, Malecot	14215	TUBEST
Tubes, Gastrostomy, Mushroom	14216	TUBEST
Tubes, Gastrostomy, Other	15264	TUBEST
Tubes, Lacrimal	17003	TUBEST
Tubes, Laryngectomy	12292	TUBEST
Tubes, Myringotomy	15323	TUBEST
Tubes, Nasogastric	14221	TUBEST
Tubes, Nephrostomy	14224	TUBEST
Tubes, Rebreathing	14226	TUBEST
Tubes, Rectal	14227	TUBEST
Tubes, Stomach Evacuator	14230	TUBEST

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Tubes, Tracheal	14085	TUBEST
Tubes, Tracheal, Laser-Resistant	17548	TUBEST
Tubes, Tracheostomy	14096	TTCEQU
Tubes, Other	15278	TUBEST
Tubing, Conductive	14240	TUBEST
Tubing, Nonconductive	15668	TUBEST
Tubing, Oxygen Connecting	12875	MGASEQ
Tubing, Polyethylene	14246	TUBEST
Tubing, Polyvinyl Chloride	14247	TUBEST
Tubing, Silicone	15277	TUBEST
Tubing, Suction	16779	ASDVAC
Tubing, Other	15265	TUBEST
Tunnelers, Other	*****	SGINAC
Tunnelers, Vascular	16985	SGINAC

**U**

Ultrasonic Surgical Cutter/Coagulators	18049	SGINAC
Ultrasound Units, Physical Therapy	11248	ELTHEQ
Ultrasound Therapy Units	*****	DIATEQ
Uncompounded Substances	*****	EMOLPR
Ureterorenoscopes	17690	ENDOSC
Ureteroscopes	15788	ENDOSC
Ureterotomes	14285	SGINAC
Urethrographs	14287	URODEQ
Urethrometers	14288	URODEQ
Urethroscopes	14289	ENDOSC
Urethrotomes	17633	SGINAC
Urinary Catheterization Kits, Foley	16781	CATKEQ
Urinary Catheterization/Collection Kits	14292	CATKEQ
Urinary Collection Bags	14298	UROCOL
Urinary Collection Bags, Infant	16782	UROCOL
Urinary Collection Units, Precision Measurement	15310	UROCOL
Urinary Drainage Units	14297	UROCOL
Urodynamic Measurement Systems	14307	URODEQ
Uroflowmeters	14308	URODEQ
Uterine Activity Monitors	17922	GYNOBE
Uterine Device	*****	GYNOBE
Uterine Packing Devices	14313	SGINAC

**V**

Vacuum Extractors, Obstetrical	14317	GYNOBE
Vacuum Mixing Devices, Orthopaedic Cement	17891	ORTHEQ
Vaginometers	14323	GYNOBE
Vaginoscopes	15630	ENDOSC

**ECRI UMDNS DEVICE TERM****ECRI IMDC****ADG CODE**

Valves	14325	RESCAQ
Valves, Automatic Shut-Off, Heart-Lung Bypass	17463	CPBYEQ
Valves, Demand	17169	VBDVAC
Valves, Ear	14328	GRFIMP
Valves, Eye	14329	GRFIMP
Valves, Hydrocephalic	17090	SHUNTS
Valves, Nonbreathing	14332	VBDVAC
Valves, Positive End-Expiratory Pressure	14337	VBDVAC
Valves, Suction Control, Ventricular Vent	17581	VBDVAC
Valves, Speech, Tracheostomy	17563	SPTHUN
Valvulotomes	14339	SGINAC
Vascular Prostheses, Human or Animal Origin, Coating	*****	PROVAO
Vectorcardiographs	14345	CARDIO
Vein Finders	14346	LOSEFI
Ventilation Monitors	12678	VBDVAC
Ventilator IMV Controllers	14359	VBDVAC
Ventilators, High-Frequency	15783	VBDVAC
Ventilators, Negative-Pressure	17877	VBDVAC
Ventilators, Pressure Cycled	14360	VBDVAC
Ventilators, Time Cycled	14361	VBDVAC
Ventilators, Volume, Mobile	17429	VBDVAC
Ventilators, Volume, Portable	17423	VBDVAC
Ventilators, Other	14355	VBDVAC
Ventriculosopes	14364	ENDOSC
Vibrators	14369	AWPSDA
Vibrometers	14371	NEURDE
Vibrators, Bone, Audiometric	15980	AUDIEQ
Vibrophonocardiographs	14372	CARDIO
Video Image Processors, Endoscopic	18034	ENDOSC
Video Image Processors, Operating Microscope	18036	MICROP
Visual Function Analyzers	14382	OPHDEQ
Vital Signs Monitors	17678	CARDMA
Vitrectomy Units	14386	OPHINS

**W**

Warmers, Anesthetic	16371	ANEQAC
Warmers, Blood, Dry Heat	10449	BLOODW
Warmers, Blood, Microwave	10450	BLOODW
Warmers, Blood, Radio-Frequency	15644	BLOODW
Warmers, Blood, Water Bath	17160	BLOODW
Warmers, Radiant, Adult	13249	THERCO
Warmers, Radiant, Infant, Mobile	17433	THERCO
Warmers, Radiant, Infant, Transport	13251	THERCO
Washers, Bone Screw	17856	GRFIMP

**ECRI UMDNS DEVICE TERM****ECRI IMDC    ADG CODE**

Washer/Sterilising	14413	STERIL
Water Purification Systems, Reverse Osmosis	14437	HDHF EQ
Wetting Solutions, Contact Lens	17431	COLECP
Wheelchairs, Powered	16214	AMBEQP
Wire Drivers	14455	SGINAC
Wires, Bone, Implantable	*****	GRFIMP
Wires, Ligature	14457	SGINAC
Wound Drainage Kits	16521	DRKTAC
Wound Drainage Kits, Closed, Craniotomy	15318	DRKTAC
Wraps, Burn	17604	BDDRST

**X**

X-ray Generators	16602	RADGEQ
X-ray Image Recorders	16512	RADGEQ
Xenon Systems	14467	NUCMEQ



## **Therapeutic Goods Administration**

PO Box 100 Woden ACT 2606 Australia

Email: [info@tga.gov.au](mailto:info@tga.gov.au) Phone: 1800 020 653 Fax: 02 6232 8605

[www.tga.gov.au](http://www.tga.gov.au)