

10. Finished Product Release and Expiry Specifications:

TEST	SPECIFICATION		Method Reference
Description	A clear, colourless or almost colourless solution free from visible particulate.		Visual
Identification	[REDACTED]		[REDACTED]
Colour & Clarity of Solution	[REDACTED]		[REDACTED]
Assay	Release: [REDACTED]	Expiry: [REDACTED]	[REDACTED]
Related substances [REDACTED]	Release: [REDACTED]	Expiry: [REDACTED]	[REDACTED]
pH	[REDACTED]		[REDACTED]
Extractable Volume	NLT nominal fill		BP/PhEur
Sterility	Complies with test		Ph. Eur.
Bacterial Endotoxins	[REDACTED]		[REDACTED]
Particulate Matter	[REDACTED]		[REDACTED]

11. Test methodology:

Copies of all In-house test methods are provided in **Appendix 1**.