

Draft of January 2001

1

u6

Draft of January 2001


2


665

Draft of January 2001

4


u3





Draft of January 2001


6





Draft of January 2001


7


uc

















DOSAGE AND ADMINISTRATION

Method of administration

Before using Trileptal oral suspension, the bottle must be shaken well and the dose prepared immediately afterwards. The prescribed amount of oral suspension should be withdrawn from the bottle using the oral dosing syringe supplied. Trileptal oral suspension may be swallowed directly from the syringe or can be mixed in a small glass of water just prior to administration. After each use, the bottle should be closed and the outside of the oral syringe should be wiped with a dry, clean tissue.

PRESENTATION

Oral suspension

- *Trileptal 60 mg/mL: off-white to slightly reddish brown oral suspension in 250 mL bottles. The bottles are supplied with a 10 mL oral dosing syringe and a press-in bottle adaptor.*